

Please refer to www.invacare.ie for Local Distributor

Account Name :	Account No :
Address :	Order No :
Tel No:	Order Qty :
Delivery address :	Order Date :

Matrx Posture Back Deep

Ireland order form

Maximum User Weight: 227 kg

WITDH

HEIGHT	cm		36	41	43	46	51
		inch	14	16	17	18	20
	31	12					
	41	16					
51	20						



Deep MatRX Posture Back Backrest

Chest Width	Backrest Width					Fits W/ch Widths		With Extended Mounting Pins EMP*	
			31cm (12")	41cm (16")	51cm (20")				
≤26cm (10")	36cm (14")		PBD1412NH	PBD1416NH		33cm - 38cm	(13" - 15")	41cm (16")	
≤31cm (12")	41cm (16")		PBD1612NH	PBD1616NH	PBD1620NH	38cm - 43cm	(15" - 17")	46cm (18")	
≤33cm (13")	43cm (17")		PBD1712NH	PBD1716NH	PBD1720NH	41cm - 46cm	(16" - 18")	48cm (19")	
≤36cm (14")	46cm (18")		PBD1812NH	PBD1816NH	PBD1820NH	43cm - 48cm	(17" - 19")	51cm (20")	
≤41cm (16")	51cm (20")		PBD2012NH	PBD2016NH	PBD2020NH	48cm - 53cm	(19" - 21")	56cm (22")	

*Not suitable for HD users (max. 136kg/300lb weight capacity)

Features

5" Contour Depth, 40° angle adjustment, 4" depth adjustment, 20° rotation, 136kg (21 stone) user weight

Code	X	Description		Price
PB Deep				
PBD	NH	Posture Back Deep		
1545001		Posture Back Deep PBD1412NH	330mm-380mm	310mm €532
1545002		Posture Back Deep PBD1416NH	330mm-380mm	410mm €532
1545003		Posture Back Deep PBD1612NH	380mm-430mm	310mm €532
1545058		Posture Back Deep PBD1616NH	380mm-430mm	410mm €532
1545059		Posture Back Deep PBD1620NH	380mm-430mm	510mm €532
1545060		Posture Back Deep PBD1712NH	410mm-460mm	310mm €532
1545061		Posture Back Deep PBD1716NH	410mm-460mm	410mm €532
1545062		Posture Back Deep PBD1720NH	410mm-460mm	510mm €796
1545063		Posture Back Deep PBD1812NH	430mm-480mm	310mm €532
1545064		Posture Back Deep PBD1816NH	430mm-480mm	410mm €796
1545065		Posture Back Deep PBD1820NH	430mm-480mm	510mm €796
1545066		Posture Back Deep PBD2012NH	480mm-530mm	310mm €532
1545067		Posture Back Deep PBD2016NH	480mm-530mm	410mm €532
1545068		Posture Back Deep PBD2020NH	480mm-530mm	510mm €796

Hardware					
1545275		Easy Set mouting hardware (pair) - plastic	Fits 19mm, 22mm and 25mm back canes		€357
1545276		Easy Set mouting hardware (pair) - plastic	Fits 29mm back canes		€357
1545277		Fixed Space Saver Hardware			€260
1545278		Split Mounting hardware			€357
1544776		Heavy Duty Mounting Hardware (pair) - metal	Fits 19mm, 22mm and 25mm back canes		€337

Harness interface				
CHIPB		Harness Interface to allow attachment of various chest supports without drilling of the back shell	Width, angle and height adjustability for precise orientation of shoulder straps	€129
Privacy flap attachment to fit all backs from 30.5cm to 50.5cm wide				
1551799		Extra Small	26cm x 41cm	€80
1551800		Small	33cm x 41cm	€80
1551903		Medium	38cm x 41cm	€80
1551904		Large	43cm x 41cm	€80
Spares				
EMP		Extended Mounting Pins	Allows Matrx backs to be installed on wider wheelchairs	€69
Headrest				
1634173		250mm soft headrest pad complete with angle and depth adjustable stem		€349
Headrest Receiver to suit backrest width:				
1544189		Small Adapter for backs 12" wide		€58
1544188		Medium Adapter for backs 13-17" wide		€58
1544187		Large Adapter for backs 18-20" wide		€58
Replacement covers				
1545069		Spare Matrx PB deep cover	330mm-380mm	€109
1545070		Spare Matrx PB deep cover	330mm-380mm	€109
1545071		Spare Matrx PB deep cover	380mm-430mm	€109
1545072		Spare Matrx PB deep cover	380mm-430mm	€109
1545073		Spare Matrx PB deep cover	380mm-430mm	€109
1545074		Spare Matrx PB deep cover	410mm-460mm	€109
1545075		Spare Matrx PB deep cover	410mm-460mm	€109
1545076		Spare Matrx PB deep cover	410mm-460mm	€109
1545077		Spare Matrx PB deep cover	430mm-480mm	€109
1545078		Spare Matrx PB deep cover	430mm-480mm	€109
1545079		Spare Matrx PB deep cover	430mm-480mm	€109
1545080		Spare Matrx PB deep cover	480mm-530mm	€109
1545081		Spare Matrx PB deep cover	480mm-530mm	€109
1545082		Spare Matrx PB deep cover	480mm-530mm	€109

Special Instructions - Please supply the following information:



Name:

Email:

Phone:

Description of special request (please include relevant option code):