

Please refer to www.invacare.ie for Local Distributor

Account Name :	Account No :
Address :	Order No :
Tel No:	Order Qty :
Delivery address :	Order Date :

Matrix Posture Back

Ireland order form
Maximum User Weight: 136kg

		WITDTH				
		36	41	43	46	51
HEIGHT	cm					
	inch	14	16	17	18	20
	31	12				
	41	16				
51	20					



Standard MatRX Posture Back Backrest (no hardware)

Chest Width	Backrest Width				Fits W/ch Widths	With Extended Mounting Pins EMP*	
		31cm (12")	41cm (16")	51cm (20")			
≤26cm (10")	36cm (14")	MPB1412NH	MPB1416NH		33cm - 38cm (13" - 15")	41cm (16")	
≤31cm (12")	41cm (16")	MPB1612NH	MPB1616NH	MPB1620NH	38cm - 43cm (15" - 17")	46cm (18")	
≤33cm (13")	43cm (17")	MPB1712NH	MPB1716NH	MPB1720NH	41cm - 46cm (16" - 18")	48cm (19")	
≤36cm (14")	46cm (18")	MPB1812NH	MPB1816NH	MPB1820NH	43cm - 48cm (17" - 19")	51cm (20")	
≤41cm (16")	51cm (20")	MPB2021NH	MPB2016NH	MPB2020NH	48cm - 53cm (19" - 21")	56cm (22")	

*not suitable for HD users (max. 136kg/300lb weight capacity)

Features

5" Contour Depth, 40° angle adjustment, 4" depth adjustment, 20° rotation, 136kg (21 stone) user weight

Matrix PB		Width	Height	Price
1544783	Matrix Posture Back MPB1412NH	330mm-380mm	310mm	€458
1544784	Matrix Posture Back MPB1416NH	330mm-380mm	410mm	€458
1544785	Matrix Posture Back MPB1612NH	380mm-430mm	310mm	€458
1515727	Matrix Posture Back MPB1616NH	380mm-430mm	410mm	€458
1515728	Matrix Posture Back MPB1620NH	380mm-430mm	510mm	€681
1544786	Matrix Posture Back MPB1712NH	410mm-460mm	310mm	€681
1515729	Matrix Posture Back MPB1716NH	410mm-460mm	410mm	€681
1515730	Matrix Posture Back MPB1720NH	410mm-460mm	510mm	€681
1544787	Matrix Posture Back MPB1812NH	430mm-480mm	310mm	€458
1515731	Matrix Posture Back MPB1816NH	430mm-480mm	410mm	€681
1515732	Matrix Posture Back MPB1820NH	430mm-480mm	510mm	€681
1544788	Matrix Posture Back MPB2012NH	480mm-510mm	310mm	€458
1518131	Matrix Posture Back MPB2016NH	480mm-510mm	410mm	€681
1518132	Matrix Posture Back MPB2020NH	480mm-510mm	510mm	€681

Hardware



Easy Set hardware



Fixed hardware



Split



Heavy Duty hardware

1545275	Easy set hardware for all Matrix backs	Fits 19mm, 22mm and 25mm back canes	€357
1545276	Easy set hardware to fit 29mm back canes	Fits 29mm back canes	€357
1545277	Fixed Space Saver Mounting Hardware		€260
1545278	Split Mounting hardware		€357
1544776	Heavy Duty Easy Set Mounting Hardware	Fits 19mm, 22mm and 25mm back canes	€337

Harness interface				
CHIPB		Harness Interface to allow attachment of various chest supports without drilling of the back shell	Width, angle and height adjustability for precise orientation of shoulder straps	€129
Privacy flap attachment to fit all backs from 30.5cm to 50.5cm wide				
1551799		Extra Small	26cm x 41cm	€80
1551800		Small	33cm x 41cm	€80
1551903		Medium	38cm x 41cm	€80
1551904		Large	43cm x 41cm	€80
Spares				
EMP		Extended Mounting Pins	Allows Matrix backs to be installed on wider wheelchairs	€69
Headrest				
1634173		250mm soft headrest pad complete with angle and depth adjustable stem (All hardware included)		€349
Headrest Receiver to suit backrest width:				
1544189		Small Adapter for backs 12" wide		€58
1544188		Medium Adapter for backs 13-17" wide		€58
1544187		Large Adapter for backs 18-20" wide		€58

Special Instructions - Please supply the following information:

Name:

Email:

