

Medley[®] Ergo

More options. More choice.



Medley Ergo profiling beds
Trusted by healthcare professionals



Customised designs

The Medley Ergo range of medical beds fits perfectly at home and long-term care settings.

With more options, each design aims to achieve comfort, rest and sleep - while prioritising the needs of user and carer.

Choose from:

Medley Ergo

Medley Ergo Low



Safety

All Invacare beds meet IEC 60601-2-52 safety standards.



Quality

Our bed designs and specifications adhere to the highest standards.



Care-led

Made for the caring environment.



Comfort

Every Invacare bed promotes rest, dignity and independence.



Homely Design

Medley Ergo embodies comfort, warmth and personal touches.



* Medley Ergo shown with optional Select bed end panels and side rails. Invacare mattress sold separately.

More options. More choice.



The Medley Ergo range is easy to use and handle as well as designed with options to maximise comfort and ensure practicality at the touch of a button.



Medley Ergo



Elevated leg section



Medley Ergo in tilt

01 Customised care
A range of options to meet various needs

11°
Supportive tilt
For positioning and pressure management.

02 Quick & easy setup
Simple tool-free assembly.

03 Adjustable positioning
Backrest, thigh & legrest adjustment. Legrest elevates manually.

04 Height adjustable
Each model adjusts for comfortable working heights and optimal patient care.

Intuitive hand control
Clear visual display for ease of use. Optional lockable tilt available.

Homely design
Optional Select Beech wood bed end covers available to fit into any home.

Easy to clean
Suitable to wipe with mild detergents making it ideal for single or multiple use environments. See User Manual.

Interchangeable options
Healthcare environments can benefit from the versatility of using the range.

Range of side rails
Choose from Bella with integrated handles, Aria, Verso HC or Scala 2.

08

07

05

06

Note the tilt function is only practised for users under medical supervision.

Medley Ergo shown with optional Select bed end covers and optional Bella side rails..

Medley® Ergo Low

Positioning the bed at a lower level increases security and reassurance for users and carers.
21 - 61 cm or 28 - 68 cm.



Medley® Ergo options

Customised care with a choice of options.



Transport trolley



Customer Support & After Sales

Ask our team. We are always here to help.

Tel: 01656 776222 & 01656 776333



Bella side rail in Beech wood



Aria side rail full length in steel



Scala 2 side rail collapsible in steel*



Lifting pole with adjustable handle



Support handle 30 x 40 cm



Support handle 50 x 40 cm



Alternative handset with lockable tilt



Side rail cover for additional safety

*available in a range of height to accommodate foam and dynamic mattress system

Technical specifications



Width inside
90 cm



Width outside
102 cm



Length inside
208 cm



Length outside
224 cm



Backrest angle
0 - 70°



Tilt angle
11°



Height adjustment
40 - 80 cm
or 33 - 73cm



Heaviest part
17 kg



Mattress support
80 - 23 - 30
- 67 cm



Max user weight
145 kg



Total weight
64 kg



Thigh angle
0 - 33°

Medley Ergo Low

All technical details as above except:

Height adjustment:
21 - 61 cm
or 28 - 68 cm

Contact your Authorised
Distributor:

Colours



Steel



Beech



Complies to the safety standard,
IEC 60601-2-52.



ISO 13485 certified

T: +44 1656 776222
uk@invacare.com
www.invacare.co.uk
www.invacare.ie

Invacare UK & Ireland
Unit 4 Pencoed Technology
Park
Bridgend CF35 5AQ

