

Invacare Medley Ergo Community Bed

Product Data Sheet



Product Range Available

- Medley Ergo
- Medley Ergo low
- Medley Ergo Select

Standard Features

- Comply to the Bed standard IEC 60601-2-52
- Four sectioned profiling mattress support
- Newly optimised mattress platform sections to reduce shear and friction
- Increased comfort when sitting
- Minimised risk of sliding downwards during profiling
- Suitable for shorter clients
- Can be positioned in 2 heights
- Wide range of side rails and side rail covers available
- All models are available on a transporter kit

Medley Low

- Low height of 210mm to maximise safety for vulnerable clients

Medley Select

- More attractive, homely design
- Covered bed ends
- Equipped with a Rastofix to enable individual positioning for the legs.

Technical Specification

Medley Ergo & Low

- Width outside: 1020mm
- Width Inside: 900mm

- Length outside: 2240mm
- Length inside: 2080mm
- Mattress support dimensions: 800 (Head section) 230 (Seat plate) x 300 (thigh section) x 670mm (lower leg section)
- Backrest Angle: 0-70°
- Thigh Angle: 0-24°
- Legrest lift 11°
- Maximum user weight 145kg (23 stone)
- Safe Working Load 180kg (28.5 stone)

Medley Ergo

- Height adjustment: 400 – 800mm or 330 – 730mm

Medley Ergo Low

- Height adjustment 280 – 680mm or 210 – 610mm

Optional accessories available including:

- Full length side rails metal or wooden
- Scala 2 range side rails
- Crash mat for the Medley Ergo Low bed
- Lifting pole
- Height extender for full length wooden side rails
- Transport kit
- 15cm Mattress support extension

Supporting Documentation

- Medley Ergo range brochure PN: 1554900
- Bed Accessory brochure*
- Beds & PAC DL leaflet PN: 1568402
- Bed rail ruler PN: 1556575

**Contact Marketing department*

Manufacturers Contact Details

Invacare Ltd, Unit 4 Pencoed Technology Park, Bridgend, CF35 5AQ

Tel: 01656 776222 **E-mail:** UK@invacare.com

Web: www.invacare.co.uk