

mo-vis

SPECIALIST CONTROLS



Specialist controls combine
with the Invacare Power Series



Introducing Mo-Vis

Mo-Vis specialist controls combine with the Invacare Power Series with either Modulite or Ultra Low Maxx seating, to provide practical, tailored mobility solutions around your needs.



All-round, All-round light and All-round HD

Developed for all-round use and is suitable for the main control input of powerchair users. It can be used as a standard, chin or attendant joystick. The All-round Joystick Light has the same specifications, but with less than half the force of the standard required, with the HD developed for users with increased force.



Scoot control

An intuitive and safe steering device intended to support an attendant in controlling or manoeuvring a power chair, both inside and outside.



Micro joystick

An ultra-sensitive proportional joystick developed for users with very limited strength, movement range, or endurance.



Multi joystick

A small proportional joystick, which requires significantly less power and movement compared to a normal joystick.



Foot control

A proportional control, movements of the foot can be converted into steering and control of the powerchair.



Proportional and Switch head control

Highly adjustable to the user's needs.



Sip and Puff head control

Providing enhanced mobility for individuals with little extremity control.



Hand warmer

Keeps the area around the wheelchair's input device warm to help improve the flexibility of cold hands and fingers.



Variety of switches and buttons

To assist in the operation of connected devices.

salesuk@dhcare.com
www.dhcare.co.uk

Invacare International GmbH is a member of the DHCare group of companies. Invacare International GmbH reserves the right to modify existing product models at any time. Product specifications may vary across regions. ©2026 Invacare International GmbH. All rights reserved. Mo-vis_EU_26082026

