


Posture belt

Rea® Spirea2  
 Rea® Spirea3  
 Rea® Spirea4NG  
 Rea® Focus  
 Invacare® SpinX



- Ⓐ = M6x55
- Ⓑ = M6x45
- Ⓒ = M6x12

<p>Rea® Spirea2 16"</p> <p>1543658</p>	<p>Rea® Spirea2 20" - 24"</p> <p>1543658</p>	<p>Rea® Spirea3</p> <p>1543659</p>
<p>Rea® Spirea4NG</p> <p>1543660</p>	<p>Invacare® SpinX</p> <p>1543662</p>	<p>Rea® Focus</p> <p>1543660</p> <p>Rea® Focus wide: Ⓐ              Rea® Focus Laguna: Ⓑ</p>

