

# Mover Range

Transfer Platform



**Active Transfer Solutions**

Smooth, confident sit-to-stand transfers in any setting



# Designed for active transfers

The Mover range provides smooth, convenient and dignified transfers between bed, wheelchair, toilet or chair. Built for active users who can bear weight on their legs, both models help patients stand up safely while supporting caregivers with ergonomic handling – every day, in every room.

Both Mover Plus and Mover Aqua share a narrow base frame, double castors for easy manoeuvring in confined spaces, and comfortable cushions. Constant eye contact between user and carer keeps the transfer calm and confident.

Choose Mover Plus for daily indoor transfers in care settings and at home. Choose Mover Aqua when transfers happen in bathrooms, showers or other wet environments where hygiene matters most.



> Mover Plus



> Mover Aqua



## Intended use

Mover Plus is built for everyday indoor transfers on flat surfaces. Mover Aqua is designed for use in wet environments – including showers – without compromising on safety or stability.




## User capacity

Mover Aqua and Mover Plus support up to **150 kg** – matching the demands of bariatric care and high-load environments.



## Materials & hygiene

Mover Aqua features water-resistant materials, PU cushions and stainless-steel components for easy cleaning in high-hygiene settings. Mover Plus uses standard padded cushions ideal for everyday indoor care.



## Technical data ▼

More information about this product, including the user manual, please visit our website: [www.dhcare.com](http://www.dhcare.com)



Max user weight ▼



Width - external ▼



Width - step ▼



Height ▼

Mover Aqua	150 kg	63.5cm	38cm	107cm
Mover Plus	150 kg	63.5cm	38cm	107cm



Length ▼



Seat height ▼



Weight ▼

Mover Aqua	96cm	80cm	28 kg
Mover Plus	96cm	80cm	28 kg

Tel: +xxx  
[xxx@dhcare.com](mailto:xxx@dhcare.com)  
[www.dhcare.com](http://www.dhcare.com)

DHCare Ltd.  
4D/E, The Ark, 201 Talgarth Road,  
Hammersmith,  
London, W6 8BJ, UK

