

Ocean Ergo

ACCESSORIES COMPARISON MATRIX

Compatibility key:

Ocean Ergo only
Ocean only
Ocean & Ocean Ergo

Standard soft seat, soft seat insert, pan and pan holder are standard features across all models

Accessories Ocean/Ocean Ergo	Code	New Aquatec Ocean Ergo				Aquatec Ocean			
		Ocean Ergo	Ocean Ergo SP	Ocean Ergo XL	Ocean Ergo XL SP	Ocean Transit	Ocean SP	Ocean XL	Ocean XL SP*
Soft seat Ocean Ergo	AP1603142	✓	✓	✓	✓	✗	✗	✗	✗
Soft seat Ocean	16341	✗	✗	✗	✗	✓	✓	✓	✓
Soft seat insert Ocean Ergo	AP1603143	✓	✓	✓	✓	✗	✗	✗	✗
Soft seat insert Ocean	16343	✗	✗	✗	✗	✓	✓	✓	✓
Small soft seat Ocean Ergo	AP1603145	✓	✓	✓	✓	✗	✗	✗	✗
Small soft seat Ocean	16342	✗	✗	✗	✗	✓	✓	✓	✓
Small soft seat insert Ocean Ergo	AP1603152	✓	✓	✓	✓	✗	✗	✗	✗
Splash guard / pommel Ocean Ergo	AP1603155	✓	✓	✓	✓	✗	✗	✗	✗
Splash guard / pommel Ocean	16511	✗	✗	✗	✗	✓	✓	✓	✓
Universal soft seat Ocean Ergo	AP1603154	✓	✓	✓	✓	✗	✗	✗	✗
Universal soft seat Ocean	1526032	✗	✗	✗	✗	✓	✓	✓	✓
XL back Ocean Ergo	AP1604015	✓	✓	✓	✓	✗	✗	✗	✗
XL back Ocean	1534979	✗	✗	✗	✗	✓	✓	✓	✓
Toilet Pan holder Ocean Ergo	AP1602078	✓	✓	✓	✓	✗	✗	✗	✗
Toilet Pan holder Ocean	1527985	✗	✗	✗	✗	✓	✓	✓	✓
Front anti-tipper Ocean Ergo Models	AP1601986	✓	✓	✓	✓	✗	✗	✗	✗
Comfort soft backrest Ocean Ergo Models	AP1604136	✓	✓	✓	✓	✗	✗	✗	✗
Comfort soft backrest Ocean Ergo XL Models	AP1604016	✓	✓	✓	✓	✗	✗	✗	✗
Ergonomic soft seat	1531259	✓	✓	✓	✓	✓	✓	✓	✓
Variable soft seat	1544469	✓	✓	✓	✓	✓	✓	✓	✓
Ocean skirt	1627829	✓	✓	✓	✓	✓	✓	✓	✓
Plastic backrest	1575299	✓	✓	✓	✓	✓	✓	✓	✓
Soft backrest cushion	1558769	✓	✓	✓	✓	✓	✓	✓	✓
Armrest locking device	1470077	✓	✓	✓	✓	✓	✓	✓	✓
Chest belt	1470081	✓	✓	✓	✓	✓	✓	✓	✓
Pelvic belt	1470082	✓	✓	✓	✓	✓	✓	✓	✓
Safety bar	10167-10	✓	✓	✗	✗	✓	✓	✗	✗
Safety bar XL	1470078	✗	✗	✓	✓	✗	✗	✓	✓
Lateral support	1535077	✓	✓	✓	✓	✓	✓	✓	✓
24" wheels (set)	16379	✓	✓	✓	✓	✓	✓	✓	✓
Calf and amputee support L/H	AP1643574	✓	✓	✓	✓	✓	✓	✓	✓
Calf and amputee support R/H	AP1643573	✓	✓	✓	✓	✓	✓	✓	✓
Calf support with foot support L/H	AP1643572	✓	✓	✓	✓	✓	✓	✓	✓
Calf support with foot support R/H	AP1643543	✓	✓	✓	✓	✓	✓	✓	✓
Calf Support	1531819	✓	✓	✓	✓	✓	✓	✓	✓
Shortened leg rests L/H	1591006	✓	✓	✓	✓	✓	✓	✓	✓
Shortened leg rests R/H	1591003	✓	✓	✓	✓	✓	✓	✓	✓
Toilet Pan	10230	✓	✓	✓	✓	✓	✓	✓	✓

New accessories for the Aquatec Ocean Ergo ▼



► **Standard soft seat**

With key shaped cut out. Ergonomically shaped.



► **Soft seat insert**

Can be used with or without the soft seat when commode function is not required.



► **Small soft seat**

With small, oval opening for smaller and young adults. Ergonomically shaped.



► **Small soft seat insert**

Can be used with the small soft seat when the commode function is not required.



► **Splash guard**

Clips on easily for use with the standard plastic seat, standard soft seat and small soft seat. May also help as a pommel.



► **Universal soft seat**

Comfortable, closed seat surface.



► **XL back**

Increases the width between the armrests by 80mm. Can be retrofitted to standard Ocean Ergo models.



► **Toilet pan**

Ergonomic handle allows removal of the pan from the holder easily

Capacity scale inside the pan

Lid holder ensures pan and lid stay together



► **Comfort backrest**

For additional comfort and easy cleaning.

Find out more at www.safepatienthandling.co.uk