



The *Invacare H440 Omega Eco* is the right choice if you want to combine function, comfort and design. With the snap button height adjustment and the innovative toilet recess insert, this toilet chair is simple and effortless to use for everyone.



- Modern and light design for the perfect fit in the living environment
- Adjustable seating height from 405 mm to 585 mm – simple to operate
- Easy handling with the recess insert
- Perfect over the toilet position
- Easy to clean cover material
- High load capacity: 140 kg

Technical data

Order numbers	H440 Omega Eco
Commode chair, 5 pack	1542903

Dimensions	
Seat width	450 mm
Seat depth	430 mm
Seat height	405 - 585 mm
Total width	580 mm
Total depth	530 mm
Total height	790 - 970 mm
Height to armrests	620 - 800 mm
Height, ground to lower seat edge:	
With holder for toilet pan	335 - 515 mm
Without holder for toilet pan	345 - 525 mm
Height floor to seat plate	405 - 585 mm
Width between the armrests	470 mm
Width, cut out in seat plate	200 mm
Depth, cut out in seat plate	255 mm

More info	
Load capacity	140 kg
Weight (kg)	7,5 kg
Colour	Lavender blue

Invacare International Sàrl
 Route de Cité-Ouest 2 - 1196 Gland
 Switzerland
 Tel: +41 (0)22 354 60 10
 Fax: +41 (0)22 354 60 11
hqeuropa@invacare.com
www.invacare.eu.com

© 2012 Invacare International Sàrl
 All rights reserved.

All information quoted is believed to be correct at time of print.
 Invacare reserves the right to alter product specifications without prior consultation.

Omega Eco - EU - 08/2012