

- Standard
 Optional field

Mandatory fields:¹⁾
 The order cannot be processed without details!

- ORDER**
 COST ESTIMATE

Invoice address

(Stamp or completely legible)

 Cust. No.

 Company

 Street

 ZIP/Town /
Shipping address (if different to the invoice address)

 Company

 Street

 ZIP/Town /
Type of order

- Mounting by customer²⁾ (Please fill out the "Wheelchair details" section completely!)
 Mounting by Alber in the factory³⁾ 1565790 100.00 EUR

 Consignment

twion M24 1565810 *

Scope of delivery:

- 2 drive wheels, 6 km/h, 24" wheel size
- Standard model: 6.3 kg per wheel
- 2 built-in battery packs in the wheels
- 1 charger

Wheelchair details

 Wheelchair manufacturer
 (e.g. Kirschall)

 Model (e.g. Compact)

 Year of construction

 Variants (e.g. Hemi)



 Drum brakes⁵⁾ Yes No
 ↳ Note: Drum brakes must be removed!


Equipped with vehicle retention system (please consult with Alber Customer Service)

Rear wheel extension available

Available wheel position

- 1 = rear
 2 = centre
 3 = front



 Seat width (seat tube measured on the outside) cm


 Seat height (rear, measured to lower edge of seat sling) cm

Important accessories

- Bracket to fit on the wheelchair 1467178 205.00 EUR
 (in a few cases price could be different)
- Anti-tippers with jack-up function 1489214 (pair) 262.00 EUR

Comments (mounting, etc.)

.....

.....

.....

* Please refer to your price list.



- Standard **Mandatory fields:¹⁾**
 Optional field The order cannot be processed without details!

Cover

- | | | | |
|--|---------|--------|-----------|
| <input checked="" type="radio"/> White cover | 1565889 | (pair) | 0.00 EUR |
| <input type="checkbox"/> Black cover | 1565888 | (pair) | 25.00 EUR |

Push-rim

- | | | | |
|--|---------|--------|----------|
| <input checked="" type="radio"/> Matt aluminium, coated black | 1565890 | (pair) | 0.00 EUR |
| <input type="checkbox"/> Stainless steel, brushed (+400 g per wheel) | 1565891 | (pair) | 0.00 EUR |

Tyres

- | | | | |
|---|---------|--------|-----------|
| <input checked="" type="radio"/> Schwalbe Marathon Plus Evolution | 1565894 | (pair) | 0.00 EUR |
| <input type="checkbox"/> Schwalbe Ultremo ZX (-300 g per wheel) | 1565893 | (pair) | 40.00 EUR |
| <input type="checkbox"/> Schwalbe Airborne (-100 g per wheel) | 1565892 | (pair) | 40.00 EUR |

Accessories

- | | | |
|--|---------|------------|
| <input type="checkbox"/> Mobility Plus package | 1565789 | 327.00 EUR |
|--|---------|------------|

Only for smartphones with the Android operating system and able to connect to the internet.

Comprising: Outer ring with exclusive design, booklet with operating instructions, licence code

Functions:

- Reading off driving parameters including battery status, speed, number of kilometres covered
- 2 additional ECO and SPORT modes
- Increase assist up to 10 km/h
- BlueDrive remote control using a smartphone
- Records journeys and trip profiles

- | | | |
|--|---------|-----------|
| <input type="checkbox"/> Brake extension | 1467780 | 22.00 EUR |
|--|---------|-----------|

Manual wheel axle retrofit^{4) 5)} (compatible with twion wheels)

- | | | | |
|---|---------|--------|------------|
| <input type="checkbox"/> Retrofit, stub axle to \varnothing 1/2" quick-release axle | 1565005 | (pair) | 85.00 EUR* |
|---|---------|--------|------------|

*in a few cases price could be different

Date / / Signature of ordering party

As at January 1, 2014. All prices are subject to the statutory VAT, and include packaging ex works. The General Terms and Conditions of Ulrich Alber GmbH shall apply, subject to changes.

- 1) In order to ensure that your order is processed smoothly and promptly, we require ALL the information in the shaded fields.
- 2) Mounting kits can only be sent to trained specialist dealers (mounting training).
- 3) Please send the wheelchair carriage with a copy of this order as soon as possible and carriage-paid to Alber.

- 4) If the manual wheels are used in combination with twion, these must be fitted with 1/2" stub axles, without drum brakes and 1" wide tyres.
- 5) If drum brakes are installed, these must be removed as they are not compatible with twion drive wheels. Standard lever brakes are used instead.

