



INVACARE France OPERATIONS SAS
Route de Saint Roch, 37230 Fondettes, France
Tel: 33 (0)247 626 466 Fax: 33 (0)247 626 488

Invacare® Action³NG CE
Local Homologation Number :
ISO - MDA , TUV Logos
Crossbar Folding Mechanism, Seat Angle Possibility : 0°-3°, 6° (Delivered Always 3°).
Rear Wheel Position : Standard, Passive, Active, Heel Straps and Tool Kit.
User Maximum Weight : 125Kg
Price In Local Currency & Euros (If Required) Of The Standard Product.

Customer Name :
Dealer :
Account No :
Address :
Tel No :
Delivery address :
Order No :
Order Qty :
Order Date :

Please make your selection with a cross in the box : Standard feature Free of Charge (NC9999) Non Conformance with ACTION3 Option 9999.

Gjeldende fra 1. juli 2018

→ Customer weight =		Kg
Code	Description	
MODEL		
CEB0000	<input checked="" type="checkbox"/> ACTION ³ NG	(NC : 0491, 0691, 1737, 0791, 1837, 1151, 1181, 1182, 1183, 1211, 1212, 1481)
SPECIAL		
N	<input type="checkbox"/> NO SPECIAL	(NC: 0991)
CEBSPEC	<input type="checkbox"/> SPECIAL	
BASE MODEL FOR PRICE LIST		
CEB0026	<input checked="" type="checkbox"/> BASE MODEL ACTION ³ NG	
ALBER		
CEB0080	<input type="checkbox"/> BRACKET FOR VIAMOBIL V25	(NC : 0110, 0140, 0150, 0160, 0170, 0300, 0310, 0315, 0461, 0462, 0463, 0464, 0466, 0908, 1170, 1201, 1212, 1331, 1332, 1345, 1346, 1365, 1370, 1375, 1430, 1440, 1441)
CEB0090	<input type="checkbox"/> BRACKET FOR VIAMOBIL ECO	(NC : 0140, 0150, 0160, 0170, 0290, 0300, 0310, 0315, 0461, 0462, 0463, 0464, 0466, 0908, 0940, 1201, 1212, 1331, 1332, 1345, 1346, 1365, 1370, 1375, 1430, 1440, 1441, 1451)
PROPULSION TYPE		
CEB0100	<input checked="" type="checkbox"/> SELF PROPULSION	(NC : 1120, 1130)
CEB0110	<input type="checkbox"/> TRANSIT PROPULSION WITH RH TIPPING AID	(NC : 0251, 0290, 0300, 0310, 0315, 0351, 0428, 0431, 0511, 0561, 0738, 0739, 1051, 1052, 1060, 1070, 1080, 1081, 1090, 1095, 1100, 1110, 1114, 1116, 1117, 1118, 1119, 1140, 1150, 1051, 1180, 1181, 1182, 1183, 1185, 1190, 1212, 1260, 1261, 1275, 1320, 1330, 1331, 1332, 1365, 1390)
CEB0140	<input type="checkbox"/> DUAL HANDRIM RIGHT HAND PROPULSION (NEW DHR SYSTEM), QR AXLE	(NC : 0300, 0310, 0315, 0351, 0361, 0511, 0561, 0738, 0739, 1051, 1070, 1081, 1090, 1095, 1100, 1110, 1114, 1116, 1117, 1118, 1119, 1130, 1140, 1150, 1151, 1170, 1181, 1182, 1183, 1185, 1212, 1260, 1261, 1280, 1330, 1331, 1332, 1365)
CEB0150	<input type="checkbox"/> DUAL HANDRIM LEFT HAND PROPULSION (NEW DHR SYSTEM), QR AXLE	(NC : 0300, 0310, 0315, 0351, 0361, 0511, 0561, 0738, 0739, 1051, 1070, 1081, 1090, 1095, 1100, 1110, 1114, 1116, 1117, 1118, 1119, 1130, 1140, 1150, 1151, 1170, 1181, 1182, 1183, 1185, 1212, 1260, 1261, 1280, 1320, 1331, 1332, 1365)
CEB0160	<input type="checkbox"/> OAD BY LEVER RH PROPULSION	(NC : 0251, 0290, 0300, 0310, 0315, 0351, 0361, 0511, 0561, 0640, 0650, 0655, 0695, 0696, 0732, 0738, 0739, 0755, 1051, 1052, 1070, 1081, 1090, 1095, 1100, 1110, 1114, 1116, 1117, 1118, 1119, 1130, 1140, 1150, 1151, 1180, 1181, 1182, 1183, 1185, 1190, 1201, 1212, 1260, 1261, 1280, 1320, 1330, 1331, 1332, 1346, 1365, 1390)
CEB0170	<input type="checkbox"/> OAD BY LEVER LH PROPULSION	(NC : 0251, 0290, 0300, 0310, 0315, 0351, 0361, 0511, 0561, 0640, 0650, 0655, 0695, 0696, 0738, 0739, 0740, 0750, 0755, 0795, 0796, 0832, 1051, 1052, 1070, 1081, 1090, 1095, 1100, 1110, 1114, 1116, 1117, 1118, 1119, 1130, 1140, 1150, 1151, 1180, 1181, 1182, 1183, 1185, 1190, 1201, 1212, 1260, 1261, 1280, 1320, 1330, 1331, 1332, 1346, 1365, 1390)
SEAT		
SEAT WIDTH		
CEB0180	<input type="checkbox"/> 380 MM SEAT WIDTH	(NC : 0251, 0361, 0372, 0421, 0441, 0442, 0851, 0852, 0855, 0881, 1481)
CEB0190	<input type="checkbox"/> 405 MM SEAT WIDTH	(NC : 0251, 0361)
CEB0200	<input type="checkbox"/> 430 MM SEAT WIDTH	(NC: 0853)
CEB0210	<input type="checkbox"/> 455 MM SEAT WIDTH	(NC: 0853)
CEB0220	<input type="checkbox"/> 480 MM SEAT WIDTH	(NC : 0361, 0372, 0421, 0441, 0442, 0881, 0851, 0852, 0853, 0854, 1481)
CEB0230	<input type="checkbox"/> 505 MM SEAT WIDTH	(NC: 0851, 0852, 0853, 0854, 0855, 0856)
SEAT DEPTH		
CEB0240	<input type="checkbox"/> 400 MM SEAT DEPTH	(NC: 0361, 0879, 0883, 1181)
CEB0250	<input type="checkbox"/> 450 MM SEAT DEPTH	(NC: 0853, 0878, 0883)
CEB0251	<input type="checkbox"/> 500 MM SEAT DEPTH (with Active Rear Wheel Support ONLY).	(NC : 0320, 0330, 0350, 0351, 0361, 0372, 0851, 0852, 0853, 0854, 0855, 0856, 0869, 0871, 0873, 0874, 0878, 0879, 0908, 1201, 1365)
CEB0252	<input type="checkbox"/> 400 TO 450 MM ADJUSTABLE SEAT DEPTH SET AT 430 MM	(NC : 0320, 0330, 0350, 0351, 0361, 0372, 0883)
SEAT HEIGHT		
CEB0260	<input type="checkbox"/> 510 MM SEAT HEIGHT	(NC : 0511, 0561, 1737, 0738, 0739, 1837, 1049, 1050, 1051, 1052, 1100, 1110, 1117)
CEB0270	<input checked="" type="checkbox"/> 485 MM SEAT HEIGHT	(NC: 0511, 0561, 1051)
CEB0280	<input type="checkbox"/> 460 MM SEAT HEIGHT	(NC : 1000, 1010, 1020, 1030, 1130)
CEB0290	<input type="checkbox"/> 435 MM SEAT HEIGHT	(NC : 0935, 1000, 1010, 1020, 1030, 1052, 1130, 1331, 1332, 1365, 1370, 1440, 1441)
CEB0300	<input type="checkbox"/> 410 MM SEAT HEIGHT	(NC : 0351, 0361, 0511, 0561, 0660, 0670, 1737, 0738, 0739, 0760, 0770, 1837, 0838, 0839, 0935, 1000, 1010, 1020, 1030, 1049, 1050, 1052, 1060, 1070, 1080, 1081, 1090, 1095, 1114, 1118, 1119, 1130, 1181, 1182, 1280, 1331, 1332, 1365, 1370, 1440, 1441)
CEB0310	<input type="checkbox"/> 385 MM SEAT HEIGHT	(NC : 0351, 0361, 0511, 0561, 0600, 0611, 0620, 0660, 0670, 1737, 0738, 0739, 0760, 0770, 1837, 0838, 0839, 0935, 1000, 1010, 1020, 1030, 1049, 1050, 1052, 1060, 1070, 1080, 1081, 1090, 1095, 1099, 1100, 1110, 1114, 1117, 1118, 1119, 1130, 1280, 1331, 1332, 1365, 1375, 1440, 1441)
CEB0315	<input type="checkbox"/> 360 MM SEAT HEIGHT, 0° SEAT ANGLE (AVAILABLE WITH LOW FORK, 5" FRONT WHEEL & 20" REAR WHEELS)	(NC : 0351, 0361, 0511, 0561, 0600, 0611, 0620, 0660, 0670, 1737, 0738, 0739, 0760, 0770, 1837, 0838, 0839, 0935, 1000, 1010, 1020, 1030, 1049, 1050, 1052, 1060, 1070, 1080, 1081, 1090, 1095, 1099, 1100, 1110, 1114, 1117, 1118, 1119, 1130, 1280, 1331, 1332, 1365, 1375, 1440, 1441)



Code	Description		
SEATING OPTIONS TYPE			
CEB0840	NO SEATING OPTION	(NC : 0849, 0851, 0852, 0853, 0854, 0855, 0856, 0869, 0871, 0873, 0874, 0875, 0876, 0878, 0879, 0881, 0882, 0883, 0905, 0906, 0907, 0908, 0930, 0935, 0936, 0940, 0945, 0952)	
SEAT CUSHION TYPE			
CEB0849	BASIC 50 MM FOAM CUSHION	NYLON & ANTI-SLIP COVER (NC: 0952)	
CEB0869	MATRIX FLO-TECH LITE CUSHION	(NC : 0883, 0952)	
CEB0871	MATRIX FLO-TECH CONTOUR VISCO CUSHION	(NC : 0883, 0952)	
CEB0873	MATRIX FLO-TECH IMAGE CUSHION	(NC : 0883, 0952)	
CEB0874	MATRIX FLO-TECH SOLUTION XTRA CUSHION	(NC : 0883, 0952)	
CEB0875	MATRIX VI CUSHION	ADDITIONAL DELIVERY TIME + 5 DAYS (NC: 0952)	
CEB0876	MATRIX LIBRA CUSHION	ADDITIONAL DELIVERY TIME + 5 DAYS (NC: 0952)	
SEAT CUSHION DEPTH			
CEB0878	FOR 400 MM SEAT DEPTH		
CEB0879	FOR 450 MM SEAT DEPTH		
CEB0883	FOR 500 MM SEAT DEPTH	(NC: 0852)	
OTHER SEATING OPTIONS			
CEB0881	TRUNK SUPPORT COMFORT DARTEX (PAIR)	ONLY FOR LAGUNA SHELL	
CEB0882	LATERAL POSITIONING PADS (PAIR)		
CEB0905	NECKREST DARTEX (INCLUDING BACKBRACE EXCEPT FOR LAGUNA)	(NC : 0906, 0907, 0908, 0940, 0945)	
CEB0906	HEADREST DARTEX (INCLUDING BACKBRACE EXCEPT FOR LAGUNA)	SIDWISE ADJUSTABLE (NC : 0907, 0908, 0940, 0945)	
CEB0907	HEADREST WITH CHEEK SUPPORT DARTEX (INCLUDING BACKBRACE EXCEPT FOR LAGUNA)	SIDWISE ADJUSTABLE (NC : 0908, 0940, 0945)	
CEB0908	MATRIX SOFT PAD HEADREST (SOFT PAD 250MM) WITH FIXED HARDWARE	(ONLY FOR MATRIX BACK) (NC : 0940, 0945)	
CEB0930	ADJUSTABLE CALF STRAP		
CEB0935	LONGER FOOTREST TUBE (UL= 450-530 MM)	(NC:0936)	
CEB0936	SHORTER FOOTREST SUPPORT (PAIR) (UL= 230-430 MM)		
CEB0940	BACKBRACE	(NC: 0945)	
CEB0945	FOLDABLE BACKBRACE	(NC: 1451)	
COLOURS			
CEB0961	SLATE GREY (RAL7024 FINE TEXTURE)		
CEB0964	PEARL WHITE		
CEB0965	ELECTRIC GREEN		
CEB0966	MATT BLACK		
CEB0967	BLUE POP		
CEB0968	HYPNOTIC PURPLE		
CEB0969	RED MATT		
CEB0980	ANTHRACITE BLACK		
CEB0991	SPECIAL COLOUR	ON DEMAND ONLY	
CASTER HOUSING TYPE			
CEB0995	STANDARD		
CEB0996	LIGHT WEIGHT		
FRONT WHEELS			
CEB1000	200 X 32 MM (8" X 1.1/4")	PNEUMATIC TYRE	
CEB1010	200 X 32 MM (8" X 1.1/4")	BLACK SOLID TYRE	
CEB1020	200 X 50 MM (8" X 2")	PNEUMATIC TYRE	
CEB1030	200 X 50 MM (8" X 2")	BLACK SOLID TYRE	
CEB1049	150 X 32 MM (6" X 1.1/4")	SOLID TYRE MEDIUM HARD	
CEB1050	150 X 32 MM (6" X 1.1/4")	BLACK SOLID TYRE	
CEB1051	125 X 25 MM (5" X 1")	SOLID TYRE	
CEB1052	140 X 45 MM (5" 1/2 X 1" 3/4)	SOLID TYRE MEDIUM HARD	
REAR WHEELS			
REAR WHEEL TYPE			
CEB1060	610 X 35 MM (24" x 1.3/8")	PNEUMATIC TYRE	SPOKED RIM (NC: 1261, 1181, 1182)
CEB1070	610 X 35 MM (24" x 1.3/8")	PNEUMATIC TYRE	PLASTIC RIM (NC : 1151, 1181, 1182, 1260, 1261, 1280, 1331, 1332)
CEB1080	600 X 35 MM (24" x 1.3/8")	SOLID TYRE	SPOKED RIM (NC: 1181, 1182, 1261, 1320, 1330, 1331, 1332)
CEB1081	600 X 35 MM (24" x 1.3/8")	BLACK SOLID TYRE	SPOKED RIM (NC: 1181, 1182, 1261, 1320, 1330, 1331, 1332)
CEB1090	600 X 35 MM (24" x 1.3/8")	SOLID TYRE	PLASTIC RIM (NC : 1151, 1181, 1182, 1260, 1261, 1280, 1320, 1330, 1331, 1332)
CEB1095	595 X 25 MM (24" x 1")	NOTCHED PNEUMATIC TYRE	SPOKED RIM (NC : 1181, 1182, 1261, 1280, 1320, 1330, 1331, 1332)
CEB1100	560 X 35 MM (22" x 1.3/8")	PNEUMATIC TYRE	SPOKED RIM (NC : 1151, 1181, 1182, 1260)
CEB1110	555 X 35 MM (22" x 1.3/8")	SOLID TYRE	SPOKED RIM (NC : 1151, 1181, 1182, 1260, 1320, 1330, 1331, 1332)
CEB1116	515 X 35 MM (20" x 1.3/8")	SOLID TYRE	SPOKED RIM (NC: 1151, 1170, 1181, 1182, 1260, 1280, 1320, 1330, 1331, 1332)
CEB1117	560 X 25 MM (22" x 1") LIGHTWEIGHT	PNEUMATIC TYRE	RADIAL SPOKED (NC: 1151, 1170, 1201, 1230, 1260, 1261, 1280, 1320, 1330, 1331, 1332)
CEB1118	595 X 25 MM (24" x 1") LIGHTWEIGHT	PNEUMATIC TYRE	RADIAL SPOKED (NC: 1170, 1201, 1230, 1260, 1261, 1280, 1320, 1330, 1331, 1332)
CEB1119	595 X 25 MM (24" x 1") LIGHTWEIGHT	MARATHON PLUS PNEUMATIC TYRE	RADIAL SPOKED (NC: 1170, 1201, 1230, 1260, 1261, 1280, 1320, 1330, 1331, 1332)
CEB1114	595 X 25 MM (24" x 1.3/8") LIGHTWEIGHT	BLACK SOLID TYRE	RADIAL SPOKED (NC: 1170, 1201, 1260, 1261, 1280, 1320, 1330, 1331, 1332)
CEB1130	305 X 45 MM (12" x 2")	SOLID TYRE	PLASTIC RIM (NC : 1140, 1150, 1151, 1180, 1181, 1182, 1260, 1280, 1320, 1330, 1331, 1332, 1380, 1390)
HANDRIM TYPE			
CEB1140	ALUMINIUM PRE-ANODIZED HANDRIM		
CEB1150	HANDRIM COATED IN BLACK ANTI-SLIP PLASTIC		
REAR WHEELS AXLE TYPE			
CEB1170	FIXED AXLE	(NC : 1181, 1182, 1185, 1365)	
CEB1180	QUICK RELEASE AXLE		
REAR WHEELS POSITION			
CEB1185	ACTIVE REAR WHEELS POSITION	37.5MM FRONTWARDS/BACKREST CANE AXLE (with Active Rear Wheel Support) (NC: 1320, 1330, 1331, 1332, 1345, 1346, 1365)	
CEB1190	STANDARD REAR WHEELS POSITION	25MM FRONTWARDS/BACKREST CANE AXLE (12.5MM with Active Rear Wheel Support and Seat Depth 500MM)	
CEB1200	PASSIVE REAR WHEELS POSITION	25MM BACKWARDS/BACKREST CANE AXLE (12.5MM with Active Rear Wheel Support and Seat Depth 500MM)	
CEB1201	DUAL AMPUTEE REAR WHEELS POSITIONS	TRANSIT=75MM BACKWARDS/BACKREST CANE AXLE, ADJUSTABLE TO 50MM BY (NC: 1212, 1275, 1280, 1320, 1330, 1331, 1332, 1365)	
PUMP			
CEB1220	WITHOUT PUMP		
CEB1230	PUMP AND SUPPORT		
CEB1235	PUMP - HIGH PRESSURE		
CHASSIS OPTIONS TYPE			
CEB1240	NO CHASSIS OPTION	(NC : 1241, 1260, 1261)	
CEB1241	FRONT WHEELS EXTENSION (+ 50 MM OUTWARD)		
CEB1260	SPOKE GUARDS PAIR 600 MM (24") BLACK	(NC: 1261)	
CEB1261	SPOKE GUARDS PAIR 560 & 515 MM (22" & 20") CLEAR		
BRAKES TYPE			
CEB1270	PUSH TO LOCK BRAKES		
CEB1275	PULL TO LOCK BRAKES		
CEB1280	ATTENDANT HUB BRAKES (INCLUDING PUSH TO LOCK BRAKES)		
CEB1320	OAD WIRE CONTROLLED BRAKE RH	(NC : 1390)	
CEB1330	OAD WIRE CONTROLLED BRAKE LH	(NC : 1390)	
SAFETY OPTIONS TYPE			
CEB1340	NO SAFETY OPTION	(NC : 1345, 1346, 1360, 1361, 1362, 1363, 1364, 1365, 1370, 1375, 1380, 1390)	
CEB1345	TRANSPORTATION KIT INCLUDING POSTURE BELT	(NC: 1360, 1363, 1364)	
CEB1346	COMFORT TRANSPORTATION KIT INCLUDING POSTURE BELT & BACK BRACE	(Headrest & Trunk Supports have to be removed) (NC: 1360, 1363, 1364)	
CEB1360	POSTURE BELT	(NC: 1361, 1362, 1363, 1364)	
CEB1361	HIP BELT, 2-PT, CTR-PULL, PB	M X 460MM PADS, FLAT-MOUNT (NC: 1362, 1363, 1364)	
CEB1362	HIP BELT, 2-PT, CTR-PULL, PB	L X 620MM PADS, FLAT-MOUNT (NC: 1363, 1364)	
CEB1363	HIP BELT, 4-PT, CTR-PULL, PB	M X 460MM PADS, FLAT-MOUNT (NC: 1363, 1364)	
CEB1364	HIP BELT, 4-PT, CTR-PULL, PB	L X 620MM PADS, FLAT-MOUNT (NC: 1363, 1364)	
CEB1365	TRANSFER KIT (100 MM (4") WHEELS)	ANTI-BARRIER	
CEB1370	REAR ANTI-TIPPERS PAIR	(NC : 1375, 1440, 1441)	
CEB1375	REAR ANTI-TIPPERS PAIR WITH TIPPING AID	(REA TYPE) (NC : 1440, 1441)	
CEB1390	REFLECTORS KIT		
CEB1390	BRAKE EXTENSION PAIR		
ATTENDANT OPTIONS TYPE			
CEB1420	NO ATTENDANT OPTION	(NC : 1430, 1440, 1441, 1451)	
CEB1430	CANE HOLDER		
CEB1440	TIPPING AID RH		
CEB1441	TIPPING AID LH		
CEB1451	DRIP STAND SUPPORT WITH DRIPSTAND (WITH BACKBRACE 0940 ONLY)		
ACCESSORIES TYPE			
CEB1480	NO ACCESSORY	(NC : 1480, 1481, 1482, 1483, 1484, 1485)	
CEB1480	ACRYLIC TRAY	(NC: 1481, 1482, 1483, 1484, 1485)	
CEB1483	HALF TABLE TRAY, CLEAR, FLIP UP, RH	(NC: 1482, 1484)	
CEB1484	HALF TABLE TRAY, CLEAR, FLIP UP, LH	(NC: 1482, 1484)	
CEB1482	LOCKING DEVICE FO COMFORT TRAY	(WITH CEA1481 ONLY)	
CEB1485	COMFORT PADDING FOR HALF TABLE TRAY	(WITH CEA1483 OR CEA1484 ONLY)	
FREIGHT & PRICE			
FREIGHT PRICE			
TOTAL PRICE			
January 2018 (3/3)			

