



# Manual Wheelchair European Prescription Form

EPF\_S3\_XLT\_Max\_180701

Created/Modified by: Jasmine Bernasconi

INVACARE® XLT Max



Local Homologation number, ISO, MDA, TUV logos

We reserve the rights to have a + or - tolerance of 1 cm on the following measurements: seat width, seat depth, SHf, SHr and UL.

We reserve the rights to make technical changes without notice !

Customer Name :

Dealer :

Account No :

Address :

Tel No :

Delivery address :

Order No :

Order Qty :

Order Date :

please make your selection with a cross in the box  Standard feature  Free of Charge (NC9999) non conformance with option 9999

→ Max Customer weight = 180 Kg

Code	Description	Non conformance	
BTA0010	Invacare XLT Max		<input checked="" type="checkbox"/>
BTA8000	Special (Full description and drawing by fax or E-mail)		<input type="checkbox"/>

**PROPULSION TYPE**

BTA0020	Standard wheel hub		<input checked="" type="checkbox"/>
BTA0030	Assistant manoeuvred drumbrakes	0570;0580;0630;0970;1050	<input type="checkbox"/>

**FRAME TYPE**

BTA0040	Active frame, knee angle 75°		<input checked="" type="checkbox"/>
---------	------------------------------	--	-------------------------------------

**SEAT DIMENSIONS****SEAT WIDTH**

BTA0050	505 mm		<input type="checkbox"/>
BTA0060	555 mm		<input type="checkbox"/>
BTA0070	605 mm		<input type="checkbox"/>

**SEAT DEPTH**

BTA0080	500 mm (adjustable from seat depth 460 to 500 mm)		<input type="checkbox"/>
BTA0090	550 mm (adjustable from seat depth 510 to 550 mm)		<input type="checkbox"/>
BTA0100	600 mm (adjustable from seat depth 560 to 600 mm)	0650;	<input type="checkbox"/>

**SEAT HEIGHT****FRONT SEAT TO FLOOR HEIGHT = SHf**

BTA0110	440 mm	0200;0260;0270;0410;0420;0430;0440;0895; 0900; 0910;0920;0930	<input type="checkbox"/>
BTA0120	450 mm	0200;0270;0410;0420;0430;0440;0895;0900; 0910; 0920;0930	<input type="checkbox"/>
BTA0130	460 mm	0200;0270;0410;0420;0430;0440;0880;0920;0930	<input type="checkbox"/>
BTA0140	470 mm	0200;0420;0430;0440;0920;0930	<input type="checkbox"/>
BTA0150	480 mm	0430;0440	<input type="checkbox"/>
BTA0160	490 mm	0430;0440	<input type="checkbox"/>
BTA0170	500 mm	0440	<input type="checkbox"/>
BTA0180	510 mm	0440	<input type="checkbox"/>
BTA0190	520 mm	0200;0210;0880	<input type="checkbox"/>

**SHf / Castor Size**

5,5" : 440-510 mm  
6" : 460-520 mm  
7" : 480-520 mm

**REAR SEAT TO FLOOR HEIGHT = SHr, SHr ≤ SHf**

BTA0200	390 mm		<input type="checkbox"/>
BTA0210	400 mm		<input type="checkbox"/>
BTA0220	410 mm		<input type="checkbox"/>
BTA0230	420 mm		<input type="checkbox"/>
BTA0240	430 mm		<input type="checkbox"/>
BTA0250	440 mm		<input type="checkbox"/>
BTA0260	450 mm		<input type="checkbox"/>
BTA0270	460 mm		<input type="checkbox"/>

**SHr / Rear wheel size**

24" : 390-460 mm

**KNEE TO HEEL LENGTH, without seat cushion, UL ≤ SHf-40 mm**

BTA0280	270 mm (Only with Footrest high mounted)	0670;0680;0690	<input type="checkbox"/>
BTA0290	280 mm (Only with Footrest high mounted)	0670;0680;0690	<input type="checkbox"/>
BTA0300	290 mm (Only with Footrest high mounted)	0670;0680;0690	<input type="checkbox"/>
BTA0310	300 mm (Only with Footrest high mounted)	0670;0680;0690	<input type="checkbox"/>
BTA0320	310 mm (Only with Footrest high mounted)	0670;0680;0690	<input type="checkbox"/>
BTA0330	320 mm (Only with Footrest high mounted)	0670;0680;0690	<input type="checkbox"/>
BTA0340	330 mm (Only with Footrest high mounted)	0670;0680;0690	<input type="checkbox"/>
BTA0350	340 mm (Only with Footrest high mounted)	0670;0680;0690	<input type="checkbox"/>
BTA0360	350 mm (Only with Footrest high mounted)	0670;0680;0690	<input type="checkbox"/>
BTA0370	360 mm (Only with Footrest high mounted)	0670;0680;0690	<input type="checkbox"/>
BTA0380	370 mm (Only with Footrest high mounted)	0670;0680;0690	<input type="checkbox"/>
BTA0390	380 mm	0700;	<input type="checkbox"/>
BTA0400	400 mm	0700;	<input type="checkbox"/>
BTA0410	420 mm	0700;	<input type="checkbox"/>
BTA0420	440 mm	0700;	<input type="checkbox"/>
BTA0430	460 mm	0700;	<input type="checkbox"/>
BTA0440	480 mm	0700;	<input type="checkbox"/>

**BACKREST****BACKREST TYPE**

BTA0450	Velcro tension adjustable, height adjustable, infinitely variable angle adjustable and foldable backrest		<input checked="" type="checkbox"/>
---------	--	--	-------------------------------------

**BACKREST HEIGHT**

BTA0510	400 mm (Height adjustable from 400 to 440 mm by 20 mm steps)		<input checked="" type="checkbox"/>
BTA0520	420 mm (Height adjustable from 400 to 440 mm by 20 mm steps)		<input type="checkbox"/>
BTA0530	440 mm (Height adjustable from 400 to 440 mm by 20 mm steps)		<input type="checkbox"/>
BTA0540	450 mm (Height adjustable from 450 to 490 mm by 20 mm steps)		<input type="checkbox"/>
BTA0550	470 mm (Height adjustable from 450 to 490 mm by 20 mm steps)		<input type="checkbox"/>
BTA0560	490 mm (Height adjustable from 450 to 490 mm by 20 mm steps)		<input type="checkbox"/>

**PUSH HANDLES**

BTA0570	Fixed push handles		<input checked="" type="checkbox"/>
BTA0580	Height adjustable push handles, external, short grip		<input type="checkbox"/>
BTA0590	Height adjustable push handles, external, long grip for drum brakes		<input type="checkbox"/>
BTA0600	Height adjustable push handles, external, including neckrest, long grip for drum brakes		<input type="checkbox"/>
BTA0610	Height adjustable push handles, external, including headrest, long grip for drum brakes		<input type="checkbox"/>
BTA0620	Height adjustable push handles, external, including headrest with cheek support, long grip for drum brakes		<input type="checkbox"/>
BTA0630	Without push handles		<input type="checkbox"/>

**SEAT****SEAT TYPE**

BTA0640	Velcro tension- and depth adjustable, infinitely variable on 40 mm.		<input checked="" type="checkbox"/>
---------	---	--	-------------------------------------

**SEAT CUSHIONS**

BTA0650	Flo-tech Contour		<input checked="" type="checkbox"/>
BTA0660	Without seat cushion		<input type="checkbox"/>

LEGRESTS / FOOTRESTS			
<b>FOOTREST INCL. CALF STRAP</b>			
BTA0670	Foot bow rigid		◆
BTA0680	Foot bow rigid with cover		□
BTA0690	Foot board, depth and angle adjustable		□
BTA0700	Footrest, high mounted, for knee to heel length 270 to 370 mm		□
<b>ARMRESTS/ SIDE GUARDS</b>			
<b>ARMRESTS/ SIDE GUARD LH TYPE</b>			
BTA0710	Side guard, removable, LH		□
BTA0720	Mudguard, removable, LH		□
BTA0730	Armrest "Rio" height adjustable, short pad (280 mm), LH		◆
BTA0740	Armrest Comfort height adjustable, wide pad, LH		□
BTA0750	Without armrest, LH	1110	□
<b>ARMRESTS/ SIDE GUARD RH TYPE</b>			
BTA0760	Side guard, removable, RH		□
BTA0770	Mudguard, removable, RH		□
BTA0780	Armrest "Rio" height adjustable, short pad (280 mm), RH		◆
BTA0790	Armrest Comfort height adjustable, wide pad, RH		□
BTA0800	Without armrest, RH	1110	□
<b>FRAME</b>			
<b>FRAME COLOURS</b>			
BTA0810	Anthracite Black		◆
BTA0820	Grey		□
BTA0830	Blue		□
BTA0840	Red		□
BTA0850	Orange		□
BTA0860	Blasted titanium		□
<b>CASTOR ATTACHMENT TYPES</b>			
BTA0870	Housing with angle adjustment		◆
<b>FRONT CASTORS,</b>			
BTA0880	140x40 soft (5,5")		□
BTA0895	150x30 medium hard (6")		◆
BTA0900	150x30 soft (6")		□
BTA0910	150x32 pneumatic (6")		□
BTA0920	180x40 soft (7")		□
BTA0930	180x40 pneumatic (7")		□
<b>FORKS</b>			
BTA0940	Fixed steps		◆
<b>REAR WHEELS SIZE/TYPE</b>			
BTA0950	24" Standard		◆
<b>TYRE TYPE</b>			
BTA0960	Tyres solid		◆
BTA0970	Tyres black marathon evolution		□
<b>HANDRIM TYPE</b>			
BTA0980	Handrim aluminium		◆
BTA0990	Handrim stainless		□
BTA1000	Handrim plastic coated		□
BTA1010	Handrim foam coated		□
BTA1020	Handrim black grip		□
<b>REAR WHEEL FIXATION</b>			
BTA1030	Standard rear wheel fixation		◆
<b>CAMBER ANGLE</b>			
BTA1040	0° camber		□
BTA1050	2° camber		◆
<b>WHEEL OPTIONS</b>			
BTA1060	Spokeguard Frost		□
<b>SAFETY</b>			
<b>BRAKE TYPE</b>			
BTA1070	Standard brake		◆
<b>SAFETY OPTIONS</b>			
BTA1080	Pelvic belt		□
BTA1090	Brake lever extension, fold away, LH		□
BTA1100	Brake lever extension, fold away, RH		□
BTA1110	Auto lock for armrest		□
BTA1120	Anti-tip devices, LH (For safety reasons we recommend the use of two anti-tip devices.)	1180	◆
BTA1130	Anti-tip devices, RH (For safety reasons we recommend the use of two anti-tip devices.)	1190	◆
BTA1140	Reflectors kit (wheels and chassis)		□
<b>ACCESSORIES</b>			
BTA1150	Pump		□
BTA1160	Cane holder, LH		□
BTA1170	Cane holder, RH		□
BTA1200	Tool kit		◆
BTA1210	Frame cover		□
<b>SPECIAL</b>			
<b>FREIGHT &amp; PRICE</b>			
<b>PACKAGING</b>			
<b>FREIGHT PRICE</b>			
<b>TOTAL PRICE</b>			

**SHf / Castor Size**

5,5" : 440-510 mm  
6" : 460-520 mm  
7" : 480-520 mm

**SHr / Rear wheel size**

24" : 390-460 mm