

Orion[®] Pro

Feel the Freedom

Now with Lithium Battery
and Emergency Brake



Comfort, safety and reliability
with a powerful drive performance



Feel the difference

Our scooters have been designed together with healthcare experts, to provide unrivalled standards of quality, safety and comfort.

Through continuous dedicated research and an unrivalled experience in designing for mobility, at Invacare we pride ourselves in scooters that improve quality of life. Being able to regain confidence, trust and independence from your scooter is important.

From ergonomic steering controls to the adaptable cushioned seats and adjustable seat risers, every detail of our scooters is considered and refined to make getting about feel great.

Quality



Market-leading build quality with exceptional attention to detail.

Safety



Compliant with the latest European Medical Device Standards.

Comfort



Adjustable design to provide lasting comfort.



Feel in control with Orion® Pro

Introducing our Orion® Pro with lithium battery.
A high-performing scooter, designed for reliable performance,
comfort and safety wherever you go.

Built for longer journeys with the optional lithium
battery, the Orion Pro delivers an exceptionally long
driving range and total peace of mind.

Drive further, for longer and without worrying about
running out of power. Giving you the confidence to go
wherever you want, whenever you want!

Advanced suspension and large wheels ensure a
smooth, stable ride over uneven ground, around
corners and over obstacles. Incredibly easy to steer and
start, it keeps you confidently in control.



Improved driving control

- Ergonomic steering for added control
- Designed for all-terrain performance
- Anti-splash guards



Lithium power for ultimate freedom

- Optional high-performance lithium battery
- Significantly extended driving range on a single charge
- Enjoy worry-free journeys without range anxiety
- Spend more time driving and less time charging
- Designed for maximum reliability and peace of mind
- Consistent performance even as charge decreases
- Maintenance free and a long service life



Designed for
**improved traction
around corners**
and safe, all-terrain
performance



Safety in mind

- Automatic speed reduction for safe cornering
- Two-step disengaging lever to prevent free-wheeling
- Brake lights and indicators
- Headlights for improved visibility
- Emergency brake button



Built around comfort

- High quality cushioned seat
- Adjustable steering column
- Adjustable seat height
- Flip up armrests and swivel seat makes getting on and off easy



Adjustable steering for greater comfort



Emergency brake for intuitive safety

- Conveniently positioned on the handlebar
- Allows immediate braking in case of danger
- Provides a more intuitive braking experience, especially for users accustomed to the traditional bicycle braking
- Gives a confident and reassuring drive with increased safety



Add to your scooter

Tailor the Orion® Pro to your taste and needs.
Turn over to explore the optional
extras available.



Swivel seat for easy transfer



High quality lighting



Emergency brake button

Add to your scooter

We have a wide range of optional extras available, so you can customise your Orion® Pro scooter. Ask your retailer to discover the full range.



Rear basket



Side & rear view mirrors



Front & rear bumper



Mobile phone holder



Crutch holder



Seatbelt



Cupholder



Interchangeable shroud panels

Technical data

For more information about this product, including the user manual, please visit our website www.invacare.eu.com



Seat width ▼

510 mm



Seat depth ▼

470 mm



Seat to floor height ▼

Min. 440 mm
Max. 465 mm



Backrest height ▼

490 mm



Driving unit width ▼

690 mm



Total length ▼

1300 mm



Empty weight ▼

126 kg



The heaviest part ▼

62 kg



Max. user weight ▼

150 kg



Battery capacity ▼

2 x 12 V / 100 Ah lithium



Turning diameter ▼

2500 mm



Max. climbable obstacle height ▼

100 mm



Ground clearance ▼

100 mm



Max. safe slope ▼

10 °



Motor capacity ▼

550W / max. 1600W



Max. driving distance ▼

70 km



Max. driving speed ▼

18km/h

Colours*



Alum Blue



Black Vinyl

* Scooter colours may vary slightly. Please check with your retailer.

Upholstery

Tel: +41 61 487 70 70
hqeuropa@invacare.com
www.invacare.eu.com

Invacare International GmbH
Neuhofweg 51
4147 Aesch BL
Switzerland

