

# INVACARE® TOP END® PRO™-2 ALL SPORT WHEELCHAIR





The Invacare® Top End® Pro™-2 All Sport Wheelchair is all about sports. If you want to get into the game, this is the chair for you. Just name it: basketball, tennis, softball, table tennis, badmitton, floor hockey or any other sport that is played on a court and the Top End Pro-2 Wheelchair can get you started. The Pro-2 All Sport Wheelchair is offered in three seat widths with lots of adjustments which is perfect for clubs and aspiring athletes. Get in the game.

### **Features**

Aluminum frame: choice of 3 seat widths

Adjustable seat depth 13"-17"

Screw on seat upholstery

Adjustable front and rear seat heights

Adjustable clamp-on tubular recessed footrest

Adjustable center of gravity

Removable single swivel anti-tip bar

Camber: 15°

25" High performance wheels with high pressure

clincher tires

Quick release axles

Non-folding adjustable back: angle and height

Adjustable tension back upholstery

Aluminum side guards

Cushion

Positioning strap

Rear frame protectors

250 lb. weight capacity

## Model #: P2AS













- 2. Adjustable seat depth: 13" 17"
- 3. Removable single swivel anti-tip bar
- 4. Adjustable seat height: front and rear
- 5. 5K Freewheel option



# Invacare Corporation

www.invacare.com

#### HEA

One Invacare Way Elyria, Ohio 44035-4190 (800) 333-6900

#### Canada

570 Matheson Blvd. E., Unit 8 Mississauga, Ontario L4Z 4G4 Canada (800) 668-5324

# Top End

www.topendwheelchair.com Sports and Recreation Products 4501 63rd Circle North Pinellas Park, Florida 33781-5914 (727) 522-8677 (800) 532-8677 Fax: (727) 522-1007

©2014

Invacare Corporation.
All rights reserved.
Trademarks are identified by the symbols ™ and ®.
All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.
Specifications are subject to change without notification.
Form No. 13-378
Rev. 12/14

# **Options**

Basketball carrying bag

Click straps

Metal O-Ring positioning belt

Extra seat cushion

The 5K Freewheel option is great for short road races too, so this chair truly can be used for ALL sports

Note: Some assembly is required.

