



# Invacare® Top End® Pro BB™ Wheelchair



## Product Description:

The Invacare Top End Pro™ BB Wheelchair is designed for aspiring basketball athletes who want durability and quick, easy adjustments, all at a great price. Adjustments that can take up to four hours on other chairs, take only 15 minutes in the Pro BB. Plus, it is so easy to order that you only need to measure seat width and depth, select a wheel size using our wheel guide and pick a color. Now that is easy!

## About New Smart Forms:

Top End is excited to introduce a new way to order Top End products. Our new Smart Forms provide an easy and electronic way to fill out order forms, save copies, and submit to them via e-mail. Please see below for the features and useful information.

### How to use:

- Use any Adobe Acrobat application to electronically fill out the order form
- Use the default selections to order a product with the most common configuration
- Use the **"Reset"** button located in the forms footer to reset selections to their defaults
- Use the **"Clear All"** button located in the forms footer to clear selections, including defaults. This feature allows for a blank order form to be printed and completed by hand
- Use the **"E-mail to Top End order box"** button located on the last page of the form to e-mail the completed form directly to the [topendorders@invacare.com](mailto:topendorders@invacare.com) e-mail inbox
- Forms can be completed electronically, printed, and e-mailed or faxed to [topendorders@invacare.com](mailto:topendorders@invacare.com) or 727-522-1007 if access to Acrobat Reader for Windows is unavailable

*(Acrobat Reader for Windows with java-script enabled is required to save and attach to e-mail)*

### Features:

- Tabulated options for an organized and standardized work-flow
- New more specific option codes which will appear on confirmations to improve customer / manufacturing communication
- Common special requests have been given option codes
- Built in exclusions to help eliminate errors
- Ability to save the form after it is filled out
- Ability to e-mail the form inside Acrobat Reader rather than saving a copy and attaching to an e-mail
- Ability to clear the entire form so it can be printed and completed by hand in the field when a computer is not available
- Ability to reset form to default selections



# Invacare® Top End® Pro BB™ Wheelchair



## Order Request Type:

☐ Quote ☐ CAD/Quote ☐ Order

PO #

## Provider / Customer Information:

Provider Information:		End User Information:	
Date:	Name:	Height:	
IVC Acct #:	Ship to Address:	Weight:	
Contact:	Address Line 2:	Leg Length:	
Phone #:	City:	Hip Width:	
Fax/Email:	State:	Zip:	Phone #:
<input type="radio"/> Ship to Provider		<input type="radio"/> Ship to End User	
		Email:	

**Complete chair starts at: \$2,395.00**

### Base Model:

Item #:	Description:	Price:
<input type="radio"/> 10000	Top End Pro BB Wheelchair	\$2395

### Sizing:

Item #:	Description:	Price:
<b>Seat Width:</b>		
<input type="radio"/> 15005	10" Seat Width	\$0
<input type="radio"/> 15007	11" Seat Width	\$0
<input type="radio"/> 15009	12" Seat Width	\$0
<input type="radio"/> 15011	13" Seat Width	\$0
<input type="radio"/> 15013	14" Seat Width	\$0
<input type="radio"/> 15015	15" Seat Width	\$0
<input type="radio"/> 15017	16" Seat Width	\$0
<input type="radio"/> 15019	17" Seat Width	\$0
<input type="radio"/> 15021	18" Seat Width	\$0
<input type="radio"/> 15023	19" Seat Width	\$0
<input type="radio"/> 15025	20" Seat Width	\$0
<b>Seat Depth:</b>		
<input type="radio"/> 15034	10" Seat Depth	\$0
<input type="radio"/> 15036	11" Seat Depth	\$0
<input type="radio"/> 15038	12" Seat Depth	\$0
<input type="radio"/> 15040	13" Seat Depth	\$0
<input type="radio"/> 15042	14" Seat Depth	\$0
<input type="radio"/> 15044	15" Seat Depth	\$0
<input type="radio"/> 15046	16" Seat Depth	\$0
<input type="radio"/> 15048	17" Seat Depth	\$0
<input type="radio"/> 15050	18" Seat Depth	\$0
<b>Back Height:</b>		
<input type="radio"/> 15126	Adjustable 9" to 11"	\$0
<input type="radio"/> 15127	Adjustable 10" to 14"	\$0
<input type="radio"/> 15128	Adjustable 14" to 18"	\$0
<b>Seat Height: (see product guide)</b>		
<input type="radio"/> 15099	Adjustable Front: 19" to 21" (Tall) Adjustable Rear: 15" to 21"	\$0
<input type="radio"/> 15100	Adjustable Front: 17" to 19" (Short) Adjustable Rear: 13" to 17"	\$0

### Wheelbase: (front to rear axle)

Item #:	Description:	Price:
<input type="radio"/> 75066	15" Wheelbase	\$0
<input type="radio"/> 75068	17" Wheelbase	\$0
<input type="radio"/> 75070	19" Wheelbase	\$0

### Wheel Camber Angle:

Item #:	Description:	Price:
<input type="radio"/> 45008	15 Degree Camber	\$0
<input type="radio"/> 45009	18 Degree Camber	\$0
<input type="radio"/> 45010	20 Degree Camber	\$0

### Center of Gravity: (See Product Diagram 1)

Item #:	Description:	Price:
<input type="radio"/> 45072	Adjustable 3.5" to 5.5"	\$0

### Back Type:

Item #:	Description:	Price:
<input type="radio"/> 15159	Non-Fold-Down, Angle Adjustable back	\$0
<input type="radio"/> 15165	Fold-Down back	\$250

### Side Guards:

Item #:	Description:	Price:
<input type="radio"/> 20004	Fixed, Aluminum, Clamp-on	\$125
<input type="radio"/> 20003	Fold-In, Aluminum - Large	\$125

### Anti-Tip:

Item #:	Description:	Price:
<input type="radio"/> 20030	Single, Adjustable Height	\$0
<input type="radio"/> 20031	Double, Adjustable Height	\$120

### Footrest:

Item #:	Description:	Price:
<input type="radio"/> 60001	Platform, Height Adjustable	\$0
<input type="radio"/> 60030	Click Strap, Single Pad	\$75



# Invacare® Top End® Pro BB™ Wheelchair



Wheel Options: (see exclusions)		
Item #:	Description:	Price:
<b>Size and Style:</b>		
<input type="radio"/> 50001	24" Sun High Performance	\$0
<input type="radio"/> 50005	24" Spinerger SLX - Black	\$895
<input type="radio"/> 50009	24" Spinerger SLX - Yellow	\$895
<input type="radio"/> 50013	24" Spinerger SLX - White	\$895
<input type="radio"/> 50017	24" Spinerger SLX - Red	\$895
<input type="radio"/> 50021	24" Spinerger SLX - Blue	\$895
<input type="radio"/> 50025	24" Spinerger SLX - Green	\$895
<input type="radio"/> 50029	24" Spinerger SLX - Pink	\$895
<input type="radio"/> 50116	24" Sun Equalizer Cross Spoke (see definitions)	\$250
<input type="radio"/> 50002	25" Sun High Performance	\$0
<input type="radio"/> 50006	25" Spinerger SLX - Black	\$895
<input type="radio"/> 50010	25" Spinerger SLX - Yellow	\$895
<input type="radio"/> 50014	25" Spinerger SLX - White	\$895
<input type="radio"/> 50018	25" Spinerger SLX - Red	\$895
<input type="radio"/> 50022	25" Spinerger SLX - Blue	\$895
<input type="radio"/> 50026	25" Spinerger SLX - Green	\$895
<input type="radio"/> 50030	25" Spinerger SLX - Pink	\$895
<input type="radio"/> 50117	25" Sun Equalizer Cross Spoke (see definitions)	\$250
<input type="radio"/> 50003	26" Sun High Performance	\$0
<input type="radio"/> 50007	26" Spinerger SLX - Black	\$895
<input type="radio"/> 50011	26" Spinerger SLX - Yellow	\$895
<input type="radio"/> 50015	26" Spinerger SLX - White	\$895
<input type="radio"/> 50019	26" Spinerger SLX - Red	\$895
<input type="radio"/> 50023	26" Spinerger SLX - Blue	\$895
<input type="radio"/> 50027	26" Spinerger SLX - Green	\$895
<input type="radio"/> 50031	26" Spinerger SLX - Pink	\$895
<input type="radio"/> 50118	26" Sun Equalizer Cross Spoke (see definitions)	\$250
<b>Handrim Size/Material:</b>		
<input type="radio"/> 50300	24" Aluminum Anodized (6 tab)	\$0
<input type="radio"/> 50301	25" Aluminum Anodized (6 tab)	\$0
<input type="radio"/> 50312	25" Aluminum 12 tab	\$150
<input type="radio"/> 50302	26" Aluminum Anodized (6 tab)	\$0
<input type="radio"/> 50313	26" Aluminum 12 tab	\$150
<b>Handrim Style/Mounting:</b>		
<input type="radio"/> 50308	Standard Mount	\$0
<input type="radio"/> 50309	Close Mount	\$0
<input type="radio"/> 50310	Close Mount, Cut-Tabs - Aluminum	\$40
<input type="radio"/> 50311	Close Mount, Cut-Tabs - Titanium	\$80
<b>Tires:</b>		
<input type="radio"/> 50200	24" Primo, Grey	\$0
<input type="radio"/> 50204	24" Primo, Red	\$75
<input type="radio"/> 50208	24" Primo, Yellow	\$75
<input type="radio"/> 50201	25" Primo, Grey	\$0
<input type="radio"/> 50205	25" Primo, Red	\$75
<input type="radio"/> 50209	25" Primo, Yellow	\$75
<input type="radio"/> 50202	26" Primo, Grey	\$0
<input type="radio"/> 50206	26" Primo, Red	\$75
<input type="radio"/> 50210	26" Primo, Yellow	\$75

Axles:		
Item #:	Description:	Price:
<input type="radio"/> 50400	Quick Release Axles	\$0

Upholstery:		
Item #:	Description:	Price:
<b>Material/Color:</b>		
<input type="radio"/> 30001	Nylon, Black	\$0
<input type="radio"/> 30002	Nylon, Woodland Camouflage	\$100
<input type="radio"/> 30003	Nylon, Digital Camouflage	\$100
<b>Seat Style:</b>		
<input type="radio"/> 30019	Screw-On Seat Upholstery	\$0
<b>Back Style:</b>		
<input type="radio"/> 30012	Adjustable Tension Back Upholstery	\$0
<b>Options:</b>		
<input type="radio"/> 30004	Leather Flaps (see definitions)	\$100
<input type="radio"/> 30005	2" Soft Cushion	\$75
<input type="radio"/> 30006	3" Soft Cushion	\$75
<input type="radio"/> 30008	2" Hard Cushion	\$125
<input type="radio"/> 30009	3" Hard Cushion	\$125
<input type="radio"/> 30010	4" Hard Cushion	\$125

Options and Accessories:		
Item #:	Description:	Price:
<b>Safety:</b>		
<input type="radio"/> 25004	Rear Frame Protector (Required for Tournament play)	\$20
<b>Positioning:</b>		
<input type="radio"/> 25066	Calf Strap	\$0
<input type="radio"/> 25001	Click Strap, 8" Small	\$75
<input type="radio"/> 25002	Click Strap, 10" Medium	\$80
<input type="radio"/> 25003	Click Strap, 12" Large	\$86
<input type="radio"/> 25012	Metal O-Ring Strap	\$60
<input type="radio"/> 25013	Hook and Loop Strap	\$40
<b>Accessories:</b>		
<input type="radio"/> 25011	Basketball Carrying Bag	\$75
<b>Apparel:</b>		
<input type="radio"/> 25009	T-Shirt, Large	\$10
<input type="radio"/> 25010	T-Shirt, Extra Large	\$10



# Invacare® Top End® Pro BB™ Wheelchair




## Color Options:

Item #:	Color:	Hue:	Price:
<b>Frame Color:</b>			
<input type="radio"/> 55002	Red	Tangerine Red	\$0
<input type="radio"/> 55003	Red	Raspberry Red	\$0
<input type="radio"/> 55005	Orange	Simple Orange	\$0
<input type="radio"/> 55007	Yellow	Safety Yellow	\$0
<input type="radio"/> 55008	Green	Grasshopper Green	\$0
<input type="radio"/> 55010	Blue	Midnight Blue	\$0
<input type="radio"/> 55011	Blue	Royal Blue	\$0
<input type="radio"/> 55012	Blue	Electric Blueberry	\$0
<input type="radio"/> 55014	Purple	Grape Madness	\$0
<input type="radio"/> 55015	Pink	Bubble-Gum Pink	\$0
<input type="radio"/> 55016	White	White Snow	\$0
<input type="radio"/> 55017	Black	Stealth Black	\$0
<input type="radio"/> 55018	Black	Wet Black	\$0
<input type="radio"/> 55019	Black	Black Prism	\$0
<input type="radio"/> 55020	Silver	Silver Metallic	\$0
<input type="radio"/> 55023	Grey	Stealth Grey	\$0

## Product Illustration:



## Product Diagrams:



## Definitions:

<b>Sun Equalizer Wheels:</b>	All Sun Equalizer wheels have a crossed spoke lacing pattern. 24" have 36 spokes and come with 6 tab aluminum anodized handrims, 25" and 26" have 48 spokes and come with 12 tab aluminum anodized handrims.
<b>Leather flaps:</b>	The leather flap upgrade includes a 3" leather flaps for reinforced heavy duty upholstery and an aluminum back frame plug to help prevent tearing.

## Exclusions:

<b>Wheels/tires:</b>	Must choose the same size tires and handrims as wheels.
----------------------	---

## Product Guides:

### Wheelbase Guide:

Seat Depth:	10" to 14"	15" to 16"	17" to 18"
Recommended Wheelbase:	15"	17"	19"

### Wheel Size Guide:

Rear Seat Height:	Short	Short or Tall	Tall
Player Specifications:	Junior's/Short arms	Adult/3.0 or lower	Adult/3.5 or higher
Recommended Wheel Size:	24"	25"	26"

Note: All Special requests are subject to approval, additional charges, and additional lead time

Special Request 1)	
Special Request 2)	
Special Request 3)	
Special Request 4)	
Special Request 5)	