^{Invacare®} | **Pronto® M41**™





New versatile mobility

The *Invacare Pronto M41* is the compact, manoeuvrable indoor chair that also has excellent outdoor driving capabilities. It combines smooth mid-wheel drive and SureStep® suspension to allow an individual to navigate small spaces and offers a great drive with a top speed of 6 km/h. With a choice of seating options the *Pronto M41* is a great versatile chair to suit a range of clients.

Its sleek, stylish, candy red shrouding ensures the **Pronto M41** looks are as engaging as its performance. The cross compatibility of a number of options and accessories means the **Pronto M41** fits straight into the expansive Invacare range of mobility products.





Pronto M41 with Cab Seat

Features and Options

Optimum Comfort

Pronto M41 has been designed with SureStep technology to ensure excellent stability and manoeuvrability in the narrowest of places. The tension adjustable backrest is a standard feature with the Standard seat option, and is angle adjustable -10° to 30°. The ability to swivel 90° in either direction allows easy access and this combined with height adjustable armrests, a one piece flip-up

footboard and a standard 50 mm cushion ensures maximum comfort for day long use. The Cab Seat Option with padded back and integral headrest also offers the benefit of flip up armrests for easy transfers. Suitable for indoor and outdoor environments, the **Pronto M41** provides a perfect solution for those clients who are looking for a chair that is extremely comfortable and versatile.



Centre wheel drive

Highly manoeuvrable with centre wheel drive





SureStep®

Maximum stability and safety with the SureStep feature.



Cab Seat with Flip Armrests

For ease of transfer the armrest of the Pronto M41 with cabseat flips up and down easily.



Technical data



Standard Seat 440 - 480 mm Cab Seat 410 - 460 mm



410, 460, 510 mm 406 - 457 mm



92 - 115°



940 - 1100 mm



480, 540 mm ¹

597 mm⁶

Pronto M41 275 - 465 mm ²



455, 480, 505, 530, 555 mm ⁴



610 mm



980 - 1030 front riggings



1020 mm







2 x 32 Ah



Pronto M41

50 mm

7° / 12.3% 3

Fixed tilt 5°



Pronto M41

Frame colour

21 km

6 km/h

1/ measured without seat cushion on Standard Seat

2/ certain legrest lengths are restricted to a special legrest type 3/ the max. safe slope depends on the wheelchair configuration.

- For further details please see the user manual 4/ manual settable
- 5/ Fixed tilt position 5° 6/ With headrest

Invacare Limited Pencoed Technology Park - Pencoed -Bridgend CF35 5AQ - United Kingdom Tel: +44 1 656 776 222 Fax: +44 1656 776 220 E-mail: uk@invacare.com Sales Order E-mail: ordersuk@invacare.com

Invacare Ireland Ltd Unit 5 - Seatown Business Campus Seatown Road - SWORDS - County Dublin - Ireland Tel. +353 1 8107084 Fax +353 1 8107085 Email: Ireland@invacare.com www.invacare.ie

www.invacare.co.uk

© 2014 Invacare International Sàrl All rights reserved.

All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior

Pronto M41 - UK - 01/2014 - 1552649



Yes, you can.





Candy red