

Tel: 01 810 7084 Fax: 01 810 7085
Email: ordersireland@invacare.com
Online Spares available at www.invacare.ie

JANUARY 2020

Customer Ref :	Account No :
Contact No :	Order No :
Delivery Address :	Order Qty :
	Order Date :



Pronto M41 Ireland

Max user weight 136kg (21.4 stone)



Manual Positioning

Standard Features



- Compact, manoeuvrable, **MWD Pronto M41 Base**
- Telescopic adjustable Seat Widths from 330mm - 530mm (13" - 21") with additional 50mm (2") by arms
- Seat Depths 410mm - 560mm Stepless Adjustment (16" - 22")
- Tension Adjustable Backrest
- Standard T-Armrests, detachable with Short armpad
- Manual settable seat Tilt (0° - 4° - 6°)
- Fixed seat heights of 450mm - 470mm - 490mm.
- Manual settable backrest Recline (0° - 30° in 3.75° steps)
- Standard 80° hangers with fixed footplate
- Shark2 REM D01 control system with standard bracket
- Candy Red Shroud Colour
- 255mm x 80mm (10" x 3") Solid Insert / Puncture Proof Grey Drive Wheels
- 6 km/h (4mph) In-line Motors and 35Ah Batteries
- 8A (200 - 250V) Charger
- 4-Point tie-down system
- Lapbelt with autolock (adjustable both sides)



* Configurations may vary from the image shown

Standard Base Price : €3,638

	X	Please make your selection with a cross in the box	◆ Standard feature	○ Free of Charge
Customer weight =		Kg		
Code		Description	Ireland	
AQE0059		Pronto M41 "Modulite fixed seating"	◆	
AQX0059		Pronto M41 "Modulite" Special Basic Model-PLEASE COMPLETE DETAILS AT THE END OF THE FORM	POA	
SEATING TYPE : MODULITE MULTI ADJUSTABLE SEATING SYSTEM				
AQE2100	X	Modulite Seating Unit	€3,638	
SEAT FRAME AND SEAT PLATE				
AQE2111	X	Telescopic Seat Frame with metal seat plate (with stepless adjustment on armrests, legrests and backrest)	◆	

MODULITE SEAT SIZE SETTING CHART:

Invacare® **Modulite™**
An evolution in seating designed around you.



			SEAT WIDTHS				
			330mm* (13")	AQE2121 380mm (15")	AQE2122 430mm (17")	AQE2123 480mm (19")	AQE2124 530mm (21")
Seat Depth AQE2129	Standard Seat Rail 410 - 510mm	AQE2132 (16") 410mm	○	○	○	○	○
		AQE2134 (18") 460mm	○	○	○	○	○
		AQE2136 (20") 510mm	○	○	○	○	○
Seat Depth AQE2130	Extended Seat Rail 510 - 560mm	AQE2136 (20") 510mm	○	○	○	○	○
		AQE2138 (22") 560mm	○	○	○	○	○

*330mm (13") seat width achieved on a 380mm seat frame (2121) with inrigged armrests (2183)

1. Seat depth 410mm: cannot combine T-Armrests (2191) with flexible lapbelt mount (2113)
2. Seat cushion only available with 560mm seat depth (2138) with seat width 530mm (2124)
3. Seat width 530mm (2124) not possible in combination with following armrests (2192/4)
4. Seat depth 560mm (2138) When in combination with Flex Tension Adjustable Backrest the seat depth is reduced
5. Modulite Teen seating available, seat width and depth from 305mm (12")

SEAT OPTION

AQESWO	"H" Frame Interface Kit	€460
--------	-------------------------	------

SEAT CUSHION (sizes according to seat frame settings)

AQE2142	Matrix Flo-tech Lite	€120
AQE2143	Matrix Flo-tech Contour Visco	Not available with seat width 530mm (21") or black spacer fabric €120
AQE2147	Matrix Libra, 90mm thick	Not available with seat width 530mm with 410mm seat depth, seat depth 560mm only available with 530mm seat width or black spacer fabric €447
AQE2140	Without seat cushion	○

BACKREST SUPPORT

AQE2162	Backposts 540mm (Tall)	Must be selected for flex3 adjustable backrest ◆
AQE2163	Backposts 480mm (Short)	Not available with flex3 adjustable backrest ○

PUSH HANDLES FOR BACKPOSTS


AQE2107	Backposts WITH push handles	◆
AQE2108	Backposts WITHOUT push handles, Push handles will be cut down on the vertical section	Not available with attendant controls, dual controls, tension adjustable backrest 2171 in conjunction with a headrest, headrest 1495 €46

BACKREST UPHOLSTERY / CUSHION

AQE2171	Basic Tension Adjustable, (540mm or 480mm upholstery height) Straps easily accessible and adjustable from the rear, Black fabric only	Not available with laterals, headrest option 1496 ◆
AQE2179	Flex3 Tension Adjustable (540mm - 480mm) Black Dartex only	Not available with flip up or following armrests, chin remote and soft armpads €173
AQE2172	Flex3 adjustable Backrest Plates (without Backrest Cushion)	€257
AQE2173	Flex3 adjustable Backrest Plates with Laguna Backrest cushion	€417
AQE2174	Matrix Posture Back, height 400mm (Mesh fabric)	Not available with laterals, 530mm seat width, headrest option 1496 €519
AQE2175	Matrix Posture Back, height 500mm (Mesh fabric)	Not available with laterals, 530mm seat width, headrest option 1496 €519
AQE2176	Matrix Elite Back, height 500mm (Mesh fabric)	Not available with laterals, 530mm seat width, headrest option 1496, seat width of 380mm 2121 in combination with flip up, following armrests €519
AQE2177	Matrix Elite TR Back, height 500mm (Mesh fabric)	Not available with laterals, 530mm seat width, headrest option 1496 €519
AQE2178	Matrix Elite Deep Back, height 500mm (Mesh fabric)	Not available with laterals or headrest option 1496 €672
AQE2170	Without Backrest Upholstery	Not available with laterals or headrest option 1496 ○

UPHOLSTERY TYPE

AQE2151	Black Nylon	◆
AQE2152	Black Dartex	○
AQE2155	Black Spacer fabric	○




ARMREST (Pair)			
AQE2191	T-Armrest Standard Range of adjustment from seat plate 245mm to 315mm		◆
AQE2181	T-Armrest Extended height Range of adjustment from seat plate 295mm to 365mm		€46
AQE2183	In-rigged Armrests (reduces seat width by 25mm each side)	Not available with hip support	€27
AQE2193	Flip-up Armrest (height and depth adjustable) (height adjustable from 230mm - 360mm)	Not available with laterals, combination of 530mm seat plate and flex3 adjustable backrest plates.	€39
AQE2192	Following Armrest (height adjustable from 230mm - 300mm)	Not available with laterals, combination of 530mm seat plate and flex3 adjustable backrest plates.	€190
AQE2194	Following Armrest (height adjustable from 300mm - 360mm)	Not available with laterals, combination of 530mm seat plate and flex3 adjustable backrest plates.	€190

ARMPAD (Pair)			
AQE2195	Short Armpad (325mm)		◆
AQE2196	Long Armpad (375mm)	Not available with seat depth 410mm (16")	○
AQE2197	Long Ergonomic Armpad (397mm)	Not available with seat depth 410mm (16")	○
AQE2199	Soft Armpad (REA Azalea Wide Pad)	Not available with seat depth 410mm (16")	€98














HIP SUPPORTS, Multi adjustable LH (for flip up and following armrests only)			
AQE2184	SMALL 140mm x 100mm	Not available with 410mm seat depth and manual centre legrests	€259
AQE2186	MEDIUM 230mm x 100mm	Not available with 410mm seat depth and manual centre legrests	€259
AQE2188	LARGE 330mm x 100mm	Not available with 410mm seat depth and manual centre legrests	€259


HIP SUPPORTS, Multi adjustable RH (for flip up and following armrests only)			
AQE2185	SMALL 140mm x 100mm	Not available with 410mm seat depth and manual centre legrests	€259
AQE2187	MEDIUM 230mm x 100mm	Not available with 410mm seat depth and manual centre legrests	€259
AQE2189	LARGE 330mm x 100mm	Not available with 410mm seat depth and manual centre legrests	€259

EXTENDED HEIGHT (+2") MOUNTING HARDWARE FOR HIP SUPPORTS			
AQE2165	Extended height mounting hardware for hip supports (8.1" / 206mm high)		○



HEADREST / NECKREST (Dartex covered) - Adequate headrest option must be selected			
AQE1350	Neckrest: height and depth adjustable		€81
AQE1432	Headrest: height, depth and sideways adjustable		€226
AQE2323	Elan Standard Pad 150mm (6") wide with standard cover complete with Elan Headrest Hardware	Not available with head controls	€369
AQE2324	Elan Standard Pad 255mm (10") wide with standard cover complete with Elan Headrest Hardware	Not available with head controls	€369
AQE2322	Elan Standard Pad 355mm (14") wide with standard cover complete with Elan Headrest Hardware	Not available with head controls	€369
AQE2325	Elan Occipital Pad 230mm x 180mm (9" x 7") with standard cover complete with Elan Headrest Hardware	Not available with head controls	€369
AQE2326	Elan Occipital Pad 305mm x 200mm (12" x 8") with standard cover complete with Elan Headrest Hardware	Not available with head controls	€369
AQE2327	Elan 4 Point Headrest Pad 280mm x 255mm (11" x 10") with standard cover, complete with Elan Headrest Hardware	Not available with head controls	€499
AQE1427	Headrest with cheek support: height and depth adjustable		€226
AQE1496	All bracketry as option 1432 without headrest cushion		€84
AQE1495	Push bar dia. 22mm in conjunction with telescopic seat frame / backposts, No Headrest - no mounting hardware.		○

		Rigid pad 1446			Multi-adjustable pad 1445
LATERAL SUPPORT (Pair)					
AQE1446		Rigid pad (Azalea type, Dartex)	Not available with flex3 tension backrest		€168
AQE1445		Multi-adjustable pad (Azalea type, Dartex)			€168
SEAT ADJUSTMENTS					
AQE1091	X	Manually settable Seat Tilt			◆
SEAT TO FLOOR HEIGHT (only to choose with option 1091)					
AQE0354		Seat height: 450mm			○
AQE0364		Seat height: 470mm			○
AQE0382		Seat height: 490mm			○
BACKREST RECLINER					
AQE0975	X	Manual Settable Recliner (0° - 30° in 3.75° steps) Foldable			◆
PLEASE CHOOSE EITHER CENTRE MOUNTED OR SIDE MOUNTED LEGRESTS					
MANUAL CENTRE MOUNT 0556		Standard 80° 0526	VARI F 0503	VARI A 0533	ADM 0532
					
CENTRE MOUNTED LEGRESTS (INCLUDING CALF PAD AND ANGLE & WIDTH ADJUSTABLE FOOTPLATE 240mm depth)					
AQE0556		Manual : Hanger, Centre mounted, angle settable (0° - 90°): knee to heel length 280mm - 385mm			€327
SIDE MOUNTED LEGRESTS					
AQE0526		Standard 80° hangers	Not available with angle adjustable footplate		◆
AQE0503		VARI F Angle adjustable with tools 0° - 70°, length 290mm - 460mm			○
AQE0533		VARI A Angle adjustable with hand wheel 0° - 70°, length 290mm - 460mm		Not available with fixed footplate	€127
AQE0532		ADM Manual Angle adjustable with quick release lever 0° - 80°, length 290mm - 460mm		Not available with fixed footplate	€221
AQE0501		Without Hanger			○
FOOTPLATE TYPE (Pair)					
AQE0512		Fixed Footplate			◆
AQE0504		Angle Adjustable Footplate			○
AQE0513		Angle and depth adjustable Footplate (Adapt)			€116
SHARK CONTROLS					
REMOTE & BRACKET					
AQE1801		Standard Remote Shark 2 REM D01			◆
AQE1207		Standard Bracket RH			○
AQE1208		Standard Bracket LH			○
AQE1707		Swing Away + height adjustable RH high end			○
AQE1708		Swing Away + height adjustable LH high end			○
REMOTE DX2 SYSTEM					
REM Type A	REM 24SD	G90A	Chin control	Dual control	IDC Dual Control
					
REMOTE & BRACKET					
AQE1638		REM Type A, (enables operation up to 4 electric adjustments)			○
AQE1601		ACS REM 24SD			€272
AQE1614		G90A			€201
AQE1207		Standard Bracket RH			○
AQE1208		Standard Bracket LH			○
AQE1707		Swing Away + height adjustable RH High End			€144
AQE1708		Swing Away + height adjustable LH High End			€144

ATTENDANT REMOTE & BRACKET			
AQE1643		REM Type A (Up to 4 electric actuators)	€18
AQE1602		ACS REM 24SD	○
AQE1616		G90A	○
AQE1207		Standard Bracket RH	○
AQE1208		Standard Bracket LH	○
CHIN CONTROL & BRACKET			
AQE1719		Chin Remote with standard Joystick (Always to be mounted on the Master Remote's side)	(Only compatible with option 1601) €2,293
AQE1289		Swing Away Bracket RH, side Manual	Only available with RH user bracket ○
AQE1288		Swing Away Bracket LH, side Manual	Only available with LH user bracket ○
AQE1287		Swing Away Bracket RH, Powered	Only available with RH user bracket €1,546
AQE1286		Swing Away Bracket LH, Powered	Only available with LH user bracket €1,546
DUAL CONTROL & BRACKET (can only be selected as an addition to user remote)			
AQE1720		Compact Dual Control	€495
AQE1207		Standard Bracket RH	Only available with LH user bracket ○
AQE1208		Standard Bracket LH	Only available with RH user bracket ○
INVACARE IDC DUAL CONTROL & BRACKET (Only available with ACS)			
AQE1645		Invacare IDC Dual control	€681
AQE1221		IDC Dual control Detachable Bracket RH	○
AQE1222		IDC Dual control Detachable Bracket LH	○
DRIVE UNIT OPTIONS			
FRAME COLOUR (COLOURED PARTS IS A DETACHABLE AND CHANGEABLE PLASTIC SHROUD)			
AQE0227	X	Candy Red	◆
CENTRE WHEELS (DRIVING WHEELS)			
AQE0781	X	Centre Wheel 255mm x 80mm (10"-3") Solid Insert / Puncture Proof Grey	◆
PERFORMANCE OPTIONS			
MOTOR TYPE			
AQE1282	X	6km/h (4mph)	◆
BATTERY			
AQE1261	X	35Ah AGM pair	◆
CHARGER			
AQE1343		OptiCharge Charger 8A (200 - 250V)	◆
AQE1313		OptiCharge Charger 10A (200 - 250V) Charging time will be 20% faster	€158
SAFETY ACCESSORIES			
TRANSPORTATION KIT			
AQE1494	X	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory)	◆
FIXATION BELT			
AQE1329		Lapbelt (adjustable both sides) fitted to base of backrest	○
AQE2113		Lapbelt (adjustable both sides) adjustable mounting on seat rail extrusion	Not available with seat depth 410mm and T-armrest ◆
BODYPOINT PELVIC BELTS (Includes mounting kit)			
AQE2701		Evoflex Pelvic Stabilizer, Rehab Latch size L (padded length 620mm)	Not available with seat width 380mm, 430mm, seat depth 410mm, hip supports €106
AQE2702		Evoflex Pelvic Stabilizer, Rehab Latch size M (padded length 460mm)	Not available seat width 480mm, 530mm, seat depth 410mm, hip supports €106
AQE2703		2-Point Hip Belt, Centre pull - push button, Cinch mount size L (padded length 620mm)	Not available with seat width 380mm, 430mm, seat depth 410mm, hip supports €98
AQE2704		2-Point Hip Belt, Centre pull - push button, Cinch mount size M (padded length 460mm)	Not available seat width 480mm, 530mm, seat depth 410mm, hip supports €98
AQE2705		2-Point Hip Belt, Dual pull, Side release button, Cinch mount size L (padded length 620mm)	Not available with seat width 380mm, 430mm, seat depth 410mm, hip supports €98
AQE2706		2-Point Hip Belt, Dual pull, Side release button, Cinch mount size M (padded length 460mm)	Not available seat width 480mm, 530mm, seat depth 410mm, hip supports €98
ACCESSORIES			
BAG HOOKS			
			
AQE2128		Bag hooks (Max load up to 6 kg = impact on max user weight)	Not available with flex3 backrest 2172/2173, attendant brackets, dual control €86
CANE HOLDER (in combination with attendant or dual control the cane holder has to be mounted on the opposite side of the remote)			
AQE1312		Cane Holder RH	€58
AQE1363		Cane Holder LH	€58
TRAY (PLEASE CHOOSE OPPOSITE HAND TO USER CONTROL)			
AQE1404		Tray Swing Away RH	€201
AQE1405		Tray Swing Away LH	€201
TOTAL PRICE			



Powerchair Alternative Options

Need a different option not on the form?

Take a look at our Powerchair Alternative Options form, available as standard

Powerchair Alternative Option required:



Special Instructions - Please supply the following information:

Name:

Email:

Phone:

Fax:

Description of special request (please include relevant option code or part number):