

Tel: 01 810 7084 Fax: 01 810 7085
Email: ordersireland@invacare.com
Online Spares available at www.invacare.ie

JANUARY 2020

Customer Ref :	Account No :
Contact No :	Order No :
Delivery Address :	Order Qty :
	Order Date :



Pronto M41 Ireland

Max user weight 136kg (21.4 stone)



Powered Positioning

Standard Features



- Compact, manoeuvrable, **MWD Pronto M41 Base**
- Telescopic adjustable Seat Widths from 330mm - 530mm (13" - 21") with additional 50mm (2") by arms
- Seat Depths 410mm - 560mm Stepless Adjustment (16" - 22")
- Tension Adjustable Backrest
- Standard T-Armrests, detachable with Short armpad
- Powered seat Tilt (0° - 30°)
- Manual settable backrest Recline (0° - 30°)
- Standard 80° hangers with fixed footplate
- REM A Joystick with standard bracket
- Candy Red Shroud Colour
- 255mm x 80mm (10" x 3") Solid Insert/Puncture Proof Grey Drive Wheels
- 6 km/h (4mph) In-line Motors and 35Ah Batteries
- 8A (200 - 250V) Charger
- 4-Point tie-down system
- Lapbelt with autolock (adjustable both sides)



* Configurations may vary from the image shown

Standard Base Price : €4,568

	X	Please make your selection with a cross in the box	◆ Standard feature	○ Free of Charge
Customer weight = Kg				
Code		Description		Ireland
AQD0059		Pronto M41 "Modulite"	◆	
AQX0059		Pronto M41 "Modulite" Special Basic Model-PLEASE COMPLETE DETAILS AT THE END OF THE FORM		POA
SEATING TYPE : MODULITE MULTI ADJUSTABLE SEATING SYSTEM				
AQD2100	X	Modulite Seating Unit		€4,568
SEAT FRAME AND SEAT PLATE				
AQD2111	X	Telescopic Seat Frame with metal seat plate (with stepless adjustment on armrests, legrests and backrest)	◆	

MODULITE SEAT SIZE SETTING CHART:

			SEAT WIDTHS				
			330mm* (13")	AQD2121 380mm (15")	AQD2122 430mm (17")	AQD2123 480mm (19")	AQD2124 530mm (21")
Seat Depth AQD2129	Standard Seat Rail 410 - 510mm	AQD2132 (16") 410mm	○	○	○	○	○
		AQD2134 (18") 460mm	○	○	○	○	○
		AQD2136 (20") 510mm	○	○	○	○	○
Seat Depth AQD2130	Extended Seat Rail 510 - 560mm	AQD2136 (20") 510mm	○	○	○	○	○
		AQD2138 (22") 560mm	○	○	○	○	○

Invacare **Modulite**[™]
An evolution in seating designed around you.



*330mm (13") seat width achieved on a 380mm seat frame (2121) with inrigged armrests (2183)

1. Seat depth 410mm: cannot combine T-Armrests (2191) with flexible lapbelt mount (2113)
2. No seat cushion available with 560mm seat depth (2138)
3. Seat width 530mm (2124) not possible in combination with following armrests (2192/4)
4. Seat depth 560mm (2138) When in combination with Flex Tension Adjustable Backrest the seat depth is reduced by 10mm
5. Modulite Teen seating available, seat width and depth from 305mm (12") - See New Teen Form

SEAT OPTION

AQDSWO	"H" Frame Interface Kit	€460
--------	-------------------------	------

SEAT CUSHION (sizes according to seat frame settings)

AQD2142	Matrx Flo-tech Lite	€116
AQD2143	Matrx Flo-tech Contour Visco	Not available with black spacer fabric €158
AQD2147	Matrx Libra, 90mm thick	Not available with seat width 530mm with 410mm seat depth, seat depth 560mm only available with 530mm seat width, black spacer fabric €483
AQD2140	Without seat cushion	○

BACKREST SUPPORT

AQD2162	Backposts 540mm (Tall)	Must be selected for flex3 adjustable backrest ◆
AQD2163	Backposts 480mm (Short)	Not available with flex3 adjustable backrest ○

PUSH HANDLES FOR BACKPOSTS


AQD2107	Backposts WITH push handles	◆
AQD2108	Backposts WITHOUT push handles Push handles will be cut down on the vertical section.	Not available with attendant controls, dual controls, tension adjustable backrest 2171 in conjunction with a headrest, headrest 1495 €46

BACKREST UPHOLSTERY / CUSHION

AQD2171	Basic Tension Adjustable, (540mm or 480mm upholstery height) Straps easily accessible and adjustable from the rear, Black fabric only	Not available with laterals, headrest option 1496 ◆
AQD2179	Flex3 Tension Adjustable (540mm - 480mm) Black Dartex only	Not available with flip up, following armrests, chin remote, soft armpads €173
AQD2172	Flex3 adjustable Backrest Plates (without Backrest Cushion)	€257
AQD2173	Flex3 adjustable Backrest Plates with Laguna Backrest cushion	€417
AQD2174	Matrx Posture Back, height 400mm (Mesh fabric)	Not available with laterals, 530mm seat width, headrest option 1496 €570
AQD2175	Matrx Posture Back, height 500mm (Mesh fabric)	Not available with laterals, 530mm seat width, headrest option 1496 €570
AQD2176	Matrx Elite Back, height 500mm (Mesh fabric)	Not available with laterals, 530mm seat width, headrest option 1496, seat width of 380mm 2121 in combination with flip up, following armrests €570
AQD2177	Matrx Elite TR Back, height 500mm (Mesh fabric)	Not available with laterals, 530mm seat width, headrest option 1496 €570
AQD2178	Matrx Elite Deep Back, height 500mm (Mesh fabric)	Not available with laterals, headrest option 1496 €632
AQD2170	Without Backrest Upholstery	Not available with laterals, headrest option 1496 ○

UPHOLSTERY TYPE

AQD2151	Black fabric	◆
AQD2152	Black Dartex	○
AQD2155	Black Spacer fabric	○



ARMREST (Pair)			
AQD2191		T-Armrest Standard Range of adjustment from seat plate 245mm to 315mm	◆
AQD2181		T-Armrest Extended height Range of adjustment from seat plate 295mm to 365mm	€46
AQD2183		In-rigged Armrests (reduces seat width by 25mm each side)	Not available with hip support €27
AQD2193		Flip-up Armrest (height and depth adjustable) (height adjustable from 230mm - 360mm)	Not available with laterals or in combination with 530mm seat plate and flex3 adjustable backrest plates. €66
AQD2192		Following Armrest (height adjustable from 230mm - 300mm)	Not available with laterals or in combination with 530mm seat plate and flex3 adjustable backrest plates. €190
AQD2194		Following Armrest (height adjustable from 300mm - 360mm)	Not available with laterals or in combination with 530mm seat plate and flex3 adjustable backrest plates. €190


ARMPADS (Pair)			
AQD2195		Short Armpad (325mm)	◆
AQD2196		Long Armpad (375mm)	○
AQD2197		Long Ergonomic Armpad (397mm)	○
AQD2199		Soft Armpad (REA Azalea Wide Pad)	Not available with seat depth 410mm (16") €98

HIP SUPPORTS, Multi adjustable LH (for flip up and following armrests only)			
AQD2184		SMALL 140mm x 100mm	Not available with 410mm seat depth and manual centre legrests €259
AQD2186		MEDIUM 230mm x 100mm	Not available with 410mm seat depth and manual centre legrests €259
AQD2188		LARGE 330mm x 100mm	Not available with 410mm seat depth and manual centre legrests €259


HIP SUPPORTS, Multi adjustable RH (for flip up and following armrests only)			
AQD2185		SMALL 140mm x 100mm	Not available with 410mm seat depth and manual centre legrests €259
AQD2187		MEDIUM 230mm x 100mm	Not available with 410mm seat depth and manual centre legrests €259
AQD2189		LARGE 330mm x 100mm	Not available with 410mm seat depth and manual centre legrests €259

EXTENDED HEIGHT (+2") MOUNTING HARDWARE FOR HIP SUPPORTS			
AQD2165		Extended height mounting hardware for hip supports (8.1" / 206mm high)	○


SEAT ACCESSORIES



HEADREST / NECKREST (Dartex covered) - Adequate headrest option must be selected			
AQD1350		Neckrest: height and depth adjustable	€81
AQD1432		Headrest: height, depth and sideways adjustable	€226
AQD1427		Headrest with cheek support: height and depth adjustable	€226
AQD2323		Elan Standard Pad 150mm (6") wide with standard cover complete with Elan Headrest Hardware	Not available with head controls €369
AQD2324		Elan Standard Pad 255mm (10") wide with standard cover complete with Elan Headrest Hardware	Not available with head controls €369
AQD2322		Elan Standard Pad 355mm (14") wide with standard cover complete with Elan Headrest Hardware	Not available with head controls €369
AQD2325		Elan Occipital Pad 230mm x 180mm (9" x 7") with standard cover complete with Elan Headrest Hardware	Not available with head controls €369
AQD2326		Elan Occipital Pad 305mm x 200mm (12" x 8") with standard cover complete with Elan Headrest Hardware	Not available with head controls €369
AQD2327		Elan 4 Point Headrest Pad 280mm x 255mm (11" x 10") with standard cover, complete with Elan Headrest Hardware	Not available with head controls €497
AQD1496		All bracketry as option 1432 without headrest cushion	€84
AQD1495		Push bar dia. 22mm in conjunction with telescopic seat frame / backposts, No Headrest - no mounting hardware.	○



**Rigid pad
1446**



**Multi-adjustable pad
1445**

LATERAL SUPPORT (Pair)

AQD1446	Rigid pad (Azalea type, Dartex)	€168
AQD1445	Lateral Support with extended arm, multi-adjustable pad	€168

SEAT LIFTER / SEAT TILT

		Tilt Range	SEAT HEIGHT 420mm (16")	Selection
AQD1093	Electric Tilt with COG Shift	0° - 30°	◆	
AQD1094	Seat Lifter with electric tilt (COG Shift), Not possible with Rem A Control	0° - 30°	€2,048	



BACKREST RECLINER

		Recline Range	Selection
AQD0975	Manually settable recline	0° - 30°	◆
AQD0977	Powered Recliner with quick release for easy folding of backrest	0° - 30°	€318



PLEASE CHOOSE EITHER CENTRE MOUNTED OR SIDE MOUNTED LEGRESTS

MANUAL CENTRE MOUNT 0556	POWERED CENTRE MOUNT 0589	80° standard 0526	VARI F 0503	VARI A 0533	ADM 0532	ADE 0531

CENTRE MOUNTED LEGRESTS (INCLUDING CALF PAD AND ANGLE ADJUSTABLE FOOTPLATE 240mm depth)

AQD0556	Manual : Hanger, centre mounted, angle settable (0° - 90°): knee to heel length 280mm - 385mm	€306
AQD0589	LNX Powered : Centre mounted legrest, length compensating	€1,548







SIDE MOUNTED LEGRESTS

FOOTREST HANGER TYPE, Swing away and detachable (Pair)

AQD0526	Standard 80° hangers	Not available with angle adjustable footplate	◆
AQD0503	VARI F Angle adjustable with tools 0° - 70°, length 290mm - 460mm		○
AQD0533	VARI A Angle adjustable with hand wheel 0° - 70°, length 290mm - 460mm	Not available with fixed footplate	€127
AQD0532	ADM Manual Angle adjustable with quick release lever 0° - 80°, length 290mm - 460mm	Not available with fixed footplate	€221
AQD0531	ADE Powered Angle adjustable powered by joystick 0° - 80°, length 290mm - 460mm	Not available with fixed footplate	€653
AQD0501	Without Hanger		○

FOOTPLATE TYPE (Pair)

AQD0512	Fixed Footplate	◆
AQD0504	Angle Adjustable Footplate	○
AQD0513	Angle and depth adjustable Footplate (Adapt)	€116

REMOTE DX2 SYSTEM					
REM Type A	REM 24SD	G90A	Chin control	Dual control	IDC Dual Control
					
REMOTE & BRACKET					
AQD1638		REM Type A (enables operation up to 4 electric adjustments)			◆
AQD1601		ACS REM 24SD			€272
AQD1614		G90A			€201
AQD1207		Standard Bracket RH			○
AQD1208		Standard Bracket LH			○
AQD1707		Swing Away + height adjustable RH High End			€144
AQD1708		Swing Away + height adjustable LH High End			€144
ATTENDANT REMOTE & BRACKET					
AQD1643		REM Type A (Up to 4 electric actuators)			€18
AQD1602		ACS REM 24SD			€215
AQD1616		G90A			€215
AQD1207		Standard Bracket RH			○
AQD1208		Standard Bracket LH			○
CHIN CONTROL & BRACKET					
AQD1719		Chin Remote with standard Joystick (Always to be mounted on the Master Remote's side)	(Only compatible with option 1601)		€2,293
AQB1289		Swing Away Bracket RH, side Manual	Only available with RH user bracket		○
AQD1288		Swing Away Bracket LH, side Manual	Only available with LH user bracket		○
AQD1287		Swing Away Bracket RH, Powered	Only available with RH user bracket		€1,546
AQD1286		Swing Away Bracket LH, Powered	Only available with LH user bracket		€1,546
DUAL CONTROL & BRACKET (can only be selected as an addition to user remote)					
AQD1720		Compact Dual Control			€495
AQD1207		Standard Bracket RH	Only available with LH user bracket		○
AQD1208		Standard Bracket LH	Only available with RH user bracket		○
INVACARE IDC DUAL CONTROL & BRACKET					
AQD1645		Invacare IDC Dual control			€681
AQD1221		IDC Dual control Detachable Bracket RH			○
AQD1222		IDC Dual control Detachable Bracket LH			○
DRIVE UNIT OPTIONS					
FRAME COLOUR (COLOURED PARTS IS A DETACHABLE AND CHANGEABLE PLASTIC SHROUD)					
AQD0227	X	Candy Red			◆
CENTRE WHEELS (DRIVING WHEELS)					
AQD0781	X	Centre Wheel 255mm x 80mm (10"- 3") Solid Insert / Puncture Proof Grey			◆
PERFORMANCE OPTIONS					
MOTOR TYPE					
AQD1282	X	6 km/h (4mph)			◆
BATTERY					
AQD1261	X	35Ah AGM pair			◆
CHARGER					
AQD1343		OptiCharge Charger 8A (200 - 250V)			◆
AQD1313		OptiCharge Charger 10A (200 - 250V) Charging time will be 20% faster			€158
SAFETY ACCESSORIES					
TRANSPORTATION KIT					
AQD1494	X	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory)			◆
FIXATION BELT					
AQD1329		Lapbelt (adjustable both sides) fitted to base of backrest			○
AQD2113		Lapbelt (adjustable both sides) adjustable mounting on seat rail extrusion	Not available with seat depth 410mm (16") and T-armrest		◆

BODYPOINT PELVIC BELTS (Includes mounting kit)

AQD2701	Evoflex Pelvic Stabilizer, Rehab Latch size L (padded length 620mm)	Not available with seat width 380mm, 430mm, seat depth 410mm, hip supports	€106
AQD2702	Evoflex Pelvic Stabilizer, Rehab Latch size M (padded length 460mm)	Not available seat width 480mm, 530mm, seat depth 410mm, hip supports	€106
AQD2703	2-Point Hip Belt, Centre pull - push button, Cinch mount size L (padded length 620mm)	Not available with seat width 380mm, 430mm, seat depth 410mm, hip supports	€98
AQD2704	2-Point Hip Belt, Centre pull - push button, Cinch mount size M (padded length 460mm)	Not available seat width 480mm, 530mm, seat depth 410mm, hip supports	€98
AQD2705	2-Point Hip Belt, Dual pull, Side release button, Cinch mount size L (padded length 620mm)	Not available with seat width 380mm, 430mm, seat depth 410mm, hip supports	€98
AQD2706	2-Point Hip Belt, Dual pull, Side release button, Cinch mount size M (padded length 460mm)	Not available seat width 480mm, 530mm, seat depth 410mm, hip supports	€98

ACCESSORIES

BAG HOOKS



AQD2128	Bag hooks (Max load up to 6 kg = impact on max user weight)	Not available with flex3 backrest 2172/2173, attendant brackets, dual control	€86
---------	---	---	-----

CANE HOLDER

AQD1312	Cane Holder RH		€39
AQD1363	Cane Holder LH		€37

TRAY (PLEASE CHOOSE OPPOSITE HAND TO USER CONTROL)

AQD1404	Tray Swing Away RH		€350
AQD1405	Tray Swing Away LH		€350

TOTAL PRICE



Need a different option not on the form?

Take a look at our Powerchair Alternative Options form, available as standard

Powerchair Alternative Option required:



Special Instructions - Please supply the following information:

Name:

Email:

Phone:

Fax:

Description of special request (please include relevant option code or part number):