



Tel: (00353) 1 8107084 Fax: (00353) 1 8107085  
Email: ordersireland@invacare.com  
Online Spares available at www.invacare.ie

JANUARY 2020

Customer Ref :	Account No :
Contact No :	Order No :
Delivery Address :	Order Qty :
	Order Date :



# Pronto M41 Modulite Teen

Max user weight 136kg (21.4 stone)

Invacare® **Modulite™**  
An evolution in seating designed around you.

## Ireland

Standard Features



## Powered Positioning

- Compact, manoeuvrable, **MWD Pronto M41 Base**
- Telescopic adjustable Seat Widths from 305mm - 380mm (12" - 15")
- Seat Depths 305mm - 560mm Stepless Adjustment (12" - 22")
- Choice of backrest (Tension adjustable, Flex 2 Laguna, Flex3 Mistral and Matrix)
- Detachable T-Armrests, height, width and depth adjustable with standard armpad
- Powered seat tilt (0° - 30°) and manually settable backrest recline (0° - 30°)
- Swing away and detachable legrest hangers with fixed footplate
- REM A Joystick with standard bracket
- Candy Red Shroud Colour
- 255mm x 80mm (10" x 3") Solid Insert/Puncture Proof Grey Drive Wheels
- 6km/h (4mph) In-line Motors and 32Ah Batteries
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)
- 8A (200 - 250V) Charger



**Standard Base Price : €3,836**

\* Configurations may vary from the image shown

	X	Please make your selection with a cross in the box	◆ Standard feature	○ Free of Charge
		Customer weight =	Kg	
Code		Description		Retail
AQX0059		Pronto M41 Modulite Teen	◆	
AQX0059		Pronto M41 Modulite Teen Special Basic Model-PLEASE COMPLETE DETAILS AT THE END OF THE FORM		POA
SWO0025	X	Invacare ONE to ONE fitting service	◆	
<b>SEATING TYPE : MODULITE MULTI ADJUSTABLE SEATING SYSTEM</b>				
AQX2100	X	Modulite Seating Unit		€3,836
<b>SEAT FRAME AND SEAT PLATE</b>				
AQD2111	X	Telescopic Seat Frame with metal seat plate (with stepless adjustment on armrests, legrests and backrest)	◆	

**MODULITE SEAT SIZE SETTING CHART:**

			SEAT WIDTHS		
			**305mm (12")	*330mm (13")	380mm (15")
Seat Depth AQD2129	Standard Seat Rail 410mm-510mm	(12") 305mm	○	○	○
		(13") 330mm	○	○	○
		(14") 355mm	○	○	○
		(15") 380mm	○	○	○
		(16") 410mm	○	○	○
		(18") 460mm	○	○	○
		(20") 510mm	○	○	○
		Seat Depth AQD2130	Extended Seat Rail 510mm-560mm	(20") 510mm	○
(22") 560mm	○			○	○



\*\*305mm (12") seat width achieved on a 380mm seat frame (2121) with inrigged armrest (2183) and additional inrigged bracket  
\*330mm (13") seat width achieved on a 380mm seat frame (2121) with inrigged armrests (2183)

**SEAT CUSHION (sizes according to seat frame settings)**

AQD2142	Matrix Flo-tech Lite	Only Available in a 15" seat width	€116
AQD2147	Matrix Libra, 90mm thick	Only Available in a 15" seat width, Black Spacer fabric	€483
<b>Matrix Mini Series Cushions</b>			
1626906	PS cushion with Kidabra print 310mm x 310mm, 12" x 12" <b>MAX User Weight 68kg</b>		€194
1626907	PS cushion with Kidabra print 310mm x 360mm, 12" x 14" <b>MAX User Weight 68kg</b>		€194
1626908	PS cushion with Kidabra print 310mm x 410mm, 12" x 16" <b>MAX User Weight 68kg</b>		€194
1627456	PS cushion with Kidabra print 330mm x 310mm, 13" x 12" <b>MAX User Weight 68kg</b>		€194
1627457	PS cushion with Kidabra print 330mm x 360mm, 13" x 14" <b>MAX User Weight 68kg</b>		€194
1627458	PS cushion with Kidabra print 330mm x 410mm, 13" x 16" <b>MAX User Weight 68kg</b>		€194
1626909	PS cushion with Kidabra print 360mm x 310mm, 14" x 12" <b>MAX User Weight 68kg</b>		€194
1626910	PS cushion with Kidabra print 360mm x 360mm, 14" x 14" <b>MAX User Weight 68kg</b>		€194
1626911	PS cushion with Kidabra print 360mm x 410mm, 14" x 16" <b>MAX User Weight 68kg</b>		€194
AQD2140	Without seat cushion		○

**BACKREST SUPPORT**

AQD2162	Backposts 540mm (Tall)	Must be selected for Flex3 adjustable backrest	◆
AQD2163	Backposts 480mm (Short)	Not available with Flex3 Adjustable backrest	○

**PUSH HANDLES FOR BACKPOSTS**












AQD2107	Backposts WITH push handles		◆
AQD2108	Backposts WITHOUT push handles Push handles will be cut down on the vertical section.	Not available with attendant controls, dual controls, tension adjustable backrest AQD2171 in conjunction with a headrest, Headrest AQD1495	€46

**BACKREST UPHOLSTERY / CUSHION**


AQD2171	Basic Tension Adjustable, (540mm or 480mm upholstery height) Straps easily accessible and adjustable from the rear, Black fabric only	Not available with laterals, headrest option 1496 Only available in 15" seat width	◆
AQD2172	Flex3 adjustable Backrest Plates (without Backrest Cushion)	Only available in a 15" seat width	€257
AQD2173	Flex3 adjustable Backrest Plates with Laguna Backrest cushion	Only available in a 15" seat width	€417
SP1562143	Mistral Flex 3 Backrest shell - shorter & narrower than the Laguna Flex 3	Only available in a 15" seat width	€132
SP1431668-TR26	Mistral Backrest Cushion, Black Dartex	Only available in a 15" seat width	€139.34
SP1431668-TR35	Mistral Backrest Cushion, Black Plush	Only available in a 15" seat width	€151.94

**MATRIX MINI SERIES BACKREST: Extended Mounting Pins Required with all Mini Backs**

EMP	Extended Mounting Pins		
MINI1207	Mini Back with Kidabra print 310mm x 180mm, 12" x 7" <b>MAX User Weight 68kg</b>		€586
MINI1210	Mini Back with Kidabra print 310mm x 260mm, 12" x 10" <b>MAX User Weight 68kg</b>		€586
MINI1212	Mini Back with Kidabra print 310mm x 310mm, 12" x 12" <b>MAX User Weight 68kg</b>		€586
MINI1214	Mini Back with Kidabra print 310mm x 360mm, 12" x 14" <b>MAX User Weight 68kg</b>		€586
MINI1307	Mini Back with Kidabra print 330mm x 180mm, 13" x 7" <b>MAX User Weight 68kg</b>		€586
MINI1310	Mini Back with Kidabra print 330mm x 260mm, 13" x 10" <b>MAX User Weight 68kg</b>		€586
MINI1312	Mini Back with Kidabra print 330mm x 310mm, 13" x 12" <b>MAX User Weight 68kg</b>		€586
MINI1314	Mini Back with Kidabra print 330mm x 360mm, 13" x 14" <b>MAX User Weight 68kg</b>		€586
MINI1407	Mini Back with Kidabra print 360mm x 180mm, 14" x 7" <b>MAX User Weight 68kg</b>		€586
MINI1410	Mini Back with Kidabra print 360mm x 260mm, 14" x 10" <b>MAX User Weight 68kg</b>		€586
MINI1412	Mini Back with Kidabra print 360mm x 310mm, 14" x 12" <b>MAX User Weight 68kg</b>		€586
MINI1414	Mini Back with Kidabra print 360mm x 360mm, 14" x 14" <b>MAX User Weight 68kg</b>		€586
AQD2170	Without Backrest Upholstery	Not available with laterals or Headrest AQD1496	○

BACKREST OPTIONS			
SP1571542 RH SP1571553 LH	Angle adjustable backrest support brackets	Not compatible with Powered Backrest recliner	€45.34
<b>UPHOLSTERY TYPE</b>			
AQD2151	Black fabric		◆
AQD2152	Black Dartex		○
<b>AQD2155</b>	<b>Black Spacer fabric</b>		○
<b>ARMREST (Pair)</b>			
<b>T-Armrest 2191</b>	<b>In-rigged armrest 2183</b>	<b>Hip supports</b>	<b>Short armpad</b>
			
<b>Long armpad</b>			
			
AQD2191	T-Armrest Standard Range of adjustment from seat plate 245mm to 315mm		◆
AQD2181	T-Armrest Extended height Range of adjustment from seat plate 295mm to 365mm		€46
AQD2183	In-rigged Armrests (reduces seat width by 25mm each side)	Must Be selected with a 12" & 13" Seat Width	€27
SP1582159	Inrigged Armrest Standard inrigged arm option MUST BE SELECTED with additional brackets.	Must Be selected with a 12" & 13" Seat Width, Not Available with hip support	€36.35
<b>ARMPADS (Pair)</b>			
AQD2195	Short Armpad (325mm)		◆
AQD2196	Long Armpad (375mm)	Not Available with Seat depth 410mm (16")	○
<b>HIP SUPPORTS, Multi adjustable LH (for flip up and following armrests only)</b>			
AQD2184	SMALL 140mm x 100mm	Not available with 410mm seat depth and manual centre legrests	€259
AQD2186	MEDIUM 230mm x 100mm	Not available with 410mm seat depth and manual centre legrests	€259
<b>HIP SUPPORTS, Multi adjustable RH (for flip up and following armrests only)</b>			
AQD2185	SMALL 140mm x 100mm	Not available with 410mm seat depth and manual centre legrests	€259
AQD2187	MEDIUM 230mm x 100mm	Not available with 410mm seat depth and manual centre legrests	€259
SEAT ACCESSORIES			
<b>HEADREST / NECKREST (Dartex covered) - Adequate headrest option must be selected</b>			
	<b>Headrest 1432</b>		<b>Matrix Elan Adjustment</b>
			<b>Elan Standard Pad, 2 sizes</b>
			<b>Elan Occipital Pad</b>
			<b>Elan 4 Point Pad</b>
			<b>Matrix Elan Vs Matrix mini Elan</b>
<b>ELAN SERIES HEADREST'S</b>			
1603255	Elan Standard Pad - Small (160mm x 90mm) Standard cover		€162
1603254	Elan Standard Pad - Medium (260mm x 130mm) Standard cover		€162
1603291	Elan Occipital Pad (230mm x 180mm) Standard cover		€162
1629822	Elan 4 Point Headrest Pad (Small) Standard cover		€300
<b>Please select Mounting hardware type for Elan headrest's</b>			
1603250	Matrix Elan 6-point Hardware Infinite adjustments		€238
1603251	Matrix Elan 6-point Mini (paediatric) 100mm shorter post		€238
1603249	Matrix standard multi-axis hardware Limited adjustability		€238
1603252	Elan adapter for Flex 3 back		€20
<b>STANDARD HEADREST OPTIONS</b>			
AQD1432	Headrest: height, depth and sideways adjustable		€ 226
AQD1496	All bracketry as option 1432 without headrest cushion		€84
AQD1495	Push bar dia. 22mm in conjunction with telescopic seat frame / backposts, No Headrest - no mounting hardware.		○

**LATERAL SUPPORT (Pair) For Flex3 Backs**



**Lateral support  
1446**



**Lateral multi-adjustable pad  
1445**

AQD1446		Rigid pad (Azalea type, Dartex)	€168
AQD1445		Lateral Support with extended arm, multi-adjustable pad	€168

**MATRX MINI LATERAL SUPORTS**

1627466		Swing-away lateral Hardware (each)	€156
EPAD-XS		Extra small lateral pad (each)	€81
1627465		Chest harness interface kit	€58

**SEAT LIFTER / SEAT TILT**

		SEAT HEIGHT			
		Tilt Range	420mm (16")	Selection	
AQD1093	Electric Tilt with COG Shift	0° - 30°	◆		
AQD1094	Seat Lifter with electric tilt (COG Shift), <b>Not possible with Rem A Control</b>	0° - 30°	€2,048		




**BACKREST RECLINER**



		Recline Range		Selection	
AQD0975	Manually setttable recline	0° - 30°	◆		
AQD0977	Powered Recliner with quick release for easy folding of backrest	0° - 30°	€318		

**PLEASE CHOOSE EITHER CENTRE MOUNTED OR SIDE MOUNTED LEGRESTS**

**MANUAL CENTRE MOUNT  
0556**



**POWERED CENTRE MOUNT  
0589**



**80° standard  
0526**



**CENTRE MOUNTED LEGRESTS (INCLUDING CALF PAD AND ANGLE ADJUSTABLE FOOTPLATE 240mm depth)**

AQD0556		Manual : Hanger, centre mounted, angle setttable (0° - 90°): knee to heel length 280mm - 385mm	€306
AQD0589		LNX Powered : Centre mounted legrest, length compensating	€1,548

**FOOTPLATE TYPE**

ARM0559		Angle Adjustable Footplate	○
---------	--	----------------------------	---

**SIDE MOUNTED LEGRESTS**

**FOOTREST HANGER TYPE, Swing away and detachable (Pair)**







AQD0526		Standard hanger 80°	Not Available with Angle adjustable footplate	◆
1606173-0601		Black short 90° Legrest Hanger x 2 - Action style.		€88.49
1606172-0601		Black short 80° Legrest Hanger x 2 - Action style.		€86.79
AQD0501		Without Hanger		○

**NB if either of the above legrest options is selected , the following items will be required:**

SP1531195		Legrest support carrier LH	€34.04
SP1531196		Legrest support carrier RH	€34.04

**FOOTPLATE TYPE (Pair), Standard legrest options.**

AQD0512		Fixed Footplate	◆
AQD0504		Angle Adjustable Footplate	○
AQD0513		Angle and depth adjustable Footplate (Adapt)	€116

FOOTPLATE OPTIONS: for use with Short 80° & 90° Legrest			
1428660 LH 1428654 RH		Angle and depth adjustable Footplate pair. Including stem	€32.41
1530620 LH 1530614 RH		Fixed Footplate pair. Including stem	€17.72
1492150 LH 1492149 RH		Upmounted Footplate	€108.79
REMOTE DX2 SYSTEM			
<b>REM Type A</b>	<b>REM 24SD</b>	<b>G90A</b>	<b>Chin control</b>
			
			<b>Dual control</b>
			
			<b>IDC Dual Control</b>
			
REMOTE & BRACKET			
AQD1638		REM Type A (enables operation up to 4 electric adjustments)	◆
AQD1601		ACS REM 24SD	€272
AQD1614		G90A	€201
AQD1207		Standard Bracket RH	○
AQD1208		Standard Bracket LH	○
AQD1707		Swing Away + height adjustable RH High End	€144
AQD1708		Swing Away + height adjustable LH High End	€144
ATTENDANT REMOTE & BRACKET			
AQD1643		REM Type A (Up to 4 electric actuators)	€18
AQD1602		ACS REM 24SD	€215
AQD1616		G90A	€215
AQD1207		Standard Bracket RH	○
AQD1208		Standard Bracket LH	○
CHIN CONTROL & BRACKET			
AQD1719		Chin Remote with standard Joystick (Always to be mounted on the Master Remote's side)	(Only compatible with Option1601 ) €2,293
AQB1289		Swing Away Bracket RH, side Manual	Only available with RH user bracket ○
AQD1288		Swing Away Bracket LH, side Manual	Only available with LH user bracket ○
AQD1287		Swing Away Bracket RH, Powered	Only available with RH user bracket €1,546
AQD1286		Swing Away Bracket LH, Powered	Only available with LH user bracket €1,546
DUAL CONTROL & BRACKET (can only be selected as an addition to user remote)			
AQD1720		Compact Dual Control	€495
AQD1207		Standard Bracket RH	Only available with LH user bracket ○
AQD1208		Standard Bracket LH	Only available with RH user bracket ○
INVACARE IDC DUAL CONTROL & BRACKET			
AQD1645		Invacare IDC Dual control	€681
AQD1221		IDC Dual control Detachable Bracket RH	○
AQD1222		IDC Dual control Detachable Bracket LH	○
DRIVE UNIT OPTIONS			
FRAME COLOUR (COLOURED PARTS IS A DETACHABLE AND CHANGEABLE PLASTIC SHROUD)			
AQD0227	X	Candy Red	◆
CENTRE WHEELS (DRIVING WHEELS)			
AQD0781	X	Centre Wheel 255mm x 80mm (10"- 3") Solid Insert / Puncture Proof Grey	◆
PERFORMANCE OPTIONS			
MOTOR TYPE			
AQD1282	X	6 km/h (4mph)	◆
BATTERY			
AQD1261	X	32Ah AGM pair	◆
CHARGER			
AQD1343		OptiCharge Charger 8A (200 - 250V)	◆
AQD1313		OptiCharge Charger 10A (200 - 250V) Charging time will be 20% faster	€158
SAFETY ACCESSORIES			
TRANSPORTATION KIT			
AQD1494	X	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory)	◆
FIXATION BELT			
AQD1329		Lapbelt (adjustable both sides) fitted to base of backrest	○
AQD2113		Lapbelt (adjustable both sides) adjustable mounting on seat rail extrusion	Not Available with seat depth 410mm (16") and T-Armrest ◆

BODYPOINT BELTS		
HB216-S38-B2	HB216 2-point belt Small x 380mm Pads, cinch mount	€84
HB216-M46-B2	HB216 2-point belt Medium x 460mm Pads, cinch mount	€84
HB219-S38-B2	HB219 2-point belt Small x 380mm Pads, cinch mount	€96
HB219-M46-B2	HB219 2-point belt Medium x 460mm Pads, cinch mount	€96
HB235-S32-B2	HB235 2-point belt Extra small x, 320mm Pads, cinch mount	€100
HB235-S38-B2	HB235 2-point belt Small x 380mm Pads, cinch mount	€100
EB235S	Evoflex Small pelvic stabilizer, push button	€120
EB225M	Evoflex Medium pelvic stabilizer, Latch	€120
EB225L	Evoflex Large pelvic stabilizer, Latch	€120

BODYPOINT HARNESES		
SH250XXS-B2	PivotFit XX-Small, front pull, cinch mount, dynamic trunk support	€143
SH250XS-B2	PivotFit X-Small, front pull, cinch mount, dynamic trunk support	€143
SH250S-B2	PivotFit Small, front pull, cinch mount, dynamic trunk support	€143
SH250M-B2	PivotFit Medium, front pull, cinch mount, dynamic trunk support	€143
SH250XXS-C2-B2	PivotFit XX-Small, shoulder harness, Plain top strap, Cinch mount	€152
SH250XS-C2-B2	PivotFit Small, shoulder harness, Plain top strap, Cinch mount	€152
SH250M-C2-B2	PivotFit Medium, shoulder harness, Plain top strap, Cinch mount	€152
SH260XXS-B2	PivotFit XX-small, rear pull, cinch mount, standard	€152
SH260XS-B2	PivotFit X-small, rear pull, cinch mount, standard	€152
SH260S-B2	PivotFit Small, rear pull, cinch mount, standard	€152
SH260M-B2	PivotFit Medium, rear pull, cinch mount, standard	€152
SH270XXS-B2	PivotFit XX-Small, Front-Pull cinch mount, dynamic trunk support	€152
SH270XS-B2	PivotFit X-Small, Front-Pull cinch mount, dynamic trunk support	€152
SH270S-B2	PivotFit Small, Front-Pull, cinch mount, dynamic trunk support	€152
SH270M-B2	PivotFit Medium, Front-Pull cinch mount, dynamic trunk support	€152
SH270XXS-C2-B2	PivotFit XXS, Front pull, Plain top strap, cinch mount	€158
SH270XS-C2-B2	PivotFit XS, Front pull, Plain top strap, cinch mount	€158
SH270S-C2-B2	PivotFit Small, Front pull, Plain top strap, cinch mount	€158
SH270M-C2-B2	PivotFit Medium, Front pull, Plain top strap, cinch mount	€158
SH280XXS-B2	PivotFit XX-Small, rear pull, cinch mount, dynamic	€158
SH280XS-B2	PivotFit X-Small, rear pull, cinch mount, dynamic	€158
SH280S-B2	PivotFit Small, rear pull, cinch mount, dynamic	€158
SH280M-B2	PivotFit Medium, rear pull, cinch mount, dynamic	€158
1627465	Chest harness interface kit	€48

### ACCESSORIES

#### BAG HOOKS *NEW*



AQD2128	Bag hooks (Max load up to 6 kg = impact on max user weight)	Not available with Flex3 Backrest 2172/2173, Attendant Brackets, Dual Control	€72
---------	---	---	-----

#### CANE HOLDER

AQD1312	Cane Holder RH	€39
AQD1363	Cane Holder LH	€37

#### TRAY (PLEASE CHOOSE OPPOSITE HAND TO USER CONTROL)

AQD1404	Tray Swing Away RH	€350
AQD1405	Tray Swing Away LH	€350

#### TOTAL PRICE



Need a different option not on the form?  
Take a look at our Powerchair Alternative Options form, available as standard.

Powerchair Alternative Option required:



Please refer to our Invacare Unique special order form, for more complex custom requirements. Available at [www.invacare.co.uk](http://www.invacare.co.uk)

**Special Instructions - Please supply the following information:**

**Name:**

**Email:**

**Phone:**

**Fax:**

**Description of special request (please include relevant option code or part number):**