

Tel: 01 810 7084 Fax: 01 810 7085  
Email: ordersireland@invacare.com  
Online Spares available at www.invacare.ie

JANUARY 2020

Customer Ref :	Account No :
Contact No :	Order No :
Delivery Address :	Order Qty :
	Order Date :



# Pronto M41 Ireland

## Standard Seating

Max user weight 136kg (21.4 stone)

### Standard Features




- Standard Seating System
- Adjustable Seat Widths from 440mm - 480mm (17" - 19")
- Adjustable Seat Depth Set at 460mm (410mm - 510mm) (16", 18" or 20")
- Backrest Height 480mm or 540mm (19" - 21.25")
- Black Nylon Tension Adjustable Backrest Upholstery
- Detachable and height adjustable armrest assembly with extended armpad
- Choice of 5 Seat Heights 455mm - 555mm (18" - 22")
- Manually settable backrest Recline
- Adjustable flip-up centre footboard
- Shark2 REM D01 control with standard bracket
- Candy Red Shrouds
- 255mm x 80mm (10" x 3") Solid Insert / Puncture Proof Grey Drive Wheels
- 6 km/h (4mph) Motors and 35Ah AGM Batteries
- 8A (200-250V) Charger
- 4-Point tie-down system



**Standard Base Price: €3,068**

\* Configurations may vary from the image shown

	X	Please make your selection with a cross in the box	◆ Standard feature	○ Free of Charge
Customer weight =      Kg				
Code		Description		Ireland
AQB0000		Pronto M41 Basic Model	◆	
AQZ0000		Pronto M41 Special Basic Model - PLEASE COMPLETE DETAILS AT THE END OF THE FORM		POA
<b>SEATING TYPE</b>				
AQB1150	x	Standard Seat Unit		€3,068
<b>SEAT WIDTH &amp; SEAT CUSHION</b>				
AQB0343		Seat Width of Seat Frame 440mm - 480mm (17" - 19")	◆	
AQB1336		Seat Cushion (width 430mm, 50mm thick)	◆	
<b>SEAT DEPTH</b>				
AQB0327	x	460mm (adjustable to 410mm or 510mm) (16", 18" or 20")	◆	
<b>BACKREST HEIGHT (without Seat Cushion)</b>				
AQB0448		480mm (19")	◆	
AQB0454		540mm (21.25")	○	
<b>BACK CUSHION</b>				
AQB1178	x	Tension Adjustable (Black Nylon)	◆	
<b>UPHOLSTERY TYPE</b>				
AQB0268	x	Upholstery colour Black (Nylon)	◆	

<b>ARMREST (Pair)</b>			
AQB0606	<b>x</b>	Detachable Armrest, height adjustable (290mm - 360mm without seat cushion)	◆
<b>ARMPADS (Pair)</b>			
AQB0649	<b>x</b>	Armpad extended (+ 50mm)	◆
<b>SEAT ACCESSORIES</b>			
<b>HEADREST / NECKREST (Dartex) - Adequate headrest option must be selected</b>			
AQB1432		Headrest : Depth, Width and Height adjustable	€226
AQB1495		Push bar dia. 22mm in conjunction with telescopic seat frame / backposts, No Headrest - no mounting hardware.	○
<b>SEAT HEIGHT</b>			
AQB0325		455mm (18")	○
AQB0339		480mm (19")	○
AQB0387		505mm (20")	○
AQB0338		530mm (21")	○
AQB0389		555mm (22")	○
<b>SEAT ADJUSTMENTS</b>			
<b>BACKREST RECLINER</b>			
AQB0959	<b>x</b>	Manual settable reclining angle by hand screws (4°/14°/21.5°/29°/36.5°/44°)	◆
<b>CENTRE MOUNTED LEGREST</b>			
<b>CENTRE MOUNTED FIXED FOOTBOARD</b>			
AQB0595	<b>x</b>	Adjustable flip-up footboard (no hangers)	◆
<b>CONTROL OPTION</b>			
<b>REMOTE &amp; BRACKET</b>			
AQB1801	<b>x</b>	Standard Remote Shark 2 REM D01	◆
AQB1207		Standard Bracket RH	○
AQB1208		Standard Bracket LH	○
AQB1707		Swing Away and Height Adjustable Bracket RH	€144
AQB1708		Swing Away and Height Adjustable Bracket LH	€144
<b>DRIVE UNIT OPTIONS</b>			
<b>SHROUD COVER COLOUR</b>			
AQB0227	<b>x</b>	Candy Red	◆
<b>WHEELS</b>			
AQB0781	<b>x</b>	Centre Wheel 255mm x 80mm (10" - 3") Solid Insert / Puncture Proof Grey	◆
<b>PERFORMANCE OPTIONS</b>			
AQB1282		Motor type 3 km/h (2mph)	◆
AQB1294		Motor type 6 km/h (4mph)	○
AQB1261	<b>x</b>	Battery 35Ah AGM pair	◆
<b>CHARGER</b>			
AQB1343		OptiCharge Charger 8A (200 - 250V)	◆
AQB1313		OptiCharge Charger 10A (200 - 250V) Charging time will be 20% faster	€158
<b>SAFETY OPTIONS</b>			
<b>TRANSPORTATION KIT</b>			
AQB1494	<b>x</b>	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) <b>(Adequate headrest is compulsory)</b>	◆
<b>FIXATION BELT</b>			
AQB1329	<b>x</b>	Metal Buckle Lap Belt (Adjustable from Both Sides)	◆
			<b>Total Price</b>
			
<p>Need a different option not on the form? Take a look at our Powerchair Alternative Options form, available as standard <b>Powerchair Alternative Option required:</b></p>			



**Special Instructions - Please supply the following information:**

**Name:**

**Email:**

**Phone:**

**Fax:**

**Description of special request (please include relevant option code or part numbers):**