

Tel: 01 810 7084 Fax: 01 810 7085  
Email: ordersireland@invacare.com  
Online Spares available at www.invacare.ie

JANUARY 2020

Customer Ref :	Account No :
Contact No :	Order No :
Delivery Address :	Order Qty :
	Order Date :



# Pronto M41 Ireland

## One piece Seat plate

Max user weight 136kg (21.4 stone)

### Standard Features



- Compact, manoeuvrable, **MWD Pronto M41 Base**
- One Piece Seat Plate width 430mm (17") (adjustable by arms 440mm - 480mm, 17" - 19") with standard foam cushion (50mm, 2")
- Seat Depths 410mm - 510mm (16" - 20")
- Tension adjustable backrest upholstery
- Standard T-Armrests, detachable with Short armpad
- Fixed seat Tilt (0° - 4° - 6°)
- Fixed seat heights of 450mm - 470mm - 490mm
- Manual settable backrest Recline (0° - 30°)
- Adjustable flip-up centre footboard
- Shark 2 REM D01 control with standard bracket
- Candy Red Shroud Colour
- 255mm x 80mm (10" x 3") Solid Insert / Puncture Proof Grey Drive Wheels
- 6 km/h (4mph) In-line Motors and 35Ah Batteries
- 8A (200 - 250V) Charger
- 4-Point tie-down system
- Lapbelt with autolock (adjustable both sides)



\* Configurations may vary from the image shown

**Standard Base Price : €3,020**

	X	Please make your selection with a cross in the box	◆ Standard feature	○ Free of Charge
Customer weight =            Kg				
Code		Description		Ireland
AQC0059		Pronto M41 "one piece seating"	◆	
AQX0059		Pronto M41 "one piece seating" Special Basic Model-PLEASE COMPLETE DETAILS AT THE END OF THE FORM		POA
<b>SEATING TYPE</b>				
AQC2100	X	One Piece Standard Seating Unit		€3,020
<b>SEAT FRAME AND SEAT PLATE</b>				
AQC2110	X	One Piece Standard Seat Plate	◆	
<b>SEAT WIDTH SETTING FOR ADJUSTABLE SEAT FRAME</b>				
AQC2122	X	430mm - 480mm (17" - 19") Seat Plate (additional 50mm, 2" adjustment on armrests)	◆	
<b>SEAT DEPTH SETTING</b>				
AQC2132		410mm (16") (adjustable 410mm - 510mm stepless adjustment)		○
AQC2134		460mm (18") (adjustable 410mm - 510mm stepless adjustment)		○
AQC2136		510mm (20") (adjustable 410mm - 510mm stepless adjustment)		○

<b>SEAT CUSHION (sizes according to seat frame settings)</b>			
AQC2141		Standard Foam Cushion 50mm, 2" thick	◆
AQC2140		without seat cushion	○
<b>BACKREST SUPPORT</b>			
AQC2161	X	Backrest frame, detachable, height 480mm 19" (one size, manual settings only)	◆
<b>BACKREST UPHOLSTERY / CUSHION</b>			
AQC2171	X	Tension Adjustable Backrest (540mm or 480mm upholstery height) Straps easily accessible and adjustable from the rear	◆
<b>UPHOLSTERY TYPE</b>			
AQC2151	X	Black fabric	◆
<b>ARMREST (Pair)</b>			
AQC2191	X	Single post 'T-Armrest' removable and height adjustable Range of adjustment from seat plate 245mm to 315mm	◆
<b>ARMPAD (Pair)</b>			
AQC2195		Short Armpad (325mm)	◆
AQC2196		Long Armpad (375mm)	○
			Not available with seat depth 410mm (16")
<b>SEAT ACCESSORIES</b>			
<b>HEADREST / NECKREST (Dartex) - Adequate headrest option must be selected</b>			
AQC1350		Neckrest: height and depth adjustable	€81
AQC1495		Push bar dia. 22mm in conjunction with telescopic seat frame / backposts, No Headrest - no mounting hardware.	◆
<b>SEAT ADJUSTMENTS</b>			
AQC1091	X	Manually settable Seat Tilt	◆
<b>SEAT TO FLOOR HEIGHT (only to choose with option 1091)</b>			
AQC0354		Seat height 450mm	○
AQC0364		Seat height 470mm	○
AQC0382		Seat height 490mm	○
<b>BACKREST RECLINER</b>			
AQC0975		Manual Settable Recliner (in 3.75° steps)	◆
<b>PLEASE CHOOSE EITHER CENTRE MOUNTED OR SIDE MOUNTED LEGRESTS</b>			
<b>CENTRE MOUNTED FIXED FOOTBOARD</b>			
AQC0595		Adjustable flip-up footboard (no hangers)	◆
<b>SIDE MOUNTED LEGRESTS</b>			
<b>FOOTREST HANGER TYPE (Pair)</b>			
AQC0526		Standard 80° hangers	○
<b>FOOTPLATE TYPE (Pair)</b>			
AQC0512		Fixed Footplate	◆
AQC0513		Angle and depth adjustable Footplate (Adapt)	€116
<b>CONTROL OPTION</b>			
<b>REMOTE &amp; BRACKET</b>			
AQC1801	X	Standard Remote Shark2 REM D01	◆
AQC1207		Standard Bracket RH	○
AQC1208		Standard Bracket LH	○
AQC1707		Swing Away + height adjustable RH high end	○
AQC1708		Swing Away + height adjustable LH high end	○
<b>DRIVE UNIT OPTIONS</b>			
<b>FRAME COLOUR (COLOURED PARTS IS A DETACHABLE AND CHANGEABLE PLASTIC SHROUD)</b>			
AQC0227	X	Candy Red	◆
<b>CENTRE WHEELS (DRIVING WHEELS)</b>			
AQC0781	X	Centre Wheel 255mm x 80mm (10" - 3") Solid Insert / Puncture Proof Grey	◆
<b>PERFORMANCE OPTIONS</b>			
<b>MOTOR TYPE</b>			
AQC1294		3 km/h (2mph)	○
AQC1282		6 km/h (4mph)	◆
<b>BATTERY</b>			
AQC1261	X	35Ah AGM pair	◆
<b>CHARGER</b>			
AQC1343		OptiCharge Charger 8A (200 - 250V)	◆
AQC1313		OptiCharge Charger 10A (200 - 250V) Charging time will be 20% faster	€158

**SAFETY ACCESSORIES****TRANSPORTATION KIT**

AQC1494	X	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) <b>(Adequate headrest is compulsory)</b>	◆
<b>FIXATION BELT</b>			
AQC1329	X	Lapbelt (adjustable both sides) fitted to base of backrest	◆

**ACCESSORIES****CANE HOLDER**

AQC1312		Cane Holder RH	€58
AQC1363		Cane Holder LH	€58
<b>TRAY (PLEASE CHOOSE OPPOSITE HAND TO USER CONTROL)</b>			
AQC1404		Tray Swing Away RH	€201
AQC1405		Tray Swing Away LH	€201

**TOTAL PRICE**

Need a different option not on the form?

Take a look at our Powerchair Alternative Options form, available as standard

**Powerchair Alternative Option required:**



**Special Instructions - Please supply the following information:**

**Name:**

**Email:**

**Phone:**

**Fax:**

**Description of special request (please include relevant option code or part number):**