

Please refer to www.invacare.ie for Local Distributor

Customer name:	Account No :
Contact no:	Order No :
Delivery address:	Order Qty :
	Order Date :

Propad Revolve Si (Silicon gel insert)
Retail order Form
Maximum User Weight 190Kg (30 Stone) *

		WIDTH									
cm		36	38	41	43	46	48	51	56	61	
inch		14	15	16	17	18	19	20	22	24	
DEPTH	36	■	■	■							
	38	■	■	■							
	41	■	■	■		■					
	43				■	■	■				
	46			■		■			■	■	
	51								■	■	


Can be used any way up in any direction



Lightweight comfortable foam that moulds to patient shape, providing excellent comfort and pressure redistributing properties



Durable high resilient foam that can be used anyway up




Utilises silicon fluid sac from Flo-tech range for pressure redistribution

Cushion protected with a multi stretch water-resistant cover

One piece moulded design excellent for side transfers

Standard features:

- Very High risk
- 3" height
- Flat - any way up
- Foam and fluid combination
- Sewn cover
- Weight limit: 30 stone



Code	X	Description	Price
Complete Standard cushion up to *159kg user weight			
1629479		Revolve Si with platilon liner 14x14" (36x36cm)	Sewn blue wipedown cover € 229.00
1629480		Revolve Si with platilon liner 14x15" (36x38cm)	Sewn blue wipedown cover € 229.00
1629481		Revolve Si with platilon liner 14x16" (36x41cm)	Sewn blue wipedown cover € 229.00
1629482		Revolve Si with platilon liner 15x14" (38x36cm)	Sewn blue wipedown cover € 229.00
1629483		Revolve Si with platilon liner 15x15" (38x38cm)	Sewn blue wipedown cover € 229.00
1629484		Revolve Si with platilon liner 15x16" (38x41cm)	Sewn blue wipedown cover € 229.00
1629485		Revolve Si with platilon liner 16x14" (41x36cm)	Sewn blue wipedown cover € 229.00
1629486		Revolve Si with platilon liner 16x15" (41x38cm)	Sewn blue wipedown cover € 229.00
1601977		Revolve Si with platilon liner 16x16" (41x41cm)	Sewn blue wipedown cover € 229.00
1601994		Revolve Si with platilon liner 16x18" (41x46cm)	Sewn blue wipedown cover € 229.00
1601993		Revolve Si with platilon liner 17x17" (43x43cm)	Sewn blue wipedown cover € 229.00
1601995		Revolve Si with platilon liner 18x16" (46x41cm)	Sewn blue wipedown cover € 229.00
1629487		Revolve Si with platilon liner 18x17" (46x43cm)	Sewn blue wipedown cover € 229.00
1601996		Revolve Si with platilon liner 18x18" (46x46cm)	Sewn blue wipedown cover € 229.00
1601997		Revolve Si with platilon liner 19x17" (48x43cm)	Sewn blue wipedown cover € 229.00
1601998		Revolve Si with platilon liner 20x18" (51x46cm)	Sewn blue wipedown cover € 229.00
1602472		Revolve Si with platilon liner 20x20" (51x51cm)	Sewn blue wipedown cover € 229.00

1629621	Revolve Si with platilon liner 14x14" (36x36cm)	Sewn black wipedown cover	€ 229.00
1629622	Revolve Si with platilon liner 14x15" (36x38cm)	Sewn black wipedown cover	€ 229.00
1629623	Revolve Si with platilon liner 14x16" (36x41cm)	Sewn black wipedown cover	€ 229.00
1629624	Revolve Si with platilon liner 15x14" (38x36cm)	Sewn black wipedown cover	€ 229.00
1629625	Revolve Si with platilon liner 15x15" (38x38cm)	Sewn black wipedown cover	€ 229.00
1629626	Revolve Si with platilon liner 15x16" (38x41cm)	Sewn black wipedown cover	€ 229.00
1629627	Revolve Si with platilon liner 16x14" (41x36cm)	Sewn black wipedown cover	€ 229.00
1629628	Revolve Si with platilon liner 16x15" (41x38cm)	Sewn black wipedown cover	€ 229.00
1627596	Revolve Si with platilon liner 16x16" (41x41cm)	Sewn black wipedown cover	€ 229.00
1627598	Revolve Si with platilon liner 16x18" (41x46cm)	Sewn black wipedown cover	€ 229.00
1627597	Revolve Si with platilon liner 17x17" (43x43cm)	Sewn black wipedown cover	€ 229.00
1627599	Revolve Si with platilon liner 18x16" (46x41cm)	Sewn black wipedown cover	€ 229.00
1629629	Revolve Si with platilon liner 18x17" (46x43cm)	Sewn black wipedown cover	€ 229.00
1627600	Revolve Si with platilon liner 18x18" (46x46cm)	Sewn black wipedown cover	€ 229.00
1627601	Revolve Si with platilon liner 19x17" (48x43cm)	Sewn black wipedown cover	€ 229.00
1627602	Revolve Si with platilon liner 20x18" (51x46cm)	Sewn black wipedown cover	€ 229.00
1627603	Revolve Si with platilon liner 20x20" (51x51cm)	Sewn black wipedown cover	€ 229.00
1629562	Revolve Si with Aquaguard zip 14x14" (36x36cm)	Welded baltic blue wipedown cover	€ 252.00
1629563	Revolve Si with Aquaguard zip 14x15" (36x38cm)	Welded baltic blue wipedown cover	€ 252.00
1629564	Revolve Si with Aquaguard zip 14x16" (36x41cm)	Welded baltic blue wipedown cover	€ 252.00
1629565	Revolve Si with Aquaguard zip 15x14" (38x36cm)	Welded baltic blue wipedown cover	€ 252.00
1629566	Revolve Si with Aquaguard zip 15x15" (38x38cm)	Welded baltic blue wipedown cover	€ 252.00
1629567	Revolve Si with Aquaguard zip 15x16" (38x41cm)	Welded baltic blue wipedown cover	€ 252.00
1629568	Revolve Si with Aquaguard zip 16x14" (41x36cm)	Welded baltic blue wipedown cover	€ 252.00
1629569	Revolve Si with Aquaguard zip 16x15" (41x38cm)	Welded baltic blue wipedown cover	€ 252.00
1602138	Revolve Si with Aquaguard zip 16x16" (41x41cm)	Welded baltic blue wipedown cover	€ 252.00
1602140	Revolve Si with Aquaguard zip 16x18" (41x46cm)	Welded baltic blue wipedown cover	€ 252.00
1602139	Revolve Si with Aquaguard zip 17x17" (43x43cm)	Welded baltic blue wipedown cover	€ 252.00
1602141	Revolve Si with Aquaguard zip 18x16" (46x41cm)	Welded baltic blue wipedown cover	€ 252.00
1629570	Revolve Si with Aquaguard zip 18x17" (46x46cm)	Welded baltic blue wipedown cover	€ 252.00
1602142	Revolve Si with Aquaguard zip 18x18" (46x46cm)	Welded baltic blue wipedown cover	€ 252.00
1602143	Revolve Si with Aquaguard zip 19x17" (48x43cm)	Welded baltic blue wipedown cover	€ 252.00
1602144	Revolve Si with Aquaguard zip 20x18" (51x46cm)	Welded baltic blue wipedown cover	€ 252.00
1602615	Revolve Si with Aquaguard zip 20x20" (51x51cm)	Welded baltic blue wipedown cover	€ 252.00

* Maximum user weight dependant on size

Complete Bariatric cushion up to 190kg User Weight

1630661	Revolve Si with platilon liner 22x18" (56x46cm)	Sewn blue wipedown cover	€ 229.00
1630662	Revolve Si with platilon liner 22x20" (56x51cm)	Sewn blue wipedown cover	€ 229.00
1630663	Revolve Si with platilon liner 24x18" (61x46cm)	Sewn blue wipedown cover	€ 229.00
1630664	Revolve Si with platilon liner 24x20" (61x51cm)	Sewn blue wipedown cover	€ 229.00
1630681	Revolve Si with platilon liner 22x18" (56x46cm)	Sewn black wipedown cover	€ 229.00
1630682	Revolve Si with platilon liner 22x20" (56x51cm)	Sewn black wipedown cover	€ 229.00
1630683	Revolve Si with platilon liner 24x18" (61x46cm)	Sewn black wipedown cover	€ 229.00
1630684	Revolve Si with platilon liner 24x20" (61x51cm)	Sewn black wipedown cover	€ 229.00
1630673	Revolve Si with Aquaguard zip 22x18" (56x46cm)	Welded baltic blue wipedown cover	€ 252.00
1630674	Revolve Si with Aquaguard zip 22x20" (56x51cm)	Welded baltic blue wipedown cover	€ 252.00
1630675	Revolve Si with Aquaguard zip 24x18" (61x46cm)	Welded baltic blue wipedown cover	€ 252.00
1630676	Revolve Si with Aquaguard zip 24x20" (61x51cm)	Welded baltic blue wipedown cover	€ 252.00

TOTAL PRICE:

Special Instructions - Extended lead times will apply

Name: _____ **Email:** _____ **Phone:** _____



Description of special request (please include relevant option code):