

Propad[®] Overlay

Pressure redistributing overlay

Placed on a conventional mattress for increased comfort and pressure redistribution.



High quality foam overlay

For enhanced pressure redistribution



Yes, you can.[®]

Transforms conventional mattresses



Pressure Redistribution

- › Larger, more supportive cutouts at the centre and smaller cutouts at each end, enables more even pressure redistribution.
- › Design conforms to the form of the body and cradles bony prominences to reduce pressure points.



Shear & Friction Reduction

- › Independent, geometrically cut surface cutouts reduce shear and friction during repositioning.
- › Multi-stretch, water-resistant and vapour-permeable cover further reduces shear and friction.



Comfort

- › Castellated foam conforms to movement and body position, providing excellent comfort.
- › The base of the cover has an anti-slip material, which prevents slippage.

Technical data ▼

For more information about this product, including the user manual, please visit our website: www.invacare.co.uk



Sizes Available ▼



Max. User Weight ▼



Max. Cleaning Temperature ▼

Height	Sizes Available	Max. User Weight	Max. Cleaning Temperature
85 mm height	1875 x 680 mm 1875 x 875 mm 1970 x 860 mm 2000 x 850 mm	108 kg	95 °C
100 mm height	1780 x 780 mm	108 kg	95 °C

Manufactured to: BS EN 597-1 Cig, BS EN 597-2 Match and BS 7177 (ignition Source 5).

Essential nursing care is pivotal in pressure injury management. This mattress may positively contribute to the outcome of a pressure injury prevention plan that includes education, clinical judgment and action planning based on vulnerability, which are fundamental factors in the prevention of pressure injuries. A range of assessment scales can be used as a formal method of assessing risk from pressure injury development, and should be used in conjunction with an informal assessment (informed nursing judgement) as part of a comprehensive and individualised plan.

T: +44 1656 776222
uk@invacare.com
www.invacare.co.uk
www.invacare.ie

Invacare UK & Ireland
 Unit 4 Pencoed Technology
 Park
 Bridgend CF35 5AQ



Yes, you can.®



© 2024 Invacare International GmbH. All rights reserved. All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation. Propad® Overlay - uk - 11/2024