



## A Refined Scooter Control Solution

The new **R-series 90** means a refined drive experience, blending innovative performance improvements with an unprecedented balance of power and size.

### Performance Features:

- 90 Amps of power when you need it.
- Improved motor matching algorithms for improved curb-climbing and hill-starting capabilities.
- Dynamic's unique Speed Reduction Wiper (SRW) provides improved safety by increasing stability and control when turning.
- Drop in replacement with industry standard connections and mounting support.
- Programmable acceleration curves, zero roll back on slopes.
- 2 drive profiles, brake and reverse lights, reversing beeper and electronic park brake release.
- Intelligent motor and battery management features provide automatic power flow optimisation.
- Advanced diagnostics to provide indepth logging and analysis of battery performance / abuse.
- A Smart Servicing Scheduler calculates and indicates when servicing is due.

# R-series 90



## Physical Specifications:

Supply Voltage	24V DC
Operating Voltage	18–32V DC
Peak Current	70 Amps
Boost Current	20 Amps
Boost Duration	1–10 seconds
Brake Voltage	24V DC
Charging Current	8 Amps RMS

Protection Rating	Electronics to IPx5
Operating Temperature	-25°C to 50°C
Storage Temperature	-40°C to 65°C
Materials	Die-cast aluminium base with polymer cover (UL94 V-0 rating)
Part Numbers	DR90-A01

Corporate Office New Zealand  
Ph: +64-3-962-2519  
Fax: +64-3-962-2966  
Email: [info@dynamiccontrols.com](mailto:info@dynamiccontrols.com)

USA  
Ph: +1-440-979-0657  
Fax: +1-440-979-1028  
Email: [usa@dynamiccontrols.com](mailto:usa@dynamiccontrols.com)

Asia  
Ph: +886-955-335-243  
Fax: +886 225-981-562  
Email: [asia@dynamiccontrols.com](mailto:asia@dynamiccontrols.com)

Europe  
Ph: +44-1905-772-321  
Fax: +44-1905-827-520  
Email: [europe@dynamiccontrols.com](mailto:europe@dynamiccontrols.com)



**dynamic**™

Dynamic Controls is the world's leading manufacturer of electronic controls for power wheelchairs and scooters. Dynamic was established in 1972 and is headquartered in New Zealand. Regional centres are located in Europe, United States, Asia, and Australasia. ISO 13485 accredited – Dynamic goes above and beyond industry standard expectations to ensure customers receive the best products possible.

For more information about Dynamic Controls visit [www.dynamiccontrols.com](http://www.dynamiccontrols.com)