

Raizer[®] M

Lifting Chair



Safe Moving & Handling

Falls Management



Feel safe and comfortable with Raizer M

Raizer M is a manually operated lifting chair for use in case of a fall. Quick assembly and ease of use allows the person to be lifted from floor to almost standing in just a few minutes.

Raizer M is most suited for use in a private home and it can be operated by anyone. It takes no special skills to operate it, and it cannot be assembled incorrectly.

Legs and backrests fit in only where they belong, and clear assembly

indicators verify that the parts have been assembled correctly.

The chair is operated by means of a crank handle – in an ergonomically correct manner and with no strain on the helper.



Always Ready

Manually operated by crank handle.



Max. 160 kg

Safe lifting capacity for a wide range of users.



Few minutes

From floor to near-standing in just a few minutes.



> **Optional headrest**
Prevents neck strain during the lift.

> **Seat belt included**
Safely stabilises the fallen person throughout.

> **Crank handle**
Insertable from both sides – ergonomic and effort-free.

Technical data ▼



Lying Dimensions ▼

132cm x 88cm x
27cm



Standing Dimensions ▼

54cm x 88cm x
127cm



SWL ▼

Max. 150kg



Weight ▼

Seat: 7.2kg
Legs and backrest:
5.2kg
Total: 12.4kg
