



**Yes, you can:**

Tel: 01656 776222  
 Email: ordersuk@invacare.com  
 Online Spares available at: www.invacare.co.uk

# Manual Invacare Ireland Prescription Form

07 JANUARY 2026  
 Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## Rea Azalea Self Propel Tall

### Standard Features

- Max User Weight 135KG/21 Stone
- Crash Tested ISO 7176-Part 19
- 24" Self Propulsion with User Manoeuvred Tyre Brakes and Solid Black Tyres
- Choice of Seat Widths 390mm, 440mm or 490mm (+50mm Width Adjustment)
- Seat Height 500mm
- Depth Adjustable Seat 480mm-550mm
- Without Seat Cushion, Supplied with Seat Board
- Flex 3 Backrest Plate, Width and Height Adjustable
- Passad 2 Backrest Cushion with Multi-stretch Black PU coated fabric Cover and Push Handles Connected
- Assistant Operated Gas Piston Seat Tilt +1° - +24°
- Backrest Recline 30° with Tilt Lock
- Angle Adjustable Legrests with Calf Support, Height and Depth Adjustable Footplates
- Height Adjustable Ducktail Pad Armrests with Autolock
- Pearl Silver Frame Colour
- Chassis Extender with Solid Black Front Castors
- Quick Release Self Propel Wheels, Solid Tyres and Aluminium Handrims
- Pelvic Belt, Adjustable both sides
- Anti Tipplers

**Standard Base price: €3,146**



Please make your selection with a cross in the box    ○ Free of charge    ◆ standard feature    POA special request

Customer weight =                      KG

Code	Description		
<b>BASE PRODUCT*</b>			
BNF9003	Rea Azalea Self Propel Tall		€3,146
BNF8000	SPECIAL REQUEST - PLEASE INCLUDE QUOTATION REFERENCE IN SPECIAL DESCRIPTION AT THE END OF THE FORM WHERE APPLICABLE		POA
<b>PROPULSION*</b>			
BNF0025	Self Propel 24" wheels with user brakes	(NC 1076)	◆
BNF0030	Self Propel 24" wheels with assistant hub brakes	(NC 1074)	€162
<b>SEAT DIMENSIONS*</b>			
BNF0036	Tall (+ 50 mm In seat height, depth and chassis length)		◆

<b>SEAT WIDTH*</b>				
BNF0040		390 mm/15" (Small chassis)	(NC 0244, 0245, 0246, 0247)	○
BNF0050		440 mm/17" (Medium chassis)	(NC 0241, 0242, 0246, 0247)	○
BNF0060		490 mm/19" (Large chassis)	(NC 0241, 0242, 0243, 0244)	○
<b>SEAT WIDTH EXTENSION*</b>				
BNF0065		+ 50 mm seat width pre-adjustment on armrests and legrests		○
BNF0085		Without seat width extension, set at 0mm adjustment but can be adjusted wider if required		○
<b>SEAT DEPTH*</b>				
BNF0105		480-550 mm / 19" -22" (Azalea Tall - Incl Chassis Extension)		◆
<b>SEAT HEIGHT*</b>				
BNF0125		500 mm / 20" (Azalea Tall)		◆
<b>SEAT ANGLE</b>				
BNF0127		Negative Seat Angle -1 To +23° (For attendant operated gas piston in negative tilt)	(NC 0360)	○
<b>SEAT TYPE*</b>				
BNF0130		Depth Adjustable		◆
<b>SEAT CUSHION*</b>				
BNF0165		Cushion Floshape seat upholstery	(NC 1111)	€72
BNF0167		Cushion Vicair Adjuster O2 10cm, complete with Black Comfair cover	(NC 0191, 0201, 0202)	€696
BNF0168		Cushion Matrx Libra, complete with Black cover	(NC 0191, 0201, 0202)	€645
BNF0169		Cushion Flo-Tech Solution Xtra, complete with Black cover	(NC 0191, 0201, 0202)	€566
BNF0179		Cushion Vicair Vector O2 100mm with 2 top covers (Additional 15 day leadtime)	(NC 0182, 0183, 0186, 0187, 0188, 0189, 0191, 0201, 0202)	€800
BNF0180		Without Seat Cushion, Seat Plate only	(NC 0191, 0201, 0202)	○
<b>SEAT DEPTH OF MATRX LIBRA AND FLO-TECH SOLUTION XTRA CUSHIONS (OPTIMISED CUSHION DEPTH TO SUIT SEAT DEPTH SETTING)</b>				
BNF0182		Short cushion depth: Libra delivered at 405mm deep, Solution Xtra delivered at 430mm deep		○
BNF0183		Long cushion depth: Libra delivered at 455mm deep, Solution Xtra delivered at 455mm deep		○
<b>MATRX LIBRA/FLOTECH SOLUTION XTRA CUSHION ACCESSORIES</b>				
BNF0186		Tapered Wedge Pair		€28
BNF0187		Abductor		€58
BNF0188		Lateral Side Wedge		€46
BNF0189		Anterior Front Wedge		€43
<b>SEAT UPHOLSTERY*</b>				
BNF0191		Black Plush TR35		○
BNF0201		Multi-stretch Black PU coated fabric		◆
BNF0202		Black Spacer TR36		○
BNF1111		No upholstery, select with options 0167 Vicair Adjuster O2 or 0180 Without seat cushion		○
<b>BACKREST WIDTH EXTENSION*</b>				
BNF0205		+ 50 mm Backrest width pre-adjustment, set at +50mm but can be adjusted if required	(NC 0221, 0223, 0225)	○
BNF0207		Without Backrest width pre-adjustment, set at 0mm but can be adjusted if required		○

<b>BACKREST*</b>			
BNF0210		Flex 3 Backrest, Fixed Plate (Width Adjustable)	(NC 0233, 0234, 0235, 0240, 0241, 0242, 0243, 0244, 0245, 0246, 0247, 0343) ◆
BNF0219		New!! Flex 3+ Backrest features tension adjustment, width adjustment, and fore/aft positioning of the vertical backpost tubes to accommodate a wider range of clinical requirements.	(NC 0233, 0235, 0234, 0241, 0242, 0243, 0244, 0245, 0246, 0247, 0277, 0343) €289
BNF0220		Flex 3 Backrest features tension adjustment and width adjustment of the vertical backpost tubes.	(NC 0233, 0235, 0234, 0241, 0242, 0243, 0244, 0245, 0246, 0247, 0277, 0343) €126
BNF0221		Matrx Elite E2 (includes hardware mounted to frame, backrest supplied in box)	(NC 0240, 0250, 0260, 0270, 0305, 0311, 0321, 0322, 0800, 0810, 0820, 0830, 0880, 0890, 0900, 0920) €668
BNF0223		Matrx Elite E2 Deep (includes hardware mounted to frame, backrest supplied in box)	(NC 0240, 0250, 0260, 0270, 0305, 0311, 0321, 0322, 0800, 0810, 0820, 0830, 0880, 0890, 0900, 0920) €739
BNF0225		Matrx Posture Back (includes hardware mounted to frame, backrest supplied in box)	(NC 0235, 0240, 0241, 0242, 0243, 0244, 0245, 0246, 0247, 0250, 0260, 0270, 0305, 0311, 0321, 0322, 0335, 0800, 0810, 0820, 0830, 0880, 0890, 0900, 0920) €1,038
BNF0230		Backrest Tubes Without Backrest	(NC 0233, 0234, 0235, 0240, 0241, 0242, 0243, 0244, 0245, 0246, 0247, 0250, 0260, 0270, 0305, 0311, 0321, 0322, 0345, 0740, 0800, 0810, 0820, 0830, 0880, 0890, 0900, 0920) ○
<b>MATRIX BACKREST HEIGHT</b>			
BNF0233		Matrx Backrest Height 410 mm / 16"	○
BNF0235		Matrx Backrest Height 460 mm / 18"	○
BNF0234		Matrx Backrest Height 510 mm / 20"	○
<b>MATRIX BACKREST WIDTH</b>			
BNF0241		Matrx Backrest Width 360mm / 14" (Chest width 265mm / 10 1/2")	○
BNF0242		Matrx Backrest Width 380mm / 15" (Chest width 290mm / 11 1/2")	○
BNF0243		Matrx Backrest Width 410mm / 16" (Chest width 315mm / 12 1/2")	○
BNF0244		Matrx Backrest Width 430mm / 17" (Chest width 340mm / 13")	○
BNF0245		Matrx Backrest Width 460mm / 18" (Chest width 370mm / 14 1/2")	○
BNF0246		Matrx Backrest Width 480mm / 19" (Chest width 390mm / 15")	○
BNF0247		Matrx Backrest Width 510mm / 20" (Chest width 420mm / 16 1/2")	○
<b>BACKREST CUSHION*</b>			
BNF0240		Cover for tension adjustable upholstery	(NC 2222) ○
BNF0250		Cushion, Laguna (Lateral Support)	(NC 2222) ○
BNF0260		Cushion, Mistral 2 (Waist Support)	(NC 0322, 2222) €179
BNF0270		Cushion, Passad 2 (Shoulder Support)	(NC 2222) ◆
BNF0276		Vicair Multifunctional O2 pressure relieving backrest cushion, supplied with Black mesh cover (Additional 5 day lead time)	(NC 0305, 0311, 0321, 0322) €1,000
BNF0277		Vicair Multifunctional O2 Flex pressure relieving backrest cushion, supplied with Black mesh cover (Additional 5 day lead time)	(NC 0305, 0311, 0321, 0322) €1,000
BNF0300		Without Backrest Cover/ Cushion	(NC 0305, 0311, 0321, 0322) ○
<b>WIDER BACKREST COVER</b>			
BNF0305		Extra large backrest cover or cushion. One size wider backrest cushion or tension cover fitted if selected	○
<b>BACKREST UPHOLSTERY*</b>			
BNF0311		Black Plush TR35	○
BNF0321		Multi-stretch Black PU coated fabric	◆
BNF0322		Black Spacer TR36	○
BNF2222		No upholstery, select with Matrx, and Without backrest options	○

<b>PUSH HANDLES*</b>			
BNF0335		Push Handles connected	◆
BNF0343		Push Handles connected,with Rea head/neckrest attachment	€96
BNF0345		Push Bar, angle adjustable	€123
<b>SEAT TILT AND BACKREST REGULATION*</b>			
BNF0350		Assistant operated gas piston seat tilt and backrest recline (NC 5555)	◆
BNF0360		Negative Seat Angle -1 To +23° (For attendant operated gas piston in negative tilt) (NC 1263)	€958
<b>TILT LOCK*</b>			
BNF1263		Tilt Lock (locking of pre-set angle)	◆
BNF5555		No Tilt Lock, Select with 0360 Electric operated seat and backrest	○
<b>CENTRAL LEGREST</b>			
BNF0370		Central Legrest, includes calfpad and footrest. Also select options 0420, 0470, 3333, 4444 (NC 0376, 0378, 0380, 0390, 0410, 0430, 0440, 0460, 0475, 0476, 0477, 0720, 0930, 0940, 0941, 0942, 0952, 0954, 0956, 1290)	€377
<b>LEGREST HANGER ATTACHMENT*</b>			
BNF0375		Standard position for legrest mounting	◆
SWO0189		Standard position with 0-30° Legrest Abduction function (NC 0370, 0475, 0476, 0477)	€122
BNF0376		High position 50mm higher for leg length compensation	○
BNF0378		90° Hanger attachment to achieve a 90° legrest position on the Elevating Legrest (NC 0390, 0410, 0440, 0460)	○
<b>LEGREST HANGER LH*</b>			
BNF0380		Elevating Legrest, angle adjustable 0 to 75° swing in/swing out with knee pad, calf support with locking clip LH (NC 0475, 0490, 0942, 0945, 6666)	◆
BNF0390		Legrest Hanger fixed 80 deg, swing in/swing out LH (NC 0928, 0929, 0930)	○
BNF0410		Legrest Amputee LH (NC 0475, 0476, 0477, 0490, 0520, 0928, 0929, 0930, 0941, 0942, 0945, 0952, 0954, 0956)	€171
BNF0420		Without Legrest, LH (NC 0475, 0476, 0477, 0490, 0520, 0928, 0929, 0930, 0941, 0942, 0945, 0952, 0954, 0956)	○
<b>LEGREST HANGER RH*</b>			
BNF0430		Elevating Legrest, angle adjustable 0 to 75°swing In/swing out with knee pad and calf support with locking clip RH (NC 0475, 0550, 0942, 0945, 6666)	◆
BNF0440		Legrest Hanger fixed 80 deg, swing in/swing out RH (NC 0931, 0932, 0940)	○
BNF0460		Legrest Amputee RH (NC 0475, 0476, 0477, 0550, 0580, 0931, 0932, 0940, 0941, 0942, 0945, 0952, 0954, 0956)	€171
BNF0470		Without Legrest, RH (NC 0475, 0476, 0477, 0550, 0580, 0931, 0932, 0940, 0941, 0942, 0945, 0952, 0954, 0956)	○
<b>FOOTREST LH*</b>			
BNF0475		Footrest one piece, depth and angle adjustable short tube (NC 0477, 0490, 0520, 0550, 0580, 0930, 0940, 0952, 0954, 0956)	€145
BNF0476		Footrest one piece, depth and angle adjustable long tube (NC 0477, 0490, 0520, 0550, 0580, 0930, 0940, 0952, 0954, 0956)	€145
BNF0477		Footbox, angle adjustable, padded with short tube (NC 0490, 0520, 0550, 0580, 0930, 0940, 0941, 0942, 0945, 0952, 0954, 0956)	€687
BNF0490		Footrest lower part, depth and angle adjustable with short tube and heel strap, LH. Not compatible with Elevating legrest.	○
BNF0520		Footrest lower part, depth and angle adjustable with long tube and heel strap, LH. All knee to heel ranges can be achieved when used with Elevating legrests.	○
BNF3333		Without footrest, only select with option 0370 Central legrest, 0410 Amputee legrest and 0420 Without legrest	○

FOOTREST RH*			
BNF0550		Footrest lower part, depth and angle adjustable with short tube and heel strap, RH. Not compatible with Elevating legrest.	○
BNF0580		Footrest lower part, depth and angle adjustable with long tube and heel strap, RH. All knee to heel ranges can be achieved when used with Elevating legrests.	○
BNF4444		Without footrest, only select with options 0370 Central legrest, 0460 Amputee legrest, 0470 Without legrest, 0475 and 0476 One piece footrests or 0477 Footbox	○
ARMREST LH*			
BNF0600		Armrest Ducktail Pad LH, height adjustable with autolock (NC 1232)	◆
BNF0607		Armrest Reinforced Wide Pad, LH, height adjustable with autolock (NC 1230, 1231, 1232, 1250, 1260)	€112
BNF0610		Armrest Padded LH, height adjustable with autolock (NC 0680, 0695, 1230, 1231, 1232, 1250, 1260, 9999)	€79
BNF0615		Armrest Hemiplegic LH, height adjustable with autolock (NC 0695, 1230, 1231, 1232, 1250, 1260)	€222
BNF0620		Without Armrest LH Select when choosing BNA1232 Half Tray LH, Narrow Armrest is included (NC 0695, 1230, 1231, 1250, 1260)	○
ARMREST RH*			
BNF0630		Armrest Ducktail Pad RH, height adjustable with autolock (NC 1233)	◆
BNF0637		Armrest Reinforced Wide Pad, RH, height adjustable with autolock (NC 1230, 1231, 1233, 1250, 1260)	€112
BNF0640		Armrest Padded RH, height adjustable with autolock. (NC 0680, 0695, 1230, 1231, 1233, 1250, 1260, 9999)	€79
BNF0645		Armrest Hemiplegic, RH, height adjustable with autolock (NC 0695, 1230, 1231, 1233, 1240, 1245, 1250, 1260)	€222
BNF0650		Without Armrest RH Select when choosing BNF1233 Half Tray RH, Narrow Armrest is included (NC 0695, 1230, 1231, 1240, 1245, 1250, 1260)	○
ARMREST OPTIONS			
BNF0680		Side Support Cover, PU fabric cover for the sideguards (Removable with tools) (NC 0695, 9999)	€63
BNF9999		Side Support Cover with foam cushion. 2cm thick cushion supplied complete with PU fabric cover. Reduces seat width by 2cm for a comfortable fit. Consists of options 0680 and 0690 (Removable with tools) (NC 0695)	€63
BNF0695		Pad for sideguard, removable with hook and loop fastening	€63
SEAT OPTIONS			
BNF0710		Incontinence Cover	€43
BNF0720		Pommel Abduction Cushion (NC 0722)	€145
BNF0722		2 cm Increased seat height plate	€89
BNF0725		Tilt Scale for seat tilt	€53
HEADREST OPTIONS*			
BNF0730		Headrest, Sideway Adjustable Fixed hardware included (NC 0742, 0744, 0746, 0750, 0770, 0782)	◆
BNF0750		Headrest with cheek support, sideway adjustable Fixed hardware included (NC 0740, 0742, 0744, 0746, 0770, 0782)	€135
BNF0731		Matrx Elan Standard Pad (254 mm / 10") (NC 0740, 0770, 0782, 0790, 8888)	€189
BNF0732		Matrx Elan Standard Pad (355 mm / 14") (NC 0740, 0770, 0782, 0790, 8888)	€189
BNF0733		Matrx Elan Occipital Pad (229 mm / 9") (NC 0740, 0770, 0782, 0790, 8888)	€189
BNF0734		Matrx Elan Occipital Pad (305 mm / 12") (NC 0740, 0770, 0782, 0790, 8888)	€189
BNF0735		Matrx Elan 4-Point Pad (229W x 177H mm / 9"W x 7"H) (NC 0740, 0770, 0782, 0790, 8888)	€357
BNF0736		Matrx Elan 4-Point Pad (280W x 254H mm / 11"W x 10"H) (NC 0740, 0770, 0782, 0790, 8888)	€357
BNF9991		Without headrest (NC 0742, 0744, 0746)	○

HEADREST HARDWARE*				
BNF0742		Matrx Elan Hardware	(NC 0740, 0770, 0790, 1206, 1207, 1300)	€285
BNF0744		Matrx Elan Extended Hardware (Extra link)	(NC 0740, 0770, 0790, 1206, 1207, 1300)	€340
BNF0746		Matrx Elan LOXX Hardware	(NC 0740, 0770, 0790, 1206, 1207, 1300)	€285
BNF8888		Without Elan headrest hardware (Only select with option 0730 & 0750 Headrests)		○
OTHER HEADREST OPTIONS				
BNF0740		Cover for headrest		€39
BNF0770		Neckrest		○
BNF9992		Without neckrest		○
BNF0782		Alternative headrest/neckrest mounting bracket		€24
BNF0790		Handles Clinic Kit for Head-/Neckrest (Levers which replace screws for easy adjustment)		€61
TRUNK SUPPORT LH				
BNF0800		Trunk Support Multi Adjustable with standard arm LH	(NC 0900)	○
BNF0820		Trunk Support Multi Adjustable with extended height arm LH	(NC 0900)	○
BNF0880		NEW!! Trunk Support Swingaway LH, multi adjustable	(NC 0900)	€128
TRUNK SUPPORT RH				
BNF0810		Trunk Support Multi Adjustable with standard arm RH	(NC 0900)	○
BNF0830		Trunk Support Multi Adjustable with extended height arm RH	(NC 0900)	○
BNF0890		NEW!! Trunk Support Swingaway RH, multi adjustable	(NC 0900)	€128
BACKREST OPTIONS				
BNF0900		Trunk Support Attachment to mount laterals later. If trunk supports are chosen, the attachment bracket is already included		○
BNF0915		Tilt Scale for backrest recline		€53
BNF0920		Wedges (One Pair)		€14
ELEVATING LEGREST OPTIONS*				
BNF0928		Calf pad cover, for LH elevating legrest		○
BNF0929		Knee pad for LH elevating legrest		◆
BNF0931		Calf pad cover, for RH elevating legrest		○
BNF0932		Knee pad for RH elevating legrest		◆
BNF6666		Without Elevating Legrest Options (Select with Fixed legrests, Amputee legrests and Without legrests only)		○
LEGREST OPTIONS				
BNF0930		Footplate Extender LH	(NC 0941)	€40
BNF0940		Footplate Extender RH	(NC 0941)	€40
BNF0941		Cover Footplate, Black Plush Fabric		€66
BNF0942		Calf Strap Padded		€88
BNF0945		Calf Strap (for One piece footrest)		○
BODYPOINT ANKLE HUGGERS				
BNF0952		Bodypoint Ankle Huggers Large (Supplied with fixing kit, unmounted)		€189
BNF0954		Bodypoint Ankle Huggers Medium (Supplied with fixing kit, unmounted)		€189
BNF0956		Bodypoint Ankle Huggers Small (Supplied with fixing kit, unmounted)		€189

<b>CHASSIS COLOUR*</b>			
BNF0944	<input type="checkbox"/>	Silver Pearl	◆
BNF0946	<input type="checkbox"/>	Azurite (Dark Blue)	○
BNF0947	<input type="checkbox"/>	Sand (Beige)	○
BNF0948	<input type="checkbox"/>	Electric Green	○
BNF0949	<input type="checkbox"/>	Happy Red	○
<b>CHASSIS SAFETY*</b>			
BNF0950	<input type="checkbox"/>	Chassis extender, castors extend 50mm forward for added stability	◆
<b>CASTOR ATTACHMENT*</b>			
BNF0960	<input type="checkbox"/>	Fixed castor attachment	◆
<b>FRONT CASTORS*</b>			
BNF1015	<input type="checkbox"/>	200x32 mm (8" x 1 1/4") Black Solid Tyre	€25
BNF1025	<input type="checkbox"/>	200x45 mm (8" X 1 3/4") Black Solid tyre	◆
<b>REAR WHEEL SIZE*</b>			
BNF1062	<input type="checkbox"/>	24" / 600 mm Self Propel rear wheel with QR attachment	◆
<b>REAR WHEEL TYPE*</b>			
BNF1074	<input type="checkbox"/>	Standard self propel with colourless hub, clear anodized rim and Stainless steel spokes	◆
BNF1076	<input type="checkbox"/>	Hub Brake 22" or 24" Self Propel with colourless hub, clear anodized rim and Stainless steel spokes	○
<b>REAR WHEEL TYRE*</b>			
BNF1082	<input type="checkbox"/>	Pneumatic Profile Black 1 3/8" tyres	○
BNF1092	<input type="checkbox"/>	Solid Black tyres	◆
<b>HANDRIM*</b>			
BNF1100	<input type="checkbox"/>	Aluminium handrim	◆
BNF1120	<input type="checkbox"/>	Plastic coated handrim	€61
<b>SPOKEGUARD</b>			
BNF1160	<input type="checkbox"/>	Frost spokeguard	€205
<b>BRAKE*</b>			
BNF1170	<input type="checkbox"/>	User manoeuvred tyre brakes, includes extended handles	◆
BNF1180	<input type="checkbox"/>	One Arm Brake LH	€154
BNF1185	<input type="checkbox"/>	One Arm Brake RH	€154
<b>PELVIC BELT*</b>			
BNF1200	<input type="checkbox"/>	Pelvic belt	◆
BNF1201	<input type="checkbox"/>	Bodypoint 2-Point Centre-pull Push Button Medium (Mounted to chair)	€171
BNF1202	<input type="checkbox"/>	Bodypoint 2-Point Centre-pull Push Button Large (Mounted to chair)	€171
BNF1203	<input type="checkbox"/>	Bodypoint 4-Point Centre-pull Push Button Medium (Mounted to chair)	€161
BNF1204	<input type="checkbox"/>	Bodypoint 4-Point Centre-pull Push Button Large (Mounted to chair)	€161
<b>BODYPOINT HARNESS</b>			
BNF1206	<input type="checkbox"/>	Bodypoint Pivotfit Harness Medium (Supplied with fixing kit, unmounted)	€171
BNF1207	<input type="checkbox"/>	Bodypoint Pivotfit Harness Large (Supplied with fixing kit, unmounted)	€171
<b>TRANSPORTATION BRACKETS*</b>			
BNF1225	<input type="checkbox"/>	Transportation Brackets, mounted to the front of the wheelchair	€16

OTHER OPTIONS		
BNF1230		Fixed Table Tray without elbow pads includes mounting kit (NC 1237, 1260) €173
BNF1231		Fixed Table Tray with elbow Pads includes mounting kit (NC 1237, 1260) €173
BNF1232		Half table tray, clear plastic, flip up, depth adjustable, left side (armrest with straight narrow pad included) (NC 1235, 1260) €138
BNF1233		Half table tray, clear, plastic, flip up, depth adjustable, right side (armrest with straight narrow pad included) (NC 1235, 1260) €138
BNF1240		Swing away table without elbow pads (NC 1237, 1250, 1260) €242
BNF1245		Swing away table with elbow pads (NC 1237, 1250, 1260) €242
BNF1235		Table tray padded cushion (To be ordered in combination with the fixed & swingaway table tray) €63
BNF1237		Half table tray padded cushion €63
BNF1250		Locking device for table tray €17
BNF1260		Mounting kit for table (if ordered separately) €10
BNF1290		Handles Clinic Kit for backrest, armrest, legrest and seat frame (Levers which replace screws for easy adjustment) €67
BNF1300		Harness attachment , must be selected in combination with Head/Neckrest €98
BNF1310		Attachment for positioning belt €39

<b>TOTAL PRICE</b>	
--------------------	--

**Non Conforming Options:**

BNF0182	Short cushion depth: Libra delivered at 405mm deep, Solution Xtra delivered at 430mm deep
BNF0183	Long cushion depth: Libra delivered at 455mm deep, Solution Xtra delivered at 455mm deep
BNF0186	Tapered Wedge Pair
BNF0187	Abductor
BNF0188	Lateral Side Wedge
BNF0189	Anterior Front Wedge
BNF0191	Black Plush TR35
BNF0201	Multi-stretch Black PU coated fabric
BNF0202	Black Spacer TR36
BNF1111	No upholstery, select with options 0167 Vicair Adjuster O2 or 0180 Without seat cushion
BNF0221	Matrx Elite E2 (includes hardware mounted to frame, backrest supplied in box)
BNF0223	Matrx Elite E2 Deep (includes hardware mounted to frame, backrest supplied in box)
BNF0225	Matrx Posture Back (includes hardware mounted to frame, backrest supplied in box)
BNF0233	Matrx Backrest Height 410 mm / 16"
BNF0235	Matrx Backrest Height 460 mm / 18"
BNF0234	Matrx Backrest Height 510 mm / 20"
BNF0241	Matrx Backrest Width 360mm / 14" (Chest width 265mm / 10 1/2")
BNF0242	Matrx Backrest Width 380mm / 15" (Chest width 290mm / 11 1/2")
BNF0243	Matrx Backrest Width 410mm / 16" (Chest width 315mm / 12 1/2")
BNF0244	Matrx Backrest Width 430mm / 17" (Chest width 340mm / 13")
BNF0245	Matrx Backrest Width 460mm / 18" (Chest width 370mm / 14 1/2")
BNF0246	Matrx Backrest Width 480mm / 19" (Chest width 390mm / 15")
BNF0247	Matrx Backrest Width 510mm / 20" (Chest width 420mm / 16 1/2")
BNF0240	Cover for tension adjustable upholstery
BNF0250	Cushion, Laguna (Lateral Support)
BNF0260	Cushion, Mistral 2 (Waist Support)
BNF0270	Cushion, Passad 2 (Shoulder Support)
BNF0277	Vicair Multifunctional O2 Flex pressure relieving backrest cushion, supplied with Black mesh cover (Additional 5 day lead time)
BNF0305	Extra large backrest cover or cushion. One size wider backrest cushion or tension cover fitted if selected
BNF0311	Black Plush TR35
BNF0321	Multi-stretch Black PU coated fabric
BNF0322	Black Spacer TR36
BNF2222	No upholstery, select with Matrx, and Without backrest options
BNF0335	Push Handles connected
BNF0343	Push Handles connected,with Rea head/neckrest attachment
BNF0345	Push Bar, angle adjustable
BNF0360	Negative Seat Angle -1 To +23° (For attendant operated gas piston in negative tilt)
BNF1263	Tilt Lock (locking of pre-set angle)
BNF5555	No Tilt Lock, Select with 0360 Electric operated seat and backrest
BNF0370	Central Legrest, includes calfpad and footrest. Also select options 0420, 0470, 3333, 4444
SWO0189	Standard position with 0-30° Legrest Abduction function
BNF0376	High position 50mm higher for leg length compensation
BNF0378	90° Hanger attachment to achieve a 90° legrest position on the Elevating Legrest

BNF0380	Elevating Legrest, angle adjustable 0 to 75° swing in/swing out with knee pad, calf support with locking clip LH
BNF0390	Legrest Hanger fixed 80 deg, swing in/swing out LH
BNF0410	Legrest Amputee LH
BNF0430	Elevating Legrest, angle adjustable 0 to 75° swing In/swing out with knee pad and calf support with locking clip RH
BNF0440	Legrest Hanger fixed 80 deg, swing in/swing out RH
BNF0460	Legrest Amputee RH
BNF0475	Footrest one piece, depth and angle adjustable short tube
BNF0476	Footrest one piece, depth and angle adjustable long tube
BNF0477	Footbox, angle adjustable, padded with short tube
BNF0490	Footrest lower part, depth and angle adjustable with short tube and heel strap, LH. Not compatible with Elevating legrest.
BNF0520	Footrest lower part, depth and angle adjustable with long tube and heel strap, LH. All knee to heel ranges can be achieved when used with Elevating legrests.
BNF0550	Footrest lower part, depth and angle adjustable with short tube and heel strap, RH. Not compatible with Elevating legrest.
BNF0580	Footrest lower part, depth and angle adjustable with long tube and heel strap, RH. All knee to heel ranges can be achieved when used with Elevating legrests.
BNF0680	Side Support Cover, PU fabric cover for the sideguards (Removable with tools)
BNF9999	Side Support Cover with foam cushion. 2cm thick cushion supplied complete with PU fabric cover. Reduces seat width by 2cm for a comfortable fit. Consists of options 0680 and 0690 (Removable with tools)
BNF0695	Pad for sideguard, removable with hook and loop fastening
BNF0720	Pommel Abduction Cushion
BNF0722	2 cm Increased seat height plate
BNF0750	Headrest with cheek support, sideways adjustable Fixed hardware included
BNF0742	Matrx Elan Hardware
BNF0744	Matrx Elan Extended Hardware (Extra link)
BNF0746	Matrx Elan LOXX Hardware
BNF8888	Without Elan headrest hardware (Only select with option 0730 & 0750 Headrests)
BNF0740	Cover for headrest
BNF0770	Neckrest
BNF0782	Alternative headrest/neckrest mounting bracket
BNF0790	Handles Clinic Kit for Head-/Neckrest (Levers which replace screws for easy adjustment)
BNF0800	Trunk Support Multi Adjustable with standard arm LH
BNF0820	Trunk Support Multi Adjustable with extended height arm LH
BNF0880	NEW!! Trunk Support Swingaway LH, multi adjustable
BNF0810	Trunk Support Multi Adjustable with standard arm RH
BNF0830	Trunk Support Multi Adjustable with extended height arm RH
BNF0890	NEW!! Trunk Support Swingaway RH, multi adjustable
BNF0900	Trunk Support Attachment to mount laterals later. If trunk supports are chosen, the attachment bracket is already included
BNF0920	Wedges (One Pair)
BNF0928	Calf pad cover, for LH elevating legrest
BNF0929	Knee pad for LH elevating legrest
BNF0931	Calf pad cover, for RH elevating legrest
BNF0932	Knee pad for RH elevating legrest
BNF6666	Without Elevating Legrest Options (Select with Fixed legrests, Amputee legrests and Without legrests only)
BNF0930	Footplate Extender LH
BNF0940	Footplate Extender RH
BNF0941	Cover Footplate, Black Plush Fabric

BNF0942	Calf Strap Padded
BNF0945	Calf Strap (for One piece footrest)
BNF0952	Bodypoint Ankle Huggers Large (Supplied with fixing kit, unmounted)
BNF0954	Bodypoint Ankle Huggers Medium (Supplied with fixing kit, unmounted)
BNF0956	Bodypoint Ankle Huggers Small (Supplied with fixing kit, unmounted)
BNF1074	Standard self propel with colourless hub, clear anodized rim and Stainless steel spokes
BNF1076	Hub Brake 22" or 24" Self Propel with colourless hub, clear anodized rim and Stainless steel spokes
BNF1206	Bodypoint Pivotfit Harness Medium (Supplied with fixing kit, unmounted)
BNF1207	Bodypoint Pivotfit Harness Large (Supplied with fixing kit, unmounted)
BNF1230	Fixed Table Tray without elbow pads includes mounting kit
BNF1231	Fixed Table Tray with elbow Pads includes mounting kit
BNF1232	Half table tray, clear plastic, flip up, depth adjustable, left side (armrest with straight narrow pad included)
BNF1233	Half table tray, clear, plastic, flip up, depth adjustable, right side (armrest with straight narrow pad included)
BNF1240	Swing away table without elbow pads
BNF1245	Swing away table with elbow pads
BNF1235	Table tray padded cushion (To be ordered in combination with the fixed & swingaway table tray)
BNF1237	Half table tray padded cushion
BNF1250	Locking device for table tray
BNF1260	Mounting kit for table (if ordered separately)
BNF1290	Handles Clinic Kit for backrest, armrest, legrest and seat frame (Levers which replace screws for easy adjustment)
BNF1300	Harness attachment , must be selected in combination with Head/Neckrest



For a more complex customer requirement,  
please email: [orders@invacare.com.au](mailto:orders@invacare.com.au)

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):