



**Yes, you can:**

Tel: 01656 776222  
Email: ordersuk@invacare.com  
Online Spares available at: www.invacare.co.uk

07 JANUARY 2026  
Customer Ref:

|                          |                    |
|--------------------------|--------------------|
| <b>Customer Ref*</b>     | <b>Account No*</b> |
| <b>Contact No*</b>       | <b>Order No*</b>   |
| <b>Delivery Address*</b> | <b>Order Qty*</b>  |
|                          | <b>Order date*</b> |

### Rea Azalea Ireland Stock Spec 3

#### Standard Features

- Max user weight: 135kg (21.5 stone)
- Crash tested ISO 7176-Part 19
- Transit Propulsion 16" Black solid tyre, attendant hub brakes
- Small, medium or large chassis, +50mm seat width adjustment
- Seat height 450mm/18"
- Depth adjustable seat 430mm-500mm/17"-20"
- Flex 3 Tension adjustable backrest
- Black Dartex Passad 2 backrest cushion
- Black Dartex Flo-shape seat cushion
- Connected push handles
- Assistant operated gas piston seat and backrest recline with tilt lock
- Height adjustable wide pad armrests
- Removable, angle adjustable legrests with calf pads and angle adjustable footplates
- Sideways adjustable headrest
- 200 X 45mm Black medium hard front castors
- User manoeuvred tyre brakes
- Pelvic belt
- Anti-tippers
- Pearl Silver grey frame



Please make your selection with a cross in the box    ☐ Free of charge    ☒ standard feature    POA special request

Customer weight =                      KG

| Code                      | Description           |        |
|---------------------------|-----------------------|--------|
| <b>Rea Azalea Transit</b> |                       |        |
| 1542419                   | W 390mm x D 430-500mm | €3,171 |
| 1542420                   | W 440mm x D 430-500mm | €3,171 |
| 1542421                   | W 490mm x D 430-500mm | €3,171 |

**TOTAL PRICE**

**Non Conforming Options:**



For a more complex customer requirement,  
please email: [orders@invacare.com.au](mailto:orders@invacare.com.au)

**Special Instructions**

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):