



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

07 JANUARY 2026
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Rea Azalea Ireland Stock Spec 4

Standard Features

- Max user weight: 135kg (21.5 stone)
- Crash tested ISO 7176-Part 19
- Transit Propulsion 16" Black solid tyre, attendant hub brakes
- Small, medium or large chassis, +100mm seat width adjustment
- Depth adjustable seat 430mm-500mm/17"-20"
- Seat height 450mm/18"
- Flex 3 Backrest Plate
- Black Dartex Passad 2 backrest cushion
- Without seat cushion
- Connected push handles
- Assistant operated gas piston seat and backrest recline with tilt lock
- Height adjustable wide pad armrests
- Removable, angle adjustable legrests with calf pads and angle adjustable footplates
- Sideways adjustable headrest
- Trunk supports, multi-adjustable extended height
- 200 X 45mm Black medium hard front castors
- User manoeuvred tyre brakes
- Pelvic belt
- Anti-tippers
- Pearl Silver grey frame



Please make your selection with a cross in the box ☐ Free of charge ☒ standard feature POA special request

Customer weight = KG

Code	Description	
Rea Azalea Transit		
1542449	W 390mm x D 430-500mm	€2,941
1542450	W 440mm x D 430-500mm	€2,941
1542451	W 490mm x D 430-500mm	€2,941

TOTAL PRICE	
--------------------	--

Non Conforming Options:



For a more complex customer requirement,
please email: orders@invacare.com.au

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):