



**Yes, you can:**

Tel: 01656 776222  
 Email: ordersuk@invacare.com  
 Online Spares available at: www.invacare.co.uk

# Manual Invacare Ireland Prescription Form

07 JANUARY 2026  
 Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## Rea Azalea Ireland Stock Spec 5

### Standard Features

- Max user weight: 135kg (21.5 stone)
- Crash tested ISO 7176-Part 19
- Transit Propulsion 16" Black solid tyre, attendant hub brakes
- Small, medium or large chassis, without seat width adjustment
- Depth adjustable seat 430mm-500mm/17"-20"
- Seat height 450mm/18"
- Flex 3 Backrest Plate
- Black Dartex Passad 2 backrest cushion
- Black Dartex Flo-shape seat cushion
- Connected push handles
- Assistant operated gas piston seat and backrest recline with tilt lock
- Height adjustable wide pad armrests
- Removable, angle adjustable legrests with calf pads and angle adjustable footplates
- Sideways adjustable headrest
- Trunk supports, multi-adjustable extended height
- 200 X 45mm Black medium hard front castors
- User manoeuvred tyre brakes
- Pelvic belt
- Anti-tippers
- Pearl Silver grey frame



Please make your selection with a cross in the box    ○ Free of charge    ◆ standard feature    POA special request

Customer weight =                      KG

Code	Description	
<b>Rea Azalea Transit</b>		
1542580	W 390mm x D 430-500mm	€3,171
1542581	W 440mm x D 430-500mm	€3,171
1542582	W 490mm x D 430-500mm	€3,171

<b>TOTAL PRICE</b>	
--------------------	--

**Non Conforming Options:**



For a more complex customer requirement,  
please email: [orders@invacare.com.au](mailto:orders@invacare.com.au)

**Special Instructions**

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):