



Yes, you can:

Tel: 01656 776222
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

07 JANUARY 2026
 Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Rea Azalea Transit

Standard Features

- Max User Weight 135KG/21 Stone
- Crash Tested ISO 7176-Part 19
- 16" Transit Propulsion with Attendant Hub Brake and Solid Black Tyres
- Choice of Seat Widths 390mm, 440mm or 490mm (+50mm or +100mm Width Adj)
- Choice of 2 Seat Heights 400mm or 450mm
- Depth Adjustable Seat 430mm-500mm
- Without Seat Cushion, Supplied with Seat Board
- Flex 3 Backrest Plate, Width and Height Adjustable
- Passad 2 Backrest Cushion with Multi-stretch Black PU coated fabric Cover and Push Handles Connected
- Assistant Operated Gas Piston Seat Tilt +1° - +24°
- Backrest Recline 30° with Tilt Lock
- Angle Adjustable Legrests with Calf Support, Height and Depth Adjustable Footplates
- Height Adjustable Wide Pad Armrests with Autolock
- Pearl Silver Frame Colour, Solid Black Front Castors
- User Manoeuvred Tyre Brakes
- Pelvic Belt, Adjustable both sides
- Anti Tippers
- Front Transportation Brackets

Standard Base price: €2,948



Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code	Description	
BASE PRODUCT*		
BNE9002	Rea Azalea Transit	€2,948
BNE8000	SPECIAL REQUEST	POA
PROPULSION*		
BNE0020	Transit 16" wheels with hub brakes	◆
SEAT WIDTH*		
BNE0040	390 mm/15" (Small chassis) (NC 0244, 0245, 0246, 0247)	○
BNE0050	440 mm/17" (Medium chassis) (NC 0241, 0242, 0246, 0247)	○
BNE0060	490 mm/19" (Large chassis) (NC 0206, 0241, 0242, 0243, 0244)	○

SEAT WIDTH EXTENSION*			
BNE0065		+ 50 mm pre-adjustment on armrests and legrests	○
BNE0075		+ 100 mm pre-adjustment on armrests and legrests (NC 0110)	○
BNE0085		Without seat width extension, set at 0mm adjustment but can be adjusted wider if required (NC 0206)	○
SEAT DEPTH*			
BNE0100		430-500 mm / 17" - 20"	◆
SEAT HEIGHT*			
BNE0110		400 mm / 16" Seat Height (NC 0370, 1015, 1025)	○
BNE0120		450 mm / 18" Seat Height	○
SEAT ANGLE			
BNE0127		Negative Seat Angle -1 To +23° (For attendant operated gas piston in negative tilt) (NC 0360)	○
SEAT TYPE*			
BNE0130		Depth Adjustable	◆
BNE0140		Sliding Seat (NC 0720)	€192
SEAT CUSHION*			
BNE0165		Cushion Flosshape seat upholstery (NC 1111, 0182, 0183, 0186, 0187, 0188, 0189)	€72
SWO0190		Cushion Flosshape seat upholstery Shorter depth 450mm deep. Only use with seat depth in shortest position (NC 1111, 0182, 0183, 0186, 0187, 0188, 0189)	€72
BNE0167		Cushion Vicair Adjuster O2 10cm, complete with Black Comfair cover (NC 0182, 0183, 0186, 0187, 0188, 0189, 0201, 0202)	€696
BNE0168		Cushion MatrX Libra, complete with Black cover (NC 0201, 0202)	€645
BNE0169		Cushion Flo-Tech Solution Xtra, complete with Black cover (NC 0201, 0202)	€566
BNE0179		Cushion Vicair Vector O2 100mm with 2 top covers (Additional 15 day leadtime) (NC 0182, 0183, 0186, 0187, 0188, 0189, 0201, 0202)	€800
BNE0180		Without Seat Cushion, Seat Plate only (NC 0182, 0183, 0201, 0202)	○
SEAT DEPTH OF MATRX LIBRA AND FLO-TECH SOLUTION XTRA CUSHIONS (OPTIMISED CUSHION DEPTH TO SUIT SEAT DEPTH SETTING)			
BNE0182		Short cushion depth: Libra delivered at 405mm deep, Solution Xtra delivered at 430mm deep	○
BNE0183		Long cushion depth: Libra delivered at 455mm deep, Solution Xtra delivered at 455mm deep	○
MATRX LIBRA/FLOTECH SOLUTION XTRA CUSHION ACCESSORIES			
BNE0186		Tapered Wedge Pair	€28
BNE0187		Abductor	€58
BNE0188		Lateral Side Wedge	€46
BNE0189		Anterior Front Wedge	€43
SEAT UPHOLSTERY*			
BNE0201		Multi-stretch Black PU coated fabric	◆
BNE0202		Black Spacer TR36	○
BNE1111		No upholstery, select with options 0167 Vicair Adjuster O2 or 0180 Without seat cushion	○
BACKREST WIDTH EXTENSION*			
BNE0205		+ 50 mm Backrest width pre-adjustment, set at +50mm but can be adjusted if required (NC 0221, 0223, 0225)	○
BNE0206		+ 100 mm Backrest width pre-adjustment, set at +100mm but can be adjusted if required (NC 0220, 0221, 0223, 0225)	○
BNE0207		Without Backrest width pre-adjustment, set at 0mm but can be adjusted if required	○

BACKREST*			
BNE0210		Flex 3 Backrest, Fixed Plate Width Adjustable. Headrest attachment built into backrest plate	(NC 0233, 0234, 0235, 0240, 0241, 0242, 0243, 0244, 0245, 0246, 0247, 0343) ◆
BNE0219		New!! Flex 3+ Backrest features tension adjustment, width adjustment, and fore/aft positioning of the vertical backpost tubes to accommodate a wider range of clinical requirements.	(NC 0233, 0235, 0234, 0241, 0242, 0243, 0244, 0245, 0246, 0247, 0277, 0343) €289
BNE0220		Flex 3 Backrest features tension adjustment and width adjustment of the vertical backpost tubes.	(NC 0233, 0235, 0234, 0241, 0242, 0243, 0244, 0245, 0246, 0247, 0277, 0343) €126
BNE0221		Matrx Elite E2 (includes hardware mounted to frame, backrest supplied in box)	(NC 0240, 0250, 0260, 0270, 0305, 0311, 0321, 0322, 0800, 0810, 0820, 0830, 0880, 0890, 0900, 0920, 7777) €668
BNE0223		Matrx Elite E2 Deep (includes hardware mounted to frame, backrest supplied in box)	(NC 0240, 0250, 0260, 0270, 0305, 0311, 0321, 0322, 0800, 0810, 0820, 0830, 0880, 0890, 0900, 0920, 7777) €739
BNE0225		Matrx Posture Back (includes hardware mounted to frame, backrest supplied in box)	(NC 0235, 0240, 0241, 0242, 0243, 0244, 0245, 0246, 0247, 0250, 0260, 0270, 0305, 0311, 0321, 0322, 0335, 0800, 0810, 0820, 0830, 0880, 0890, 0900, 0920) €1,038
BNE0230		Backrest Tubes Without Backrest	(NC 0233, 0234, 0235, 0240, 0241, 0242, 0243, 0244, 0245, 0246, 0247, 0250, 0260, 0270, 0305, 0311, 0321, 0322, 0345, 0740, 0800, 0810, 0820, 0830, 0880, 0890, 0900, 0920) ○
MATRIX BACKREST HEIGHT			
BNE0233		Matrx Backrest Height 410 mm / 16"	○
BNE0235		Matrx Backrest Height 460 mm / 18"	○
BNE0234		Matrx Backrest Height 510 mm / 20"	○
MATRIX BACKREST WIDTH			
BNE0241		Matrx Backrest Width 360mm / 14" (Chest width 265mm / 10 1/2")	○
BNE0242		Matrx Backrest Width 380mm / 15" (Chest width 290mm / 11 1/2")	○
BNE0243		Matrx Backrest Width 410mm / 16" (Chest width 315mm / 12 1/2")	○
BNE0244		Matrx Backrest Width 430mm / 17" (Chest width 340mm / 13")	○
BNE0245		Matrx Backrest Width 460mm / 18" (Chest width 370mm / 14 1/2")	○
BNE0246		Matrx Backrest Width 480mm / 19" (Chest width 390mm / 15")	○
BNE0247		Matrx Backrest Width 510mm / 20" (Chest width 420mm / 16 1/2")	○
BNE7777		No backrest width, select with Matrx Posture back only	○
BACKREST CUSHION*			
BNE0240		Cover for tension adjustable upholstery	(NC 2222) ○
BNE0250		Cushion, Laguna (Lateral Support)	(NC 2222) ○
BNE0260		Cushion, Mistral 2 (Waist Support)	(NC 0322, 2222) €179
BNE0270		Cushion, Passad 2 (Shoulder Support)	(NC 2222) ◆
BNE0276		Vicair Multifunctional O2 pressure relieving backrest cushion, supplied with Black mesh cover (Additional 5 day lead time)	(NC 0305, 0311, 0321, 0322) €1,000
BNE0277		Vicair Multifunctional O2 Flex pressure relieving backrest cushion, supplied with Black mesh cover (Additional 5 day lead time)	(NC 0305, 0311, 0321, 0322) €1,000
BNE0300		Without Backrest Cover/ Cushion	(NC 0305, 0311, 0321, 0322) ○
WIDER BACKREST COVER			
BNE0305		Extra large backrest cover or cushion. One size wider backrest cushion or tension cover fitted if selected	○
BACKREST UPHOLSTERY*			
BNE0311		Black Plush TR35	○
BNE0321		Multi-stretch Black PU coated fabric	◆
BNE0322		Black Spacer TR36	○
BNE2222		No upholstery, select with Matrx, and Without backrest options	○

PUSH HANDLES*			
BNE0335		Push Handles connected	◆
BNE0343		Push Handles connected,with Rea head/neckrest attachment	€96
BNE0345		Push Bar, Angle Adjustable	€123
SEAT TILT AND BACKREST REGULATION*			
BNE0350		Assistant operated gas piston seat tilt and backrest recline (NC 5555)	◆
BNE0360		Electrical operated seat tilt and backrest recline (Tilt range is -1° To +23°) (NC 1263)	€958
TILT LOCK*			
BNE1263		Tilt Lock (locking of pre-set angle)	◆
BNE5555		No Tilt Lock, Select with 0360 Electric operated seat and backrest	○
CENTRAL LEGREST			
BNE0370		Central Legrest, includes calfpad and footrest. Also select options 0420, 0470, 3333, 4444 (NC 0376, 0378, 0380, 0390, 0400, 0410, 0430, 0440, 0450, 0460, 0475, 0476, 0477, 0720, 0930, 0940, 0941, 0942, 0952, 0954, 0956, 1290)	€377
LEGREST HANGER ATTACHMENT*			
BNE0375		Standard position for legrest mounting	◆
SWO0189		Standard position with 0-30° Legrest Abduction function (NC 0370, 0475, 0476, 0477)	€116
BNE0376		High position 50mm higher for leg length compensation	○
BNE0378		90° Hanger attachment to achieve a 90° legrest position on the Elevating Legrest only (NC 0390, 0400, 0410, 0440, 0450, 0460)	○
LEGREST HANGER LH*			
BNE0380		Elevating Legrest, angle adjustable 0 to 75° swing in/swing out with knee pad, calf support with locking clip LH (NC 0475, 0490, 0942, 0945, 6666)	◆
BNE0390		Legrest Hanger fixed 80 deg, swing in/swing out LH (NC 0928, 0929, 0930)	○
BNE0400		Legrest Hanger fixed 90 deg, swing in/swing out LH (NC 0520, 0928, 0929, 0930)	○
BNE0410		Legrest Amputee LH (NC 0475, 0476, 0477, 0490, 0520, 0928, 0929, 0930, 0941, 0942, 0945, 0952, 0954, 0956)	€171
BNE0420		Without Legrest, LH (NC 0475, 0476, 0477, 0490, 0520, 0928, 0929, 0930, 0941, 0942, 0945, 0952, 0954, 0956)	○
LEGREST HANGER RH*			
BNE0430		Elevating Legrest, angle adjustable 0 to 75°swing In/swing out with knee pad and calf support with locking clip RH (NC 0475, 0550, 0942, 0945, 6666)	◆
BNE0440		Legrest Hanger fixed 80 deg, swing in/swing out RH (NC 0931, 0932, 0940)	○
BNE0450		Legrest Hanger fixed 90 deg, swing in/swing out RH (NC 0580, 0931, 0932, 0940)	○
BNE0460		Legrest Amputee RH (NC 0475, 0476, 0477, 0550, 0580, 0931, 0932, 0940, 0941, 0942, 0945, 0952, 0954, 0956)	€171
BNE0470		Without Legrest, RH (NC 0475, 0476, 0477, 0550, 0580, 0931, 0932, 0940, 0941, 0942, 0945, 0952, 0954, 0956)	○
FOOTREST LH*			
BNE0475		Footrest one piece, depth and angle adjustable short tube (NC 0477, 0490, 0520, 0550, 0580, 0930, 0940, 0952, 0954, 0956)	€145
BNE0476		Footrest one piece, depth and angle adjustable long tube (NC 0477, 0490, 0520, 0550, 0580, 0930, 0940, 0952, 0954, 0956)	€145
BNE0477		Footbox, angle adjustable, padded with short tube (NC 0490, 0520, 0550, 0580, 0930, 0940, 0941, 0942, 0945, 0952, 0954, 0956)	€687
BNE0490		Footrest lower part, depth and angle adjustable with short tube and heel strap, LH. Not compatible with Elevating legrest.	○
BNE0520		Footrest lower part, depth and angle adjustable with long tube and heel strap, LH. All knee to heel ranges can be achieved when used with Elevating legrests.	○
BNE3333		Without footrest, only select with option 0370 Central legrest, 0410 Amputee legrest and 0420 Without legrest	○

FOOTREST RH*			
BNE0550		Footrest lower part, depth and angle adjustable with short tube and heel strap, RH. Not compatible with Elevating legrest.	○
BNE0580		Footrest lower part, depth and angle adjustable with long tube and heel strap, RH. All knee to heel ranges can be achieved when used with Elevating legrests.	○
BNE4444		Without footrest, only select with options 0370 Central legrest, 0460 Amputee legrest, 0470 Without legrest, 0475 and 0476 One piece footrests or 0477 Footbox	○
ARMREST LH*			
BNE0600		Armrest Ducktail Pad LH, height adjustable with autolock (NC 1232)	○
BNE0605		Armrest Wide Pad LH, height adjustable with autolock. Compatible with 16" wheels only (NC 1232)	◆
BNE0607		Armrests Reinforced Wide Pad, LH, height adjustable with autolock (NC 1230, 1231, 1232, 1250, 1260)	€112
BNE0610		Armrest Padded LH, height adjustable with autolock (NC 0680, 0695, 1230, 1231, 1232, 1250, 1260, 9999)	€79
BNE0615		Armrest Hemiplegic LH, height adjustable with autolock (NC 0695, 1230, 1231, 1232, 1240, 1245, 1250, 1260)	€222
BNE0620		Without Armrest LH Select when choosing BNE1232 Half Tray LH, Narrow Armrest is included (NC 0695, 1230, 1231, 1240, 1245, 1250, 1260)	○
ARMREST RH*			
BNE0630		Armrest Ducktail Pad RH, height adjustable with autolock (NC 1233)	○
BNE0635		Armrest Wide Pad RH, height adjustable with autolock. Compatible with 16" wheels only (NC 1233)	◆
BNE0637		Armrest Reinforced Wide Pad, RH, height adjustable with autolock (NC 1230, 1231, 1233, 1250, 1260)	€112
BNE0640		Armrest Padded RH, height adjustable with autolock. (NC 0680, 0695, 1230, 1231, 1233, 1250, 1260, 9999)	€79
BNE0645		Armrest Hemiplegic, RH, height adjustable with autolock (NC 0695, 1230, 1231, 1233, 1240, 1245, 1250, 1260)	€222
BNE0650		Without Armrest RH Select when choosing BNE1233 Half Tray RH, Narrow Armrest is included (NC 0695, 1230, 1231, 1240, 1245, 1250, 1260)	○
ARMREST OPTIONS			
BNE0680		Side Support Cover, PU fabric cover for the sideguards (Removable with tools) (NC 0695, 9999)	€63
BNE9999		Side Support Cover with foam cushion. 2cm thick cushion supplied complete with PU fabric cover. Reduces seat width by 2cm for a comfortable fit. Consists of options 0680 and 0690 (Removable with tools) (NC 0695)	€63
BNE0695		Pad for sideguard, removable with hook and loop fastening	€63
SEAT OPTIONS			
BNE0710		Incontinence Cover	€43
BNE0720		Pommel Abduction Cushion	€145
BNE0725		Tilt Scale for seat tilt	€53

HEADREST OPTIONS*			
BNE0730		Headrest, Sideway Adjustable Fixed hardware included (NC 0742, 0744, 0746, 0750, 0770, 0782)	◆
BNE0750		Headrest with cheek support, sideway adjustable Fixed hardware included (NC 0740, 0742, 0744, 0746, 0770, 0782)	€128
BNE0731		Matrx Elan Standard Pad (254 mm / 10") (NC 0740, 0770, 0782, 0790)	€189
BNE0732		Matrx Elan Standard Pad (355 mm / 14") (NC 0740, 0770, 0782, 0790)	€189
BNE0733		Matrx Elan Occipital Pad (229 mm / 9") (NC 0740, 0770, 0782, 0790)	€189
BNE0734		Matrx Elan Occipital Pad (305 mm / 12") (NC 0740, 0770, 0782, 0790)	€189
BNE0735		Matrx Elan 4-Point Pad (229W x 177H mm / 9"W x 7"H) (NC 0740, 0770, 0782, 0790)	€357
BNE0736		Matrx Elan 4-Point Pad (280W x 254H mm / 11"W x 10"H) (NC 0740, 0770, 0782, 0790)	€357
BNE9991		Without headrest (NC 0742, 0744, 0746)	○
HEADREST HARDWARE*			
BNE0742		Matrx Elan Hardware (NC 0740, 0770, 0790, 1206, 1207, 1300)	€285
BNE0744		Matrx Elan Extended Hardware (Extra link) (NC 0740, 0770, 0790, 1206, 1207, 1300)	€340
BNE0746		Matrx Elan LOXX Hardware (NC 0740, 0770, 0790, 1206, 1207, 1300)	€285
BNE8888		Without Elan headrest hardware (Only select with option 0730 & 0750 Headrests) (NC 0732, 0733, 0734, 0735, 0736)	○
OTHER HEADREST OPTIONS			
BNE0740		Cover for headrest	€39
BNE0770		Neckrest	○
BNE9992		Without neckrest	○
BNE0782		Alternative headrest/neckrest mounting bracket	€24
BNE0790		Handles Clinic Kit for Head/Neckrest (Levers which replace screws for easy adjustment)	€61
TRUNK SUPPORT LH			
BNE0800		Trunk Support Multi Adjustable with standard arm LH (NC 0900)	○
BNE0820		Trunk Support Multi Adjustable with extended height arm LH (NC 0900)	○
BNE0880		NEW!! Trunk Support Swingaway LH, multi adjustable (NC 0900)	€128
TRUNK SUPPORT RH			
BNE0810		Trunk Support Multi Adjustable with standard arm RH (NC 0900)	○
BNE0830		Trunk Support Multi Adjustable with extended height arm RH (NC 0900)	○
BNE0890		NEW!! Trunk Support Swingaway RH, multi adjustable (NC 0900)	€128
BACKREST OPTIONS			
BNE0900		Trunk Support Attachment to mount laterals later. If trunk supports are chosen, the attachment bracket is already included	○
BNE0915		Tilt Scale for backrest recline	€48
BNE0920		Wedges (One Pair)	€14
ELEVATING LEGREST OPTIONS*			
BNE0928		Calf pad cover, for LH elevating legrest	○
BNE0929		Knee pad for LH elevating legrest	◆
BNE0931		Calf pad cover, for RH elevating legrest	○
BNE0932		Knee pad for RH elevating legrest	◆
BNE6666		Without Elevating Legrest Options (Select with Fixed legrests, Amputee legrests and Without legrests only)	○

LEGREST OPTIONS			
BNE0930	<input type="checkbox"/>	Footplate Extender LH (NC 0941)	€40
BNE0940	<input type="checkbox"/>	Footplate Extender RH (NC 0941)	€40
BNE0941	<input type="checkbox"/>	Cover Footplate, Black Plush Fabric	€66
BNE0942	<input type="checkbox"/>	Calf Strap Padded	€88
BNE0945	<input type="checkbox"/>	Calf Strap (for One piece footrest)	○
BODYPOINT ANKLE HUGGERS			
BNE0952	<input type="checkbox"/>	Bodypoint Ankle Huggers Large (Supplied with fixing kit, unmounted)	€189
BNE0954	<input type="checkbox"/>	Bodypoint Ankle Huggers Medium (Supplied with fixing kit, unmounted)	€189
BNE0956	<input type="checkbox"/>	Bodypoint Ankle Huggers Small (Supplied with fixing kit, unmounted)	€189
CHASSIS COLOUR*			
BNE0944	<input type="checkbox"/>	Silver Pearl	◆
BNE0946	<input type="checkbox"/>	Azurite (Dark Blue)	○
BNE0947	<input type="checkbox"/>	Sand (Beige)	○
BNE0948	<input type="checkbox"/>	Electric Green	○
BNE0949	<input type="checkbox"/>	Happy Red	○
CHASSIS SAFETY			
BNE0950	<input type="checkbox"/>	Chassis extender, castors extend 50mm forward for added stability	○
CASTOR ATTACHMENT*			
BNE0960	<input type="checkbox"/>	Fixed castor attachment	◆
FRONT CASTORS*			
BNE0992	<input type="checkbox"/>	150x30 mm (6" x 1 1/4") Medium Hard Black tyre	○
BNE1015	<input type="checkbox"/>	200x32 mm (8" x 1 1/4") Black Solid Tyre	○
BNE1025	<input type="checkbox"/>	200x45 mm (8" X 1 3/4") Black Solid tyre	◆
REAR WHEEL SIZE*			
BNE1047	<input type="checkbox"/>	16" / 405 mm Transit rear wheel with fixed attachment	◆
REAR WHEEL TYPE*			
BNE1076	<input type="checkbox"/>	Hub Brake 16" Transit with Black Plastic Hub and Rim	◆
REAR WHEEL TYRE*			
BNE1092	<input type="checkbox"/>	Solid Black tyres	◆
BRAKE*			
BNE1170	<input type="checkbox"/>	User manoeuvred tyre brakes, includes extended handles	◆
PELVIC BELT*			
BNE1200	<input type="checkbox"/>	Pelvic belt	◆
BNE1201	<input type="checkbox"/>	Bodypoint 2-Point Centre-pull Push Button Medium (Mounted to chair)	€171
BNE1202	<input type="checkbox"/>	Bodypoint 2-Point Centre-pull Push Button Large (Mounted to chair)	€171
BNE1203	<input type="checkbox"/>	Bodypoint 4-Point Centre-pull Push Button Medium (Mounted to chair)	€161
BNE1204	<input type="checkbox"/>	Bodypoint 4-Point Centre-pull Push Button Large (Mounted to chair)	€161
BODYPOINT HARNESS			
BNE1206	<input type="checkbox"/>	Bodypoint Pivotfit Harness Medium (Supplied with fixing kit, unmounted)	€171
BNE1207	<input type="checkbox"/>	Bodypoint Pivotfit Harness Large (Supplied with fixing kit, unmounted)	€171

TRANSPORTATION BRACKETS*			
BNE1225		Transportation Brackets, mounted to the front of the wheelchair	€16
OTHER OPTIONS			
BNE1230		Fixed Table Tray without elbow pads includes mounting kit (NC 1237, 1260)	€173
BNE1231		Fixed Table Tray with elbow Pads includes mounting kit (NC 1237, 1260)	€173
BNE1232		Half table tray, clear plastic, flip up, depth adjustable, left side (armrest with straight narrow pad included) (NC 1235, 1260)	€138
BNE1233		Half table tray, clear, plastic, flip up, depth adjustable, right side (armrest with straight narrow pad included) (NC 1235, 1260)	€138
BNE1240		Swing away table without elbow pads (NC 1237, 1250, 1260)	€242
BNE1245		Swing away table with elbow pads (NC 1237, 1250, 1260)	€242
BNE1235		Table tray padded cushion (To be ordered in combination with the fixed & swingaway table tray)	€63
BNE1237		Half table tray padded cushion	€63
BNE1250		Locking device for table tray	€17
BNE1260		Mounting kit for table (if ordered separately)	€10
BNE1290		Handles Clinic Kit for backrest, armrest, legrest and seat frame (Levers which replace screws for easy adjustment)	€67
BNE1300		Harness attachment , must be selected in combination with Head/Neckrest	€98
BNE1305		Front castor locks, includes front tiedown brackets. Allows castors to remain parallel for easy access into confined spaces such as lifts or vehicles	€240
BNE1310		Attachment for positioning belt	€39

TOTAL PRICE	
--------------------	--

Non Conforming Options:

BNE0110	400 mm / 16" Seat Height
BNE0182	Short cushion depth: Libra delivered at 405mm deep, Solution Xtra delivered at 430mm deep
BNE0183	Long cushion depth: Libra delivered at 455mm deep, Solution Xtra delivered at 455mm deep
BNE0186	Tapered Wedge Pair
BNE0187	Abductor
BNE0188	Lateral Side Wedge
BNE0189	Anterior Front Wedge
BNE0201	Multi-stretch Black PU coated fabric
BNE0202	Black Spacer TR36
BNE1111	No upholstery, select with options 0167 Vicair Adjuster O2 or 0180 Without seat cushion
BNE0206	+ 100 mm Backrest width pre-adjustment, set at +100mm but can be adjusted if required
BNE0220	Flex 3 Backrest features tension adjustment and width adjustment of the vertical backpost tubes.
BNE0221	Matrx Elite E2 (includes hardware mounted to frame, backrest supplied in box)
BNE0223	Matrx Elite E2 Deep (includes hardware mounted to frame, backrest supplied in box)
BNE0225	Matrx Posture Back (includes hardware mounted to frame, backrest supplied in box)
BNE0233	Matrx Backrest Height 410 mm / 16"
BNE0235	Matrx Backrest Height 460 mm / 18"
BNE0234	Matrx Backrest Height 510 mm / 20"
BNE0241	Matrx Backrest Width 360mm / 14" (Chest width 265mm / 10 1/2")
BNE0242	Matrx Backrest Width 380mm / 15" (Chest width 290mm / 11 1/2")
BNE0243	Matrx Backrest Width 410mm / 16" (Chest width 315mm / 12 1/2")
BNE0244	Matrx Backrest Width 430mm / 17" (Chest width 340mm / 13")
BNE0245	Matrx Backrest Width 460mm / 18" (Chest width 370mm / 14 1/2")
BNE0246	Matrx Backrest Width 480mm / 19" (Chest width 390mm / 15")
BNE0247	Matrx Backrest Width 510mm / 20" (Chest width 420mm / 16 1/2")
BNE7777	No backrest width, select with Matrx Posture back only
BNE0240	Cover for tension adjustable upholstery
BNE0250	Cushion, Laguna (Lateral Support)
BNE0260	Cushion, Mistral 2 (Waist Support)
BNE0270	Cushion, Passad 2 (Shoulder Support)
BNE0277	Vicair Multifunctional O2 Flex pressure relieving backrest cushion, supplied with Black mesh cover (Additional 5 day lead time)
BNE0305	Extra large backrest cover or cushion. One size wider backrest cushion or tension cover fitted if selected
BNE0311	Black Plush TR35
BNE0321	Multi-stretch Black PU coated fabric
BNE0322	Black Spacer TR36
BNE2222	No upholstery, select with Matrx, and Without backrest options
BNE0335	Push Handles connected
BNE0343	Push Handles connected,with Rea head/neckrest attachment
BNE0345	Push Bar, Angle Adjustable
BNE0360	Electrical operated seat tilt and backrest recline (Tilt range is -1° To +23°)
BNE1263	Tilt Lock (locking of pre-set angle)
BNE5555	No Tilt Lock, Select with 0360 Electric operated seat and backrest
BNE0370	Central Legrest, includes calfpad and footrest. Also select options 0420, 0470, 3333, 4444

SWO0189	Standard position with 0-30° Legrest Abduction function
BNE0376	High position 50mm higher for leg length compensation
BNE0378	90° Hanger attachment to achieve a 90° legrest position on the Elevating Legrest only
BNE0380	Elevating Legrest, angle adjustable 0 to 75° swing in/swing out with knee pad, calf support with locking clip LH
BNE0390	Legrest Hanger fixed 80 deg, swing in/swing out LH
BNE0400	Legrest Hanger fixed 90 deg, swing in/swing out LH
BNE0410	Legrest Amputee LH
BNE0430	Elevating Legrest, angle adjustable 0 to 75° swing In/swing out with knee pad and calf support with locking clip RH
BNE0440	Legrest Hanger fixed 80 deg, swing in/swing out RH
BNE0450	Legrest Hanger fixed 90 deg, swing in/swing out RH
BNE0460	Legrest Amputee RH
BNE0475	Footrest one piece, depth and angle adjustable short tube
BNE0476	Footrest one piece, depth and angle adjustable long tube
BNE0477	Footbox, angle adjustable, padded with short tube
BNE0490	Footrest lower part, depth and angle adjustable with short tube and heel strap, LH. Not compatible with Elevating legrest.
BNE0520	Footrest lower part, depth and angle adjustable with long tube and heel strap, LH. All knee to heel ranges can be achieved when used with Elevating legrests.
BNE0550	Footrest lower part, depth and angle adjustable with short tube and heel strap, RH. Not compatible with Elevating legrest.
BNE0580	Footrest lower part, depth and angle adjustable with long tube and heel strap, RH. All knee to heel ranges can be achieved when used with Elevating legrests.
BNE0680	Side Support Cover, PU fabric cover for the sideguards (Removable with tools)
BNE9999	Side Support Cover with foam cushion. 2cm thick cushion supplied complete with PU fabric cover. Reduces seat width by 2cm for a comfortable fit. Consists of options 0680 and 0690 (Removable with tools)
BNE0695	Pad for sideguard, removable with hook and loop fastening
BNE0720	Pommel Abduction Cushion
BNE0750	Headrest with cheek support, sideways adjustable Fixed hardware included
BNE0732	Matrx Elan Standard Pad (355 mm / 14")
BNE0733	Matrx Elan Occipital Pad (229 mm / 9")
BNE0734	Matrx Elan Occipital Pad (305 mm / 12")
BNE0735	Matrx Elan 4-Point Pad (229W x 177H mm / 9"W x 7"H)
BNE0736	Matrx Elan 4-Point Pad (280W x 254H mm / 11"W x 10"H)
BNE0742	Matrx Elan Hardware
BNE0744	Matrx Elan Extended Hardware (Extra link)
BNE0746	Matrx Elan LOXX Hardware
BNE0740	Cover for headrest
BNE0770	Neckrest
BNE0782	Alternative headrest/neckrest mounting bracket
BNE0790	Handles Clinic Kit for Head/Neckrest (Levers which replace screws for easy adjustment)
BNE0800	Trunk Support Multi Adjustable with standard arm LH
BNE0820	Trunk Support Multi Adjustable with extended height arm LH
BNE0880	NEW!! Trunk Support Swingaway LH, multi adjustable
BNE0810	Trunk Support Multi Adjustable with standard arm RH
BNE0830	Trunk Support Multi Adjustable with extended height arm RH
BNE0890	NEW!! Trunk Support Swingaway RH, multi adjustable
BNE0900	Trunk Support Attachment to mount laterals later. If trunk supports are chosen, the attachment bracket is already included
BNE0920	Wedges (One Pair)

BNE0928	Calf pad cover, for LH elevating legrest
BNE0929	Knee pad for LH elevating legrest
BNE0931	Calf pad cover, for RH elevating legrest
BNE0932	Knee pad for RH elevating legrest
BNE6666	Without Elevating Legrest Options (Select with Fixed legrests, Amputee legrests and Without legrests only)
BNE0930	Footplate Extender LH
BNE0940	Footplate Extender RH
BNE0941	Cover Footplate, Black Plush Fabric
BNE0942	Calf Strap Padded
BNE0945	Calf Strap (for One piece footrest)
BNE0952	Bodypoint Ankle Huggers Large (Supplied with fixing kit, unmounted)
BNE0954	Bodypoint Ankle Huggers Medium (Supplied with fixing kit, unmounted)
BNE0956	Bodypoint Ankle Huggers Small (Supplied with fixing kit, unmounted)
BNE1015	200x32 mm (8" x 1 1/4") Black Solid Tyre
BNE1025	200x45 mm (8" X 1 3/4") Black Solid tyre
BNE1206	Bodypoint Pivotfit Harness Medium (Supplied with fixing kit, unmounted)
BNE1207	Bodypoint Pivotfit Harness Large (Supplied with fixing kit, unmounted)
BNE1230	Fixed Table Tray without elbow pads includes mounting kit
BNE1231	Fixed Table Tray with elbow Pads includes mounting kit
BNE1232	Half table tray, clear plastic, flip up, depth adjustable, left side (armrest with straight narrow pad included)
BNE1233	Half table tray, clear, plastic, flip up, depth adjustable, right side (armrest with straight narrow pad included)
BNE1240	Swing away table without elbow pads
BNE1245	Swing away table with elbow pads
BNE1235	Table tray padded cushion (To be ordered in combination with the fixed & swingaway table tray)
BNE1237	Half table tray padded cushion
BNE1250	Locking device for table tray
BNE1260	Mounting kit for table (if ordered separately)
BNE1290	Handles Clinic Kit for backrest, armrest, legrest and seat frame (Levers which replace screws for easy adjustment)
BNE1300	Harness attachment , must be selected in combination with Head/Neckrest



For a more complex customer requirement,
please email: orders@invacare.com.au

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):