



Yes, you can:

# Manual Ireland Prescription Form

Please refer to [www.invacare.ie](http://www.invacare.ie) for Local Distributor

NOVEMBER 2023  
Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## Rea Azalea Transit One Solution

### Standard Features

- Max User Weight 135KG/21 Stone
- Crash Tested ISO 7176-Part 19
- 16" Transit Propulsion with Attendant Hub Brake and Solid Black Tyres
- Choice of Seat Widths 390mm, 440mm or 490mm (+50mm or +100mm Width Adj)
- Choice of 2 Seat Heights 400mm or 450mm
- Depth Adjustable Seat 430mm-500mm
- Without Seat Cushion, Supplied with Seat Board
- Flex 3 Backrest Plate, Width and Height Adjustable
- Passad 2 Backrest Cushion with Multi-stretch Black PU coated fabric Cover and Push Handles Connected
- Assistant Operated Gas Piston Seat Tilt +1° - +25°
- Backrest Recline 30° with Tilt Lock
- Angle Adjustable Legrests with Calf Support, Height and Depth Adjustable Footplates
- Height Adjustable Wide Pad Armrests with Autolock
- Pearl Silver Frame Colour, Solid Black Front Castors
- User Manoeuvred Tyre Brakes
- Pelvic Belt, Adjustable both sides
- Anti Tippers
- Front Transportation Brackets

**Standard Base price: €2,668**



Please make your selection with a cross in the box     Free of charge     standard feature    POA special request

Customer weight =                      KG

Code	Description	
<b>PROPULSION*</b>		
BNE0020	Transit 16" wheels with hub brakes	<input checked="" type="checkbox"/>
<b>SEAT WIDTH*</b>		
BNE0040	390 mm/15" (Small chassis) (NC 0244, 0245, 0246, 0247)	<input type="checkbox"/>
BNE0050	440 mm/17" (Medium chassis) (NC 0241, 0242, 0246, 0247)	<input type="checkbox"/>
BNE0060	490 mm/19" (Large chassis) (NC 0241, 0242, 0243, 0244)	<input type="checkbox"/>
<b>SEAT WIDTH EXTENSION*</b>		
BNE0065	+ 50 mm pre-adjustment on armrests and legrests	<input type="checkbox"/>
BNE0075	+ 100 mm pre-adjustment on armrests and legrests (NC 0110)	<input type="checkbox"/>
BNE0085	Without seat width extension, set at 0mm adjustment but can be adjusted wider if required	<input type="checkbox"/>

<b>SEAT DEPTH*</b>			
BNE0100	<input type="checkbox"/>	430-500 mm / 17" - 20"	◆
<b>SEAT HEIGHT*</b>			
BNE0110	<input type="checkbox"/>	400 mm / 16" Seat Height	(NC 0370, 1015, 1025) ○
BNE0120	<input type="checkbox"/>	450 mm / 18" Seat Height	○
<b>SEAT ANGLE*</b>			
BNE0127	<input type="checkbox"/>	Negative Seat Angle -1 To +23°	○
BNE0000	<input type="checkbox"/>	Standard Seat Angle +1° to +25°	◆
<b>SEAT TYPE*</b>			
BNE0130	<input type="checkbox"/>	Depth Adjustable	◆
<b>SEAT CUSHION*</b>			
BNE0165	<input type="checkbox"/>	Cushion Flosshape seat upholstery	(NC 1111) €65
BNE0167	<input type="checkbox"/>	Cushion Vicair Adjuster O2 10cm (Black)	(NC 0201, 0202) €663
SWO0190	<input type="checkbox"/>	Cushion Flosshape seat upholstery Shorter depth 450mm deep. Only use with seat depth in shortest position	(NC 0202, 1111) €65
BNE0180	<input type="checkbox"/>	Without Seat Cushion, Seat Plate only	(NC 0201, 0202) ○
<b>SEAT UPHOLSTERY*</b>			
BNE0201	<input type="checkbox"/>	Multi-stretch Black PU coated fabric	◆
BNE0202	<input type="checkbox"/>	Black Spacer TR36	○
BNE1111	<input type="checkbox"/>	No upholstery, select with options 0167 Vicair Adjuster O2 or 0180 Without seat cushion	○
<b>BACKREST WIDTH EXTENSION*</b>			
BNE0207	<input type="checkbox"/>	Without Backrest width pre-adjustment. Backrest width adjustment is not applicable to MatrX backrests	○
<b>BACKREST*</b>			
BNE0221	<input type="checkbox"/>	Matrx Elite (includes hardware mounted to frame, backrest supplied in box)	(NC 0335, 1206, 1207, 7777) €636
BNE0223	<input type="checkbox"/>	Matrx PB Elite Deep (includes hardware mounted to frame, backrest supplied in box)	(NC 0335, 1206, 1207, 7777) €704
BNE0225	<input type="checkbox"/>	Matrx Posture Back (includes hardware mounted to frame, backrest supplied in box)	(NC 0235, 0241, 0242, 0243, 0244, 0245, 0246, 0247, 0335) €989
<b>MATRIX BACKREST HEIGHT*</b>			
BNE0233	<input type="checkbox"/>	Matrx Backrest Height 410 mm / 16"	○
BNE0235	<input type="checkbox"/>	Matrx Backrest Height 460 mm / 18"	○
BNE0234	<input type="checkbox"/>	Matrx Backrest Height 510 mm / 20"	○
<b>MATRIX BACKREST WIDTH CHOICE*</b>			
BNE0241	<input type="checkbox"/>	360mm Backrest width, to fit chest width 265mm	○
BNE0242	<input type="checkbox"/>	380mm Backrest width, to fit chest width 290mm	○
BNE0243	<input type="checkbox"/>	410mm Backrest width, to fit chest width 315mm	○
BNE0244	<input type="checkbox"/>	430mm Backrest width, to fit chest width 340mm	○
BNE0245	<input type="checkbox"/>	460mm Backrest width, to fit chest width 370mm	○
BNE0246	<input type="checkbox"/>	480mm Backrest width, to fit chest width 390mm	○
BNE0247	<input type="checkbox"/>	510mm Backrest width, to fit chest width 420mm	○
BNE7777	<input type="checkbox"/>	No Backrest width, select with MatrX Posture back, option 0225	○
<b>BACKREST CUSHION*</b>			
BNE0300	<input type="checkbox"/>	Without Backrest Cover/ Cushion	○
<b>BACKREST UPHOLSTERY*</b>			
BNE2222	<input type="checkbox"/>	No upholstery, select with MatrX	○

<b>PUSH HANDLES*</b>			
BNE0335		Push Handles connected	◆
BNE0343		Push Handles connected,with Rea head/neckrest attachment	€87
BNE0345		Push Bar, Angle Adjustable	€117
<b>SEAT TILT AND BACKREST REGULATION*</b>			
BNE0350		Assistant operated gas piston seat tilt and backrest recline	◆
<b>TILT LOCK*</b>			
BNE1263		Tilt Lock (locking of pre-set angle)	◆
<b>CENTRAL LEGREST</b>			
BNE0370		Central Legrest, includes calfpad and footrest. Also select options 0420, 0470, 3333, 4444	(NC 0376, 0378, 0380, 0390, 0400, 0410, 0430, 0440, 0450, 0460, 0475, 0476, 0477, 0720, 0930, 0940, 0941, 0942, 0952, 0954, 0956, 1290) €341
<b>LEGREST HANGER ATTACHMENT*</b>			
BNE0375		Standard position for legrest mounting	◆
SWO0189		Standard position with 0-30° Legrest Abduction function	€105
BNE0376		High position 50mm higher for leg length compensation	○
BNE0378		90° Hanger attachment to achieve a 90° legrest position on the Elevating Legrest only	(NC 0390, 0400, 0410, 0440, 0450, 0460) ○
<b>LEGREST HANGER LH*</b>			
BNE0380		Elevating Legrest, angle adjustable 0 to 75° swing in/swing out with knee pad, calf support with locking clip LH	(NC 0475, 0480, 0490, 0942) ◆
BNE0390		Legrest Hanger fixed 80 deg, swing in/swing out LH	(NC 0928, 0929, 0930) ○
BNE0400		Legrest Hanger fixed 90 deg, swing in/swing out LH	(NC 0480, 0520, 0928, 0929, 0930) ○
BNE0410		Legrest Amputee LH	(NC 0475, 0476, 0477, 0480, 0490, 0520, 0928, 0929, 0930, 0941, 0942, 0945) €154
BNE0420		Without Legrest, LH	(NC 0475, 0476, 0477, 0480, 0490, 0520, 0928, 0929, 0930, 0941, 0942, 0945) ○
<b>LEGREST HANGER RH*</b>			
BNE0430		Elevating Legrest, angle adjustable 0 to 75°swing In/swing out with knee pad and calf support with locking clip RH	(NC 0475, 0540, 0550, 0942, 0945, 6666) ◆
BNE0440		Legrest Hanger fixed 80 deg, swing in/swing out RH	(NC 0931, 0932, 0940) ○
BNE0450		Legrest Hanger fixed 90 deg, swing in/swing out RH	(NC 0540, 0580, 0931, 0932, 0940) ○
BNE0460		Legrest Amputee RH	(NC 0475, 0476, 0477, 0540, 0550, 0580, 0931, 0932, 0940, 0941, 0942, 0945) €154
BNE0470		Without Legrest, RH	(NC 0475, 0476, 0477, 0540, 0550, 0580, 0931, 0932, 0940, 0941, 0942, 0945) ○
<b>FOOTREST LH*</b>			
BNE0475		Footrest one piece, depth and angle adjustable short tube	(NC 0477, 0480, 0490, 0520, 0540, 0550, 0580, 0930, 0940, 0952, 0954, 0956) €131
BNE0476		Footrest one piece, depth and angle adjustable long tube	(NC 0477, 0480, 0490, 0520, 0540, 0550, 0580, 0930, 0940, 0952, 0954, 0956) €131
BNE0477		Footbox, angle adjustable, padded with short tube	(NC 0480, 0490, 0520, 0540, 0550, 0580, 0930, 0940, 0941, 0942, 0945, 0952, 0954, 0956) €622
BNE0480		Footrest lower part, fixed with short tube and heel strap, LH	(NC 0930, 0941, 0952, 0954, 0956) ○
BNE0490		Footrest lower part, depth and angle adjustable with short tube and heel strap, LH. Not compatible with Elevating legrest.	○
BNE0520		Footrest lower part, depth and angle adjustable with long tube and heel strap, LH. All knee to heel ranges can be achieved when used with Elevating legrests.	○
BNE3333		Without footrest, only select with option 0370 Central legrest, 0410 Amputee legrest and 0420 Without legrest	○

FOOTREST RH*			
BNE0540		Footrest lower part, fixed with short tube and heel strap, RH (NC 0940, 0941, 0952, 0954, 0956)	○
BNE0550		Footrest lower part, depth and angle adjustable with short tube and heel strap, RH. Not compatible with Elevating legrest.	○
BNE0580		Footrest lower part, depth and angle adjustable with long tube and heel strap, RH. All knee to heel ranges can be achieved when used with Elevating legrests.	○
BNE4444		Without footrest, only select with options 0370 Central legrest, 0460 Amputee legrest, 0470 Without legrest, 0475 and 0476 One piece footrests or 0477 Footbox	○
ARMREST LH*			
BNE0600		Armrest Ducktail Pad LH, height adjustable with autolock (NC 1232)	○
BNE0605		Armrest Wide Pad LH, height adjustable with autolock. Compatible with 16" wheels only (NC 1232)	◆
BNE0607		Armrests Reinforced Wide Pad, LH, height adjustable with autolock (NC 1230, 1231, 1232, 1250, 1260)	€101
BNE0610		Armrest Padded LH, height adjustable with autolock (NC 1230, 1231, 1232, 1250, 1260)	€71
BNE0615		Armrest Hemiplegic LH, height adjustable with autolock (NC 0695, 1230, 1231, 1232, 1250, 1260)	€201
BNE0620		Without Armrest LH Select when choosing BNE1232 Half Tray LH, Narrow Armrest is included (NC 0695, 1230, 1231, 1250, 1260)	○
ARMREST RH*			
BNE0630		Armrest Ducktail Pad RH, height adjustable with autolock (NC 1233)	○
BNE0635		Armrest Wide Pad RH, height adjustable with autolock. Compatible with 16" wheels only (NC 1233)	◆
BNE0637		Armrest Reinforced Wide Pad, RH, height adjustable with autolock (NC 1230, 1231, 1233, 1250, 1260)	€101
BNE0640		Armrest Padded RH, height adjustable with autolock. (NC 1230, 1231, 1233, 1250, 1260)	€71
BNE0645		Armrest Hemiplegic, RH, height adjustable with autolock (NC 0695, 1230, 1231, 1233, 1240, 1245, 1250, 1260)	€201
BNE0650		Without Armrest RH Select when choosing BNE1233 Half Tray RH, Narrow Armrest is included (NC 0695, 1230, 1231, 1240, 1245, 1250, 1260)	○
ARMREST OPTIONS. Side support cover 0680, must be selected if 0690 Foam for side support is required			
BNE0680		Side Support Cover, removable with tools (Wipeable cover for the sideguard) (NC 0695)	€57
BNE0690		Foam for side support cover (Internal cushioned pad, which is only possible in combination with 0680) (NC 0695)	○
BNE0695		Pad for sideguard, removable with hook and loop fastening	€57
SEAT OPTIONS			
BNE0710		Incontinence Cover	€39
BNE0720		Pommel Abduction Cushion (NC 0722)	€131
BNE0722		2 cm Increased seat height plate	€81
BNE0725		Tilt Scale for seat tilt	€48
HEADREST OPTIONS			
BNE0730		Headrest, Sideway Adjustable Fixed hardware included (NC 0750, 0770)	◆
BNE0750		Headrest with cheek support, sideway adjustable Fixed hardware included (NC 0740, 0770)	€122
BNE0731		Matrx Elan Standard Pad (254 mm / 10") (NC 0740, 0770, 0790, 1206, 1207)	€180
BNE0732		Matrx Elan Standard Pad (355 mm / 14") (NC 0740, 0770, 0790, 1206, 1207)	€180
BNE0733		Matrx Elan Occipital Pad (229 mm / 9") (NC 0740, 0770, 0790, 1206, 1207)	€180
BNE0734		Matrx Elan Occipital Pad (305 mm / 12") (NC 0740, 0770, 0790, 1206, 1207)	€180
BNE0735		Matrx Elan 4-Point Pad (229W x 177H mm / 9"W x 7"H) (NC 0740, 0770, 0790, 1206, 1207)	€340
BNE0736		Matrx Elan 4-Point Pad (280W x 254H mm / 11"W x 10"H) (NC 0740, 0770, 0790, 1206, 1207)	€340

HEADREST HARDWARE			
BNE0742	<input type="checkbox"/>	Matrx Elan Hardware (NC 0740, 0770, 0782, 0790, 1206, 1207)	€271
BNE0744	<input type="checkbox"/>	Matrx Elan Extended Hardware (Extra link) (NC 0740, 0770, 0782, 0790, 1206, 1207)	€324
BNE0746	<input type="checkbox"/>	Matrx Elan LOXX Hardware (NC 0740, 0770, 0782, 0790, 1206, 1207)	€271
OTHER HEADREST OPTIONS			
BNE0740	<input type="checkbox"/>	Cover for headrest	€35
BNE0770	<input type="checkbox"/>	Neckrest	○
BNE0782	<input type="checkbox"/>	Alternative headrest/neckrest mounting bracket	€22
BNE0790	<input type="checkbox"/>	Handles Clinic Kit for Head/Neckrest (Levers which replace screws for easy adjustment)	€55
BACKREST OPTIONS			
BNE0915	<input type="checkbox"/>	Tilt Scale for backrest recline	€44
ELEVATING LEGREST OPTIONS*			
BNE0928	<input type="checkbox"/>	Calf pad cover, for LH elevating legrest	○
BNE0929	<input type="checkbox"/>	Knee pad for LH elevating legrest	◆
BNE0931	<input type="checkbox"/>	Calf pad cover, for RH elevating legrest	○
BNE0932	<input type="checkbox"/>	Knee pad for RH elevating legrest	◆
BNE6666	<input type="checkbox"/>	Without Elevating Legrest Options (Select with Fixed legrests, Amputee legrests and Without legrests only)	○
LEGREST OPTIONS			
BNE0930	<input type="checkbox"/>	Footplate Extender LH (NC 0941)	€36
BNE0940	<input type="checkbox"/>	Footplate Extender RH (NC 0941)	€36
BNE0941	<input type="checkbox"/>	Cover Footplate, Black Plush Fabric	€60
BNE0942	<input type="checkbox"/>	Calf Strap Padded	€80
BNE0945	<input type="checkbox"/>	Calf Strap (for One piece footrest)	○
BODYPOINT ANKLE HUGGERS			
BNE0952	<input type="checkbox"/>	Bodypoint Ankle Huggers Large (Supplied with fixing kit, unmounted)	€180
BNE0954	<input type="checkbox"/>	Bodypoint Ankle Huggers Medium (Supplied with fixing kit, unmounted)	€180
BNE0956	<input type="checkbox"/>	Bodypoint Ankle Huggers Small (Supplied with fixing kit, unmounted)	€180
CHASSIS COLOUR*			
BNE0944	<input type="checkbox"/>	Silver Pearl	◆
BNE0946	<input type="checkbox"/>	Azurite (Dark Blue)	○
BNE0947	<input type="checkbox"/>	Sand (Beige)	○
BNE0948	<input type="checkbox"/>	Electric Green	○
BNE0949	<input type="checkbox"/>	Happy Red	○
CHASSIS SAFETY			
BNE0950	<input type="checkbox"/>	Chassis extender, castors extend 50mm forward for added stability	○
CASTOR ATTACHMENT*			
BNE0960	<input type="checkbox"/>	Fixed castor attachment	◆
FRONT CASTORS*			
BNE0992	<input type="checkbox"/>	150x30 mm (6" x 1 1/4") Medium Hard Black tyre	○
BNE1015	<input type="checkbox"/>	200x32 mm (8" x 1 1/4") Black Solid Tyre	○
BNE1025	<input type="checkbox"/>	200x45 mm (8" X 1 3/4") Black Solid tyre	◆
REAR WHEEL SIZE*			
BNE1047	<input type="checkbox"/>	16" / 405 mm Transit rear wheel with fixed attachment	◆

<b>REAR WHEEL TYPE*</b>			
BNE1076	<input type="checkbox"/>	Hub Brake 16" Transit with Black Plastic Hub and Rim	◆
<b>REAR WHEEL TYRE*</b>			
BNE1092	<input type="checkbox"/>	Solid Black tyres	◆
<b>BRAKE*</b>			
BNE1170	<input type="checkbox"/>	User manoeuvred tyre brakes	◆
<b>PELVIC BELT*</b>			
BNE1200	<input type="checkbox"/>	Pelvic belt	◆
BNE1201	<input type="checkbox"/>	Bodypoint 2-Point Centre-pull Push Button Medium (Mounted to chair)	€163
BNE1202	<input type="checkbox"/>	Bodypoint 2-Point Centre-pull Push Button Large (Mounted to chair)	€163
BNE1203	<input type="checkbox"/>	Bodypoint 4-Point Centre-pull Push Button Medium (Mounted to chair)	€153
BNE1204	<input type="checkbox"/>	Bodypoint 4-Point Centre-pull Push Button Large (Mounted to chair)	€153
<b>BODYPOINT HARNESS</b>			
BNE1206	<input type="checkbox"/>	Bodypoint Pivotfit Harness Medium (Supplied with fixing kit, unmounted)	€163
BNE1207	<input type="checkbox"/>	Bodypoint Pivotfit Harness Large (Supplied with fixing kit, unmounted)	€163
<b>TRANSPORTATION BRACKETS*</b>			
BNE1225	<input type="checkbox"/>	Transportation Brackets, mounted to the front of the wheelchair	€14
<b>OTHER OPTIONS</b>			
BNE1230	<input type="checkbox"/>	Fixed Table Tray without elbow pads includes mounting kit (NC 1237, 1260)	€156
BNE1231	<input type="checkbox"/>	Fixed Table Tray with elbow Pads includes mounting kit (NC 1237, 1260)	€156
BNE1232	<input type="checkbox"/>	Half table tray, clear plastic, flip up, depth adjustable, left side (armrest with straight narrow pad included) (NC 1235, 1260)	€125
BNE1233	<input type="checkbox"/>	Half table tray, clear, plastic, flip up, depth adjustable, right side (armrest with straight narrow pad included) (NC 1235, 1260)	€125
BNE1240	<input type="checkbox"/>	Swing away table without elbow pads (NC 1237, 1250, 1260)	€219
BNE1245	<input type="checkbox"/>	Swing away table with elbow pads (NC 1237, 1260)	€219
BNE1235	<input type="checkbox"/>	Table tray padded cushion (To be ordered in combination with the fixed & swingaway table tray)	€57
BNE1237	<input type="checkbox"/>	Half table tray padded cushion	€57
BNE1250	<input type="checkbox"/>	Locking device for table tray	€15
BNE1260	<input type="checkbox"/>	Mounting kit for table (if ordered separately)	€8
BNE1290	<input type="checkbox"/>	Handles Clinic Kit for backrest, armrest, legrest and seat frame (Levers which replace screws for easy adjustment)	€61
BNE1300	<input type="checkbox"/>	Harness attachment , must be selected in combination with Head/Neckrest	€89
BNE1305	<input type="checkbox"/>	Front castor locks, includes front tiedown brackets. Allows castors to remain parallel for easy access into confined spaces such as lifts or vehicles	€217
BNE1310	<input type="checkbox"/>	Attachment for positioning belt	€35

	<b>TOTAL PRICE</b>	
--	--------------------	--

**Non Conforming Options:**

BNE0110	400 mm / 16" Seat Height
BNE0201	Multi-stretch Black PU coated fabric
BNE0202	Black Spacer TR36
BNE1111	No upholstery, select with options 0167 Vicair Adjuster O2 or 0180 Without seat cushion
BNE0235	Matrix Backrest Height 460 mm / 18"
BNE0241	360mm Backrest width, to fit chest width 265mm
BNE0242	380mm Backrest width, to fit chest width 290mm
BNE0243	410mm Backrest width, to fit chest width 315mm
BNE0244	430mm Backrest width, to fit chest width 340mm
BNE0245	460mm Backrest width, to fit chest width 370mm
BNE0246	480mm Backrest width, to fit chest width 390mm
BNE0247	510mm Backrest width, to fit chest width 420mm
BNE7777	No Backrest width, select with Matrix Posture back, option 0225
BNE0335	Push Handles connected
BNE0370	Central Legrest, includes calfpad and footrest. Also select options 0420, 0470, 3333, 4444
BNE0376	High position 50mm higher for leg length compensation
BNE0378	90° Hanger attachment to achieve a 90° legrest position on the Elevating Legrest only
BNE0380	Elevating Legrest, angle adjustable 0 to 75° swing in/swing out with knee pad, calf support with locking clip LH
BNE0390	Legrest Hanger fixed 80 deg, swing in/swing out LH
BNE0400	Legrest Hanger fixed 90 deg, swing in/swing out LH
BNE0410	Legrest Amputee LH
BNE0430	Elevating Legrest, angle adjustable 0 to 75° swing In/swing out with knee pad and calf support with locking clip RH
BNE0440	Legrest Hanger fixed 80 deg, swing in/swing out RH
BNE0450	Legrest Hanger fixed 90 deg, swing in/swing out RH
BNE0460	Legrest Amputee RH
BNE0475	Footrest one piece, depth and angle adjustable short tube
BNE0476	Footrest one piece, depth and angle adjustable long tube
BNE0477	Footbox, angle adjustable, padded with short tube
BNE0480	Footrest lower part, fixed with short tube and heel strap, LH
BNE0490	Footrest lower part, depth and angle adjustable with short tube and heel strap, LH. Not compatible with Elevating legrest.
BNE0520	Footrest lower part, depth and angle adjustable with long tube and heel strap, LH. All knee to heel ranges can be achieved when used with Elevating legrests.
BNE0540	Footrest lower part, fixed with short tube and heel strap, RH
BNE0550	Footrest lower part, depth and angle adjustable with short tube and heel strap, RH. Not compatible with Elevating legrest.
BNE0580	Footrest lower part, depth and angle adjustable with long tube and heel strap, RH. All knee to heel ranges can be achieved when used with Elevating legrests.
BNE0695	Pad for sideguard, removable with hook and loop fastening
BNE0720	Pommel Abduction Cushion
BNE0722	2 cm Increased seat height plate
BNE0750	Headrest with cheek support, sideways adjustable Fixed hardware included
BNE0740	Cover for headrest
BNE0770	Neckrest
BNE0782	Alternative headrest/neckrest mounting bracket
BNE0790	Handles Clinic Kit for Head/Neckrest (Levers which replace screws for easy adjustment)

BNE0928	Calf pad cover, for LH elevating legrest
BNE0929	Knee pad for LH elevating legrest
BNE0931	Calf pad cover, for RH elevating legrest
BNE0932	Knee pad for RH elevating legrest
BNE6666	Without Elevating Legrest Options (Select with Fixed legrests, Amputee legrests and Without legrests only)
BNE0930	Footplate Extender LH
BNE0940	Footplate Extender RH
BNE0941	Cover Footplate, Black Plush Fabric
BNE0942	Calf Strap Padded
BNE0945	Calf Strap (for One piece footrest)
BNE0952	Bodypoint Ankle Huggers Large (Supplied with fixing kit, unmounted)
BNE0954	Bodypoint Ankle Huggers Medium (Supplied with fixing kit, unmounted)
BNE0956	Bodypoint Ankle Huggers Small (Supplied with fixing kit, unmounted)
BNE1015	200x32 mm (8" x 1 1/4") Black Solid Tyre
BNE1025	200x45 mm (8" X 1 3/4") Black Solid tyre
BNE1206	Bodypoint Pivotfit Harness Medium (Supplied with fixing kit, unmounted)
BNE1207	Bodypoint Pivotfit Harness Large (Supplied with fixing kit, unmounted)
BNE1230	Fixed Table Tray without elbow pads includes mounting kit
BNE1231	Fixed Table Tray with elbow Pads includes mounting kit
BNE1232	Half table tray, clear plastic, flip up, depth adjustable, left side (armrest with straight narrow pad included)
BNE1233	Half table tray, clear, plastic, flip up, depth adjustable, right side (armrest with straight narrow pad included)
BNE1240	Swing away table without elbow pads
BNE1245	Swing away table with elbow pads
BNE1235	Table tray padded cushion (To be ordered in combination with the fixed & swingaway table tray)
BNE1237	Half table tray padded cushion
BNE1250	Locking device for table tray
BNE1260	Mounting kit for table (if ordered separately)
BNE1290	Handles Clinic Kit for backrest, armrest, legrest and seat frame (Levers which replace screws for easy adjustment)



For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):