



**Yes, you can:**

Tel: 01656 776222  
 Email: ordersuk@invacare.com  
 Online Spares available at: www.invacare.co.uk

# Manual Retail Prescription Form

JANUARY 2025  
 Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## Rea Azalea Transit Stock

### Standard Features

- Max user weight: 135kg (21.5 stone)
- Crash tested ISO 7176-Part 19
- Transit Propulsion 16" Black solid tyre, attendant hub brakes with +50mm pre-width adjustment
- Small, medium or large chassis
- Seat height 450mm/18"
- Depth adjustable seat 430mm-500mm/17"-20"
- Flex 3 Backrest Plate
- Multi-stretch Black PU coated Laguna backrest cushion
- Multi-stretch Black PU coated Flo-shape seat cushion
- Connected push handles
- Assistant operated gas piston seat and backrest recline with tilt lock
- Removable, angle adjustable legrests with calf pads and angle adjustable footplates
- Height adjustable wide pad armrest (including autolock)
- Sideways adjustable headrest
- Trunk supports, multi-adjustable extended height
- 200 X 45mm Black medium hard front castors
- User manoeuvred tyre brakes
- Pelvic belt, seat belt attachment bracket
- Front transportation brackets
- Anti-tippers
- Silver grey frame



Please make your selection with a cross in the box     Free of charge     standard feature    POA special request

Customer weight =                      KG

Code	Description	
<b>Rea Azalea Transit with +5cm width adjustment</b>		
1587805	W 39cm x D 43-50cm	£2,393
1587806	W 44cm x D 43-50cm	£2,393
1587807	W 49cm x D 43-50cm	£2,393

<b>TOTAL PRICE</b>	
--------------------	--

**Non Conforming Options:**



For a more complex customer requirement  
please refer to our Invacare Unique special order form.  
email [specialsuk@invacare.com](mailto:specialsuk@invacare.com)

**Special Instructions**

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_

Description of special request (please include relevant option code or part number):  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_