



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

07 JANUARY 2026
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Rea Azalea Transit Tall

Standard Features

- Max User Weight 135KG/21 Stone
- Crash Tested ISO 7176-Part 19
- 16" Transit Propulsion with Attendant Hub Brake and Solid Black Tyres
- Choice of Seat Widths 390mm, 440mm or 490mm (+50mm or +100mm Width Adj)
- Seat Height 500mm
- Depth Adjustable Seat 430mm-500mm
- Without Seat Cushion, Supplied with Seat Board
- Flex 3 Backrest Plate, Width and Height Adjustable
- Passad 2 Backrest Cushion with Multi-stretch PU coated fabric Cover and Push Handles Connected
- Assistant Operated Gas Piston Seat Tilt +1° - +24°
- Backrest Recline 30° with Tilt Lock
- Angle Adjustable Legrests with Calf Support, Height and Depth Adjustable Footplates
- Height Adjustable Wide Pad Armrests with Autolock
- Headrest, Sideway Adjustable
- Pearl Silver Frame Colour
- Chassis Extender with Solid Black Front Castors
- User Manoeuvred Tyre Brakes
- Pelvic Belt, Adjustable both sides
- Anti Tippers



Standard Base price: €3,146

Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

Code	Description	
BASE PRODUCT*		
BNG9004	Rea Azalea Transit Tall	€3,146
BNG8000	SPECIAL REQUEST - PLEASE INCLUDE QUOTATION REFERENCE IN SPECIAL DESCRIPTION AT THE END OF THE FORM WHERE APPLICABLE	POA
PROPULSION*		
BNG0020	Transit 16" wheels with hub brakes	◆
SEAT DIMENSIONS*		
BNG0036	Tall	◆

SEAT WIDTH*				
BNG0040		390 mm/15" (Small chassis)	(NC 0244, 0245, 0246, 0247)	○
BNG0050		440 mm/17" (Medium chassis)	(NC 0241, 0242, 0246, 0247)	○
BNG0060		490 mm/19" (Large chassis)	(NC 0206, 0241, 0242, 0243, 0244)	○
SEAT WIDTH EXTENSION*				
BNG0065		+ 50 mm pre-adjustment on armrests and legrests		○
BNG0075		+ 100 mm pre-adjustment on armrests and legrests		○
BNG0085		Without seat width extension	(NC 0206)	○
SEAT DEPTH*				
BNG0105		480-550 mm / 19" -22" (Azalea Tall - Incl Chassis Extension)		◆
SEAT HEIGHT*				
BNG0125		500 mm / 20" (Azalea Tall)		◆
SEAT ANGLE				
BNG0127		Negative Seat Angle -1 To +23° (For attendant operated gas piston in negative tilt)	(NC 0360)	○
SEAT TYPE*				
BNG0130		Depth Adjustable		◆
SEAT CUSHION*				
BNG0165		Cushion Floshape seat upholstery	(NC 1111, 0182, 0183, 0186, 0187, 0188, 0189)	€72
BNG0167		Cushion Vicair Adjuster O2 10cm, complete with Black Comfair cover	(NC 0182, 0183, 0186, 0187, 0188, 0189, 0201, 0202)	€696
BNG0168		Cushion Matrx Libra, complete with Black cover	(NC 0201, 0202)	€645
BNG0169		Cushion Flo-Tech Solution Xtra, complete with Black cover	(NC 0201, 0202)	€566
BNG0179		Cushion Vicair Vector O2 100mm with 2 top covers (Additional 15 day leadtime)	(NC 0182, 0183, 0186, 0187, 0188, 0189, 0201, 0202)	€800
BNG0180		Without Seat Cushion, Seat Plate only	(NC 0201, 0202)	○
SEAT DEPTH OF MATRX LIBRA AND FLO-TECH SOLUTION XTRA CUSHIONS (OPTIMISED CUSHION DEPTH TO SUIT SEAT DEPTH SETTING)				
BNG0182		Short cushion depth: Libra delivered at 405mm deep, Solution Xtra delivered at 430mm deep		○
BNG0183		Long cushion depth: Libra delivered at 455mm deep, Solution Xtra delivered at 455mm deep		○
MATRX LIBRA/FLOTECH SOLUTION XTRA CUSHION ACCESSORIES				
BNG0186		Tapered Wedge Pair		€28
BNG0187		Abductor		€58
BNG0188		Lateral Side Wedge		€46
BNG0189		Anterior Front Wedge		€43
SEAT UPHOLSTERY*				
BNG0201		Multi-stretch Black PU coated fabric		◆
BNG0202		Black Spacer TR36		○
BNG1111		No upholstery, select with options 0167 Vicair Adjuster O2 or 0180 Without seat cushion		○
BACKREST WIDTH EXTENSION*				
BNG0205		+ 50 mm Backrest width pre-adjustment, set at +50mm but can be adjusted if required	(NC 0221, 0223, 0225)	○
BNG0206		+ 100 mm Backrest width pre-adjustment, set at +100mm but can be adjusted if required	(NC 0220, 0221, 0223, 0225)	○
BNG0207		Without Backrest width pre-adjustment, set at 0mm but can be adjusted if required		○

BACKREST*			
BNG0210		Flex 3 Backrest, Fixed Plate (Width Adjustable)	(NC 0233, 0234, 0235, 0240, 0241, 0242, 0243, 0244, 0245, 0246, 0247, 0343) ◆
BNG0219		New!! Flex 3+ Backrest features tension adjustment, width adjustment, and fore/aft positioning of the vertical backpost tubes to accommodate a wider range of clinical requirements.	(NC 0233, 0235, 0234, 0241, 0242, 0243, 0244, 0245, 0246, 0247, 0277, 0343) €289
BNG0220		Flex 3 Backrest features tension adjustment and width adjustment of the vertical backpost tubes.	(NC 0233, 0235, 0234, 0241, 0242, 0243, 0244, 0245, 0246, 0247, 0277, 0343) €126
BNG0221		Matrx Elite E2 (includes hardware mounted to frame, backrest supplied in box)	(NC 0240, 0250, 0260, 0270, 0305, 0321, 0322, 0800, 0810, 0820, 0830, 0880, 0890, 0900, 0920) €668
BNG0223		Matrx Elite E2 Deep (includes hardware mounted to frame, backrest supplied in box)	(NC 0240, 0250, 0260, 0270, 0305, 0321, 0322, 0800, 0810, 0820, 0830, 0880, 0890, 0900, 0920) €739
BNG0225		Matrx Posture Back (includes hardware mounted to frame, backrest supplied in box)	(NC 0235, 0240, 0241, 0242, 0243, 0244, 0245, 0246, 0247, 0250, 0260, 0270, 0305, 0321, 0322, 0335, 0800, 0810, 0820, 0830, 0880, 0890, 0900, 0920) €1,038
BNG0230		Backrest Tubes Without Backrest	(NC 0233, 0234, 0235, 0240, 0241, 0242, 0243, 0244, 0245, 0246, 0247, 0250, 0260, 0270, 0305, 0321, 0322, 0345, 0740, 0800, 0810, 0820, 0830, 0880, 0890, 0900, 0920) ○
MATRIX BACKREST HEIGHT			
BNG0233		Matrx Backrest Height 410 mm / 16"	○
BNG0235		Matrx Backrest Height 460 mm / 18"	○
BNG0234		Matrx Backrest Height 510 mm / 20"	○
MATRIX BACKREST WIDTH			
BNG0241		Matrx Backrest Width 360mm / 14" (Chest width 265mm / 10 1/2")	○
BNG0242		Matrx Backrest Width 380mm / 15" (Chest width 290mm / 11 1/2")	○
BNG0243		Matrx Backrest Width 410mm / 16" (Chest width 315mm / 12 1/2")	○
BNG0244		Matrx Backrest Width 430mm / 17" (Chest width 340mm / 13")	○
BNG0245		Matrx Backrest Width 460mm / 18" (Chest width 370mm / 14 1/2")	○
BNG0246		Matrx Backrest Width 480mm / 19" (Chest width 390mm / 15")	○
BNG0247		Matrx Backrest Width 510mm / 20" (Chest width 420mm / 16 1/2")	○
BACKREST CUSHION*			
BNG0240		Cover for tension adjustable upholstery	(NC 2222) ○
BNG0250		Cushion, Laguna (Lateral Support)	(NC 2222) ○
BNG0260		Cushion, Mistral 2 (Waist Support)	(NC 0322, 2222) €179
BNG0270		Cushion, Passad 2 (Shoulder Support)	(NC 2222) ◆
BNG0276		Vicair Multifunctional O2 pressure relieving backrest cushion, supplied with Black mesh cover (Additional 5 day lead time)	(NC 0305, 0321, 0322) €1,000
BNG0277		Vicair Multifunctional O2 Flex pressure relieving backrest cushion, supplied with Black mesh cover (Additional 5 day lead time)	(NC 0305, 0321, 0322) €1,000
BNG0300		Without Backrest Cover/ Cushion	(NC 0305, 0321, 0322) ○
WIDER BACKREST COVER			
BNG0305		Extra large backrest cover or cushion. One size wider backrest cushion or tension cover fitted if selected	○
BACKREST UPHOLSTERY*			
BNG0321		Multi-stretch PU coated fabric	○
BNG0322		Black Spacer TR36	○
BNG2222		No upholstery, select with Matrx, and Without backrest options	○

PUSH HANDLES*			
BNG0335		Push Handles connected	◆
BNG0343		Push Handles connected,with Rea head/neckrest attachment	€96
BNG0345		Push Bar, angle adjustable	€123
SEAT TILT AND BACKREST REGULATION*			
BNG0350		Assistant operated gas piston seat tilt and backrest recline (NC 5555)	◆
BNG0360		Electrical operated seat tilt and backrest recline (Tilt range is -1° To +23°) (NC 1263)	€958
TILT LOCK*			
BNG1263		Tilt Lock (locking of pre-set angle)	◆
BNG5555		No Tilt Lock, Select with 0360 Electric operated seat and backrest	○
CENTRAL LEGREST			
BNG0370		Central Legrest, includes calfpad and footrest. Also select options 0420, 0470, 3333, 4444 (NC 0189, 0376, 0378, 0380, 0390, 0410, 0430, 0440, 0460, 0475, 0476, 0477, 0720, 0930, 0940, 0941, 0942, 0952, 0954, 0946, 1290)	€377
LEGREST HANGER ATTACHMENT*			
BNG0375		Standard position for legrest mounting	◆
SWO0189		Standard position with 0-30° Legrest Abduction function (NC 0370, 0475, 0476, 0477)	€122
BNG0376		High position 50mm higher for leg length compensation	○
BNG0378		90° Hanger attachment to achieve a 90° legrest position on the Elevating Legrest (NC 0390, 0410, 0440, 0460)	○
LEGREST HANGER LH*			
BNG0380		Elevating Legrest, angle adjustable 0 to 75° swing in/swing out with knee pad, calf support with locking clip LH (NC 0475, 0490, 0942, 0945, 6666)	◆
BNG0390		Legrest Hanger fixed 80 deg, swing in/swing out LH (NC 0928, 0929, 0930)	○
BNG0410		Legrest Amputee LH (NC 0475, 0476, 0477, 0490, 0520, 0928, 0929, 0930, 0941, 0942, 0945, 0952, 0954, 0956)	€171
BNG0420		Without Legrest, LH (NC 0475, 0476, 0477, 0490, 0520, 0928, 0929, 0930, 0941, 0942, 0945, 0952, 0954, 0956)	○
LEGREST HANGER RH*			
BNG0430		Elevating Legrest, angle adjustable 0 to 75°swing In/swing out with knee pad and calf support with locking clip RH (NC 0475, 0550, 0942, 0945, 6666)	◆
BNG0440		Legrest Hanger fixed 80 deg, swing in/swing out RH (NC 0931, 0932, 0940)	○
BNG0460		Legrest Amputee RH (NC 0475, 0476, 0477, 0550, 0580, 0931, 0932, 0940, 0941, 0942, 0945, 0952, 0954, 0956)	€171
BNG0470		Without Legrest, RH (NC 0475, 0476, 0477, 0550, 0580, 0931, 0932, 0940, 0941, 0942, 0945, 0952, 0954, 0956)	○
FOOTREST LH*			
BNG0475		Footrest one piece, depth and angle adjustable short tube (NC 0477, 0490, 0520, 0550, 0580, 0930, 0940, 0952, 0954, 0956)	€145
BNG0476		Footrest one piece, depth and angle adjustable long tube (NC 0477, 0490, 0520, 0550, 0580, 0930, 0940, 0952, 0954, 0956)	€145
BNG0477		Footbox, angle adjustable, padded with short tube (NC 0490, 0520, 0550, 0580, 0930, 0940, 0941, 0942, 0945, 0952, 0954, 0956)	€688
BNG0490		Footrest lower part, depth and angle adjustable with short tube and heel strap, LH. Not compatible with Elevating legrest.	○
BNG0520		Footrest lower part, depth and angle adjustable with long tube and heel strap, LH. All knee to heel ranges can be achieved when used with Elevating legrests.	○
BNG3333		Without footrest, only select with option 0370 Central legrest, 0410 Amputee legrest and 0420 Without legrest	○

FOOTREST RH*			
BNG0550		Footrest lower part, depth and angle adjustable with short tube and heel strap, RH. Not compatible with Elevating legrest.	○
BNG0580		Footrest lower part, depth and angle adjustable with long tube and heel strap, RH. All knee to heel ranges can be achieved when used with Elevating legrests.	○
BNG4444		Without footrest, only select with options 0370 Central legrest, 0460 Amputee legrest, 0470 Without legrest, 0475 and 0476 One piece footrests or 0477 Footbox	○
ARMREST LH*			
BNG0600		Armrest Ducktail Pad LH, height adjustable with autolock (NC 1232)	○
BNG0605		Armrest Wide Pad LH, height adjustable with autolock. Compatible with 16" wheels only (NC 1232)	◆
BNG0607		Armrests Reinforced Wide Pad, LH, height adjustable with autolock (NC 1230, 1231, 1232, 1250, 1260)	€112
BNG0610		Armrest Padded LH, height adjustable with autolock (NC 0680, 0695, 1230, 1231, 1232, 1250, 1260, 9999)	€79
BNG0615		Armrest Hemiplegic LH, height adjustable with autolock (NC 0695, 1230, 1231, 1232, 1240, 1245, 1250, 1260)	€222
BNG0620		Without Armrest LH Select when choosing BNG1232 Half Tray LH, Narrow Armrest is included (NC 0695, 1230, 1231, 1240, 1245, 1250, 1260)	○
ARMREST RH*			
BNG0630		Armrest Ducktail Pad RH, height adjustable with autolock (NC 1233)	○
BNG0635		Armrest Wide Pad RH, height adjustable with autolock. Compatible with 16" wheels only (NC 1233)	◆
BNG0637		Armrest Reinforced Wide Pad, RH, height adjustable with autolock (NC 1230, 1231, 1233, 1250, 1260)	€112
BNG0640		Armrest Padded RH, height adjustable with autolock. (NC 0680, 0695, 1230, 1231, 1233, 1250, 1260, 9999)	€79
BNG0645		Armrest Hemiplegic, RH, height adjustable with autolock (NC 0695, 1230, 1231, 1233, 1240, 1245, 1250, 1260)	€222
BNG0650		Without Armrest RH Select when choosing BNG1233 Half Tray RH, Narrow Armrest is included (NC 0695, 1230, 1231, 1240, 1245, 1250, 1260)	○
ARMREST OPTIONS			
BNG0680		Side Support Cover, removable with tools (wipeable cover for the sideguard) (NC 0695, 9999)	€63
BNG9999		Side Support Cover with foam cushion. 2cm thick cushion supplied complete with PU fabric cover. Reduces seat width by 2cm for a comfortable fit. Consists of options 0680 and 0690 (Removable with tools) (NC 0695)	€63
BNG0695		Pad for sideguard, removable with hook and loop fastening	€63
SEAT OPTIONS			
BNG0710		Incontinence Cover	€43
BNG0720		Pommel Abduction Cushion	€145
BNG0725		Tilt Scale for seat tilt	€53

HEADREST OPTIONS*			
BNG0730		Headrest, Sideway Adjustable (NC 0742, 0744, 0746, 0750, 0770, 0782)	◆
BNG0750		Headrest with cheek support, sideway adjustable Fixed hardware included (NC 0740, 0742, 0744, 0746, 0740, 0770, 0782)	€135
BNG0731		Matrx Elan Standard Pad (254 mm / 10") (NC 0740, 0770, 0782, 0790, 8888)	€189
BNG0732		Matrx Elan Standard Pad (355 mm / 14") (NC 0740, 0770, 0782, 0790, 8888)	€189
BNG0733		Matrx Elan Occipital Pad (229 mm / 9") (NC 0740, 0770, 0782, 0790, 8888)	€189
BNG0734		Matrx Elan Occipital Pad (305 mm / 12") (NC 0740, 0770, 0782, 0790, 8888)	€189
BNG0735		Matrx Elan 4-Point Pad (229W x 177H mm / 9"W x 7"H) (NC 0740, 0770, 0782, 0790, 8888)	€357
BNG0736		Matrx Elan 4-Point Pad (280W x 254H mm / 11"W x 10"H) (NC 0740, 0770, 0782, 0790, 8888)	€357
BNG9991		Without headrest (NC 0742, 0744, 0746)	○
HEADREST HARDWARE*			
BNG0742		Matrx Elan Hardware (NC 0740, 0770, 0790, 1206, 1207, 1300)	€285
BNG0744		Matrx Elan Extended Hardware (Extra link) (NC 0740, 0770, 0790, 1206, 1207, 1300)	€340
BNG0746		Matrx Elan LOXX Hardware (NC 0740, 0770, 0790, 1206, 1207, 1300)	€285
BNG8888		Without Elan headrest hardware (Only select with option 0730 & 0750 Headrests)	○
OTHER HEADREST OPTIONS			
BNG0740		Cover for headrest	€37
BNG0770		Neckrest	○
BNG9992		Without neckrest	○
BNG0782		Alternative headrest/neckrest mounting bracket	€24
BNG0790		Handles Clinic Kit for Head-/Neckrest (Levers which replace screws for easy adjustment)	€61
TRUNK SUPPORT LH			
BNG0800		Trunk Support Multi Adjustable with standard arm LH (NC 0900)	○
BNG0820		Trunk Support Multi Adjustable with extended height arm LH (NC 0900)	○
BNG0880		NEW!! Trunk Support Swingaway LH, multi adjustable (NC 0900)	€128
TRUNK SUPPORT RH			
BNG0810		Trunk Support Multi Adjustable with standard arm RH (NC 0900)	○
BNG0830		Trunk Support Multi Adjustable with extended height arm RH (NC 0900)	○
BNG0890		NEW!! Trunk Support Swingaway RH, multi adjustable (NC 0900)	€128
BACKREST OPTIONS			
BNG0900		Trunk Support Attachment to mount laterals later. If trunk supports are chosen, the attachment bracket is already included	◆
BNG0915		Tilt Scale for backrest recline	€53
BNG0920		Wedges (One Pair)	€14
ELEVATING LEGREST OPTIONS*			
BNG0928		Calf pad cover, for LH elevating legrest	○
BNG0929		Knee pad for LH elevating legrest	◆
BNG0931		Calf pad cover, for RH elevating legrest	○
BNG0932		Knee pad for RH elevating legrest	◆
BNG6666		Without Elevating Legrest Options (Select with Fixed legrests, Amputee legrests and Without legrests only)	○

LEGREST OPTIONS			
BNG0930	<input type="checkbox"/>	Footplate Extender LH (NC 0941)	€40
BNG0940	<input type="checkbox"/>	Footplate Extender RH (NC 0941)	€40
BNG0941	<input type="checkbox"/>	Cover Footplate, Black Plush Fabric	€66
BNG0942	<input type="checkbox"/>	Calf Strap Padded	€88
BNG0945	<input type="checkbox"/>	Calf Strap (for One piece footrest)	○
BODYPOINT ANKLE HUGGERS			
BNG0952	<input type="checkbox"/>	Bodypoint Ankle Huggers Large (Supplied with fixing kit, unmounted)	€189
BNG0954	<input type="checkbox"/>	Bodypoint Ankle Huggers Medium (Supplied with fixing kit, unmounted)	€189
BNG0956	<input type="checkbox"/>	Bodypoint Ankle Huggers Small (Supplied with fixing kit, unmounted)	€189
CHASSIS COLOUR*			
BNG0944	<input type="checkbox"/>	Silver Pearl	◆
BNG0946	<input type="checkbox"/>	Azurite (Dark Blue)	○
BNG0947	<input type="checkbox"/>	Sand (Beige)	○
BNG0948	<input type="checkbox"/>	Electric Green	○
BNG0949	<input type="checkbox"/>	Happy Red	○
CHASSIS SAFETY*			
BNG0950	<input type="checkbox"/>	Chassis extender, castors extend 50mm forward for added stability	◆
CASTOR ATTACHMENT*			
BNG0960	<input type="checkbox"/>	Fixed castor attachment	◆
FRONT CASTORS*			
BNG1015	<input type="checkbox"/>	200x32 mm (8" x 1 1/4") Black Solid Tyre	€25
BNG1025	<input type="checkbox"/>	200x45 mm (8" X 1 3/4") Black Solid tyre	◆
REAR WHEEL SIZE*			
BNG1047	<input type="checkbox"/>	16" / 405 mm Transit rear wheel with fixed attachment	◆
REAR WHEEL TYPE*			
BNG1076	<input type="checkbox"/>	Hub Brake 16" Transit with Black Plastic Hub and Rim	○
REAR WHEEL TYRE*			
BNG1092	<input type="checkbox"/>	Solid Black tyres	◆
BRAKE*			
BNG1170	<input type="checkbox"/>	User manoeuvred tyre brakes, includes extended handles	◆
PELVIC BELT*			
BNG1200	<input type="checkbox"/>	Pelvic belt	◆
BNG1201	<input type="checkbox"/>	Bodypoint 2-Point Centre-pull Push Button Medium (Mounted to chair)	€171
BNG1202	<input type="checkbox"/>	Bodypoint 2-Point Centre-pull Push Button Large (Mounted to chair)	€171
BNG1203	<input type="checkbox"/>	Bodypoint 4-Point Centre-pull Push Button Medium (Mounted to chair)	€161
BNG1204	<input type="checkbox"/>	Bodypoint 4-Point Centre-pull Push Button Large (Mounted to chair)	€161
BODYPOINT HARNESS			
BNG1206	<input type="checkbox"/>	Bodypoint Pivotfit Harness Medium (Supplied with fixing kit, unmounted)	€171
BNG1207	<input type="checkbox"/>	Bodypoint Pivotfit Harness Large (Supplied with fixing kit, unmounted)	€171
TRANSPORTATION BRACKETS*			
BNG1225	<input type="checkbox"/>	Transportation Brackets, mounted to the front of the wheelchair	€16

OTHER OPTIONS		
BNG1230	Fixed Table Tray without elbow pads includes mounting kit	(NC 1237, 1260) €173
BNG1231	Fixed Table Tray with elbow Pads includes mounting kit	(NC 1237, 1260) €173
BNG1232	Half table tray, clear plastic, flip up, depth adjustable, left side (armrest with straight narrow pad included)	(NC 1235, 1260) €138
BNG1233	Half table tray, clear, plastic, flip up, depth adjustable, right side (armrest with straight narrow pad included)	(NC 1235, 1260) €138
BNG1240	Swing away table without elbow pads	(NC 1237, 1250, 1260) €242
BNG1245	Swing away table with elbow pads	(NC 1237, 1250, 1260) €242
BNG1235	Table tray padded cushion (To be ordered in combination with the fixed & swingaway table tray)	€63
BNG1237	Half table tray padded cushion	€63
BNG1250	Locking device for table tray	€17
BNG1260	Mounting kit for table (if ordered separately)	€10
BNG1290	Handles Clinic Kit for backrest, armrest, legrest and seat frame (Levers which replace screws for easy adjustment)	€67
BNG1300	Harness attachment , must be selected in combination with Head/Neckrest	€92
BNG1310	Attachment for positioning belt	€39

TOTAL PRICE	
--------------------	--

Non Conforming Options:

BNG0182	Short cushion depth: Libra delivered at 405mm deep, Solution Xtra delivered at 430mm deep
BNG0183	Long cushion depth: Libra delivered at 455mm deep, Solution Xtra delivered at 455mm deep
BNG0186	Tapered Wedge Pair
BNG0187	Abductor
BNG0188	Lateral Side Wedge
BNG0189	Anterior Front Wedge
BNG0201	Multi-stretch Black PU coated fabric
BNG0202	Black Spacer TR36
BNG1111	No upholstery, select with options 0167 Vicair Adjuster O2 or 0180 Without seat cushion
BNG0206	+ 100 mm Backrest width pre-adjustment, set at +100mm but can be adjusted if required
BNG0220	Flex 3 Backrest features tension adjustment and width adjustment of the vertical backpost tubes.
BNG0221	Matrx Elite E2 (includes hardware mounted to frame, backrest supplied in box)
BNG0223	Matrx Elite E2 Deep (includes hardware mounted to frame, backrest supplied in box)
BNG0225	Matrx Posture Back (includes hardware mounted to frame, backrest supplied in box)
BNG0233	Matrx Backrest Height 410 mm / 16"
BNG0235	Matrx Backrest Height 460 mm / 18"
BNG0234	Matrx Backrest Height 510 mm / 20"
BNG0241	Matrx Backrest Width 360mm / 14" (Chest width 265mm / 10 1/2")
BNG0242	Matrx Backrest Width 380mm / 15" (Chest width 290mm / 11 1/2")
BNG0243	Matrx Backrest Width 410mm / 16" (Chest width 315mm / 12 1/2")
BNG0244	Matrx Backrest Width 430mm / 17" (Chest width 340mm / 13")
BNG0245	Matrx Backrest Width 460mm / 18" (Chest width 370mm / 14 1/2")
BNG0246	Matrx Backrest Width 480mm / 19" (Chest width 390mm / 15")
BNG0247	Matrx Backrest Width 510mm / 20" (Chest width 420mm / 16 1/2")
BNG0240	Cover for tension adjustable upholstery
BNG0250	Cushion, Laguna (Lateral Support)
BNG0260	Cushion, Mistral 2 (Waist Support)
BNG0270	Cushion, Passad 2 (Shoulder Support)
BNG0277	Vicair Multifunctional O2 Flex pressure relieving backrest cushion, supplied with Black mesh cover (Additional 5 day lead time)
BNG0305	Extra large backrest cover or cushion. One size wider backrest cushion or tension cover fitted if selected
BNG0321	Multi-stretch PU coated fabric
BNG0322	Black Spacer TR36
BNG2222	No upholstery, select with Matrx, and Without backrest options
BNG0335	Push Handles connected
BNG0343	Push Handles connected,with Rea head/neckrest attachment
BNG0345	Push Bar, angle adjustable
BNG0360	Electrical operated seat tilt and backrest recline (Tilt range is -1° To +23°)
BNG1263	Tilt Lock (locking of pre-set angle)
BNG5555	No Tilt Lock, Select with 0360 Electric operated seat and backrest
BNG0370	Central Legrest, includes calfpad and footrest. Also select options 0420, 0470, 3333, 4444
SWO0189	Standard position with 0-30° Legrest Abduction function
BNG0376	High position 50mm higher for leg length compensation
BNG0378	90° Hanger attachment to achieve a 90° legrest position on the Elevating Legrest

BNG0380	Elevating Legrest, angle adjustable 0 to 75° swing in/swing out with knee pad, calf support with locking clip LH
BNG0390	Legrest Hanger fixed 80 deg, swing in/swing out LH
BNG0410	Legrest Amputee LH
BNG0430	Elevating Legrest, angle adjustable 0 to 75° swing In/swing out with knee pad and calf support with locking clip RH
BNG0440	Legrest Hanger fixed 80 deg, swing in/swing out RH
BNG0460	Legrest Amputee RH
BNG0475	Footrest one piece, depth and angle adjustable short tube
BNG0476	Footrest one piece, depth and angle adjustable long tube
BNG0477	Footbox, angle adjustable, padded with short tube
BNG0490	Footrest lower part, depth and angle adjustable with short tube and heel strap, LH. Not compatible with Elevating legrest.
BNG0520	Footrest lower part, depth and angle adjustable with long tube and heel strap, LH. All knee to heel ranges can be achieved when used with Elevating legrests.
BNG0550	Footrest lower part, depth and angle adjustable with short tube and heel strap, RH. Not compatible with Elevating legrest.
BNG0580	Footrest lower part, depth and angle adjustable with long tube and heel strap, RH. All knee to heel ranges can be achieved when used with Elevating legrests.
BNG0680	Side Support Cover, removable with tools (wipeable cover for the sideguard)
BNG9999	Side Support Cover with foam cushion. 2cm thick cushion supplied complete with PU fabric cover. Reduces seat width by 2cm for a comfortable fit. Consists of options 0680 and 0690 (Removable with tools)
BNG0695	Pad for sideguard, removable with hook and loop fastening
BNG0720	Pommel Abduction Cushion
BNG0750	Headrest with cheek support, sideway adjustable Fixed hardware included
BNG0742	Matrx Elan Hardware
BNG0744	Matrx Elan Extended Hardware (Extra link)
BNG0746	Matrx Elan LOXX Hardware
BNG8888	Without Elan headrest hardware (Only select with option 0730 & 0750 Headrests)
BNG0740	Cover for headrest
BNG0770	Neckrest
BNG0782	Alternative headrest/neckrest mounting bracket
BNG0790	Handles Clinic Kit for Head-/Neckrest (Levers which replace screws for easy adjustment)
BNG0800	Trunk Support Multi Adjustable with standard arm LH
BNG0820	Trunk Support Multi Adjustable with extended height arm LH
BNG0880	NEW!! Trunk Support Swingaway LH, multi adjustable
BNG0810	Trunk Support Multi Adjustable with standard arm RH
BNG0830	Trunk Support Multi Adjustable with extended height arm RH
BNG0890	NEW!! Trunk Support Swingaway RH, multi adjustable
BNG0900	Trunk Support Attachment to mount laterals later. If trunk supports are chosen, the attachment bracket is already included
BNG0920	Wedges (One Pair)
BNG0928	Calf pad cover, for LH elevating legrest
BNG0929	Knee pad for LH elevating legrest
BNG0931	Calf pad cover, for RH elevating legrest
BNG0932	Knee pad for RH elevating legrest
BNG6666	Without Elevating Legrest Options (Select with Fixed legrests, Amputee legrests and Without legrests only)
BNG0930	Footplate Extender LH
BNG0940	Footplate Extender RH
BNG0941	Cover Footplate, Black Plush Fabric
BNG0942	Calf Strap Padded

BNG0945	Calf Strap (for One piece footrest)
BNG0952	Bodypoint Ankle Huggers Large (Supplied with fixing kit, unmounted)
BNG0954	Bodypoint Ankle Huggers Medium (Supplied with fixing kit, unmounted)
BNG0956	Bodypoint Ankle Huggers Small (Supplied with fixing kit, unmounted)
BNG0946	Azurite (Dark Blue)
BNG1206	Bodypoint Pivotfit Harness Medium (Supplied with fixing kit, unmounted)
BNG1207	Bodypoint Pivotfit Harness Large (Supplied with fixing kit, unmounted)
BNG1230	Fixed Table Tray without elbow pads includes mounting kit
BNG1231	Fixed Table Tray with elbow Pads includes mounting kit
BNG1232	Half table tray, clear plastic, flip up, depth adjustable, left side (armrest with straight narrow pad included)
BNG1233	Half table tray, clear, plastic, flip up, depth adjustable, right side (armrest with straight narrow pad included)
BNG1240	Swing away table without elbow pads
BNG1245	Swing away table with elbow pads
BNG1235	Table tray padded cushion (To be ordered in combination with the fixed & swingaway table tray)
BNG1237	Half table tray padded cushion
BNG1250	Locking device for table tray
BNG1260	Mounting kit for table (if ordered separately)
BNG1290	Handles Clinic Kit for backrest, armrest, legrest and seat frame (Levers which replace screws for easy adjustment)
BNG1300	Harness attachment , must be selected in combination with Head/Neckrest



For a more complex customer requirement,
please email: orders@invacare.com.au

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):