



Yes, you can:

Tel: 01656 776222
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

07 JANUARY 2026
 Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Rea Dahlia 30

Standard Features

- Max User Weight 135KG/21 Stone
- Crash Tested ISO 7176-Part 19
- 16" Transit Propulsion with Attendant Hub Brake and Solid Black Tyres
- Or Self Propel 20", 22" or 24" Self Propel Wheels
- Choice of Seat Widths 390mm, 440mm or 490mm (+50mm or +100mm Width Adj)
- Choice of 8 Seat Heights 330mm - 500mm (13" - 20")
- Depth Adjustable Seat 430mm - 500mm (17" - 20")
- Without Seat Cushion, Supplied with Seat Board
- Flex 3 Backrest Plate, Width and Height Adjustable, with Neckrest
- Passad 2 Backrest Cushion with Multi-stretch Black PU coated fabric Cover and Push Handles Connected
- Assistant Operated Gas Piston Seat Tilt +30° and Backrest Recline 30° with Tilt Lock
- Fixed 80° Legrests Swing In Swing Out. Height, Depth and Angle Adjustable Footplates
- Height Adjustable Ducktail Pad Armrests with Autolock
- Black Frame Colour, Solid Black Front Castors, User Manoeuvred Tyre Brakes
- Pelvic Belt, Adj both sides, Anti Tippers, Front Transportation Brackets

Standard Base price: €3,392



Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code	Description		
BASE PRODUCT*			
BVA0010	DAHLIA 30 °		€3,392
BVA8000	SPECIAL REQUEST - PLEASE INCLUDE QUOTATION REFERENCE IN SPECIAL DESCRIPTION AT THE END OF THE FORM WHERE APPLICABLE		POA
PROPULSION*			
BVA0020	Transit 16" wheels with hub brakes	(NC 0140, 0150, 0160, 0190, 1555, 1565, 1575, 1584, 1588)	○
BVA0030	Self Propel 24"/22"/20" wheels with user brakes	(NC 0090, 0180, 1545, 1588)	○
BVA0042	One Arm Drive Right with Quick Release Rear Wheels	(NC 0130, 0140, 0150, 0180, 0190, 1545, 1555, 1584, 1586)	€464
BVA0044	One Arm Drive Left with Quick Release Rear Wheels	(NC 0130, 0140, 0150, 0180, 0190, 1545, 1555, 1584, 1586)	€464

SEAT WIDTH*				
BVA0050		390 mm/15" (Small chassis)	(NC 0463, 0464, 0465, 0466)	○
BVA0060		440 mm/17" (Medium chassis)	(NC 0461, 0465, 0466)	○
BVA0070		490 mm/19" (Large chassis)	(NC 0461, 0462, 0463)	○
SEAT WIDTH EXTENSION*				
BVA0080		+ 50 mm pre-adjustment on armrests and legrests		○
BVA0090		+ 100 mm pre-adjustment on armrests and legrests		○
BVA0100		Without seat width extension, set at 0mm adjustment but can be adjusted wider if required	(NC 0460, 0470)	○
SEAT DEPTH*				
BVA0110		430-500 mm / 17" - 20"		◆
SEAT HEIGHT*				
BVA0130		350 mm/ 14"	(NC 0225, 0637, 0670, 1435, 1455, 1475, 1495, 1515, 1575)	○
BVA0140		380 mm/ 15"	(NC 0200, 0220, 0227, 0230, 0637, 0670, 1435, 1440, 1495, 1515, 1545, 1565, 1575)	○
BVA0150		400 mm/ 16"	(NC 0200, 0230, 0670, 1435, 1440, 1495, 1515, 1545, 1575)	○
BVA0160		430 mm/ 17"	(NC 0200, 0230, 0670, 1435, 1440, 1495, 1515, 1545, 1555)	○
BVA0170		450 mm/ 18"	(NC 0200, 0210, 0225, 0720, 0770, 1435, 1440, 1555, 1565)	○
BVA0180		480 mm/ 19"	(NC 0200, 0210, 0225, 0227, 0720, 0770, 1435, 1440, 1455, 1475, 1555, 1565, 1575, 1725, 1727)	○
BVA0190		500 mm/ 20"	(NC 0200, 0210, 0225, 0227, 0720, 0770, 1435, 1440, 1455, 1475, 1545, 1565, 1575)	○
SEAT TILT ANGLE*				
BVA0200		Reduced Tilt: Positive Seat Angle 10°		○
BVA0210		Reduced Tilt: Positive Seat Angle 20°		○
BVA0220		Full Tilt: Positive Seat Angle 30°		○
BACKREST RECLINE ANGLE*				
BVA0225		Reduced Backrest Recline: Angle 10°		○
BVA0227		Reduced Backrest Recline: Angle 20°		○
BVA0230		Full Backrest Recline: Angle 30°		○
NEGATIVE SEAT ANGLE				
BVA0240		Negative Seat Angle -1 to aid standing transfers		○
SEAT CUSHION*				
BVA0270		Cushion Floshape Black Multi-stretch PU coated fabric	(NC 0324, 0325, 0333, 0334, 0335, 0336, 0333, 0334, 0335, 0336)	€72
SWO0271		Cushion Floshape Black Multi-stretch PU coated fabric seat upholstery Shorter depth 450mm deep. Only use with seat depth in shortest position	(NC 0324, 0325, 0333, 0334, 0335, 0336, 0333, 0334, 0335, 0336)	€72
BVA0272		Cushion Floshape Black Plush	(NC 0324, 0325, 0333, 0334, 0335, 0336, 0333, 0334, 0335, 0336)	○
BVA0274		Cushion Floshape Black Spacer	(NC 0324, 0325, 0333, 0334, 0335, 0336, 0333, 0334, 0335, 0336)	€72
BVA0292		Cushion Matrix Libra, complete with Black cover		€645
BVA0294		Cushion Flo-Tech Solution Xtra, complete with Black cover		€566
BVA0300		Cushion Vicair Adjuster O2 10cm, complete with Black Comfair cover (Additional 5 day lead time)	(NC 0324, 0325, 0333, 0334, 0335, 0336, 0333, 0334, 0335, 0336)	€696
BVA0310		Cushion Vicair Vector O2 100mm with 2 top covers (Additional 15 day leadtime)	(NC 0324, 0325, 0333, 0334, 0335, 0336, 0333, 0334, 0335, 0336)	€800
BVA0320		Without Seat Cushion, Seat Plate only		○
SEAT DEPTH OF MATRIX LIBRA AND FLO-TECH SOLUTION XTRA CUSHIONS (OPTIMISED CUSHION DEPTH TO SUIT SEAT DEPTH SETTING)				
BVA0324		Short cushion depth: Libra delivered at 405mm deep, Solution Xtra delivered at 430mm deep		○
BVA0325		Long cushion depth: Libra delivered at 455mm deep, Solution Xtra delivered at 455mm deep		○

MATRIX LIBRA/FLOTECH SOLUTION XTRA CUSHION ACCESSORIES			
BVA0333		Tapered Wedge Pair	€28
BVA0334		Abductor	€58
BVA0335		Lateral Side Wedge	€46
BVA0336		Anterior Front Wedge	€43
BACKREST TUBES*			
BVA0360		Straight Back Tubes (Flex 3 and Matrix)	◆
BACKREST RECLINE PIVOT POINT*			
BVA0370		Plus 50mm width extension, outrigs backrest pivot points by 5cm overall to allow a wider backrest to be fitted to a smaller chassis	(NC 0390, 0400, 0470) €119
BVA1111		Without Backrest Pivot Point Width Extension	◆
BACKREST*			
BVA0390		Matrix Elite E2 (includes hardware mounted to frame, backrest supplied in box)	(NC 0460, 0470, 0490, 0500, 0510, 0520, 0544, 0546, 0560, 0570, 1140, 1160, 1220, 1150, 1170, 1230, 1260, 1280, 1860) €668
BVA0400		Matrix Elite E2 Deep (includes hardware mounted to frame, backrest supplied in box)	(NC 0460, 0470, 0490, 0500, 0510, 0520, 0544, 0546, 0560, 0570, 1140, 1160, 1220, 1150, 1170, 1230, 1260, 1280, 1860) €739
BVA0430		Flex 3 Backrest, Fixed Plate (Width Adjustable)	(NC 0455, 0456, 0457, 0461, 0462, 0463, 0464, 0465, 0466, 0490) ◆
BVA0440		Flex 3 Backrest features tension adjustment and width adjustment of the vertical backpost tubes.	(NC 0455, 0456, 0457, 0461, 0462, 0463, 0464, 0465, 0466, 0470, 0546) €126
BVA0445		New!! Flex 3+ Backrest features tension adjustment, width adjustment, and fore/aft positioning of the vertical backpost tubes to accommodate a wider range of clinical requirements.	(NC 0455, 0456, 0457, 0461, 0462, 0463, 0464, 0465, 0466, 0470, 0546) €289
BVA0450		Backrest Tubes Without Backrest	(NC 0455, 0456, 0457, 0461, 0462, 0463, 0464, 0465, 0466, 0460, 0470, 0490, 0500, 0510, 0520, 0544, 0546, 1140, 1160, 1220, 1150, 1170, 1230, 1260, 1280) ○
MATRIX BACKREST HEIGHT			
BVA0455		Matrix Backrest Height 410 mm / 16"	○
BVA0456		Matrix Backrest Height 460 mm / 18"	○
BVA0457		Matrix Backrest Height 510 mm / 20"	○
MATRIX BACKREST WIDTH			
BVA0461		Matrix Backrest Width 380mm / 15" (Chest width 290mm / 11 1/2")	○
BVA0462		Matrix Backrest Width 410mm / 16" (Chest width 315mm / 12 1/2")	○
BVA0463		Matrix Backrest Width 430mm / 17" (Chest width 340mm / 13")	○
BVA0464		Matrix Backrest Width 460mm / 18" (Chest width 370mm / 14 1/2")	○
BVA0465		Matrix Backrest Width 480mm / 19" (Chest width 390mm / 15")	○
BVA0466		Matrix Backrest Width 510mm / 20" (Chest width 420mm / 16 1/2")	○
BACKREST WIDTH EXTENSION*			
BVA0460		+ 50 mm Backrest width pre-adjustment, set at +50mm but can be adjusted if required	○
BVA0470		+ 100 mm Backrest width pre-adjustment, set at +100mm but can be adjusted if required	(NC 1725, 1727) ○
BVA0480		Without Backrest width pre-adjustment, set at 0mm but can be adjusted if required	○

BACKREST CUSHION*			
BVA0490		Cover for tension adjustable upholstery	○
BVA0500		Cushion, Laguna (Lateral Support) (NC 2222)	○
BVA0510		Cushion, Mistral 2 (Waist Support) (NC 0575, 2222)	€179
BVA0520		Cushion, Passad 2 (Shoulder Support) (NC 2222)	◆
BVA0544		Vicair Multifunctional O2 pressure relieving backrest cushion, supplied with Black mesh cover (Additional 5 day lead time) (NC 0560, 0570, 0575)	€1,000
BVA0546		Vicair Multifunctional O2 Flex pressure relieving backrest cushion, supplied with Black mesh cover (Additional 5 day lead time) (NC 0560, 0570, 0575)	€1,000
BVA0550		Without Backrest Cover/ Cushion (NC 0560, 0570, 0575)	○
BACKREST UPHOLSTERY*			
BVA0560		Black Plush TR35	○
BVA0570		Black Multi-stretch Black PU coated fabric	◆
BVA0575		Black Spacer TR36	○
BVA2222		No upholstery, select with Matrx and Without backrest options	○
PUSH HANDLES*			
BVA0580		Push Bar, Angle Adjustable	€123
BVA0600		Push Handles connected	◆
SEAT TILT*			
BVA0620		Assistant operated gas piston seat tilt (NC 3333)	◆
BVA0632		User operated seat tilt, with height adjustable armrest. Limited to 10 degrees of tilt (NC 0915, 0920, 0930, 0975, 0980, 0990, 1750, 1755, 1756, 1757, 1760, 1765, 1770, 1775, 1780, 1790, 1850, 3333)	€119
BVA0637		Electrical operated seat and backrest tilt (NC 0640, 0670)	€958
BACKREST REGULATION*			
BVA0640		Assistant Operated gas piston for backrest recline	◆
BVA3333		Only select if option 0637 Electrical operated is required	○
CENTRAL LEGREST			
BVA0670		Central Legrest, includes calfpad and footrest. Also select options 0740, 0790, 4444, 5555 (NC 0681, 0690, 0692, 0700, 0710, 0720, 0730, 0750, 0760, 0770, 0780, 0805, 0855, 0865, 1040, 1310, 1320, 1330, 1340, 1350, 1351, 1352, 1353)	€377
LEGREST HANGER ATTACHMENT*			
BVA0680		Standard position for legrest mounting	◆
SWO0681		Standard position with 0-30° Legrest Abduction function (NC 0670, 0805)	€110
BVA0690		High position 50mm higher for leg length compensation	○
BVA0692		90° Hanger attachment to achieve a 90° legrest position on the Elevating Legrest (NC 0710, 0720, 0730, 0760, 0770, 0780)	○
LEGREST HANGER LH*			
BVA0700		Elevating Legrest, angle adjustable 0 to 75° swing in/swing out with knee pad, calf support with locking clip LH (NC 0770, 0855, 1330, 8888)	€139
BVA0710		Legrest Hanger fixed 80 deg, swing in/swing out LH (NC 1302, 1304, 1310)	◆
BVA0720		Legrest Hanger fixed 90 deg, swing in/swing out LH (NC 0750, 0865, 1302, 1304, 1310, 1495, 1515)	○
BVA0730		Legrest Amputee LH (NC 0805, 0825, 0855, 0865, 1302, 1304, 1310, 1330, 1340, 1350)	€171
BVA0740		Without Legrest, LH (NC 0805, 0825, 0855, 0865, 1302, 1304, 1310, 1330, 1340, 1350)	○
LEGREST HANGER RH*			
BVA0750		Elevating Legrest, angle adjustable 0 to 75° swing In/swing out with knee pad and calf support with locking clip RH (NC 0855, 1330, 8888)	€139
BVA0760		Legrest Hanger fixed 80 deg, swing in/swing out RH (NC 1312, 1314, 1320)	◆
BVA0770		Legrest Hanger fixed 90 deg, swing in/swing out RH (NC 0865, 1312, 1314, 1320, 1495, 1515)	○
BVA0780		Legrest Amputee RH (NC 0805, 0845, 0855, 0865, 1312, 1314, 1320, 1330, 1340, 1350)	€171
BVA0790		Without Legrest, RH (NC 0805, 0845, 0855, 0865, 1312, 1314, 1320, 1330, 1340, 1350)	○

FOOTREST LH*			
BVA0805		Footrest one piece, depth and angle adjustable short tube	(NC 0845, 1310, 1320, 1351, 1352, 1353) €148
BVA0825		Footrest lower part, depth and angle adjustable with heel strap, LH	(NC 1330) ○
BVA4444		Without footrest, select with options Central legrest, Amputee legrest or Without legrest	○
FOOTREST RH*			
BVA0845		Footrest lower part, depth and angle adjustable with heel strap, RH	(NC 1330) ○
BVA5555		Without footrest, select with options Central legrest, Amputee legrest, Without legrest or One piece footrest	○
FOOTREST TUBE LENGTH*			
BVA0855		Short Tube 24CM	○
BVA0865		Long Tube 36CM	○
BVA6666		No Tube length. Select with Central Legrest, Amputee legrest or Without legrest	○
ARMREST LH*			
BVA0900		Armrest Black PVC Moulded Narrow Pad LH, Height Adjustable from 200 mm - 320 mm with Autolock. Adjustment taken from seat plate to top of armpad	(NC 1756) ○
BVA0910		Armrest Ducktail Pad LH, Height Adjustable from 200 mm - 320mm with Autolock	(NC 1756) ◆
BVA0915		Armrest Extended Height, Ducktail pad LH, Adjusts from 255 mm - 355 mm with Autolock	(NC 1756) ○
BVA0920		Armrest Wide Pad LH, Height Adjustable from 200 mm - 320mm with autolock.	(NC 1756) ○
BVA0930		Armrest Hemiplegic LH, Height Adjustable from 200 mm - 320 mm with autolock	(NC 1750, 1755, 1756, 1760, 1765, 1780, 1790) €222
ARMREST RH*			
BVA0960		Armrest Black PVC Moulded Narrow Pad RH, Height Adjustable from 200 mm - 320 mm with Autolock. Adjustment taken from seat plate to top of armpad	(NC 1757) ○
BVA0970		Armrest Ducktail Pad RH, Height Adjustable from 200 mm - 320mm with Autolock	(NC 1757) ◆
BVA0975		Armrest Extended Height, Ducktail pad RH, Adjusts from 255 mm - 355 mm with Autolock	(NC 1757) ○
BVA0980		Armrest Wide Pad RH, Height Adjustable from 200 mm - 320mm with autolock.	(NC 1757) ○
BVA0990		Armrest Hemiplegic RH, Height Adjustable from 200 mm - 320 mm with autolock	(NC 1020, 1750, 1755, 1757, 1760, 1765, 1770, 1775, 1780, 1790) €222
ARMREST OPTIONS			
BVA1020		Pad for sideguard, removable with hook and loop fastening	€63
SEAT OPTIONS			
BVA1030		Incontinence Cover	€43
BVA1040		Pommel Abduction Cushion	(NC 1060) €145
BVA1050		Seat Plate Extension Plates. Width, Forward, Rear or Assymetric (Supplied as a kit of 2 metal plates unmounted)	(NC 1060) €60
BVA1060		2 cm Increased seat height plate. Allows an increased 20mm of seat height by adding a foam plate to the seat board with hook and loop	€89
BVA1070		Tilt Scale for seat tilt	€53

HEADREST OPTIONS*			
BVA1080		Headrest, Sideway Adjustable	(NC 1087, 1088, 1089, 1110) ○
BVA1100		Headrest with cheek support, sideway adjustable	(NC 1087, 1088, 1089, 1090, 1110) €135
BVA1081		Matrx Elan with Standard pad (254mm / 10")	(NC 1090, 1110, 9999) €189
BVA1082		Matrx Elan with Standard pad (355mm / 14")	(NC 1090, 1110, 9999) €189
BVA1083		Matrx Elan with Occipital pad (229mm / 9")	(NC 1090, 1110, 9999) €189
BVA1084		Matrx Elan with Occipital pad (305mm / 12")	(NC 1090, 1110, 9999) €189
BVA1085		Matrx Elan with 4-Point pad (229mm W x 177mm H / 9" W x 7" H)	(NC 1090, 1110, 9999) €357
BVA1086		Matrx Elan with 4-Point pad (280mm W x 254mm H / 11" W x 10" H)	(NC 1090, 1110, 9999) €357
BVA9991		Without headrest	(NC 1087, 1088, 1089) ○
MATRX HEADREST HARDWARE*			
BVA1087		Matrx Elan hardware	(NC 1090, 1110, 1130, 1736, 1737, 1860) €285
BVA1088		Matrx Elan extended hardware, (Extra link)	(NC 1090, 1110, 1130, 1736, 1737, 1860) €340
BVA1089		Matrx Elan LOXX Hardware	(NC 1090, 1110, 1130, 1860) €285
BVA9999		Without Elan headrest hardware (Only select with option 1080, 1100 Headrest and without headrest)	○
OTHER HEADREST OPTIONS			
BVA1090		Cover for headrest	(NC 1110) €39
BVA1110		Neckrest	◆
BVA1130		Handles Clinic Kit for Head-/Neckrest (Levers which replace screws for easy adjustment)	€58
BVA9992		Without neckrest	○
TRUNK SUPPORTS LH			
BVA1140		Trunk Support Multi Adjustable with standard arm LH	(NC 1260) ○
BVA1160		Trunk Support Multi Adjustable with extended height arm LH	(NC 1260) ○
BVA1220		NEW!! Trunk Support Swingaway LH, multi adjustable	(NC 1260) €128
TRUNK SUPPORTS RH			
BVA1150		Trunk Support Multi Adjustable with standard arm RH	(NC 1260) ○
BVA1170		Trunk Support Multi Adjustable with extended height arm RH	(NC 1260) ○
BVA1230		NEW!! Trunk Support Swingaway RH, multi adjustable	(NC 1260) €128
BACKREST OPTIONS			
BVA1260		Trunk Support Attachment. If trunk supports are chosen, the attachment bracket is already included	○
BVA1280		Wedges (One Pair)	€14
BVA1300		Tilt Scale for backrest recline	€53
ELEVATING LEGREST OPTIONS*			
BVA1302		Calf pad cover, for LH elevating legrest	○
BVA1304		Knee pad for LH elevating legrest	◆
BVA1312		Calf pad cover, for RH elevating legrest	○
BVA1314		Knee pad for RH elevating legrest	◆
BVA8888		Without Elevating Legrest Options (Select with Fixed legrests, Amputee legrests and Without legrests only)	○
LEGREST OPTIONS			
BVA1310		Footplate Extender LH	(NC 1340) €40
BVA1320		Footplate Extender RH	(NC 1340) €40
BVA1330		Calf Strap (for One piece footrest)	○
BVA1340		Cover Footplate, Black Plush Fabric	€66
BVA1350		Calf Strap Padded	€88

BODYPOINT ANKLE HUGGERS			
BVA1351	<input type="checkbox"/>	Bodypoint Ankle Huggers Large (Supplied with fixing kit, unmounted)	€189
BVA1352	<input type="checkbox"/>	Bodypoint Ankle Huggers Medium (Supplied with fixing kit, unmounted)	€189
BVA1353	<input type="checkbox"/>	Bodypoint Ankle Huggers Small (Supplied with fixing kit, unmounted)	€189
CHASSIS COLOUR*			
BVA1405	<input type="checkbox"/>	Black Frame	◆
CASTOR ATTACHMENT*			
BVA1410	<input type="checkbox"/>	Fixed	◆
FRONT CASTORS*			
BVA1435	<input type="checkbox"/>	100x27 mm (4" x 1") Softroll Black Tyre For Seat Height 330mm	○
BVA1440	<input type="checkbox"/>	125x27 mm (5" x 1") Medium Hard Black Tyre For Seat Height 350 mm	○
BVA1455	<input type="checkbox"/>	140x45 mm (5 1/2" x 1 3/4") Soft Black Tyre For Seat Height 380-450 mm	€25
BVA1475	<input type="checkbox"/>	150X30 mm (6" x 1 1/4") Medium Hard Black Tyre For Seat Height 380-450 mm	○
BVA1495	<input type="checkbox"/>	200x32 mm (8" x 1 1/4") Black Solid Tyre For Seat Height 450 - 500 mm	○
BVA1515	<input type="checkbox"/>	200x45 mm (8" x 1 3/4") Black Solid Tyre For Seat Height 450-500 mm	◆
REAR WHEEL POSITION*			
BVA1534	<input type="checkbox"/>	Standard	◆
REAR WHEEL SIZE*			
BVA1545	<input type="checkbox"/>	16" / 405 mm Transit rear wheel with fixed attachment	(NC 1584, 1588, 1605, 1640, 1650, 1660, 1710, 1725, 1727, 7777) ○
BVA1555	<input type="checkbox"/>	20" / 515 mm Self Propel rear wheel with QR attachment	(NC 1586, 1588, 1625) ○
BVA1565	<input type="checkbox"/>	22" / 560 mm Self Propel rear wheel with QR attachment	(NC 7777) ○
BVA1575	<input type="checkbox"/>	24" / 600 mm Self Propel rear wheel with QR attachment	(NC 7777) ○
REAR WHEEL TYPE*			
BVA1584	<input type="checkbox"/>	Standard self propel with colourless hub, clear anodized rim and Stainless Steel spokes	(NC 1728) ○
BVA1586	<input type="checkbox"/>	Hub Brake Self Propel or Transit 16" Transit with Black Plastic Hub and Rim 22" and 24" / 560 mm and 600 mm Self Propel with colourless hub, clear anodized rim and Stainless Steel spokes	○
BVA1588	<input type="checkbox"/>	One Arm Drive Black hub and Zinc Plated Spokes with OAD Handrim System	(NC 1605, 1650, 1660, 1710, 1728) ○
REAR WHEEL TYRE*			
BVA1605	<input type="checkbox"/>	Pneumatic Profile Black 1" 3/8 tyres	○
BVA1625	<input type="checkbox"/>	Solid Black Tyres 1"3/8	◆
BVA1640	<input type="checkbox"/>	Pneumatic low profile Black 1" tyres	○
HANDRIM			
BVA1650	<input type="checkbox"/>	Aluminium handrim	○
BVA1660	<input type="checkbox"/>	Plastic coated handrim	€61
SPOKEGUARD			
BVA1710	<input type="checkbox"/>	Frost spokeguard	€186

BRAKE*			
BVA1720		User manoeuvred tyre brakes, includes extended handles	◆
BVA1725		One Arm Brake LH	€154
BVA1727		One Arm Brake RH	€154
ATTENDANT HUB BRAKE			
BVA1728		Attendant Hub Brake	○
PELVIC BELT*			
BVA1730		Pelvic belt	◆
BVA1731		Bodypoint 2-Point Centre-pull Push Button Medium (Mounted to chair)	€171
BVA1732		Bodypoint 2-Point Centre-pull Push Button Large (Mounted to chair)	€171
BVA1733		Bodypoint 4-Point Centre-pull Push Button Medium (Mounted to chair)	€161
BVA1734		Bodypoint 4-Point Centre-pull Push Button Large (Mounted to chair)	€161
BODYPOINT HARNESS			
BVA1736		Bodypoint Pivotfit Harness Medium (Supplied with fixing kit, unmounted)	€171
BVA1737		Bodypoint Pivotfit Harness Large (Supplied with fixing kit, unmounted)	€171
OTHER OPTIONS			
BVA1750		Table Tray without elbow pads includes mounting kit (NC 1765, 1790)	€173
BVA1755		Table Tray with elbow Pads includes mounting kit (NC 1765, 1790)	€173
BVA1756		Half table tray, clear plastic, flip up, depth adjustable, left side (armrest with straight narrow pad included) (NC 1760, 1790)	€138
BVA1757		Half table tray, clear, plastic, flip up, depth adjustable, right side (armrest with straight narrow pad included) (NC 1760, 1790)	€138
BVA1770		Swing away table without elbow pads (NC 1765, 1780, 1790)	€242
BVA1775		Swing away table with elbow pads (NC 1765, 1780, 1790)	€242
BVA1760		Table tray padded cushion (To be ordered in combination with the fixed & swingaway table tray) (NC 1790)	€63
BVA1765		Half table tray padded cushion (NC 1790)	€63
BVA1780		Locking device for table tray (NC 1790)	€17
BVA1790		Mounting kit for table (if ordered separately)	€10
BVA1820		Pump	€25
BVA1830		Cane Holder	€62
BVA1840		Clinic Kit	€67
BVA1850		Harness Attachment Bracket Seat (Pair) on seat	€29
BVA1860		Harness Attachment Bracket Upper	€29
BVA1870		Kerb Lever (NC 1900)	€59
BVA1900		Accessories Tray to suit Ventilators etc. Positioned at the rear of the frame between the rear wheels	€372
BVA1930		Restrictor Clip kit, Plastic clips which attach to the gas struts to limit the seat tilt or backrest recline	€8
TRANSPORTATION BRACKETS*			
BVA1890		Transportation Brackets	◆
BVA7777		Without Transportation Brackets, only select with 20" Self Propel wheels	○

TOTAL PRICE	
--------------------	--

Non Conforming Options:

BVA0090	+ 100 mm pre-adjustment on armrests and legrests
BVA0130	350 mm/ 14"
BVA0140	380 mm/ 15"
BVA0150	400 mm/ 16"
BVA0160	430 mm/ 17"
BVA0180	480 mm/ 19"
BVA0190	500 mm/ 20"
BVA0200	Reduced Tilt: Positive Seat Angle 10°
BVA0210	Reduced Tilt: Positive Seat Angle 20°
BVA0220	Full Tilt: Positive Seat Angle 30°
BVA0225	Reduced Backrest Recline: Angle 10°
BVA0227	Reduced Backrest Recline: Angle 20°
BVA0230	Full Backrest Recline: Angle 30°
BVA0324	Short cushion depth: Libra delivered at 405mm deep, Solution Xtra delivered at 430mm deep
BVA0325	Long cushion depth: Libra delivered at 455mm deep, Solution Xtra delivered at 455mm deep
BVA0333	Tapered Wedge Pair
BVA0334	Abductor
BVA0335	Lateral Side Wedge
BVA0336	Anterior Front Wedge
BVA0390	Matrx Elite E2 (includes hardware mounted to frame, backrest supplied in box)
BVA0400	Matrx Elite E2 Deep (includes hardware mounted to frame, backrest supplied in box)
BVA0455	Matrx Backrest Height 410 mm / 16"
BVA0456	Matrx Backrest Height 460 mm / 18"
BVA0457	Matrx Backrest Height 510 mm / 20"
BVA0461	Matrx Backrest Width 380mm / 15" (Chest width 290mm / 11 1/2")
BVA0462	Matrx Backrest Width 410mm / 16" (Chest width 315mm / 12 1/2")
BVA0463	Matrx Backrest Width 430mm / 17" (Chest width 340mm / 13")
BVA0464	Matrx Backrest Width 460mm / 18" (Chest width 370mm / 14 1/2")
BVA0465	Matrx Backrest Width 480mm / 19" (Chest width 390mm / 15")
BVA0466	Matrx Backrest Width 510mm / 20" (Chest width 420mm / 16 1/2")
BVA0460	+ 50 mm Backrest width pre-adjustment, set at +50mm but can be adjusted if required
BVA0470	+ 100 mm Backrest width pre-adjustment, set at +100mm but can be adjusted if required
BVA0490	Cover for tension adjustable upholstery
BVA0500	Cushion, Laguna (Lateral Support)
BVA0510	Cushion, Mistral 2 (Waist Support)
BVA0520	Cushion, Passad 2 (Shoulder Support)
BVA0544	Vicair Multifunctional O2 pressure relieving backrest cushion, supplied with Black mesh cover (Additional 5 day lead time)
BVA0546	Vicair Multifunctional O2 Flex pressure relieving backrest cushion, supplied with Black mesh cover (Additional 5 day lead time)
BVA0560	Black Plush TR35
BVA0570	Black Multi-stretch Black PU coated fabric
BVA0575	Black Spacer TR36
BVA2222	No upholstery, select with Matrx and Without backrest options
BVA0637	Electrical operated seat and backrest tilt

BVA0640	Assistant Operated gas piston for backrest recline
BVA3333	Only select if option 0637 Electrical operated is required
BVA0670	Central Legrest, includes calfpad and footrest. Also select options 0740, 0790, 4444, 5555
SWO0681	Standard position with 0-30° Legrest Abduction function
BVA0690	High position 50mm higher for leg length compensation
BVA0692	90° Hanger attachment to achieve a 90° legrest position on the Elevating Legrest
BVA0700	Elevating Legrest, angle adjustable 0 to 75° swing in/swing out with knee pad, calf support with locking clip LH
BVA0710	Legrest Hanger fixed 80 deg, swing in/swing out LH
BVA0720	Legrest Hanger fixed 90 deg, swing in/swing out LH
BVA0730	Legrest Amputee LH
BVA0750	Elevating Legrest, angle adjustable 0 to 75° swing In/swing out with knee pad and calf support with locking clip RH
BVA0760	Legrest Hanger fixed 80 deg, swing in/swing out RH
BVA0770	Legrest Hanger fixed 90 deg, swing in/swing out RH
BVA0780	Legrest Amputee RH
BVA0805	Footrest one piece, depth and angle adjustable short tube
BVA0825	Footrest lower part, depth and angle adjustable with heel strap, LH
BVA0845	Footrest lower part, depth and angle adjustable with heel strap, RH
BVA0855	Short Tube 24CM
BVA0865	Long Tube 36CM
BVA0915	Armrest Extended Height, Ducktail pad LH, Adjusts from 255 mm - 355 mm with Autolock
BVA0920	Armrest Wide Pad LH, Height Adjustable from 200 mm - 320mm with autolock.
BVA0930	Armrest Hemiplegic LH, Height Adjustable from 200 mm - 320 mm with autolock
BVA0975	Armrest Extended Height, Ducktail pad RH, Adjusts from 255 mm - 355 mm with Autolock
BVA0980	Armrest Wide Pad RH, Height Adjustable from 200 mm - 320mm with autolock.
BVA0990	Armrest Hemiplegic RH, Height Adjustable from 200 mm - 320 mm with autolock
BVA1020	Pad for sideguard, removable with hook and loop fastening
BVA1040	Pommel Abduction Cushion
BVA1060	2 cm Increased seat height plate. Allows an increased 20mm of seat height by adding a foam plate to the seat board with hook and loop
BVA1087	Matrx Elan hardware
BVA1088	Matrx Elan extended hardware, (Extra link)
BVA1089	Matrx Elan LOXX Hardware
BVA9999	Without Elan headrest hardware (Only select with option 1080, 1100 Headrest and without headrest)
BVA1090	Cover for headrest
BVA1110	Neckrest
BVA1130	Handles Clinic Kit for Head-/Neckrest (Levers which replace screws for easy adjustment)
BVA1140	Trunk Support Multi Adjustable with standard arm LH
BVA1160	Trunk Support Multi Adjustable with extended height arm LH
BVA1220	NEW!! Trunk Support Swingaway LH, multi adjustable
BVA1150	Trunk Support Multi Adjustable with standard arm RH
BVA1170	Trunk Support Multi Adjustable with extended height arm RH
BVA1230	NEW!! Trunk Support Swingaway RH, multi adjustable
BVA1260	Trunk Support Attachment. If trunk supports are chosen, the attachment bracket is already included
BVA1280	Wedges (One Pair)
BVA1302	Calf pad cover, for LH elevating legrest

BVA1304	Knee pad for LH elevating legrest
BVA1312	Calf pad cover, for RH elevating legrest
BVA1314	Knee pad for RH elevating legrest
BVA8888	Without Elevating Legrest Options (Select with Fixed legrests, Amputee legrests and Without legrests only)
BVA1310	Footplate Extender LH
BVA1320	Footplate Extender RH
BVA1330	Calf Strap (for One piece footrest)
BVA1340	Cover Footplate, Black Plush Fabric
BVA1350	Calf Strap Padded
BVA1351	Bodypoint Ankle Huggers Large (Supplied with fixing kit, unmounted)
BVA1352	Bodypoint Ankle Huggers Medium (Supplied with fixing kit, unmounted)
BVA1353	Bodypoint Ankle Huggers Small (Supplied with fixing kit, unmounted)
BVA1435	100x27 mm (4" x 1") Softroll Black Tyre For Seat Height 330mm
BVA1440	125x27 mm (5" x 1") Medium Hard Black Tyre For Seat Height 350 mm
BVA1455	140x45 mm (5 1/2" x 1 3/4") Soft Black Tyre For Seat Height 380-450 mm
BVA1475	150X30 mm (6" x 1 1/4") Medium Hard Black Tyre For Seat Height 380-450 mm
BVA1495	200x32 mm (8" x 1 1/4") Black Solid Tyre For Seat Height 450 - 500 mm
BVA1515	200x45 mm (8" x 1 3/4") Black Solid Tyre For Seat Height 450-500 mm
BVA1545	16" / 405 mm Transit rear wheel with fixed attachment
BVA1555	20" / 515 mm Self Propel rear wheel with QR attachment
BVA1565	22" / 560 mm Self Propel rear wheel with QR attachment
BVA1575	24" / 600 mm Self Propel rear wheel with QR attachment
BVA1584	Standard self propel with colourless hub, clear anodized rim and Stainless Steel spokes
BVA1586	Hub Brake Self Propel or Transit 16" Transit with Black Plastic Hub and Rim 22" and 24" / 560 mm and 600 mm Self Propel with colourless hub, clear anodized rim and Stainless Steel spokes
BVA1588	One Arm Drive Black hub and Zinc Plated Spokes with OAD Handrim System
BVA1605	Pneumatic Profile Black 1" 3/8 tyres
BVA1625	Solid Black Tyres 1"3/8
BVA1640	Pneumatic low profile Black 1" tyres
BVA1650	Aluminium handrim
BVA1660	Plastic coated handrim
BVA1710	Frost spokeguard
BVA1725	One Arm Brake LH
BVA1727	One Arm Brake RH
BVA1728	Attendant Hub Brake
BVA1736	Bodypoint Pivotfit Harness Medium (Supplied with fixing kit, unmounted)
BVA1737	Bodypoint Pivotfit Harness Large (Supplied with fixing kit, unmounted)
BVA1750	Table Tray without elbow pads includes mounting kit
BVA1755	Table Tray with elbow Pads includes mounting kit
BVA1756	Half table tray, clear plastic, flip up, depth adjustable, left side (armrest with straight narrow pad included)
BVA1757	Half table tray, clear, plastic, flip up, depth adjustable, right side (armrest with straight narrow pad included)
BVA1770	Swing away table without elbow pads
BVA1775	Swing away table with elbow pads
BVA1760	Table tray padded cushion (To be ordered in combination with the fixed & swingaway table tray)

BVA1765	Half table tray padded cushion
BVA1780	Locking device for table tray
BVA1790	Mounting kit for table (if ordered separately)
BVA1850	Harness Attachment Bracket Seat (Pair) on seat
BVA1860	Harness Attachment Bracket Upper
BVA1900	Accessories Tray to suit Ventilators etc. Positioned at the rear of the frame between the rear wheels
BVA7777	Without Transportation Brackets, only select with 20" Self Propel wheels



For a more complex customer requirement,
please email: orders@invacare.com.au

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):