



Yes, you can:

Tel: 01 810 7084 Fax: 01 810 7085
Email: ordersireland@invacare.com
www.invacare.ie

AUGUST 2022

Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Rea Dahlia 45

Standard Features

- Max User Weight 135KG/21 Stone
- Crash Tested ISO 7176-Part 19
- 16" Transit Propulsion with Attendant Hub Brake and Solid Black Tyres
- Or Self Propel 20", 22" or 24" Self Propel Wheels
- Choice of Seat Widths 390mm, 440mm or 490mm (+50mm or +100mm Width Adj)
- Choice of 7 Seat Heights 350mm - 500mm (14" - 20")
- Depth Adjustable Seat 430mm - 500mm (17" - 20")
- Without Seat Cushion, Supplied with Seat Board
- Flex 3 Backrest Plate, Width and Height Adjustable
- Passad 2 Backrest Cushion with Black Multi-stretch PU coated fabric Cover and Push Handles Connected
- Assistant Operated Gas Piston Seat Tilt +30° and Backrest Recline 30° with Tilt Lock
- Fixed 80° Legrests Swing In Swing Out. Height, Depth and Angle Adjustable Footplates
- Height Adjustable Ducktail Pad Armrests with Autolock
- Black Frame Colour, Solid Black Front Castors, User Manoeuvred Tyre Brakes
- Pelvic Belt, Adj both sides, Anti Tippers, Front Transportation Brackets



Standard Base price: €3,503

Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

Code	Description		
BASE PRODUCT*			
BXA0015	INVACARE DAHLIA 45 °		€3,503
BXA8000	SPECIAL REQUEST - PLEASE INCLUDE QUOTATION REFERENCE IN SPECIAL DESCRIPTION AT THE END OF THE FORM WHERE APPLICABLE		POA
PROPULSION*			
BXA0020	Transit 16" wheels with hub brakes	(NC 0140, 0150, 0160, 1515, 1525, 1535, 1544, 1548, 1565, 1595, 1610, 1620, 1670, 1685, 1687)	<input type="radio"/>
BXA0035	Self Propel 24"/22"/20" Wheels	(NC 0090, 0180, 0430, 1505, 1546, 1548)	<input type="radio"/>
BXA0042	One Arm Drive RH with Quick Release Rear Wheels	(NC 0120, 0140, 0180, 0430, 0595, 1505, 1515, 1544, 1546)	€441
BXA0044	One Arm Drive LH with Quick Release Rear Wheels	(NC 0120, 0140, 0180, 0430, 0595, 1505, 1515, 1544, 1546)	€441

SEAT WIDTH*			
BXA0050		390 mm/15" (Small chassis)	○
BXA0060		440 mm/17" (Medium chassis)	○
BXA0070		490 mm/19" (Large chassis)	○
SEAT WIDTH EXTENSION*			
BXA0080		+ 50 mm seat width pre-adjustment on armrests and legrests	○
BXA0090		+ 100 mm seat width pre-adjustment on armrests and legrests	(NC 1685, 1687) ○
BXA0100		Without seat width extension, set at 0mm adjustment but can be adjusted wider if required	(NC 0420, 0430) ○
SEAT DEPTH*			
BXA0110		430-500 mm / 17" - 20"	◆
SEAT HEIGHT*			
BXA0120		350 mm / 14"	(NC 0630, 1400, 1415, 1435, 1455, 1475, 1525, 1535) ○
BXA0130		380 mm / 15"	(NC 0630, 1395, 1415, 1435, 1455, 1475, 1535) ○
BXA0140		400 mm / 16"	(NC 0630, 1395, 1400, 1455, 1475, 1505, 1525, 1535) ○
BXA0150		430 mm / 17"	(NC 0630, 1395, 1400, 1455, 1475, 1505, 1535) ○
BXA0160		450 mm / 18"	(NC 1395, 1400, 1455, 1475, 1505, 1515) ○
BXA0170		480 mm / 19"	(NC 1395, 1400, 1515, 1525) ○
BXA0180		500 mm / 20"	(NC 0680, 0730, 1395, 1400, 1415, 1435, 1515, 1525, 1535, 1685, 1687) ○
SEAT TILT ANGLE*			
BXA0190		Reduced Tilt: Positive Seat Angle 30°	○
BXA0200		Full Tilt: Positive Seat Angle 45°	◆
BACKREST RECLINE ANGLE*			
BXA0203		Reduced Backrest Recline: Angle 10°	○
BXA0205		Reduced Backrest Recline: Angle 20°	○
BXA0207		Full Backrest Recline: Angle 30°	◆
SEAT CUSHION*			
BXA0230		Cushion Flosshape	(NC 1111) ○
SWO0190		Cushion Flosshape seat upholstery Shorter depth 450mm deep. Only use with seat depth in shortest position	(NC 0305, 1111) €68
BXA0260		Cushion Vicair Adjuster O2 10cm (Black)	(NC 0300, 0305) €696
BXA0280		Without Seat Cushion, Seat Plate only	(NC 0300, 0305) ◆
SEAT UPHOLSTERY*			
BXA0300		Black Multi-stretch PU coated fabric	◆
BXA0305		Black Spacer TR36	○
BXA1111		No upholstery, select with Matrix, Vicair Adjuster O2 and Without backrest options	○
BACKREST TUBES*			
BXA0310		Bent Back Tubes (20° from 380 mm above the seatplate - Matrix only)	(NC 0330, 0390, 0400) ○
BXA0320		Straight Back Tubes (Flex 3 and Matrix)	◆
BACKREST RECLINE PIVOT POINT*			
BXA0330		Plus 50mm width extension, outrigs backrest pivot points by 5cm overall to allow a wider backrest to be fitted to a smaller chassis	(NC 0430, 0595) €113
BXA0340		Without Backrest Pivot Point Width Extension	◆

BACKREST*			
BXA0350		Matrx Elite E2 (includes hardware mounted to frame, backrest supplied in box)	(NC 0420, 0430, 0450, 0460, 0470, 0480, 0530, 0535, 1100, 1110, 1120, 1130, 1200, 1210, 1220, 1240) €636
BXA0360		Matrx PB Elite Deep E2 (includes hardware mounted to frame, backrest supplied in box)	(NC 0420, 0430, 0450, 0460, 0470, 0480, 0530, 0535, 1100, 1110, 1120, 1130, 1200, 1210, 1220, 1240) €704
BXA0390		Flex 3 Backrest, Fixed Plate (Width Adjustable)	(NC 0413, 0417, 0450, 0570) ◆
BXA0400		Flex 3 Backrest Tension (Width Adjustable)	(NC 0413, 0417, 0430, 0570) €120
BXA0410		Backrest Tubes Without Backrest	(NC 0413, 0417, 0420, 0430, 0450, 0460, 0470, 0480, 0530, 0535, 1100, 1110, 1120, 1130, 1200, 1210, 1220, 1240) ○
MATRIX BACKREST HEIGHT			
BXA0413		Matrx Backrest Height 410 mm / 16"	○
BXA0417		Matrx Backrest Height 510 mm / 20"	○
BACKREST WIDTH EXTENSION*			
BXA0420		+ 50 mm Backrest width pre-adjustment, set at +50mm but can be adjusted if required	○
BXA0430		+ 100 mm Backrest width pre-adjustment, set at +100mm but can be adjusted if required	(NC 1685, 1687) ○
BXA0440		Without Backrest width pre-adjustment, set at 0mm but can be adjusted if required	○
BACKREST CUSHION*			
BXA0450		Cover for tension adjustable upholstery	(NC 2222) ○
BXA0460		Cushion, Laguna (Lateral Support)	(NC 2222) ○
BXA0470		Cushion, Mistral 2 (Waist Support)	(NC 0535, 2222) €182
BXA0480		Cushion, Passad 2 (Shoulder Support)	(NC 2222) ◆
BXA0510		Without Backrest Cover/ Cushion	(NC 0530, 0535) ○
BACKREST UPHOLSTERY*			
BXA0530		Black Multi-stretch PU coated fabric	◆
BXA0535		Black Spacer TR36	○
BXA2222		No upholstery, select with Matrx, Vicair Adjuster O2 and Without backrest options	○
PUSH HANDLES*			
BXA0540		Push Bar, Angle Adjustable	€117
BXA0560		Push Handles connected	◆
BXA0570		Push Handles connected, with head/neckrest attachment	€91
SEAT TILT*			
BXA0580		Assistant Operated gas piston seat tilt	(NC 5555) ◆
BACKREST REGULATION*			
BXA0600		Assistant Operated gas piston for backrest recline	◆
CENTRAL LEGREST			
BXA0630		Central Legrest, includes calfpad and footrest. Also select options 0700, 0750, 3333, 4444	(NC 0641, 0650, 0655, 0660, 0670, 0680, 0690, 0710, 0720, 0730, 0740, 0760, 0770, 1000, 1270, 1280, 1300, 1310, 1800) €358
LEGREST HANGER ATTACHMENT*			
BXA0640		Standard position for legrest mounting	◆
SWO0641		Standard position with 0-30° Legrest Abduction function	€110
BXA0650		High position 50mm higher for leg length compensation	○
BXA0655		90° Hanger attachment to achieve a 90° legrest position on the Elevating Legrest	(NC 0670, 0680, 0690, 0720, 0730, 0740) ○

LEGREST HANGER LH*			
BXA0660		Elevating Legrest, angle adjustable 0 to 75° swing in/swing out with knee pad, calf support with locking clip LH (NC 0760, 0790, 1290, 6666)	€132
BXA0670		Legrest Hanger fixed 80 deg, swing in/swing out LH (NC 1266, 1267, 1270)	◆
BXA0680		Legrest Hanger fixed 90 deg, swing in/swing out LH (NC 0810, 1266, 1267, 1270, 1455, 1475)	○
BXA0690		Legrest Amputee LH (NC 0760, 0770, 0790, 0810, 1266, 1267, 1270, 1300, 1310)	€162
BXA0700		Without Legrest, LH (NC 0760, 0770, 0790, 0810, 1266, 1267, 1270, 1300, 1310)	○
LEGREST HANGER RH*			
BXA0710		Elevating Legrest, angle adjustable 0 to 75° swing In/swing out with knee pad and calf support with locking clip RH (NC 0760, 0830, 1290, 6666)	€132
BXA0720		Legrest Hanger fixed 80 deg, swing in/swing out RH (NC 1276, 1277, 1280)	◆
BXA0730		Legrest Hanger fixed 90 deg, swing in/swing out RH (NC 0850, 1276, 1277, 1280, 1455, 1475)	○
BXA0740		Legrest Amputee RH (NC 0760, 0770, 0830, 0850, 1276, 1277, 1280, 1300, 1310)	€162
BXA0750		Without Legrest, RH (NC 0760, 0770, 0830, 0850, 1276, 1277, 1280, 1300, 1310)	○
FOOTREST LH*			
BXA0760		Footrest one piece, depth and angle adjustable short tube (NC 0790, 0810, 0830, 0850, 1270, 1280)	€141
BXA0770		Footrest one piece, depth and angle adjustable long tube (NC 0790, 0810, 0830, 0850, 1270, 1280)	€141
BXA0790		Footrest lower part, depth and angle adjustable with short tube and heel strap, LH	○
BXA0810		Footrest lower part, depth and angle adjustable with long tube and heel strap, LH	○
BXA3333		Without footrest, select with options Central legrest, Amputee legrest or Without legrest	○
FOOTREST RH*			
BXA0830		Footrest lower part, depth and angle adjustable with short tube and heel strap, RH	○
BXA0850		Footrest lower part, depth and angle adjustable with long tube and heel strap, RH	○
BXA4444		Without footrest, select with options Central legrest, Amputee legrest, One Piece or Without legrest	○
ARMREST LH*			
BXA0870		Armrest Ducktail Pad LH, Height Adjustable from 200 mm - 320mm with Autolock	◆
BXA0875		Armrest Extended Height, Ducktail pad LH, Adjusts from 255 mm - 355 mm with Autolock	○
BXA0880		Armrest Wide Pad LH, Height Adjustable from 200 mm - 320mm with autolock.	○
BXA0890		Armrest Hemiplegic LH, Height Adjustable from 200 mm - 320 mm with autolock (NC 0980, 1710, 1715, 1720, 1730, 1735, 1740, 1750)	€211
ARMREST RH*			
BXA0930		Armrest Ducktail Pad RH, Height Adjustable from 200 mm - 320mm with Autolock	◆
BXA0935		Armrest Extended Height, Ducktail pad RH, Adjusts from 255 mm - 355 mm with Autolock	○
BXA0940		Armrest Wide Pad RH, Height Adjustable from 200 mm - 320mm with autolock.	○
BXA0950		Armrest Hemiplegic RH, Height Adjustable from 200 mm - 320 mm with autolock (NC 0980, 1710, 1715, 1720, 1730, 1735, 1740, 1750)	€211
ARMREST OPTIONS			
BXA0980		Pad for sideguard, removable with hook and loop fastening	€60

SEAT OPTIONS			
BXA0990		Incontinence Cover	€41
BXA1000		Pommel Abduction Cushion (NC 1020)	€138
BXA1010		Seat Plate Extension Plates. Width, Forward, Rear or Assymmetric (Supplied as a kit of 2 metal plates unmounted)	€57
BXA1020		2 cm Increased seat height plate. Allows an increased 20mm of seat height by adding a foam plate to the seat board with hook and loop	€85
BXA1030		Tilt Scale for seat tilt	€50
BACKREST OPTIONS			
BXA1040		Headrest, Sideway Adjustable (NC 1070)	○
BXA1060		Headrest with cheek support, sideway adjustable (NC 1050, 1070)	€128
BXA1050		Cover for headrest (NC 1070)	€37
BXA1070		Neckrest	◆
BXA1090		Handles Clinic Kit for Head-/Neckrest (Levers which replace screws for easy adjustment)	€58
BXA1220		Trunk Support Attachment. If trunk supports are chosen, the attachment bracket is already included	○
BXA1240		Wedges (One Pair)	€13
BXA1260		Tilt Scale for backrest recline	€50
TRUNK SUPPORTS LH			
BXA1100		Trunk Support Multi Adjustable with standard arm LH (NC 1220)	○
BXA1120		Trunk Support Multi Adjustable with extended height arm LH (NC 1220)	○
BXA1200		Trunk Support Multi Adjustable Swingaway LH (NC 1220)	€132
TRUNK SUPPORTS RH			
BXA1110		Trunk Support Multi Adjustable with standard arm RH (NC 1220)	○
BXA1130		Trunk Support Multi Adjustable with extended height arm RH (NC 1220)	○
BXA1210		Trunk Support Multi Adjustable Swingaway RH (NC 1220)	€132
ELEVATING LEGREST OPTIONS*			
BXA1266		Calf pad cover, for LH elevating legrest	◆
BXA1267		Knee pad for LH elevating legrest	◆
BXA1276		Calf pad cover, for RH elevating legrest	◆
BXA1277		Knee pad for RH elevating legrest	◆
BXA6666		Without Elevating Legrest Options (Select with Fixed legrests, Amputee legrests and Without legrests only)	○
LEGREST OPTIONS			
BXA1270		Footplate Extender LH (NC 1300)	€38
BXA1280		Footplate Extender RH (NC 1300)	€38
BXA1290		Calf Strap (for One piece footrest)	○
BXA1300		Cover Footplate, Black Plush Fabric	€63
BXA1310		Calf Strap Padded	€84
CHASSIS COLOUR*			
BXA1365		Black Frame	◆
CASTOR ATTACHMENT*			
BXA1370		Fixed	◆

FRONT CASTORS*			
BXA1395		100x27 mm (4" x 1") Softroll Black Tyre	○
BXA1400		125x27 mm (5" x 1") Medium Hard Black Tyre	○
BXA1415		140x45 mm (5 1/2" x 1 3/4") Soft Black Tyre	€24
BXA1435		150X30 mm (6" x 1 1/4") Medium Hard Black Tyre	○
BXA1455		200x32 mm (8" x 1 1/4") Black Solid Tyre	○
BXA1475		200x45 mm (8" x 1 3/4") Black Solid Tyre	◆
REAR WHEEL SIZE*			
BXA1505		16" / 405 mm Transit rear wheel with fixed attachment (NC 1544, 1548, 1565, 1595, 1610, 1620, 1670, 1685, 1687)	○
BXA1515		20" / 515 mm Self Propel rear wheel with QR attachment (NC 1546, 1548, 1585)	○
BXA1525		22" / 560 mm Self Propel rear wheel with QR attachment	○
BXA1535		24" / 600 mm Self Propel rear wheel with QR attachment	○
REAR WHEEL TYPE*			
BXA1544		Standard self propel with colourless hub, clear anodized rim and Stainless Steel spokes	○
BXA1546		Hub Brake (NC 1595) 16" Transit with Black Plastic Hub and Rim 22" and 24" / 560 mm and 600 mm Self Propel with colourless hub, clear anodized rim and Stainless Steel spokes	○
BXA1548		One Arm Drive Black hub and Zinc Plated Spokes with OAD Handrim System (NC 1565, 1610, 1620, 1670)	○
REAR WHEEL TYRE*			
BXA1565		Pneumatic Profile Black 1" 3/8 tyres	○
BXA1585		Solid Black Tyres	◆
BXA1595		Pneumatic low profile Black 1" tyres	○
HANDRIM			
BXA1610		Aluminium handrim	◆
BXA1620		Plastic coated handrim	€58
SPOKEGUARD			
BXA1670		Frost spokeguard	€195
BRAKE*			
BXA1680		User manoeuvred tyre brakes	◆
BXA1685		One Arm Brake LH	€148
BXA1687		One Arm Brake RH	€148
PELVIC BELT*			
BXA1690		Pelvic belt	◆

OTHER OPTIONS			
BXA1710		Table Tray without elbow pads includes mounting kit (NC 1750)	€164
BXA1715		Table Tray with elbow Pads includes mounting kit (NC 1750)	€164
BXA1730		Swing away table without elbow pads (NC 1740, 1750)	€230
BXA1735		Swing away table with elbow pads (NC 1740, 1750)	€230
BXA1720		Table tray padded cushion (To be ordered in combination with the fixed & swingaway table tray) (NC 1750)	€60
BXA1740		Locking device for table tray (NC 1750)	€16
BXA1750		Mounting kit for table (if ordered separately)	€9
BXA1800		Clinic Kit	€64
BXA1810		Harness Attachment Bracket Seat (Pair) on seat	€28
BXA1820		Harness Attachment Bracket Upper	€28
BXA1830		Kerb Lever (NC 1860)	€57
BXA1860		Accessories Tray to suit Ventilators etc. Positioned at the rear of the frame between the rear wheels	€353
BXA1890		Restrictor Clip kit, Plastic clips which attach to the gas struts to limit the seat tilt or backrest recline	€7
TRANSPORTATION BRACKETS*			
BXA1850		Transportation Brackets	◆

TOTAL PRICE	
--------------------	--

Non Conforming Options:

BXA0090	+ 100 mm seat width pre-adjustment on armrests and legrests
BXA0120	350 mm / 14"
BXA0140	400 mm / 16"
BXA0150	430 mm / 17"
BXA0160	450 mm / 18"
BXA0180	500 mm / 20"
BXA0300	Black Multi-stretch PU coated fabric
BXA0305	Black Spacer TR36
BXA1111	No upholstery, select with Matrix, Vicair Adjuster O2 and Without backrest options
BXA0330	Plus 50mm width extension, outrigs backrest pivot points by 5cm overall to allow a wider backrest to be fitted to a smaller chassis
BXA0390	Flex 3 Backrest, Fixed Plate (Width Adjustable)
BXA0400	Flex 3 Backrest Tension (Width Adjustable)
BXA0413	Matrix Backrest Height 410 mm / 16"
BXA0417	Matrix Backrest Height 510 mm / 20"
BXA0420	+ 50 mm Backrest width pre-adjustment, set at +50mm but can be adjusted if required
BXA0430	+ 100 mm Backrest width pre-adjustment, set at +100mm but can be adjusted if required
BXA0450	Cover for tension adjustable upholstery
BXA0460	Cushion, Laguna (Lateral Support)
BXA0470	Cushion, Mistral 2 (Waist Support)
BXA0480	Cushion, Passad 2 (Shoulder Support)
BXA0530	Black Multi-stretch PU coated fabric
BXA0535	Black Spacer TR36
BXA2222	No upholstery, select with Matrix, Vicair Adjuster O2 and Without backrest options
BXA0570	Push Handles connected, with head/neckrest attachment
BXA0595	Electrical Operated seat and backrest tilt
BXA0600	Assistant Operated gas piston for backrest recline
BXA5555	Electric backrest regulation, select with option 0595
BXA0630	Central Legrest, includes calfpad and footrest. Also select options 0700, 0750, 3333, 4444
SWO0641	Standard position with 0-30° Legrest Abduction function
BXA0650	High position 50mm higher for leg length compensation
BXA0655	90° Hanger attachment to achieve a 90° legrest position on the Elevating Legrest
BXA0660	Elevating Legrest, angle adjustable 0 to 75° swing in/swing out with knee pad, calf support with locking clip LH
BXA0670	Legrest Hanger fixed 80 deg, swing in/swing out LH
BXA0680	Legrest Hanger fixed 90 deg, swing in/swing out LH
BXA0690	Legrest Amputee LH
BXA0710	Elevating Legrest, angle adjustable 0 to 75° swing In/swing out with knee pad and calf support with locking clip RH
BXA0720	Legrest Hanger fixed 80 deg, swing in/swing out RH
BXA0730	Legrest Hanger fixed 90 deg, swing in/swing out RH
BXA0740	Legrest Amputee RH
BXA0760	Footrest one piece, depth and angle adjustable short tube
BXA0770	Footrest one piece, depth and angle adjustable long tube
BXA0790	Footrest lower part, depth and angle adjustable with short tube and heel strap, LH
BXA0810	Footrest lower part, depth and angle adjustable with long tube and heel strap, LH

BXA0830	Footrest lower part, depth and angle adjustable with short tube and heel strap, RH
BXA0850	Footrest lower part, depth and angle adjustable with long tube and heel strap, RH
BXA0980	Pad for sideguard, removable with hook and loop fastening
BXA1000	Pommel Abduction Cushion
BXA1020	2 cm Increased seat height plate. Allows an increased 20mm of seat height by adding a foam plate to the seat board with hook and loop
BXA1050	Cover for headrest
BXA1070	Neckrest
BXA1220	Trunk Support Attachment. If trunk supports are chosen, the attachment bracket is already included
BXA1240	Wedges (One Pair)
BXA1100	Trunk Support Multi Adjustable with standard arm LH
BXA1120	Trunk Support Multi Adjustable with extended height arm LH
BXA1200	Trunk Support Multi Adjustable Swingaway LH
BXA1110	Trunk Support Multi Adjustable with standard arm RH
BXA1130	Trunk Support Multi Adjustable with extended height arm RH
BXA1210	Trunk Support Multi Adjustable Swingaway RH
BXA1266	Calf pad cover, for LH elevating legrest
BXA1267	Knee pad for LH elevating legrest
BXA1276	Calf pad cover, for RH elevating legrest
BXA1277	Knee pad for RH elevating legrest
BXA6666	Without Elevating Legrest Options (Select with Fixed legrests, Amputee legrests and Without legrests only)
BXA1270	Footplate Extender LH
BXA1280	Footplate Extender RH
BXA1290	Calf Strap (for One piece footrest)
BXA1300	Cover Footplate, Black Plush Fabric
BXA1310	Calf Strap Padded
BXA1395	100x27 mm (4" x 1") Softroll Black Tyre
BXA1400	125x27 mm (5" x 1") Medium Hard Black Tyre
BXA1415	140x45 mm (5 1/2" x 1 3/4") Soft Black Tyre
BXA1435	150X30 mm (6" x 1 1/4") Medium Hard Black Tyre
BXA1455	200x32 mm (8" x 1 1/4") Black Solid Tyre
BXA1475	200x45 mm (8" x 1 3/4") Black Solid Tyre
BXA1505	16" / 405 mm Transit rear wheel with fixed attachment
BXA1515	20" / 515 mm Self Propel rear wheel with QR attachment
BXA1525	22" / 560 mm Self Propel rear wheel with QR attachment
BXA1535	24" / 600 mm Self Propel rear wheel with QR attachment
BXA1544	Standard self propel with colourless hub, clear anodized rim and Stainless Steel spokes
BXA1546	Hub Brake 16" Transit with Black Plastic Hub and Rim 22" and 24" / 560 mm and 600 mm Self Propel with colourless hub, clear anodized rim and Stainless Steel spokes
BXA1548	One Arm Drive Black hub and Zinc Plated Spokes with OAD Handrim System
BXA1565	Pneumatic Profile Black 1" 3/8 tyres
BXA1585	Solid Black Tyres
BXA1595	Pneumatic low profile Black 1" tyres
BXA1610	Aluminium handrim
BXA1620	Plastic coated handrim

BXA1670	Frost spokeguard
BXA1685	One Arm Brake LH
BXA1687	One Arm Brake RH
BXA1710	Table Tray without elbow pads includes mounting kit
BXA1715	Table Tray with elbow Pads includes mounting kit
BXA1730	Swing away table without elbow pads
BXA1735	Swing away table with elbow pads
BXA1720	Table tray padded cushion (To be ordered in combination with the fixed & swingaway table tray)
BXA1740	Locking device for table tray
BXA1750	Mounting kit for table (if ordered separately)
BXA1800	Clinic Kit
BXA1860	Accessories Tray to suit Ventilators etc. Positioned at the rear of the frame between the rear wheels



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):