



Yes, you can:

Tel: 01 810 7084 Fax: 01 810 7085

Email: ordersireland@invacare.com

Online Spares available at: www.invacare.ie

JANUARY 2022

Customer Ref:

Customer Ref :	Account No :
Contact no:	Order no:
Delivery Address:	Order Qty:
	Order Date :

Rea Azalea Transit

Standard Features

Max User Weight 135KG/21 Stone
 Crash Tested ISO 7176-Part 19
 16" Transit Propulsion with Attendant Hub Brake and Solid Black Tyres
 Choice of Seat Widths 390mm, 440mm or 490mm (+50mm or +100mm Width Adj)
 Choice of 2 Seat Heights 400mm or 450mm
 Depth Adjustable Seat 430mm-500mm
 Without Seat Cushion, Supplied with Seat Board
 Flex 3 Backrest Plate, Width and Height Adjustable
 Passad 2 Backrest Cushion with Dartex Cover and Push Handles Connected
 Assistant Operated Gas Piston Seat Tilt +1° - +25°
 Backrest Recline 30° with Tilt Lock
 Angle Adjustable Legrests with Calf Support, Height and Depth Adjustable Footplates
 Height Adjustable Wide Pad Armrests with Autolock
 Pearl Silver Frame Colour, Solid Black Front Castors
 User Manoeuvred Tyre Brakes
 Pelvic Belt, Adjustable both sides
 Anti Tippers
 Front Transportation Brackets



Standard Base Price : €2,668

Please make your selection with a cross in the box ☐ Free Of Charge ☒ Standard Feature ☐ POA ☐ Special Request

Customer weight= **KG**

Code **Description**

BASE PRODUCT *

BNE9002	☐	Rea Azalea Transit (NC)	€ 2,668
BNE8000	☐	SPECIAL REQUEST	POA

PROPULSION *

BNE0020	☐	Transit 16" wheels with hub brakes (NC)	☒
---------	--------------------------	---	----------------------------------

SEAT WIDTH *

BNE0040	☐	390 mm/15" (Small chassis) (NC)	☐
---------	--------------------------	---------------------------------	-----------------------

BNE0050		440 mm/17" (Medium chassis)	(NC)	○
BNE0060		490 mm/19" (Large chassis)	(NC 0206)	○
SEAT WIDTH EXTENSION *				
BNE0065		+ 50 mm pre-adjustment on armrests and legrests		○
BNE0075		+ 100 mm pre-adjustment on armrests and legrests	(NC 0110)	○
BNE0085		Without seat width extension, set at 0mm adjustment but can be adjusted wider if required	(NC 0206)	○
SEAT DEPTH *				
BNE0100		430-500 mm / 17" - 20"	(NC)	◆
SEAT HEIGHT *				
BNE0110		400 mm / 16" Seat Height	(NC 0370, 1015, 1025)	○
BNE0120		450 mm / 18" Seat Height	(NC)	○
SEAT ANGLE *				
BNE0127		Negative Seat Angle -1 To +23°	(NC 0360)	○
BNE0000		Standard Seat Angle +1° to +25°		◆
SEAT TYPE *				
BNE0130		Depth Adjustable		◆
BNE0140		Sliding Seat	(NC 0720)	€173
SEAT CUSHION *				
BNE0165		Cushion Flosshape seat upholstery	(NC 1111)	€65
BNE0170		Cushion Vicair Multifunctional (Black)	(NC 0201, 0202)	€593
BNE0180		Without Seat Cushion, Seat Plate only	(NC 0201, 0202)	○
SEAT UPHOLSTERY *				
BNE0201		Black Dartex TR26		◆
BNE0202		Black Spacer TR36		○
BNE1111		No upholstery, select with options 0170 Multifunctional and 0180 Without seat cushion		○
BACKREST WIDTH EXTENSION *				
BNE0205		+ 50 mm Backrest width pre-adjustment, set at +50mm but can be adjusted if required	(NC 0221, 0223, 0225, 0228)	○
BNE0206		+ 100 mm Backrest width pre-adjustment, set at +100mm but can be adjusted if required	(NC 0220, 0221, 0223, 0225, 0228)	○
BNE0207		Without Backrest width pre-adjustment, set at 0mm but can be adjusted if required		○
BACKREST *				
BNE0210		Flex 3 Backrest, Fixed Plate (Width Adjustable)	(NC 0233, 0234, 0240, 0343, 0771, 0905)	◆
BNE0220		Flex 3 Backrest Tension (Width Adjustable)	(NC 0233, 0234, 0343, 0771, 0925)	€114
BNE0221		Matrx Elite (includes hardware mounted to frame, backrest supplied in box)	(NC 0240, 0250, 0260, 0270, 0280, 0290, 0321, 0322, 0335, 0336, 0337, 0338, 0339, 0340, 0341, 0342, 0344, 0345, 0346, 0347, 0348, 0349, 0350, 0351, 0352, 0353, 0354, 0355, 0356, 0357, 0358, 0359, 0360, 0361, 0362, 0363, 0364, 0365, 0366, 0367, 0368, 0369, 0370, 0371, 0372, 0373, 0374, 0375, 0376, 0377, 0378, 0379, 0380, 0381, 0382, 0383, 0384, 0385, 0386, 0387, 0388, 0389, 0390, 0391, 0392, 0393, 0394, 0395, 0396, 0397, 0398, 0399, 0400, 0401, 0402, 0403, 0404, 0405, 0406, 0407, 0408, 0409, 0410, 0411, 0412, 0413, 0414, 0415, 0416, 0417, 0418, 0419, 0420, 0421, 0422, 0423, 0424, 0425, 0426, 0427, 0428, 0429, 0430, 0431, 0432, 0433, 0434, 0435, 0436, 0437, 0438, 0439, 0440, 0441, 0442, 0443, 0444, 0445, 0446, 0447, 0448, 0449, 0450, 0451, 0452, 0453, 0454, 0455, 0456, 0457, 0458, 0459, 0460, 0461, 0462, 0463, 0464, 0465, 0466, 0467, 0468, 0469, 0470, 0471, 0472, 0473, 0474, 0475, 0476, 0477, 0478, 0479, 0480, 0481, 0482, 0483, 0484, 0485, 0486, 0487, 0488, 0489, 0490, 0491, 0492, 0493, 0494, 0495, 0496, 0497, 0498, 0499, 0500, 0501, 0502, 0503, 0504, 0505, 0506, 0507, 0508, 0509, 0510, 0511, 0512, 0513, 0514, 0515, 0516, 0517, 0518, 0519, 0520, 0521, 0522, 0523, 0524, 0525, 0526, 0527, 0528, 0529, 0530, 0531, 0532, 0533, 0534, 0535, 0536, 0537, 0538, 0539, 0540, 0541, 0542, 0543, 0544, 0545, 0546, 0547, 0548, 0549, 0550, 0551, 0552, 0553, 0554, 0555, 0556, 0557, 0558, 0559, 0560, 0561, 0562, 0563, 0564, 0565, 0566, 0567, 0568, 0569, 0570, 0571, 0572, 0573, 0574, 0575, 0576, 0577, 0578, 0579, 0580, 0581, 0582, 0583, 0584, 0585, 0586, 0587, 0588, 0589, 0590, 0591, 0592, 0593, 0594, 0595, 0596, 0597, 0598, 0599, 0600, 0601, 0602, 0603, 0604, 0605, 0606, 0607, 0608, 0609, 0610, 0611, 0612, 0613, 0614, 0615, 0616, 0617, 0618, 0619, 0620, 0621, 0622, 0623, 0624, 0625, 0626, 0627, 0628, 0629, 0630, 0631, 0632, 0633, 0634, 0635, 0636, 0637, 0638, 0639, 0640, 0641, 0642, 0643, 0644, 0645, 0646, 0647, 0648, 0649, 0650, 0651, 0652, 0653, 0654, 0655, 0656, 0657, 0658, 0659, 0660, 0661, 0662, 0663, 0664, 0665, 0666, 0667, 0668, 0669, 0670, 0671, 0672, 0673, 0674, 0675, 0676, 0677, 0678, 0679, 0680, 0681, 0682, 0683, 0684, 0685, 0686, 0687, 0688, 0689, 0690, 0691, 0692, 0693, 0694, 0695, 0696, 0697, 0698, 0699, 0700, 0701, 0702, 0703, 0704, 0705, 0706, 0707, 0708, 0709, 0710, 0711, 0712, 0713, 0714, 0715, 0716, 0717, 0718, 0719, 0720, 0721, 0722, 0723, 0724, 0725, 0726, 0727, 0728, 0729, 0730, 0731, 0732, 0733, 0734, 0735, 0736, 0737, 0738, 0739, 0740, 0741, 0742, 0743, 0744, 0745, 0746, 0747, 0748, 0749, 0750, 0751, 0752, 0753, 0754, 0755, 0756, 0757, 0758, 0759, 0760, 0761, 0762, 0763, 0764, 0765, 0766, 0767, 0768, 0769, 0770, 0771, 0772, 0773, 0774, 0775, 0776, 0777, 0778, 0779, 0780, 0781, 0782, 0783, 0784, 0785, 0786, 0787, 0788, 0789, 0790, 0791, 0792, 0793, 0794, 0795, 0796, 0797, 0798, 0799, 0800, 0801, 0802, 0803, 0804, 0805, 0806, 0807, 0808, 0809, 0810, 0811, 0812, 0813, 0814, 0815, 0816, 0817, 0818, 0819, 0820, 0821, 0822, 0823, 0824, 0825, 0826, 0827, 0828, 0829, 0830, 0831, 0832, 0833, 0834, 0835, 0836, 0837, 0838, 0839, 0840, 0841, 0842, 0843, 0844, 0845, 0846, 0847, 0848, 0849, 0850, 0851, 0852, 0853, 0854, 0855, 0856, 0857, 0858, 0859, 0860, 0861, 0862, 0863, 0864, 0865, 0866, 0867, 0868, 0869, 0870, 0871, 0872, 0873, 0874, 0875, 0876, 0877, 0878, 0879, 0880, 0881, 0882, 0883, 0884, 0885, 0886, 0887, 0888, 0889, 0890, 0891, 0892, 0893, 0894, 0895, 0896, 0897, 0898, 0899, 0900, 0901, 0902, 0903, 0904, 0905, 0906, 0907, 0908, 0909, 0910, 0911, 0912, 0913, 0914, 0915, 0916, 0917, 0918, 0919, 0920, 0921, 0922, 0923, 0924, 0925, 0926, 0927, 0928, 0929, 0930, 0931, 0932, 0933, 0934, 0935, 0936, 0937, 0938, 0939, 0940, 0941, 0942, 0943, 0944, 0945, 0946, 0947, 0948, 0949, 0950, 0951, 0952, 0953, 0954, 0955, 0956, 0957, 0958, 0959, 0960, 0961, 0962, 0963, 0964, 0965, 0966, 0967, 0968, 0969, 0970, 0971, 0972, 0973, 0974, 0975, 0976, 0977, 0978, 0979, 0980, 0981, 0982, 0983, 0984, 0985, 0986, 0987, 0988, 0989, 0990, 0991, 0992, 0993, 0994, 0995, 0996, 0997, 0998, 0999, 1000)	€570

BNE0375		Standard position for legrest mounting	◆
SWO0189		Standard position with 0-30° Legrest Abduction function	€ 110
BNE0376		High position 50mm higher for leg length compensation	○
BNE0378		90° Hanger attachment to achieve a 90° legrest position on the Elevating Legrest (NC 0390, 0400, 0410, 0440, 0450, 0460)	○
LEGREST HANGER LH *			
BNE0380		Elevating Legrest, angle adjustable 0 to 75° swing in/swing out with knee pad and calf support LH (NC 0475, 0480, 0490, 0942, 0945)	◆
BNE0390		Legrest Hanger fixed 80 deg, swing in/swing out LH (NC)	○
BNE0400		Legrest Hanger fixed 90 deg, swing in/swing out LH (NC 0480, 0520, 0930)	○
BNE0410		Legrest Amputee LH (NC 0475, 0476, 0477, 0480, 0490, 0520, 0930, 0941, 0942, 0945)	€154
BNE0420		Without Legrest, LH (NC 0475, 0476, 0477, 0480, 0490, 0520, 0930, 0941, 0942, 0945)	○
LEGREST HANGER RH *			
BNE0430		Elevating Legrest, angle adjustable 0 to 75° swing in/swing out with knee pad and calf support RH (NC 0475, 0540, 0550, 0942, 0945)	◆
BNE0440		Legrest Hanger fixed 80 deg, swing in/swing out RH (NC)	○
BNE0450		Legrest Hanger fixed 90 deg, swing in/swing out RH (NC 0540, 0580, 0940)	○
BNE0460		Legrest Amputee RH (NC 0475, 0476, 0477, 0540, 0550, 0580, 0940, 0941, 0942, 0945)	€154
BNE0470		Without Legrest, RH (NC 0475, 0476, 0477, 0540, 0550, 0580, 0940, 0941, 0942, 0945)	○
FOOTREST LH *			
BNE0475		Footrest one piece, depth and angle adjustable short tube (NC 0477, 0480, 0490, 0520, 0540, 0550, 0580, 0930, 0940)	€131
BNE0476		Footrest one piece, depth and angle adjustable long tube (NC 0477, 0480, 0490, 0520, 0540, 0550, 0580, 0930, 0940)	€131
BNE0477		Footbox, angle adjustable, padded with short tube (NC 0480, 0490, 0520, 0540, 0550, 0580, 0930, 0940, 0941, 0945)	€622
BNE0480		Footrest lower part, fixed with short tube and heel strap, LH (NC 0930, 0941)	○
BNE0490		Footrest lower part, depth and angle adjustable with short tube and heel strap, LH	○
BNE0520		Footrest lower part, depth and angle adjustable with long tube and heel strap, LH	○
BNE3333		Without footrest, only select with option 0370 Central legrest, 0410 Amputee legrest and 0420 Without legrest	○
FOOTREST RH *			
BNE0540		Footrest lower part, fixed with short tube and heel strap, RH (NC 0940, 0941)	○
BNE0550		Footrest lower part, depth and angle adjustable with short tube and heel strap, RH	○
BNE0580		Footrest lower part, depth and angle adjustable with long tube and heel strap, RH	○

BNE4444		Without footrest, only select with options 0370 Central legrest, 0460 Amputee legrest, 0470 Without legrest, 0475 and 0476 One piece footrests or 0477 Footbox	○
ARMREST LH *			
BNE0600		Armrest Ducktail Pad LH, height adjustable with autolock (NC 1232)	○
BNE0605		Armrest Wide Pad LH, height adjustable with autolock. Compatible with 16" wheels only (NC 1232)	◆
BNE0607		Armrests Reinforced Wide Pad, LH, height adjustable with autolock (NC 1230, 1231, 1232, 1250, 1260)	€101
BNE0610		Armrest Padded LH, height adjustable with autolock (NC 1230, 1231, 1232, 1250, 1260)	€71
BNE0615		Armrest Hemiplegic LH, height adjustable with autolock (NC 0695, 1230, 1231, 1232, 1250, 1260)	€201
BNE0620		Without Armrest LH Select when choosing BNE1232 Half Tray LH, (NC 0695, 1230, 1231, 1250, 1260) Narrow Armrest is included	○
ARMREST RH *			
BNE0630		Armrest Ducktail Pad RH, height adjustable with autolock (NC 1233)	○
BNE0635		Armrest Wide Pad RH, height adjustable with autolock. Compatible with 16" wheels only (NC 1233)	◆
BNE0637		Armrest Reinforced Wide Pad, RH, height adjustable with autolock (NC 1230, 1231, 1233, 1250, 1260)	€101
BNE0640		Armrest Padded RH, height adjustable with autolock. (NC 1230, 1231, 1233, 1250, 1260)	€71
BNE0645		Armrest Hemiplegic, RH, height adjustable with autolock (NC 0695, 1230, 1231, 1233, 1240, 1245, 1250, 1260)	€201
BNE0650		Without Armrest RH Select when choosing BNE1233 Half Tray RH, Narrow Armrest is included (NC 0695, 1230, 1231, 1240, 1245, 1250, 1260)	○
ARMREST OPTIONS. Side support cover 0680, must be selected if 0690 Cushion to side support is required			
BNE0680		Side Support Cover, removable with tools (Wipeable cover for the sideguard) (NC 0695)	€57
BNE0690		Cushion to side support cover (Internal cushioned pad, which is only possible in combination with 0680) (NC 0695)	○
BNE0695		Pad for sideguard, removable with hook and loop fastening	€57
SEAT OPTIONS			
BNE0710		Incontinence Cover	€39
BNE0720		Pommel Abduction Cushion (NC 0722)	€131
BNE0722		2 cm Increased seat height plate	€81
BNE0725		Tilt Scale for seat tilt	€48
BACKREST OPTIONS			
BNE0730		Headrest, Sideway Adjustable (NC 0750, 0770)	◆

BNE0750		Headrest with cheek support, sideways adjustable (NC 0740, 0770, 0780)	€122
BNE0740		Cover for headrest	€35
BNE0770		Neckrest (NC 0771)	○
BNE0771		Headrest Matrix, 250 mm Soft Pad (NC 0780, 0790)	€345
BNE0780		Cover For Neckrest	€23
BNE0782		Alternative headrest/neckrest mounting bracket	€22
BNE0790		Handles Clinic Kit for Head-/Neckrest (Levers which replace screws for easy adjustment)	€55
BNE0900		Trunk Support Attachment to mount laterals later. If trunk supports are chosen, the attachment bracket is already included	○
BNE0905		Lateral Positioning Pads. Attaches with hook and loop between the strap and the cover of the tension adjustable upholstery	€63
BNE0910		Lumbar Support Cushion	€34
BNE0915		Tilt Scale for backrest recline	€44
BNE0920		Wedges (One Pair)	€12
BNE0925		Privacy Attachment (Material connecting backrest and seat)	€34
TRUNK SUPPORT LH			
BNE0800		Trunk Support Multi Adjustable with standard arm LH (NC 0900)	○
BNE0820		Trunk Support Multi Adjustable with extended height arm LH (NC 0900)	○
BNE0880		Trunk Support Swingaway LH (NC 0900)	€116
BNE0895		Trunk Support Multi Adjustable Swingaway LH (NC 0900)	€126
TRUNK SUPPORT RH			
BNE0810		Trunk Support Multi Adjustable with standard arm RH (NC 0900)	○
BNE0830		Trunk Support Multi Adjustable with extended height arm RH (NC 0900)	○
BNE0890		Trunk Support Swingaway RH (NC 0900)	€116
BNE0896		Trunk Support Multi Adjustable Swingaway RH (NC 0900)	€126
LEGREST OPTIONS			
BNE0930		Footplate Extender LH (NC 0941)	€36
BNE0940		Footplate Extender RH (NC 0941)	€36
BNE0941		Cover Footplate, Black Plush Fabric	€60
BNE0942		Calf Strap Padded	€80
BNE0945		Calf Strap (for One piece footrest)	○
CHASSIS COLOUR *			
BNE0944		Silver Pearl	◆
BNE0946		Azurite (Dark Blue)	○
BNE0947		Sand (Beige)	○
BNE0948		Electric Green	○
BNE0949		Happy Red	○
CHASSIS SAFETY			
BNE0950		Chassis extender, castors extend 50mm forward for added stability	○
CASTOR ATTACHMENT *			
BNE0960		Fixed castor attachment	◆

FRONT CASTORS *			
BNE0992	☐	150x30 mm (6" x 1 1/4") Medium Hard Black tyre	○
BNE1015	☐	200x32 mm (8" x 1 1/4") Black Solid Tyre	○
BNE1025	☐	200x45 mm (8" X 1 3/4") Black Solid tyre	◆
REAR WHEEL SIZE *			
BNE1047	☐	16" / 405 mm Transit rear wheel with fixed attachment (NC)	◆
REAR WHEEL TYPE *			
BNE1076	☐	Hub Brake 16" Transit with Black Plastic Hub and Rim (NC)	◆
REAR WHEEL TYRE *			
BNE1092	☐	Solid Black tyres	◆
BRAKE *			
BNE1170	☐	User manoeuvred tyre brakes	◆
SAFETY OPTIONS			
BNE1220	☐	Reflector Kit (wheels and chassis)	€23
PELVIC BELT *			
BNE1200	☐	Pelvic belt	◆
TRANSPORTATION BRACKETS *			
BNE1225	☐	Transportation Brackets, mounted to the front of the wheelchair	€14
OTHER OPTIONS			
BNE1230	☐	Fixed Table Tray without elbow pads includes mounting kit (NC 1237, 1260)	€156
BNE1231	☐	Fixed Table Tray with elbow Pads includes mounting kit (NC 1237, 1260)	€156
BNE1232	☐	Half table tray, clear plastic, flip up, depth adjustable, left side (armrest with straight narrow pad included) (NC 1235, 1260)	€125
BNE1233	☐	Half table tray, clear, plastic, flip up, depth adjustable, right side (armrest with straight narrow pad included) (NC 1235, 1260)	€125
BNE1240	☐	Swing away table without elbow pads (NC 1235, 1250, 1260)	€219
BNE1245	☐	Swing away table with elbow pads (NC 1235, 1260)	€219
BNE1235	☐	Table tray padded cushion	€57
BNE1237	☐	Half table tray padded cushion	€57
BNE1250	☐	Locking device for table tray	€15
BNE1260	☐	Mounting kit for table (if ordered separately)	€8
BNE1270	☐	Pump	€58
BNE1280	☐	Cane Holder	€56
BNE1290	☐	Handles Clinic Kit for Head-/Neckrest (Levers which replace screws for easy adjustment) (NC)	€61
BNE1300	☐	Harness attachment , must be selected in combination with Head/Neckrest	€89
BNE1305	☐	Front castor locks, includes front tiedown brackets. Allows castors to remain parallel for easy access into confined spaces such as lifts or vehicles	€217
BNE1310	☐	Attachment for positioning belt	€35

TOTAL PRICE

Non Conforming Options:

BNE0110	400 mm / 16" Seat Height
BNE0201	Black Dartex TR26
BNE0202	Black Spacer TR36
BNE1111	No upholstery, select with options 0170 Multifunctional and 0180 Without seat cushion
BNE0206	+ 100 mm Backrest width pre-adjustment, set at +100mm but can be adjusted if required
BNE0220	Flex 3 Backrest Tension (Width Adjustable)
BNE0221	Matrx Elite (includes hardware mounted to frame, backrest supplied in box)
BNE0223	Matrx PB Elite Deep (includes hardware mounted to frame, backrest supplied in box)
BNE0225	Matrx Posture Back (includes hardware mounted to frame, backrest supplied in box)
BNE0228	Matrx Elite TR (includes hardware mounted to frame, backrest supplied in box)
BNE0233	Matrx Backrest Height 410 mm / 16"
BNE0234	Matrx Backrest Height 510 mm / 20"
BNE0240	Cover for tension adjustable upholstery
BNE0250	Cushion, Laguna (Lateral Support)
BNE0260	Cushion, Mistral 2 (Waist Support)
BNE0270	Cushion, Passad 2 (Shoulder Support)
BNE0280	Cushion, Shoulder High. Higher Cushion to suit taller clients
BNE0290	Cushion, Vicair Multifunctional (Black)
BNE0321	Black Dartex TR26
BNE0322	Black Spacer TR36
BNE2222	No upholstery, select with Matrx, Multifunctional and Without backrest options
BNE0335	Push Handles connected,no head/neckrest attachment
BNE0343	Push Handles connected,with Rea head/neckrest attachment
BNE0360	Electrical operated seat tilt and backrest recline
BNE1263	Tilt Lock (locking of pre-set angle)
BNE5555	No Tilt Lock, Select with 0360 Electric operated seat and backrest
BNE0370	Central Legrest, includes calfpad and footrest. Also select options 0420, 0470, 3333, 4444
SWO0189	Standard position with 0-30° Legrest Abduction function
BNE0376	High position 50mm higher for leg length compensation
BNE0378	90° Hanger attachment to achieve a 90° legrest position on the Elevating Legrest
BNE0380	Elevating Legrest, angle adjustable 0 to 75° swing in/swing out with knee pad and calf support LH
BNE0390	Legrest Hanger fixed 80 deg, swing in/swing out LH
BNE0400	Legrest Hanger fixed 90 deg, swing in/swing out LH
BNE0410	Legrest Amputee LH

BNE0430	Elevating Legrest, angle adjustable 0 to 75° swing In/swing out with knee pad and calf support RH
BNE0440	Legrest Hanger fixed 80 deg, swing in/swing out RH
BNE0450	Legrest Hanger fixed 90 deg, swing in/swing out RH
BNE0460	Legrest Amputee RH
BNE0475	Footrest one piece, depth and angle adjustable short tube
BNE0476	Footrest one piece, depth and angle adjustable long tube
BNE0477	Footbox, angle adjustable, padded with short tube
BNE0480	Footrest lower part, fixed with short tube and heel strap, LH
BNE0490	Footrest lower part, depth and angle adjustable with short tube and heel strap, LH
BNE0520	Footrest lower part, depth and angle adjustable with long tube and heel strap, LH
BNE0540	Footrest lower part, fixed with short tube and heel strap, RH
BNE0550	Footrest lower part, depth and angle adjustable with short tube and heel strap, RH
BNE0580	Footrest lower part, depth and angle adjustable with long tube and heel strap, RH
BNE0695	Pad for sideguard, removable with hook and loop fastening
BNE0720	Pommel Abduction Cushion
BNE0722	2 cm Increased seat height plate
BNE0750	Headrest with cheek support, sideways adjustable
BNE0740	Cover for headrest
BNE0770	Neckrest
BNE0771	Headrest Matrix, 250 mm Soft Pad
BNE0780	Cover For Neckrest
BNE0790	Handles Clinic Kit for Head-/Neckrest (Levers which replace screws for easy adjustment)
BNE0900	Trunk Support Attachment to mount laterals later. If trunk supports are chosen, the attachment bracket is already included
BNE0905	Lateral Positioning Pads. Attaches with hook and loop between the strap and the cover of the tension adjustable upholstery
BNE0910	Lumbar Support Cushion
BNE0920	Wedges (One Pair)
BNE0925	Privacy Attachment (Material connecting backrest and seat)
BNE0800	Trunk Support Multi Adjustable with standard arm LH
BNE0820	Trunk Support Multi Adjustable with extended height arm LH
BNE0880	Trunk Support Swingaway LH
BNE0895	Trunk Support Multi Adjustable Swingaway LH
BNE0810	Trunk Support Multi Adjustable with standard arm RH
BNE0830	Trunk Support Multi Adjustable with extended height arm RH
BNE0890	Trunk Support Swingaway RH
BNE0896	Trunk Support Multi Adjustable Swingaway RH
BNE0930	Footplate Extender LH
BNE0940	Footplate Extender RH
BNE0941	Cover Footplate, Black Plush Fabric
BNE0942	Calf Strap Padded
BNE0945	Calf Strap (for One piece footrest)
BNE1015	200x32 mm (8" x 1 1/4") Black Solid Tyre

BNE1025	200x45 mm (8" X 1 3/4") Black Solid tyre
BNE1230	Fixed Table Tray without elbow pads includes mounting kit
BNE1231	Fixed Table Tray with elbow Pads includes mounting kit
BNE1232	Half table tray, clear plastic, flip up, depth adjustable, left side (armrest with straight narrow pad included)
BNE1233	Half table tray, clear, plastic, flip up, depth adjustable, right side (armrest with straight narrow pad included)
BNE1240	Swing away table without elbow pads
BNE1245	Swing away table with elbow pads
BNE1235	Table tray padded cushion
BNE1237	Half table tray padded cushion
BNE1250	Locking device for table tray
BNE1260	Mounting kit for table (if ordered separately)
BNE1290	Handles Clinic Kit for Head-/Neckrest (Levers which replace screws for easy adjustment)



For a more complex customer requirement
please refer to our Invacare Unique special order form.
email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions - Please supply the following information:

Name :

Email :

Phone :

Fax :

Description of special request (please include relevant option code or part number):