# **Invacare SB755 Community Bed**

# **Product Data Sheet**



## **Product Range Available**

- SB755
- SB755 Detachable



#### **Standard Features**

- Complies to the Bed standard IEC 60601-2-52 **Third Edition**
- Four sectioned profiling mattress support
- Auto regression as standard to increase comfort when activating backrest, minimising shear/friction forces
- Optimised dimensions of mattress platform for increased comfort
- Integrated extendable mattress platform, additional 100mm at head and foot end
- Wide range of bed ends, side rails and side rail covers available
- Suitable for Bariatric use 31.5 Stone Maximum User Weight

#### **Standard Features Detachable**

- Unique dismountable system is designed for minimum effort
- Designed with parts that weight no more that 25kg each
- Bed parts, mattress and accessories are stored directly on the base



SB755 Detachable

#### **Technical Specification**

#### SB755 and SB755 Detachable

Width outside: 950/1010mm
Width Inside: 850/900mm
Length outside: 2130mm
Length inside: 2080mm

- Mattress support dimensions: 800mm (Head section) 260mm (Seat plate) x 260mm (thigh section) x 680mm (lower leg section)
- Backrest Angle: 0-73°
- Thigh Angle: 0-28°
- Full Legtest lift 0-18°
- Legrest angle below horizontal 0 -16°
- Maximum user weight 200kg (31.5 stone)
- Safe Working Load 235kg (37 stone)

#### **Optional accessories available including:**

- 3 bed ends available for the SB755 and Detachable
- Full length Line and Line extendable side rails
- Scala 2 range side rails
- Lifting pole
- Support handles
- Standard or soft hand control



Anita bed end



Susanne bed end



#### Victoria bed end

## **Supporting Documentation**

- SB755 range brochure PN: 1523060

- Bed Accessory brochure\*

- Beds & PAC DL leaflet PN: 1568402

- Bed rail ruler PN: 1556575

#### **Manufacturers Contact Details**

Invacare Ltd, Unit 4 Pencoed Technology Park, Bridgend, CF35 5AQ

Tel: 01656 776222 E-mail: <u>UK@invacare.com</u>

Web: www.invacare.co.uk



<sup>\*</sup>Contact Marketing department