Invacare SB755 Wide Community Bed

Product Data Sheet



Product Range Available

- SB755 Wide



Standard Features

- Complies to the Bed standard IEC 60601-2-52 **Third Edition**
- Four sectioned profiling mattress support
- Auto regression as standard to increase comfort when activating backrest, minimising shear/friction forces
- Optimised dimensions of mattress platform for increased comfort
- Integrated extendable mattress platform, additional 100mm at head and foot end
- Wide range of side rails and side rail covers available
- Suitable for Bariatric use 31.5 Maximum user weight
- Mattress platform width of 1200mm
- Detachable format standard

Technical Specification

Width outside: 1300mm
Width Inside: 1200mm
Length outside: 2130mm
Length inside: 2080mm

- Mattress support dimensions: 800mm (Head section) 260mm (Seat plate) x 260mm (thigh section) x 680mm (lower leg section)



- Backrest Angle: 0-73°

- Thigh Angle: 0-28°

- Full Legtest lift 0-18°

- Legrest angle below horizontal 0 -16°

- Maximum user weight 200kg (31.5 stone)

- Safe Working Load 235kg (37 stone)

Optional accessories available including:

- Full length Line and Line extendable side rails
- Brit V wooden side rails
- Scala 2 range side rails
- Lifting pole
- Trolley for transportation
- Support handles



Susanne bed end

Supporting Documentation

- SB755 range brochure PN: 1523060

Bed Accessory brochure*

- Beds & PAC leaflet PN: 1568402

- Bed rail ruler PN: 1556575

Manufacturers Contact Details

Invacare Ltd, Unit 4 Pencoed Technology Park, Bridgend, CF35 5AQ

Tel: 01656 776222 E-mail: <u>UK@invacare.com</u>

Web: www.invacare.co.uk



^{*}Contact Marketing department