

Sherpa 3

by United Care



A versatile ceiling lift solution

for all care environments



A solution for all lifting needs

The compact, lightweight ceiling lift with a new Lithium-Ion battery

Sherpa 3 is a safe, simple ceiling lift for seated, horizontal and floor-to-bed transfers to and from the bed, wheelchair or toilet. It supports both low and high lifting heights, making it suitable for lifting from the floor, standing training and walking rehabilitation.

Available in 200 kg and 400 kg versions with a choice of charging, traversing and spreader bar options, Sherpa 3 can be individually adapted to nearly any care setting. It works with United Care and SystemRoMedic rails and with freestanding systems such as the Castor or Helena.



1

Flexible by design

Choose 200 kg model for manual or powered traversing, end-point or in-rail charging and 2 or 4 point spreader bars.

Select the 400 kg model for manual traversing, end-point charging and a choice of 4 or 6 point spreader bars.*

2

Safe and reassuring

Manual emergency stop, electrical and mechanical emergency lowering, overload and over-speed protection, a slack strap switch, and soft start and stop.

> Optional quick release trolley

for tool-free disconnection of the motor from the rail

> Energy-efficient Gravity Powered system

that recharges the battery by gravity

> Several spreader bar options:

2-point, 4-point and 4-point electric

> Lightweight ceiling lift

with power traverse option

> Clear, intuitive interface

with battery and fault indicator lights

> Emergency stop

with electric and mechanical emergency lowering

> Easy-to-operate, ergonomic hand control



*6-point only available with Sherpa 3 400 kg.

Technical data ▼

	 Safe working load ▼	 Width of motor ▼	 Height of motor ▼	 Depth of motor ▼
Sherpa 3, 200kg	200 kg	300 mm	190 mm	240 mm
Sherpa 3, 400kg	400 kg	300 mm	190 mm	240 mm

	 Lifting range ▼	 Sound level ▼	 Total weight ▼
Sherpa 3, 200kg	2.3 m	44 dB (A)	8.5 kg
Sherpa 3, 400kg	2.3 m	44 dB (A)	11 kg

All measurements are taken with a tolerance of + / - 5 mm

Invacare International GmbH
Neuhofweg 51, 4147 Aesch BL, Switzerland
globalcomms@dhcare.com
www.dhcare.com



Invacare International GmbH is a member of the DHCare group of companies. Invacare International GmbH reserves the right to modify existing product models at any time. Product specifications may vary across regions. ©2026 Invacare International GmbH. All rights reserved. Sherpa_EU_07082026