

# Aquatec® | **Shower stools, shower chairs**



*Invacare Alize (with arms)*

- Stable shower stool
- Large ergonomic seat surface with hygiene recess makes personal hygiene easier
- Fixed armrests
- Height-adjustable, corrosion-protected legs



*Invacare Alize shower seat*

- Ergonomic seat surface with hygiene recess
- Height-adjustable, corrosion-protected legs
- Height adjustable 500 - 560 mm

**AQUATEC®**



**Yes, you can.®**



**Aquatec Galaxy**


**3 year warranty**

- Folding, height-adjustable shower chair, perfect for travelling and compact for storage
- Ergonomic seat surface, can be supplied with or without hygiene recess
- Soft backrest for additional support


**Technical data**

Order numbers	Invacare® <i>Alize stool c/w arms</i>	Invacare <i>Alize chair</i>	Aquatec <i>Galaxy</i>
Shower stool with armrests and hygiene recess	<b>H2240/1</b>		
Shower chair with backrest and armrests, with hygiene recess		<b>H2480/1</b>	<b>1470805</b>
<b>Dimensions</b>			
Seat width	410 mm (A)	410 mm (A)	400 mm (A)
Seat depth	400 mm (B)	400 mm (B)	345 mm (B)
Seat height	470 - 550 mm (C)	500 - 560 mm (C)	470 - 620 mm (C)
Total width	590 mm (D)	590 mm (D)	580 mm (D)
Total depth	480 mm (E)	480 mm (E)	460 mm (E)
Total height	550 mm (F)	940 mm (F)	840 - 990 mm (F)
Graduation, height adjustment			
Height to armrests			630 - 780 mm (G)
Width between the armrests	450 mm (G)	450 mm (G)	510 mm (H)
Width, hygiene recess	140 mm	140 mm	130 mm
Depth, hygiene recess	170 mm	170 mm	140 mm
<b>More info</b>			
Weight	5.7 kg	7.5kg	7.5 kg
Load capacity	120 kg (19 stone)	120kg (19 stone)	100 kg (16 stone)
Colour	White	White	White
Warranty	2 years	2 years	3 years

Invacare *Alize with arms*



Invacare *Alize shower seat* / Aquatec *Galaxy*





Aquatec **Pico**

- Modular shower chair system
- Four configurations
- Textured seat surface for safe sitting
- Large hygiene recess, either at the front or back, due to two mounting options for the backrest
- Height-adjustable legs
- Rubber feet provide a secure stand
- Simple and easy assembly, no tools required
- Load capacity of up to 160 kg (25 stone)

### Options and features



**Pico**



**Pico** with armrests



**Pico** with backrest

**3 year warranty**

## Technical data

Order numbers	Aquatec <b>Pico</b>
Shower stool	<b>9.20.001</b>
Shower stool with armrests	<b>9.21.001</b>
Shower chair with backrest	<b>9.23.001</b>
Shower chair complete	<b>9.22.001</b>
Backrest	<b>16469</b>
Armrests (pair)	<b>16468</b>

Dimensions	
Seat width	430 mm (A)
Seat depth	420 mm (B)
Seat height	425 - 575 mm (C)
Total width	575 mm (D)
Total depth	520 mm (E)
Total height	825 - 975 mm (F)
Height to armrests	605 - 755 mm (G)
Width between the armrests	430 mm (H)
Width, hygiene recess	180 mm
Depth, hygiene recess	145 mm
Graduation, height adjustment	25 mm

More info	
Weight, shower stool	3,2 kg
Weight, shower stool with armrests	4,1 kg
Weight, shower chair with backrest	4,0 kg
Weight, shower chair complete	4,8 kg
Load capacity	160 kg! (25 stone)
Colour	White
Warranty	3 years

Aquatec **Pico**

