

Sliding & Transfer Aids

SystemRoMedic™



Safe and easy

Reduces friction for gentle, effortless transfers



Trusted since 1984

EasySlide is the original SystemRoMedic transfer aid – used in hospitals and nursing homes worldwide for over 40 years.

Its unique Ultra Low Friction (ULF) material on the inner surface enables caregivers to transfer users on virtually any surface with minimal effort and no lifting. Available in several models, materials and sizes, EasySlide suits

a wide range of transfer situations. Classic models are made of microfibre with polyester padding and ULF lining. Most sizes can be fitted with a replaceable or disposable protective cover for improved hygiene.

EasySlide 190 x 60

Designed for supine transfers to and from beds, x-ray tables, operating tables, shower trolleys and stretchers – and equally effective for turning and repositioning in bed.

> **ULF sliding material**

Ultra Low Friction inner layer reduces resistance on any surface.

> **Seat belt included**

Safely helps to stabilise the user during transfer.

> **Soft padded surface**

Gentle on skin – comfortable for users and caregivers alike.

> **Protective cover available**

Replaceable polyester cover with carbon fibre reduces static electricity.

> **Disposable covers available**

Single-use plastic covers for maximum hygiene and less laundry.



Ultra Low Friction

Unique ULF material reduces friction on any surface – enabling gentle transfers with minimal effort and no lifting.



Multiple models & sizes

Available in microfibre and disposable versions to suit any transfer situation or hygiene requirement.

EasySlide disposable

EasySlide disposable is a thin, pliable and versatile sliding mat on a roll – up to 60 individual mats per pack. Suitable for turning, repositioning, supine and seated transfers, leg lifting, dressing assistance and sling application. Each roll comes in a hygienic wall-mounted dispenser for easy access.

> **Up to 60 mats per roll**

Always fresh, clean and ready – supplied in practical two-packs.

> **Perforated every 50 cm**

Simply unroll and tear off the required length.

> **Tubular or open**

Use as a tube or separate into a flat sheet.

> **Eco-friendly material**

Tough, durable and not water-soluble – ideal for showering situations.

> **Instructions on every mat**

User guide printed directly on each sheet.

> **Wall dispenser available**

Hygienic storage within reach – also available as a steel tabletop dispenser.

EasySlide protective cover

Fits most EasySlide models. Made of polyester with carbon fibre to reduce static electricity. Easily washable and reusable.



> **Carbon fibre lining**
Reduces static electricity for smoother, safer transfers.

OneWaySlide


OneWaySlide is a sliding mat with many functions. The main one improves sitting positions in chairs or wheelchairs. The sliding mat effectively retards forward movement yet slides backwards, making it easy for the carer to help the user move into a more comfortable position in the chair. It is simply a case of grabbing the handles and pulling – without great use of force.



> **Ultra Low Friction**
OneWaySlide is made of polyester, polyethylene and ULF (Ultra Low Friction), Handicare's own unique sliding material.

EasySlide disposable cover

Disposable plastic covers for EasySlide reduce laundry and improve infection control. Made of bio-degradable polyethylene.



> **Biodegradable polyethylene**
Hygienic, eco-conscious and low-maintenance.

ReadySlide

ReadySlide is a thin, pliable and very easy-to-use sliding mat that facilitates various types of lying transfers of users, positioning in chairs and application of lifting slings in a seated position.



> **Ultra Low Friction**
ReadySlide is made of ULF (Ultra Low Friction), Handicare's own unique sliding material.

ReadySheet

ReadySheet is a simple and versatile assistive device for application of lifting slings and for repositioning in one location, for example, for turning or moving the user higher up in bed. ReadySheet is available in two models; one open and one tubular, and in many sizes. All variants have handles for easier positioning.



> **Ultra Low Friction**
ReadySheet is made of ULF (Ultra Low Friction), Handicare's own unique sliding material.

RollerSlide

RollerSlide is used for supine transfers of inert users between beds and an x-ray or operating tables, shower trolleys, ambulance stretchers or similar. RollerSlide bridges gaps and moves with the user during the transfer.



> **Replaceable sliding cover**
The replaceable sliding cover can easily be cleaned and disposable covers, which are available as accessories, provide for even further improved hygiene and less laundry.

FlexiMove

FlexiMove is primarily designed for moving higher up in bed, turning and sideways transfers. In these situations, pulling is easier than lifting as the underside is smooth and slides very well in the bed.



> **Three sizes available**
FlexiMove is available in three sizes with four or six handles, all for users weighing max 150 kg.

Technical data ▼



	Description ▼	Length ▼	Width ▼
EasySlide 190 x 60	190 x 60	1900 mm	600 mm
	190 x 60, complete with Protective cover	1900 mm	600 mm
OneWaySlide	50 x 40, tubular	500 mm	400 mm
	50 x 45, tubular	500 mm	450 mm
	110 x 30	1100 mm	300 mm
	130 x 40	1300 mm	400 mm
	130 x 45	1300 mm	450 mm
	130 x 55	1300 mm	550 mm
ReadySlide	195 x 70 cm, polyamid	1950 mm	700 mm
	195 x 70 cm, polyamid, with handles	1950 mm	700 mm
	110 x 70 cm, polyamid	1100 mm	700 mm
	110 x 90 cm, polyamid, extra wide	1100 mm	900 mm
	110 x 100 cm	1100 mm	1000 mm
	50 x 60 cm	500 mm	600 mm
	130 x 90 cm	1300 mm	900 mm
ReadySheet	70 x 70 cm, polyamid	700 mm	700 mm
	180 x 110, polyamid	1800 mm	1100 mm
	195 x 75, polyamid	1950 mm	750 mm
FlexiMove	102 x 43 cm	1020 mm	430 mm
	112 x 53 cm	1120 mm	530 mm
	152 x 53 cm	1520 mm	530 mm
RollerSlide	light, Plastic	850 mm	500 mm
	light, Plastic, Foldable	1750 mm	500 mm
	Plywood, Foldable	1750 mm	500 mm
	flexi, Plastic, Foldable	1750 mm	500 mm
	Plywood	850 mm	500 mm
	flexi, Plastic	850 mm	500 mm
	flexi, Plastic, Foldable, HD	1750 mm	600 mm
	flexi, Plastic, HD	850 mm	600 mm

Invacare International GmbH
 Neuhofweg 51, 4147 Aesch BL, Switzerland
globalcomms@dhcare.com
www.dhcare.com



Invacare International GmbH is a member of the DHCare group of companies. Invacare International GmbH reserves the right to modify existing product models at any time. Product specifications may vary across regions. ©2026 Invacare International GmbH. All rights reserved. EasySlide_EU_14072026