


Softform[®]

Premier Visco

Pressure-redistributing viscoelastic mattress

An exceptionally comfortable support surface that conforms to the person's body shape.




Stability and pressure redistribution

For personalised pressure injury management



Yes, you can.[®]

Effective pressure injury prevention



Pressure Redistribution

- › High quality mattress facilitates pressure reduction by minimising peak pressures.
- › Deep viscoelastic insert pad is high-density and temperature-sensitive to the body, allowing it to conform to a person's shape.
- › Modular construction enables parts to be replaced easily and increases the longevity.



Shear & Friction Reduction

- › Two-way stretch, water-resistant, polyurethane cover is durable and helps to reduce shear and friction.



Comfort

- › High quality mattress delivering an exceptionally comfortable support surface.
- › Viscoelastic foam conforms around the body for custom support and comfort.
- › Four-sided concealed welded zip reduces potential fluid ingress and allows easy detachment for inspection and laundering.

Technical data ▼

For more information about this product, including the user manual, please visit our website: www.invacare.co.uk

| | | Sizes available ▼ | | | | |
|---------------|---------|-------------------|--------|--------|---------|---------|
| | | 850 mm | 880 mm | 900 mm | 1000 mm | 1050 mm |
| 150 mm height | 1970 mm | - | • | - | - | • |
| | 2000 mm | • | • | • | • | - |
| | 2080 mm | - | - | • | - | - |



Max. User Weight ▼

247.6 kg



Max. Cleaning Temperature ▼

95 °C

Manufactured to BS EN 597-1 Cig, BS EN 597-2 Match. Foam Interior: BS 5852:1990 (Ignition Source 5); Cover BS 7175:1989 (Ignition Source 5).

Essential nursing care is pivotal in pressure injury management. This mattress may positively contribute to the outcome of a pressure injury prevention plan that includes education, clinical judgment and action planning based on vulnerability, which are fundamental factors in the prevention of pressure injuries. A range of assessment scales can be used as a formal method of assessing risk from pressure injury development, and should be used in conjunction with an informal assessment (informed nursing judgement) as part of a comprehensive and individualised plan.

T: +44 1656 776222
 uk@invacare.com
 www.invacare.co.uk
 www.invacare.ie

Invacare UK & Ireland
 Unit 4 Pencoed Technology
 Park
 Bridgend CF35 5AQ



Yes, you can.®



© 2024 Invacare International GmbH. All rights reserved. All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation. Softform Premier Visco - UK - 11/2024