



Yes, you can:

Power Retail Prescription Form

Tel: 01656 776222
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

JANUARY 2025
 Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Storm 4 Max

Standard Features

- Max User Weight 220KG (34.5 Stone)
- Max Seating
- Seat Width from 490mm - 610mm
- Seat Depths 560mm - 630mm
- multi stretch PU coated fabric Upholstery
- Seat Height 480mm - 510mm
- Manual Settable Tilt with Manual Settable Recline
- Reinforced Armrests
- Reinforced Legrests with Angle adjustable flip-back footboard
- Headrest, sideways and depth adjustable
- LiNX Control REM110 with standard bracket
- Choice of 6 Colours
- Puncture Proof Castors
- Puncture Proof Drive Wheels
- 6 km/h (4mph) Motors
- 8A (200 - 250V) Charger
- 4 Point tie down system
- 73.5 Ah (Gel) batteries
- Bolted on anti tippers

Standard Base price: £7,789



Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

Code	Description	
CHAIR TYPE*		
AWA0000	Storm 4 Max	£7,789
AWZ0000	Special Storm 4 Max (Please provide details in special section)	POA
SEATING TYPE*		
AWA1145	MAX SEAT UNIT	<input checked="" type="checkbox"/>
BACK CUSHION*		
AWA1106	Tension adjustable 630mm - 620mm high with SW55 / 610mm 570mm High with SW490mm	<input checked="" type="checkbox"/>

SEAT WIDTH*			
AWA0348	<input type="checkbox"/>	Seat Width 490mm	○
AWA0378	<input type="checkbox"/>	Seat Width 550mm	○
AWA0380	<input type="checkbox"/>	Seat Width 610mm (NC 0503)	○
SEAT CUSHION*			
AWA1457	<input type="checkbox"/>	Seat cushion Max (505mm x 535mm)	○
AWA1459	<input type="checkbox"/>	Seat cushion Max (550mm x 635mm)	○
AWA1460	<input type="checkbox"/>	Seat cushion Max (610mm x 635mm)	○
SEAT DEPTH SETTABLE IN 3 STEPS*			
AWA0396	<input type="checkbox"/>	Seat Depth 560mm	○
AWA0397	<input type="checkbox"/>	Seat Depth 595mm	○
AWA0398	<input type="checkbox"/>	Special Seat depth 630mm	POA
ADAPTER TO REDUCE SEAT DEPTH			
AWA0399	<input type="checkbox"/>	Adapter to shorten the seat depth by 130mm	£49
UPHOLSTERY TYPE*			
AWA0278	<input type="checkbox"/>	multi stretch PU coated fabric	◆
ARMRESTS INCLUDING REINFORCED DEPTH ADJUSTABLE BRACKET*			
AWA0635	<input type="checkbox"/>	Fixed armrest 250mm - 340mm height adjustable removable	○
AWA0636	<input type="checkbox"/>	Fixed armrest 290mm - 380mm height adjustable removable	○
ARMPAD*			
AWA0647	<input type="checkbox"/>	Standard Armpad 325mm	○
AWA0649	<input type="checkbox"/>	Armpad extended plus 50mm	○
PUSH BAR*			
AWA1307	<input type="checkbox"/>	Push bar - height adjustable	○
HEADREST / NECKREST - Adequate headrest option must be selected*			
AWA1432	<input type="checkbox"/>	Headrest sideways and depth adjustable	◆
AWA1427	<input type="checkbox"/>	Headrest with cheek support sideways adjustable	£195
AWA1495	<input type="checkbox"/>	Equivalent non Invacare Headrest supporting the requirements of ISO7176-19:2001 Provided independently from this order	○
SEAT HEIGHT*			
AWA1025	<input type="checkbox"/>	Fixed seat height 480mm	○
AWA1006	<input type="checkbox"/>	Fixed seat height 510mm	○
SEAT TILT*			
AWA1002	<input type="checkbox"/>	Manual Seat Tilt by spindle	○
AWA1014	<input type="checkbox"/>	Electrical Tilt 0° to 15°	£514
BACKREST RECLINER*			
AWA0958	<input type="checkbox"/>	Settable recliner by spindle (0° - 25°)	○
AWA0947	<input type="checkbox"/>	Electrical recliner (with high pivot, shearing compensation)	£454
FOOTREST HANGER TYPE (PAIR)*			
AWA0564	<input type="checkbox"/>	Reinforced hanger fixed angle (not detachable without tools) (NC 0522)	○
AWA0503	<input type="checkbox"/>	Vari F : Angle adjustable with tools 0° - 70°, length 290mm - 460mm (NC 0599)	£102
CALF PADS			
AWA0583	<input type="checkbox"/>	Calf Pads adjustable	○

FOOTBOARD*			
AWA0599		Angle adjustable flip back heavy duty footboard (NC 0503)	◆
AWA0522		Angle adjustable footboard (for max user weight 200kg) (NC 0564)	○
CONTROLS*			
AWA2201		REM110 Drive only (NC 1437)	◆
AWA2203		REM211 Drive & Actuators (NC 1437)	○
AWA2204		REM216 Drive, Actuators & Lights	£130
CONTROL BRACKETRY*			
AWA1207		Standard Bracket RH	○
AWA1208		Standard Bracket LH	○
AWA1707		Swing Away RH Bracket	£113
AWA1708		Swing Away LH Bracket	£113
AWA1709		Matrx Resolve Swing in/out bracket RH	£219
AWA1710		Matrx Resolve Swing in/out bracket LH	£219
AWA2227		Bracket RH mounted for Attendant (NC 2219, 2220)	£54
AWA2228		Bracket LH mounted for Attendant (NC 2219, 2220)	£54
DUAL CONTROL			
AWA2220		REM050 Dual Control (Please also select user control and bracket) (Please select dual control bracket) (NC 2227, 2228, 2245)	£338
AWA2219		Compact Dual Control (ACU200) (Please select dual control bracket) (NC 2227, 2228, 2245, 2201)	£338
DUAL CONTROL BRACKETEY			
AWA1207		Dual Control Standard Bracket RH (NC 2227, 2228, 2245)	○
AWA1208		Dual Control Standard Bracket LH	○
INVACARE LiNX IDC DUAL CONTROL			
AWA2245		Invacare LiNX IDC Dual control (Please select IDC dual control bracket below) (NC 2227, 2228, 2220, 2219, 2201)	£641
INVACARE LiNX IDC DUAL CONTROL BRACKETRY			
AWA1221		IDC Dual control Detachable Bracket RH (NC 2227, 2228, 2220, 2219)	○
AWA1222		LiNX IDC Dual control Detachable Bracket LH (NC 2227, 2228, 2220, 2219)	○
LINX 10-WAY SWITCH			
AWA2251		10-Way Switch DLX-SW10-A for direct Powered Seating Access (included loom to connect + additional jack for mode switch)	£431
LINX 10-WAY SWITCH BRACKET			
AWA2217		10-Way Switch Bracket RH, Mounted for User	○
AWA2218		10-Way Switch Bracket LH, Mounted for User	○
AWA2227		10-Way Switch Bracket RH, Mounted for Attendant	£54
AWA2228		10-Way Switch Bracket LH, Mounted for Attendant	£54
LINX ACCESSORY			
AWA2252		Loom providing jack sockets for Drive and Mode (includes connection for 10-way switch)	£86
FURTHER LiNX COMPONENTS / ACCESSORIES			
AWA2300		To be enclosed only (Accessory loom 2252 required if no REM400 or REM500 is selected)	○
BUDDY BUTTONS			
AWA1685		Red Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	£75
AWA1686		Yellow Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	£75
AWA1687		Green Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	£75
AWA1690		Blue Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	£75

EGG SWITCH		
AWA2254		Egg Switch, black (to be enclosed only (Accessory loom 2252 required if no REM400 or REM500 is selected) £75
PIKO BUTTONS		
AWA1711		Red Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £135
AWA1712		Yellow Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £135
AWA1713		Green Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £135
AWA1714		Orange Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £135
AWA1715		Blue Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £135
AWA1716		Black Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £135
TOGGLES / BUTTONS		
AWZ2261		Stereo Toggle - direct control for Tilt only £230
AWZ2262		Stereo Button - direct control for Tilt only £219
AWZ2263		Four-Way Toggle - direct control for Tilt & Recline £286
AWZ2264		Four-Way Button - direct control for Tilt & Recline £306
AWZ2207		Toggle / Button RH Mount ○
AWZ2208		Toggle / Button LH Mount ○
AWZ2225		Without Mounting ○
BODYPOINT JOYSTICK KNOBS		
AWA2721		3" U-Shaped, Joystick Handle, Bodypoint PC101A £62
AWA2722		4" U-Shaped, Joystick Handle, Bodypoint PC102A £62
AWA2723		Domed Rubber, Joystick Handle, Bodypoint PC107A £62
AWA2724		Mushroom Grooved, Joystick Handle, Bodypoint PC110A £62
REMOTE STOP SWITCH		
AWZ2258		Remote Stop Switch (ASL504) £897
SHROUD COLOUR*		
AWA0202		Blue ○
AWA0228		Black ○
AWA0204		Transparent Cherry Red ○
AWA0222		White ○
AWA0216		Olive Green ○
AWZ0299		Special Shroud Colour (RAL Code must be supplied) POA
FRONT WHEEL*		
AWA0840		10" Grey Pneumatic ○
AWA0858		10" Grey Puncture Proof Sealed ○
AWA0842		10" Grey Puncture Proof Solid Inlay ◆
REAR WHEEL*		
AWA0771		14" Grey Pneumatic ○
AWA0773		14" Grey Puncture Proof Sealed ○
AWA0750		14" Grey Puncture Proof Solid Inlay ◆
MOTOR TYPE*		
AWA1277		6 km/h ◆

BATTERY*			
AWA1364	<input type="checkbox"/>	73.5 AH Gel	◆
SUSPENSION*			
AWA1203	<input type="checkbox"/>	Heavy Duty Suspension	◆
PERFORMANCE ACCESSORIES			
AWA1215	<input type="checkbox"/>	G-Trac (ensures the powerchair takes the intended direction as indicated by the user control, without the need for continual correction by the client on uneven or sloping terrain)	£703
CHARGER*			
AWA1343	<input type="checkbox"/>	OptiCharge Charger 8A (200 - 250V)	◆
AWA1313	<input type="checkbox"/>	OptiCharge Charger 10A (200 - 250V) Charging time will be 20% faster	£130
LIGHTS*			
AWA1254	<input type="checkbox"/>	No Lights	◆
AWA1437	<input type="checkbox"/>	LED Light system with indicators (30 LUX)	£410
TRANSPORTATION KIT*			
AWA1494	<input type="checkbox"/>	4 Point Tie Down (Tested version for using wheelchair as a vehicle seat according to ISO7176-19:2001) Adequate headrest is compulsory	◆
FIXATION BELT*			
AWA1455	<input type="checkbox"/>	Fixation Belt Extra long with Auto Lock	◆
REAR VIEW MIRROR			
AWA1415	<input type="checkbox"/>	Rear view Mirror RH - Must be selected with 10 km/h Motor (NC 2217)	£36
AWA1339	<input type="checkbox"/>	Rear view Mirror LH - Must be selected with 10 km/h Motor (NC 2218)	£36
LUGGAGE RACK			
AWA1394	<input type="checkbox"/>	Luggage Rack	£281
ACCESSORIES			
AWA1312	<input type="checkbox"/>	Cane holder RH (In combination with Attendant bracket or Dual control the cane holder has to be mounted on the opposite side of the remote)	£36
AWA1363	<input type="checkbox"/>	Cane holder LH (In combination with Attendant bracket or Dual control the cane holder has to be mounted on the opposite side of the remote)	£36
AWA1417	<input type="checkbox"/>	Fenders for rear wheels	£104

TOTAL PRICE	
--------------------	--

Non Conforming Options:

AWA0564	Reinforced hanger fixed angle (not detachable without tools)
AWA0503	Vari F : Angle adjustable with tools 0° - 70°, length 290mm - 460mm
AWA0599	Angle adjustable flip back heavy duty footboard
AWA0522	Angle adjustable footboard (for max user weight 200kg)
AWA2201	REM110 Drive only
AWA2227	Bracket RH mounted for Attendant
AWA2228	Bracket LH mounted for Attendant
AWA2220	REM050 Dual Control (Please also select user control and bracket) (Please select dual control bracket)
AWA2219	Compact Dual Control (ACU200) (Please select dual control bracket)
AWA2245	Invacare LiNX IDC Dual control (Please select IDC dual control bracket below)
AWA2217	10-Way Switch Bracket RH, Mounted for User
AWA2218	10-Way Switch Bracket LH, Mounted for User
AWA2227	10-Way Switch Bracket RH, Mounted for Attendant
AWA2228	10-Way Switch Bracket LH, Mounted for Attendant
AWA1437	LED Light system with indicators (30 LUX)



For a more complex customer requirement
please refer to our Invacare Unique special order form.
email specialsuk@invacare.com

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):