



Yes, you can:

Powerchairs Ireland Prescription Form

NOVEMBER 2023
Customer Ref:

Please refer to www.invacare.ie for Local Distributor

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Storm 4 Max

Standard Features

- Max User Weight 220KG (34.5 Stone)
- Max Seating
- Seat Width from 490mm - 610mm
- Seat Depths 560mm - 630mm
- multi stretch PU coated fabric Upholstery
- Seat Height 480mm - 510mm
- Manual Settable Tilt with Manual Settable Recline
- Reinforced Armrests
- Reinforced Legrests with Angle adjustable flip-back footboard
- Headrest, sideways and depth adjustable
- LiNX Control REM110 with standard bracket
- Choice of 6 Colours
- Puncture Proof Castors
- Puncture Proof Drive Wheels
- 6 km/h (4mph) Motors
- 8A (200 - 250V) Charger
- 4 Point tie down system
- 73.5 Ah (Gel) batteries
- Bolted on anti tippers

Standard Base price: €11,358



Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code	Description	
CHAIR TYPE*		
AWA0000	Storm 4 Max	€11,358
AWZ0000	Special Storm 4 Max (Please provide details in special section)	POA
SEATING TYPE*		
AWA1145	MAX SEAT UNIT	◆
BACK CUSHION*		
AWA1106	Tension adjustable 630mm - 620mm high with SW55 / 610mm 570mm High with SW490mm	◆

SEAT WIDTH*			
AWA0378	<input type="checkbox"/>	Seat Width 550mm	○
AWA0380	<input type="checkbox"/>	Seat Width 610mm	○
SEAT CUSHION*			
AWA1459	<input type="checkbox"/>	Seat cushion Max (550mm x 635mm)	○
AWA1460	<input type="checkbox"/>	Seat cushion Max (610mm x 635mm)	○
SEAT DEPTH SETTABLE IN 3 STEPS*			
AWA0396	<input type="checkbox"/>	Seat Depth 560mm	○
AWA0397	<input type="checkbox"/>	Seat Depth 595mm	○
ADAPTER TO REDUCE SEAT DEPTH			
AWA0399	<input type="checkbox"/>	Adapter to shorten the seat depth by 130mm	€62
UPHOLSTERY TYPE*			
AWA0278	<input type="checkbox"/>	multi stretch PU coated fabric	◆
ARMRESTS INCLUDING REINFORCED DEPTH ADJUSTABLE BRACKET*			
AWA0635	<input type="checkbox"/>	Fixed armrest 250mm - 340mm height adjustable removable	○
ARMPAD*			
AWA0647	<input type="checkbox"/>	Standard Armpad 325mm	○
AWA0649	<input type="checkbox"/>	Armpad extended plus 50mm	○
PUSH BAR*			
AWA1307	<input type="checkbox"/>	Push bar - height adjustable	○
HEADREST / NECKREST - Adequate headrest option must be selected*			
AWA1432	<input type="checkbox"/>	Headrest sideways and depth adjustable	◆
AWA1427	<input type="checkbox"/>	Headrest with cheek support sideways adjustable	€254
AWA1495	<input type="checkbox"/>	Equivalent non Invacare Headrest supporting the requirements of ISO7176-19:2001 Provided independently from this order	○
SEAT HEIGHT*			
AWA1025	<input type="checkbox"/>	Fixed seat height 480mm	○
AWA1006	<input type="checkbox"/>	Fixed seat height 510mm	○
SEAT TILT*			
AWA1002	<input type="checkbox"/>	Manual Seat Tilt by spindle	○
AWA1014	<input type="checkbox"/>	Electrical Tilt 0° to 15°	€664
BACKREST RECLINER*			
AWA0958	<input type="checkbox"/>	Settable recliner by spindle (0° - 25°)	○
AWA0947	<input type="checkbox"/>	Electrical recliner (with high pivot, shearing compensation)	€585
FOOTREST HANGER TYPE (PAIR)*			
AWA0564	<input type="checkbox"/>	Reinforced hanger fixed angle (not detachable without tools) (NC 0522)	◆
CALF PADS			
AWA0583	<input type="checkbox"/>	Calf Pads adjustable	○
FOOTBOARD*			
AWA0599	<input type="checkbox"/>	Angle adjustable flip back heavy duty footboard	◆
AWA0522	<input type="checkbox"/>	Angle adjustable footboard (for max user weight 200kg) (NC 0564)	○
CONTROLS*			
AWA2201	<input type="checkbox"/>	REM110 Drive only (NC 1437)	◆
AWA2203	<input type="checkbox"/>	REM211 Drive & Actuators (NC 1437)	○
AWA2204	<input type="checkbox"/>	REM216 Drive, Actuators & Lights	€162

CONTROL BRACKETRY*			
AWA1207		Standard Bracket RH	○
AWA1208		Standard Bracket LH	○
AWA1707		Swing Away RH Bracket	€159
AWA1708		Swing Away LH Bracket	€159
AWA2227		Bracket RH mounted for Attendant (NC 2219, 2220)	€71
AWA2228		Bracket LH mounted for Attendant (NC 2219, 2220)	€71
DUAL CONTROL			
AWA2220		REM050 Dual Control (Please also select user control and bracket) (Please select dual control bracket) (NC 2227, 2228, 2245)	€489
AWA2219		Compact Dual Control (ACU200) (Please select dual control bracket) (NC 2227, 2228, 2245, 2201)	€466
DUAL CONTROL BRACKETEY			
AWA1207		Dual Control Standard Bracket RH (NC 2227, 2228, 2245)	○
AWA1208		Dual Control Standard Bracket LH	○
INVACARE LiNX IDC DUAL CONTROL			
AWA2245		Invacare LiNX IDC Dual control (Please select IDC dual control bracket below) (NC 2227, 2228, 2220, 2219, 2201)	€748
INVACARE LiNX IDC DUAL CONTROL BRACKETRY			
AWA1221		IDC Dual control Detachable Bracket RH (NC 2227, 2228, 2220, 2219)	○
AWA1222		LiNX IDC Dual control Detachable Bracket LH (NC 2227, 2228, 2220, 2219)	○
LINX 10-WAY SWITCH			
AWA2251		10-Way Switch DLX-SW10-A for direct Powered Seating Access (included loom to connect + additional jack for mode switch)	€539
LINX 10-WAY SWITCH BRACKET			
AWA2217		10-Way Switch Bracket RH, Mounted for User	○
AWA2218		10-Way Switch Bracket LH, Mounted for User	○
AWA2227		10-Way Switch Bracket RH, Mounted for Attendant	€32
AWA2228		10-Way Switch Bracket LH, Mounted for Attendant	€32
LINX ACCESSORY			
AWA2252		Loom providing jack sockets for Drive and Mode (includes connection for 10-way switch)	€110
FURTHER LiNX COMPONENTS / ACCESSORIES			
AWA2300		To be enclosed only (Accessory loom 2252 required if no REM400 or REM500 is selected)	○
BUDDY BUTTONS			
AWA1685		Red Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	€90
AWA1686		Yellow Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	€90
AWA1687		Green Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	€90
AWA1690		Blue Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	€90
EGG SWITCH			
AWA2254		Egg Switch, black (to be enclosed only (Accessory loom 2252 required if no REM400 or REM500 is selected)	€90

PIKO BUTTONS			
AWA1711	<input type="checkbox"/>	Red Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	€161
AWA1712	<input type="checkbox"/>	Yellow Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	€161
AWA1713	<input type="checkbox"/>	Green Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	€161
AWA1714	<input type="checkbox"/>	Orange Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	€161
AWA1715	<input type="checkbox"/>	Blue Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	€161
AWA1716	<input type="checkbox"/>	Black Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	€161
TOGGLES / BUTTONS			
AWZ2261	<input type="checkbox"/>	Stereo Toggle - direct control for Tilt only	€263
AWZ2262	<input type="checkbox"/>	Stereo Button - direct control for Tilt only	€263
AWZ2263	<input type="checkbox"/>	Four-Way Toggle - direct control for Tilt & Recline	€342
AWZ2264	<input type="checkbox"/>	Four-Way Button - direct control for Tilt & Recline	€365
AWZ2207	<input type="checkbox"/>	Toggle / Button RH Mount	○
AWZ2208	<input type="checkbox"/>	Toggle / Button LH Mount	○
AWZ2225	<input type="checkbox"/>	Without Mounting	○
BODYPOINT JOYSTICK KNOBS			
AWA2721	<input type="checkbox"/>	3" U-Shaped, Joystick Handle, Bodypoint PC101A	€75
AWA2722	<input type="checkbox"/>	4" U-Shaped, Joystick Handle, Bodypoint PC102A	€75
AWA2723	<input type="checkbox"/>	Domed Rubber, Joystick Handle, Bodypoint PC107A	€75
AWA2724	<input type="checkbox"/>	Mushroom Grooved, Joystick Handle, Bodypoint PC110A	€75
REMOTE STOP SWITCH			
AWZ2258	<input type="checkbox"/>	Remote Stop Switch (ASL504)	€1,070
SHROUD COLOUR*			
AWA0202	<input type="checkbox"/>	Blue	○
AWA0228	<input type="checkbox"/>	Black	○
AWA0220	<input type="checkbox"/>	Champagne	○
AWA0204	<input type="checkbox"/>	Transparent Cherry Red	○
AWZ0299	<input type="checkbox"/>	Special Shroud Colour (RAL Code must be supplied)	POA
FRONT WHEEL*			
AWA0840	<input type="checkbox"/>	10" Grey Pneumatic	○
AWA0858	<input type="checkbox"/>	10" Grey Puncture Proof Sealed	○
AWA0842	<input type="checkbox"/>	10" Grey Puncture Proof Solid Inlay	◆
REAR WHEEL*			
AWA0771	<input type="checkbox"/>	14" Grey Pneumatic	○
AWA0773	<input type="checkbox"/>	14" Grey Puncture Proof Sealed	○
AWA0750	<input type="checkbox"/>	14" Grey Puncture Proof Solid Inlay	◆
MOTOR TYPE*			
AWA1277	<input type="checkbox"/>	6 km/h	◆
BATTERY*			
AWA1364	<input type="checkbox"/>	73.5 AH Gel	◆
SUSPENSION*			
AWA1203	<input type="checkbox"/>	Heavy Duty Suspension	◆

PERFORMANCE ACCESSORIES

AWA1215		G-Trac (ensures the powerchair takes the intended direction as indicated by the user control, without the need for continual correction by the client on uneven or sloping terrain)	€908
---------	--	---	------

CHARGER*

AWA1343		OptiCharge Charger 8A (200 - 250V)	◆
AWA1313		OptiCharge Charger 10A (200 - 250V) Charging time will be 20% faster	€174

LIGHTS*

AWA1254		No Lights	◆
AWA1437		LED Light system with indicators	€503

TRANSPORTATION KIT*

AWA1494		4 Point Tie Down (Tested version for using wheelchair as a vehicle seat according to ISO7176-19:2001) Adequate headrest is compulsory	◆
---------	--	--	---

FIXATION BELT*

AWA1455		Fixation Belt Extra long with Auto Lock	◆
---------	--	---	---

ACCESSORIES

AWA1312		Cane holder RH (In combination with Attendant bracket or Dual control the cane holder has to be mounted on the opposite side of the remote)	€43
AWA1363		Cane holder LH (In combination with Attendant bracket or Dual control the cane holder has to be mounted on the opposite side of the remote)	€43
AWA1417		Fenders for rear wheels	€134

TOTAL PRICE

Non Conforming Options:

AWA0564	Reinforced hanger fixed angle (not detachable without tools)
AWA0522	Angle adjustable footboard (for max user weight 200kg)
AWA2201	REM110 Drive only
AWA2227	Bracket RH mounted for Attendant
AWA2228	Bracket LH mounted for Attendant
AWA2220	REM050 Dual Control (Please also select user control and bracket) (Please select dual control bracket)
AWA2219	Compact Dual Control (ACU200) (Please select dual control bracket)
AWA2245	Invacare LiNX IDC Dual control (Please select IDC dual control bracket below)
AWA2227	10-Way Switch Bracket RH, Mounted for Attendant
AWA2228	10-Way Switch Bracket LH, Mounted for Attendant
AWA1437	LED Light system with indicators



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):