

# INVACARE® TOP END® T-5 7000 Series Tennis



Yes, you can.®



# INVACARE® TOP END® T-5 7000 Series Tennis



Michael Jeremiaz - France

Invacare **Top End T-5 7000 Series Tennis** steals the game for tennis lovers. Choose the tennis chair to suit you, fully welded to your precise requirements in 7000 Series aluminium: with the **Top End T-5 7000 Series Tennis** the ball's in your court!

## New material

Comparable in weight with titanium frames, the new **Invacare Top End T-5 7000 Series Tennis** in 7000 Series aluminium is stronger, lighter and more rigid than ever. The strength of 7000 aluminium alloy is even further enhanced by **Invacare Top End's** sophisticated heat treating and powder coating process.

## Game, Set and Match

Professional players know that foot placement is key in wheelchair tennis. Thanks to the range of footrest choices, the **Top End T-5 7000 Series Tennis** is specially designed for high-level players, helping them move with minimal effort. Athletes can also choose to maximise their performance through the use of side-guards, click straps and toe clips.

For professional players who don't want to compromise!

## OPTIONS AND FEATURES



Fixed swivel anti-tip, height adjustable.



Sleek curved frame members add a sharp look, for more rigidity and durability.



Custom front frame angle, foot space width, axle position and wheelbase.



Adjustable COG option.



Robin Ammerlaan - NL

## TECHNICAL DATA



12" - 21"  
(approx. 305 - 533 mm)



12" - 19"  
(approx. 305 - 483 mm)



Front: 14.5" - 21"  
(approx. 368 - 533 mm)  
Rear: 14.5" - 21"  
(approx. 368 - 533 mm)



Fixed: 8" - 15.5"  
(approx. 203 - 398 mm)  
Adj.: 9" - 18"  
(approx. 229 - 457 mm)



Max. 250 lbs  
(approx. 113 kg)

Please refer to order form for all technical specifications.



Invacare Limited  
Pencoed Technology Park - Pencoed  
Bridgend CF35 5AQ - United Kingdom  
Tel: +44 1 656 776 222  
Fax: +44 1656 776 220  
E-mail: [uk@invacare.com](mailto:uk@invacare.com)  
Sales Order E-mail:  
[ordersuk@invacare.com](mailto:ordersuk@invacare.com)  
[www.invacare.co.uk](http://www.invacare.co.uk)

Invacare Ireland Ltd  
Unit 5 - Seatown Business Campus  
Seatown Road - SWORDS  
County Dublin - Ireland  
Tel: +353 1 8107084  
Fax: +353 1 8107085  
Email: [Ireland@invacare.com](mailto:Ireland@invacare.com)  
[www.invacare.ie](http://www.invacare.ie)

© 2012 Invacare International Sàrl  
All rights reserved.

All information quoted is believed  
to be correct at time of print.  
Invacare reserves the right to alter product  
specifications without prior consultation.

TopEnd T-5 Tennis - UK - 03/2012 - 1524494



**Yes, you can.®**