



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

Power Retail Prescription Form

SEPTEMBER 2024
 Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

TDX SP2 NB Teen

Standard Features

- Max User Weight 23.5 Stone
- Compact, manoeuvrable, CWD TDX SP2 Narrow Base
- Modulite highly adjustable seating system
- Available in Seat Widths from 305mm - 380mm (12" - 15")
- Seat Depths 305mm - 560mm Stepless Adjustment (12" - 22")
- Tension adjustable backrest
- Standard T-Armrests, detachable with Short armpad
- Manual settable seat Tilt 0° - 20°
- Manual settable backrest Recline 0° - 30°, in 3.75° steps
- Standard 80 hangers with fixed footplate
- LiNX Control REM110 with standard bracket
- Selection of 10 Stunning Shroud and Rim colours
- Centre (Driving Wheels) 355mm x 76mm (14" x 3") Black Tyres
- 6 m/h (4mph) 4 Pole Motors and 50Ah Batteries
- 8A (200 - 250V) Charger
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

Standard Base price: £5,370



Q'STRAIT
QLK-150
 Compatible

Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code	Description	
CHAIR TYPE*		
ARZ0000	TDX-SP2 NB Modulite Teen	£5,370
ARZ0000	Special TDX SP2 NB Modulite Teen (Please provide details in special section)	◆
FITTING*		
SWO0025	Invacare ONE to ONE fitting service	POA
SEATING TYPE*		
ARM2100	Modulite Seating Unit (Flex 3 seating unit)	◆

SEATING FRAME*			
ARM2111		Telescopic Seat Frame with metal seat plate (with stepless adjustment on armrests, legrests and backrest)	◆
SEAT WIDTH*			
ARM0001		305mm (12")	○
ARM0002		330mm (13")	○
ARM2121		380mm (15")	○
SEAT DEPTH RAIL LENGTH*			
ARM2129		Standard Depth Package (16" - 20", 410mm - 510mm simple stepless adjustment from the side)	○
ARM2190		Extended Depth Package (20" - 22", 510mm - 560mm simple stepless adjustment from the side) (NC 2134)	○
SEAT DEPTH*			
ARM0003		305mm (12")	○
ARM0004		330mm (13")	○
ARM0005		355mm (14")	○
ARM0006		380mm (15")	○
ARM2132		410mm (16")	○
ARM2134		460mm (18")	○
ARM2136		510mm (20")	○
ARM2138		560mm (22")	○
SEAT OPTION			
ARMSWO		"H" Frame Interface Kit (NC 0002)	£446
SEAT CUSHION*			
ARM2142		Matrx Flo-tech Lite (75mm up to 127kg) (NC 0001, 0002)	£112
ARM2147		Matrx Libra 90mm thick (NC 0001, 0002)	£391
1626906		Matrx Mini series - PS cushion with Kidabra print 310mm x 310mm, 12" x 12" (max user weight 68kg)	£232
1626907		Matrx Mini series - PS cushion with Kidabra print 310mm x 360mm, 12" x 14" (max user weight 68kg)	£232
1626908		Matrx Mini series - PS cushion with Kidabra print 310mm x 410mm, 12" x 16" (max user weight 68kg)	£232
1627456		Matrx Mini series - PS cushion with Kidabra print 330mm x 310mm, 13" x 12" (max user weight 68kg)	£232
1627457		Matrx Mini series - PS cushion with Kidabra print 330mm x 360mm, 13" x 14" (max user weight 68kg)	£232
1627458		Matrx Mini series - PS cushion with Kidabra print 330mm x 410mm, 13" x 16" (max user weight 68kg)	£232
1626909		Matrx Mini series - PS cushion with Kidabra print 360mm x 310mm, 14" x 12" (max user weight 68kg)	£232
1626910		Matrx Mini series - PS cushion with Kidabra print 360mm x 360mm, 14" x 14" (max user weight 68kg)	£232
1626911		Matrx Mini series - PS cushion with Kidabra print 360mm x 410mm, 14" x 16" (max user weight 68kg)	£232
ARM2140		Without seat cushion	◆
BACKREST SUPPORT*			
ARM2162		Backposts 540mm (Tall) (must be selected for Flex3 Laguna backrest)	◆
ARM2163		Backposts Short (480mm) (NC 2173)	○
PUSH HANDLES FOR BACKPOSTS*			
ARM2107		Backposts WITH push handles	◆
ARM2108		Backposts WITHOUT push handles, Push handles will be cut down on the vertical section. (NC 2219, 2227, 2228, 2245, 2171, 1495)	£36

BACKREST UPHOLSTERY / CUSHION*			
ARM2171		Basic Tension Adjustable Backrest (540mm or 480mm upholstery height) Straps easily accessible and adjustable from the rear, Black fabric only	(NC 0002, 2155, 1496, 2229, 1445, 1288, 1289) ◆
ARM2172		Flex3 adjustable Backrest Plates (without Backrest Cushion)	£224
ARM2173		Flex3 adjustable Backrest Plates with Laguna Backrest cushion (maximum seat depth achievable 510mm (20"))	£363
SP1562143		Mistral Flex3 Backrest shell - shorter & narrower than the Laguna Flex4	(NC 0002) £159
SP1431668-TR26		Mistral Backrest, multi stretch PU coated fabric Black Cushion	(NC 0002) £145
SP1431668-TR35		Mistral Backrest, Black Plush Cushion	(NC 0002) £163
EMP		Extended Mounting Pins	£519
MINI1207		Matrx Mini series - Mini Back with Kidabra print 310mm x 180mm, 12" x 7" (max user weight 68kg) (please also order the extending Mounting Pins)	£519
MINI1210		Matrx Mini series - Mini Back with Kidabra print 310mm x 260mm, 12" x 10" (max user weight 68kg) (please also order the extending Mounting Pins)	£519
MINI1212		Matrx Mini series - Mini Back with Kidabra print 310mm x 310mm, 12" x 12" (max user weight 68kg) (please also order the extending Mounting Pins)	£519
MINI1214		Matrx Mini series - Mini Back with Kidabra print 310mm x 360mm, 12" x 14" (max user weight 68kg) (please also order the extending Mounting Pins)	£519
MINI1307		Matrx Mini series - Mini Back with Kidabra print 330mm x 180mm, 13" x 7" (max user weight 68kg) (please also order the extending Mounting Pins)	£519
MINI1310		Matrx Mini series - Mini Back with Kidabra print 330mm x 260mm, 13" x 10" (max user weight 68kg) (please also order the extending Mounting Pins)	£519
MINI1312		Matrx Mini series - Mini Back with Kidabra print 330mm x 310mm, 13" x 12" (max user weight 68kg) (please also order the extending Mounting Pins)	£519
MINI1314		Matrx Mini series - Mini Back with Kidabra print 330mm x 360mm, 13" x 14" (max user weight 68kg) (please also order the extending Mounting Pins)	£519
MINI1407		Matrx Mini series - Mini Back with Kidabra print 360mm x 180mm, 13" x 14" (max user weight 68kg) (please also order the extending Mounting Pins)	£519
MINI1410		Matrx Mini series - Mini Back with Kidabra print 360mm x 260mm, 13" x 14" (max user weight 68kg) (please also order the extending Mounting Pins)	£519
MINI1412		Matrx Mini series - Mini Back with Kidabra print 360mm x 310mm, 13" x 14" (max user weight 68kg) (please also order the extending Mounting Pins)	£519
MINI1414		Matrx Mini series - Mini Back with Kidabra print 360mm x 360mm, 14" x 14" (max user weight 68kg) (please also order the extending Mounting Pins)	£44
ARM2170		Without Backrest Upholstery	(NC 2151, 2152, 2155) ○
BACKREST OPTIONS			
SP1571542		Angle adjustable backrest support brackets RH	£44
SP1571533		Angle adjustable backrest support brackets LH	£36
UPHOLSTERY COLOUR*			
ARM2151		Black fabric	◆
ARM2152		Black multi stretch PU coated fabric	(NC 2170) ○
ARM2155		Black Spacer fabric	(NC 2147, 2171) ○
ARZ1005		Matrx Mini In Kidabra fabric	(NC 2172, 2173) ○
ARZ1006		Matrx Minin In black fabric	(NC 2172, 2173) ○
ARZ1007		Select if No upholstery is required	(NC 2172, 2173) ○

ARMREST*			
ARM2191		T-Armrest Standard, Range of adjustment from seat plate 245mm to 315mm	◆
ARM2181		T-Armrest Extended height, Range of adjustment from seat plate 295mm - 365mm.	£36
ARM2183		T-Armrest In-rigged (reduces seat width by 25mm each side) (Must be selected with 13" & 12" seat width)	£25
SP1582159		Inrigged Armrests for 12" seat width, Standard inrigged arm option MUST BE SELECTED with additional brackets (must be selected with 12" seat width)	£37
ARMPAD*			
ARM2195		Short Armpad (325mm)	◆
ARM2196		Long Armpad (375mm)	○
HIP SUPPORTS LH (Not available with 410mm seat depth and Manual Centre legrests)			
ARM2184		SMALL 140mm x 100mm, for flip up and following armrests only (NC 0556, 0589)	£226
ARM2186		MEDIUM 230mm x 100mm, for flip up and following armrests only (NC 0556, 0589)	£226
HIP SUPPORTS RH (Not available with 410mm seat depth and Manual Centre legrests)			
ARM2185		SMALL 140mm x 100mm, for flip up and following armrests only (NC 0556, 0589)	£226
ARM2187		MEDIUM 230mm x 100mm, for flip up and following armrests only (NC 0556, 0589)	£226
MOUNTING HARDWARE FOR HIP SUPPORTS			
ARM2165		Mounting hardware for hip supports (8.1" / 206mm high)	○
HEADREST / NECKREST - Adequate headrest option must be selected*			
1603255		Elan Standard Pad - Small (160mm x 90mm) Standard cover (please select mounting Hardware type for Elan Headrest)	£124
1603254		Elan Standard Pad - Medium (260mm x 130mm) Standard cover (please select mounting Hardware type for Elan Headrest)	£124
1603291		Elan Occipital Pad (230mm x 180mm) Standard cover (please select mounting Hardware type for Elan Headrest)	£232
1629822		Elan 4 Point Headrest Pad (Small) Standard cover (please select mounting Hardware type for Elan Headrest)	£184
ARM1432		Headrest, height, depth and sideways adjustable (NC 2247)	£203
ARM1496		Non Invacare Headrest but profile mount fitting Flex3 hard back	○
ARM1495		Push bar dia. 22mm in conjunction with telescopic seat frame / backposts, No Headrest - no mounting hardware.	○
ELAN MOUNTING HARDWARE			
1603250		Matrx Elan 6-point Hardware Infinite adjustments	£184
1603251		Matrx Elan 6-point Mini (paediatric 100mm shorter post)	£166
1603249		Matrx Standard multi-axis hardware Limited adjustability	£166
1603252		Elan adapter for Flex3 back	£140
LATERAL SUPPORT			
ARM1445		Multi-adjustable pad (Azalea type, multi stretch PU coated fabric) (not available with Matrx Backrests) (NC 2171)	£150
MATRX MINI LATERAL SUPPORTS			
1627466		Swing-away lateral Hardware (each)	£86
1544292		Extra small lateral Pad (each)	£59
1627465		Chest harness interface Kit	£70

SEAT LIFTER / SEAT TILT*			
ARM1091		Manually settable Seat Tilt (0° - 20°) Seat Height 450mm	◆
ARM1092		Electric Tilt (0° - 20°) Seat Height 450mm	£362
ARM1093		Electric Tilt (0° - 30°) with CoG shift Seat Height 450mm	£961
ARM1094		Seat lifter with Electrical Tilt (300mm, 30°) with CoG shift Seat Height 450mm - 480mm (Also choose option 1265)	£1,696
SEAT HEIGHT SETTING FOR 20 DEGREES SEAT TILT*			
ARM0354		Seat height set at 450mm (17.75") (from floor to front of seat plate at 0° tilt) (NC 1094)	○
ARZ1003		Seat height set at 460mm (18") (from floor to front of seat plate at 0° tilt) (NC 1092)	○
ARM0364		Seat height set at 475mm (18.75") (from floor to front of seat plate at 0° tilt) (NC 1093, 1094)	○
SEAT TILT OPTION			
ARM2054		5° Fixed Anterior Tilt Assist mount (will reduce max tilt by 5°)	£290
BACKREST RECLINER*			
ARM0975		Manual Settable Recliner (0° - 30°, in 3.75° steps) Foldable	◆
ARM0977		Powered Recliner (0° - 30°) Foldable	£453
FOOTREST HANGER TYPE*			
ARM0556		Centre Mounted: Manual : Hanger, centre mounted, angle settable (0° - 90°) : knee to heel length 390 - 510mm, (including calf pad and angle adjustable footplate) (NC 0504, 0513)	£334
ARM0589		Centre Mounted : LNX Powered: Centre mounted legrest, length compensating (including calf pad and angle adjustable footplate) (NC 0512, 0504, 0513)	£1,470
ARM0526		Side Mounted : Standard hanger 80° (NC 0559, 0504)	◆
1606173-0601		Side Mounted : Black short 90° Legrest Hanger x 2 - Action style (please select legrest support carriers below)	£136
1606172-0601		Side Mounted : Black short 80° Legrest Hanger x 2 - Action style (please select legrest support carriers below)	£136
ARM0501		Without Hanger (NC 0559, 0504, 0513)	○
LEGREST SUPPORT			
SP1531195		Legrest support carrier LH (to be selected for the Action style hangers, 1606173-0601 and 1606172-0601)	£39
SP1531196		Legrest support carrier RH (to be selected for the Action style hangers, 1606173-0601 and 1606172-0601)	£39
AMPUTEE SUPPORT			
ARM0555		Amputee Support LH (this replaces the LH Hanger)	£174
ARM0554		Amputee Support RH (this replaces the RH Hanger)	£174

FOOTPLATE TYPE*				
ARM0559		Angle Adjustable Footplate Medium 140mm x 240mm (5.5" x 9.5") for Centre Mounted Hanger	(NC 0501)	○
ARM0512		Fixed Footplate	(NC 0556, 0589, 0501)	◆
ARM0504		Angle Adjustable Footplate	(NC 0589, 0501, 0526)	£64
ARM0513		Angle and depth adjustable Footplate (Adapt)	(NC 0589, 0501)	○
1428660		Angle and depth adjustable Footplate LH, including stem (for use with short 80° & 90° Action style Legrest)		£37
1428654		Angle and depth adjustable Footplate RH, including stem (for use with short 80° & 90° Action style Legrest)		£37
1530620		Fixed Footplate LH, including stem (for use with short 80° & 90° Action style Legrest)		£19
1530614		Fixed Footplate RH, including stem (for use with short 80° & 90° Action style Legrest)		£19
1492150		Upmounted Footplate LH (for use with short 80° & 90° Action style Legrest)		£74
1492149		Upmounted Footplate RH (for use with short 80° & 90° Action style Legrest)		£63
ARMXXXX		Without footplates for option 0501 without hanger	(NC 0589, 0526)	○
CONTROLS*				
ARM2201		REM110 Drive only	(NC 1093, 1094, 0977, 0589, 1437, 2209, 2210, 2230, 2247, 2233, 2240, 2219, 2245, 2270)	◆
ARM2203		REM211 Drive & Actuators	(NC 2209, 2210, 2230, 2247, 2233, 2240, 2270)	○
ARM2204		REM216 Drive, Actuators & Lights	(NC 2210, 2230, 2247, 2233, 2240, 2270)	£137
ARM2205		REM400, touch-screen Interface, built in Bluetooth HID	(NC 2210)	£305
ARM2206		REM500, without joystick for alternative access, built in Bluetooth HID (Secondary User input required)	(NC 2228, 2270)	£695
CONTROL BRACKETRY*				
ARM1207		Standard Bracket RH	(NC 2206)	○
ARM1208		Standard Bracket LH	(NC 2206)	○
ARM1707		Swing Away RH Bracket	(NC 2206)	£119
ARM1708		Swing Away LH Bracket	(NC 2206)	£119
ARM1709		Matrx Resolve Swing in/out bracket RH	(NC 2213, 2235, 2213, 2241, 1404)	£219
ARM1710		Matrx Resolve Swing in/out bracket LH	(NC 2214, 2236, 2214, 2242, 1405)	£219
ARM2213		Centre mount swing away to RH Nucleus Bracket	(NC 1405)	£446
ARM2214		Centre mount swing away to LH Nucleus Bracket	(NC 1405)	£446
ARM2227		Bracket RH mounted for Attendant		£57
ARM2228		Bracket LH mounted for Attendant		£57
ARZ2239		Individual mounting option per (special) - Requires additional information / specification		POA
ARM2209		REM500 Swing away RH multi-adjustable Bracket		○
ARM2210		REM500 Swing away LH multi-adjustable Bracket		○
REM400 TOGGLE OPTION				
ARM2270		Assignable Toggle Switch Module for REM400, (as default set to left ON / OFF, right speed dial)	(NC 2203, 2204)	£284
CHIN CONTROL				
ARM2229		Proportional Chin Joystick (ASL138) (Requires Primary Remote REM400 or REM500) (Please select chin control bracket below)	(NC 2203, 2204)	£1,787
CHIN CONTROL BRACKETRY				
ARM1289		Chin Control Swing Away Bracket RH Manual	(NC 2203, 2204, 1287, 1286)	○
ARM1288		Chin Control Swing Away Bracket LH Manual	(NC 2203, 2204, 1287, 1286)	○
ARM1287		Chin Control Swing Away Bracket RH Powered	(NC 2203, 2204, 1287, 1286, 2227, 1404, 1312)	£2,251
ARM1286		Chin Control Swing Away Bracket LH Powered	(NC 2203, 2204, 1287, 1286, 2228, 1405, 1363)	£2,251

DUAL CONTROL				
ARM2220		REM050 Dual Control (Please select user control and bracket) (Please select dual control bracket below)	(NC 1222, 2227, 2228)	£355
ARH2219		Compact Dual Control (DLX-ACU200) (Please select user control and bracket) (Please select dual control bracket below)	(NC 1222, 2227, 2228, 2201)	£355
DUAL CONTROL BRACKETRY				
ARM1207		Dual Control Standard Bracket RH		○
ARM1208		Dual Control Standard Bracket LH		○
INVACARE LINX IDC DUAL CONTROL				
ARM2245		Invacare IDC Dual control (Please select IDC Dual control bracket below)	(NC 2220, 2227, 2228)	£673
INVACARE LINX IDC DUAL CONTROL BRACKETRY				
ARM1221		IDC Dual control Detachable Bracket RH	(NC 2220, 2227, 2228)	○
ARM1222		IDC Dual control Detachable Bracket LH	(NC 2220, 2227, 2228)	○
LINX COMPACT REMOTE				
ARM2231		Proportional Compact Remote DLX-CR400 (Requires Remote REM400 or REM500) (Please select Proportional Bracket)	(NC 2203, 2204)	£740
ARM2232		Proportional Compact Remote Low Force DLX-CR400LF (Requires Primary Remote REM400 or REM500) (Please select Proportional Bracketry)	(NC 2203, 2204)	£990
LINX COMPACT REMOTE BRACKET				
ARM2207		Proportional Compact RH Mount Bracket	(NC 2203, 2204)	○
ARM2208		Proportional Compact LH Mount Bracket	(NC 2203, 2204)	○
ARM2213		Proportional Compact Centre Mount swing away to RH (Nucleus)	(NC 2203, 2204)	£446
ARM2214		Proportional Compact Centre Mount swing away to LH (Nucleus)	(NC 2203, 2204)	£446
ARZ2215		Individual mounting option per request (Special) Requires additional information / specification	(NC 2203, 2204)	POA
LINX HEAD CONTROL				
ARM2230		Head Control with ultrasonic proximity sensors (ATOM , ASL104) with standard wings M3 (Requires Primary Remote REM400 or REM500)	(NC 2203, 2204)	£3,835
ARM2247		Head Control with ultrasonic proximity sensors (PROTON , ASL104P) with PROTON wings (Requires Primary Remote REM400 or REM500)	(NC 2203, 2204)	£4,491
LINX MICRO EXTREMITY CONTROL				
ARZ2233		MEC Proportional Joystick (Micro Extremity Control) (Please select MEC Bracketry) (Requires Primary Remote REM400 or REM500)	(NC 2203, 2204)	£3,723
LINX MICRO EXTREMITY CONTROL BRACKET				
ARZ2235		MEC Tray Mount RH (ASL603, 100mm x 150mm) (Please select MEC Bracketry)	(NC 2203, 2204)	○
ARZ2236		MEC Tray Mount LH (ASL603, 100mm x 150mm)	(NC 2203, 2204)	○
ARZ2213		Centre Mount, swing away to RH, Nucleus	(NC 2203, 2204)	£446
ARZ2214		Centre Mount, swing away to LH, Nucleus	(NC 2203, 2204)	£446
ARZ2241		Centre Mount, swing away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 2203, 2204)	£446
ARZ2242		Centre Mount, swing away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 2203, 2204)	£446
ARZ2215		Individual mounting option per request (Special) Requires additional information / specification	(NC 2203, 2204)	POA
LINX HEAVY DUTY REMOTE				
ARZ2240		Proportional Heavy Duty Remote Joystick with mode switch (ASL133) (Please select Heavy Duty Bracketry) (Requires Primary Remote REM400 or REM500)	(NC 2203, 2204)	£2,076

LINX HEAVY DUTY REMOTE BRACKET			
ARZ2235		Tray Mount RH, (ASL603, 100mm x 150mm)	(NC 2203, 2204) ○
ARZ2236		Tray Mount LH, (ASL603, 100mm x 150mm)	(NC 2203, 2204) ○
ARZ2213		Centre Mount, swing away to RH, Nucleus	(NC 2203, 2204) £446
ARZ2214		Centre Mount, swing away to LH, Nucleus	(NC 2203, 2204) £446
ARZ2241		Centre Mount, swing away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 2203, 2204) £446
ARZ2242		Centre Mount, swing away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 2203, 2204) £446
ARZ2215		Individual mounting option per request (Special) Requires additional information / specification	POA
LINX 10-WAY SWITCH			
ARM2251		10-Way Switch DLX-SW10-A for direct Powered Seating Access (included loom to connect + additional jack for mode switch) (Please select 10 way switch bracket below)	£453
LINX 10-WAY SWITCH BRACKET			
ARM2217		10-Way Switch Bracket RH, Mounted for User	(NC 2213, 2235, 2241) ○
ARM2218		10-Way Switch Bracket LH, Mounted for User	(NC 2214, 2236, 2242) ○
ARM2227		10-Way Switch Bracket RH, Mounted for Attendant	£57
ARM2228		10-Way Switch Bracket LH, Mounted for Attendant	£57
LINX ACCESSORY			
ARM2252		Loom providing jack sockets for Drive and Mode (includes connection for 10-way switch)	£90
LINX USB CHARGER PORT			
ARM2253		USB Charger Port DLX-USB02-A	£143
LINX USB CHARGER PORT BRACKET			
ARM2207		USB Charger Port RH Mount	○
ARM2208		USB Charger Port LH Mount	○
FURTHER LINX COMPONENTS / ACCESSORIES			
ARM2300		To be enclosed only (Accessory loom 2252 required if no REM400 or REM500 is selected)	○
ARZ2301		To be connected and programmed (Special)	POA
BUDDY BUTTONS			
ARM1685		Red Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	£79
ARM1686		Yellow Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	£79
ARM1687		Green Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	£79
ARM1690		Blue Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	£79
EGG SWITCH			
ARM2254		Egg Switch, black (to be enclosed only (Accessory loom 2252 required if no REM400 or REM500 is selected)	£79

PIKO BUTTONS		
ARM1711	<input type="checkbox"/>	Red Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £142
ARM1712	<input type="checkbox"/>	Yellow Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £142
ARM1713	<input type="checkbox"/>	Green Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £142
ARM1714	<input type="checkbox"/>	Orange Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £142
ARM1715	<input type="checkbox"/>	Blue Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £142
ARM1716	<input type="checkbox"/>	Black Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £142
TOGGLES / BUTTONS		
ARZ2261	<input type="checkbox"/>	Stereo Toggle - direct control for Tilt only £242
ARZ2262	<input type="checkbox"/>	Stereo Button - direct control for Tilt only £230
ARZ2263	<input type="checkbox"/>	Four-Way Toggle - direct control for Tilt & Recline £300
ARZ2264	<input type="checkbox"/>	Four-Way Button - direct control for Tilt & Recline £321
ARZ2207	<input type="checkbox"/>	Toggle / Button RH Mount ○
ARZ2208	<input type="checkbox"/>	Toggle / Button LH Mount ○
ARZ2225	<input type="checkbox"/>	Without Mounting ○
BODYPOINT JOYSTICK KNOBS		
ARM2721	<input type="checkbox"/>	3" U-Shaped, Joystick Handle, Bodypoint PC101A £65
ARM2722	<input type="checkbox"/>	4" U-Shaped, Joystick Handle, Bodypoint PC102A £65
ARM2723	<input type="checkbox"/>	Domed Rubber, Joystick Handle, Bodypoint PC107A £65
ARM2724	<input type="checkbox"/>	Mushroom Grooved, Joystick Handle, Bodypoint PC110A £65
REMOTE STOP SWITCH		
ARZ2258	<input type="checkbox"/>	Remote Stop Switch (ASL504) £942
SHROUD COLOUR*		
ARM2640	<input type="checkbox"/>	Obsidian Black Frame Shrouds ○
ARM2641	<input type="checkbox"/>	Ruby Red Frame Shrouds ○
ARM2642	<input type="checkbox"/>	Purple Passion Frame Shrouds ○
ARM2643	<input type="checkbox"/>	Peacock Blue Frame Shrouds ○
ARM2644	<input type="checkbox"/>	Arctic Blue Frame Shrouds ○
RIM COLOUR*		
ARM2650	<input type="checkbox"/>	Obsidian Black Rim Inserts ○
ARM2651	<input type="checkbox"/>	Ruby Red Rim Inserts ○
ARM2652	<input type="checkbox"/>	Purple Passion Rim Inserts ○
ARM2653	<input type="checkbox"/>	Peacock Blue Rim Inserts ○
ARM2654	<input type="checkbox"/>	Arctic Blue Rim Inserts ○
FORKS*		
ARM2632	<input type="checkbox"/>	Double sided fork ◆
CENTRE WHEELS (DRIVING WHEELS)*		
ARM0705	<input type="checkbox"/>	DW 14" (355 mm) x 3" (75 mm), Solid inlay, dual density foam insert ○
MOTOR TYPE*		
ARM1294	<input type="checkbox"/>	3 km/h (2mph) ○
ARM1282	<input type="checkbox"/>	6 km/h (4mph) ◆
ARM1908	<input type="checkbox"/>	10 km/h (6mph) (includes LED light system and indicators) £464

BATTERY*			
ARM1257		50 Ah Gel	◆
PERFORMANCE ACCESSORIES			
ARM1215		G-Trac (ensures the powerchair takes the intended direction as indicated by the user control, without the need for continual correction by the client on uneven or sloping terrain)	£738
CHARGER*			
ARM1343		OptiCharge Charger 8A (200 - 250V)	◆
ARM1313		OptiCharge Charger 10A (200 - 250V) Charging time will be 20% faster	£137
LIGHTS*			
ARM1437		Light system with indicators & LED Headlights (included with 10 km/h motor)	£407
ARM1254		Without Light system	◆
TRANSPORTATION KIT*			
ARM1450		Preparation for DAHL Docking Station: Docking station for car implementation, please order docking station separately. (Adequate headrest is compulsory.)	£524
ARM1494		4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory, see options below)	◆
FIXATION BELT			
ARM1329		Lap Belt with Auto Lock	◆
ARM2113		Lap belt with auto lock (flexible mount in side rails) (not available in combination with 410mm seat depth (2132) and T-armrest (2191)	○
BODYPOINT BELTS			
HB216-S38-B2		HB216 2-point belt Small x 380mm Pads, cinch mount	£63
HB216-M46-B2		HB216 2-point belt Medium x 460mm Pads, cinch mount	£75
HB219-S38-B2		HB219 2-point belt Small x 380mm Pads, cinch mount	£75
HB219-M46-B2		HB219 2-point belt Medium x 460mm Pads, cinch mount	£75
HB235-S32-B2		HB235 2-point belt Extra small x, 320mm Pads, cinch mount	£75
HB235-S38-B2		HB235 2-point belt Small x 380mm Pads, cinch mount	£81
EB235S		Evoflex Small pelvic stabilizer, push button	£81
EB225M		Evoflex Medium pelvic stabilizer, Latch	£81
EB225L		Evoflex Large pelvic stabilizer, Latch	£110
BODYPOINT PivotFit			
SH290XXS-B6		X-Extra small standard PivotFit with multi-directional design and cinch mount	£110
SH290XS-B6		Extra small standard PivotFit with multi-directional design and cinch mount	£110
SH290S-B6		Small standard PivotFit with multi-directional design and cinch mount	£110
SH290M-B6		Medium standard PivotFit with multi-directional design and cinch mount	£114
SH295XXS-B6		X-Extra small dynamic PivotFit with multi-directional design and cinch mount	£114
SH295XS-B6		Extra small dynamic PivotFit with multi-directional design and cinch mount	£114
SH295S-B6		Small dynamic PivotFit with multi-directional design and cinch mount	£114
SH295M-B6		Medium dynamic PivotFit with multi-directional design and cinch mount	£92

BODYPOINT Trimline			
SH235S-B6		Small Trimline shoulder harness with multi-directional straps and cinch mount	£92
SH235M-B6		Medium Trimline shoulder harness with multi-directional straps and cinch mount	£133
BODYPOINT Stayflex			
SH351XS-B6		Extra small standard Stayflex without zip, with multidirectional straps and plastic cinch mounts	£133
SH351S-B6		Small standard Stayflex without zip, with multidirectional straps and plastic cinch mounts	£133
SH360M-B6		Medium narrow Stayflex with zip, with multidirectional straps and cinch mounts	£133
SH361M-B6		Medium narrow Stayflex without zip, with multidirectional straps and cinch mount	£133
SH350M-B6		Medium standard Stayflex with zip, with multifunctional straps and plastic cinch mounts	£76
BODYPOINT Monoflex			
SH120XS-B3		Extra Small, front pull, grommet strap	£76
SH120S-B3		Small, front pull, grommet Strap	£76
SH120M-B3		Medium, front pull, grommet strap	£76
SH120XS-B1		Extra Small, Monoflex underarm release, flat mount	£76
SH120S-B1		Small, Monoflex underarm release, flat mount	£76
SH120M-B1		Medium, Monoflex underarm release, flat mount	£76
SH122XS-B1		Extra-small Monflex centre release, flat mount	£76
SH122S-B1		Small Monflex centre release, flat mount	£76
SH122M-B1		Medium Monflex centre release, flat mount	£59
BODYPOINT HARNESES			
1627465		Chest harness interface kit	£70
CANE HOLDER			
ARM1312		Cane holder RH (In combination with Attendant bracket or Dual control the cane holder has to be mounted on the opposite side of the remote)	£38
ARM1363		Cane holder LH (In combination with Attendant bracket or Dual control the cane holder has to be mounted on the opposite side of the remote)	£38
TRAY			
ARM1404		Tray Swing Away RH	(NC 1707, 2213, 2214, 2209, 1287) £174
ARM1405		Tray Swing Away LH	(NC 1708, 2213, 2214, 2210, 1286) £174

TOTAL PRICE	
--------------------	--

Non Conforming Options:

ARM0001	305mm (12")
ARM0002	330mm (13")
ARM2134	460mm (18")
ARM2147	Matrx Libra 90mm thick
ARM2171	Basic Tension Adjustable Backrest (540mm or 480mm upholstery height) Straps easily accessible and adjustable from the rear, Black fabric only
ARM2172	Flex3 adjustable Backrest Plates (without Backrest Cushion)
ARM2173	Flex3 adjustable Backrest Plates with Laguna Backrest cushion (maximum seat depth achievable 510mm (20"))
MINI1312	Matrx Mini series - Mini Back with Kidabra print 330mm x 310mm, 13" x 12" (max user weight 68kg) (please also order the extending Mounting Pins)
ARM2170	Without Backrest Upholstery
ARM2151	Black fabric
ARM2152	Black multi stretch PU coated fabric
ARM2155	Black Spacer fabric
ARM1496	Non Invacare Headrest but profile mount fitting Flex3 hard back
ARM1495	Push bar dia. 22mm in conjunction with telescopic seat frame / backposts, No Headrest - no mounting hardware.
ARM1445	Multi-adjustable pad (Azalea type, multi stretch PU coated fabric) (not available with Matrx Backrests)
ARM1092	Electric Tilt (0° - 20°) Seat Height 450mm
ARM1093	Electric Tilt (0° - 30°) with CoG shift Seat Height 450mm
ARM1094	Seat lifter with Electrical Tilt (300mm, 30°) with CoG shift Seat Height 450mm - 480mm (Also choose option 1265)
ARM0977	Powered Recliner (0° - 30°) Foldable
ARM0556	Centre Mounted: Manual : Hanger, centre mounted, angle settable (0° - 90°) : knee to heel length 390 - 510mm, (including calf pad and angle adjustable footplate)
ARM0589	Centre Mounted : LNX Powered: Centre mounted legrest, length compensating (including calf pad and angle adjustable footplate)
ARM0526	Side Mounted : Standard hanger 80°
ARM0501	Without Hanger
ARM0559	Angle Adjustable Footplate Medium 140mm x 240mm (5.5" x 9.5") for Centre Mounted Hanger
ARM0512	Fixed Footplate
ARM0504	Angle Adjustable Footplate
ARM0513	Angle and depth adjustable Footplate (Adapt)
ARM2201	REM110 Drive only
ARM2203	REM211 Drive & Actuators
ARM2204	REM216 Drive, Actuators & Lights
ARM2206	REM500, without joystick for alternative access, built in Bluetooth HID (Secondary User input required)
ARM1707	Swing Away RH Bracket
ARM1708	Swing Away LH Bracket
ARM2213	Centre mount swing away to RH Nucleus Bracket
ARM2214	Centre mount swing away to LH Nucleus Bracket
ARM2227	Bracket RH mounted for Attendant
ARM2228	Bracket LH mounted for Attendant
ARM2209	REM500 Swing away RH multi-adjustable Bracket
ARM2210	REM500 Swing away LH multi-adjustable Bracket
ARM2270	Assignable Toggle Switch Module for REM400, (as default set to left ON / OFF, right speed dial)
ARM2229	Proportional Chin Joystick (ASL138) (Requires Primary Remote REM400 or REM500) (Please select chin control bracket below)

ARM1289	Chin Control Swing Away Bracket RH Manual
ARM1288	Chin Control Swing Away Bracket LH Manual
ARM1287	Chin Control Swing Away Bracket RH Powered
ARM1286	Chin Control Swing Away Bracket LH Powered
ARM2220	REM050 Dual Control (Please select user control and bracket) (Please select dual control bracket below)
ARH2219	Compact Dual Control (DLX-ACU200) (Please select user control and bracket) (Please select dual control bracket below)
ARM2245	Invacare IDC Dual control (Please select IDC Dual control bracket below)
ARM1222	IDC Dual control Detachable Bracket LH
ARM2213	Proportional Compact Centre Mount swing away to RH (Nucleus)
ARM2214	Proportional Compact Centre Mount swing away to LH (Nucleus)
ARM2230	Head Control with ultrasonic proximity sensors (ATOM , ASL104) with standard wings M3 (Requires Primary Remote REM400 or REM500)
ARM2247	Head Control with ultrasonic proximity sensors (PROTON , ASL104P) with PROTON wings (Requires Primary Remote REM400 or REM500)
ARZ2233	MEC Proportional Joystick (Micro Extremity Control) (Please select MEC Bracketry) (Requires Primary Remote REM400 or REM500)
ARZ2235	MEC Tray Mount RH (ASL603, 100mm x 150mm) (Please select MEC Bracketry)
ARZ2236	MEC Tray Mount LH (ASL603, 100mm x 150mm)
ARZ2213	Centre Mount, swing away to RH, Nucleus
ARZ2214	Centre Mount, swing away to LH, Nucleus
ARZ2241	Centre Mount, swing away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
ARZ2242	Centre Mount, swing away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
ARZ2240	Proportional Heavy Duty Remote Joystick with mode switch (ASL133) (Please select Heavy Duty Bracketry) (Requires Primary Remote REM400 or REM500)
ARZ2235	Tray Mount RH, (ASL603, 100mm x 150mm)
ARZ2236	Tray Mount LH, (ASL603, 100mm x 150mm)
ARZ2213	Centre Mount, swing away to RH, Nucleus
ARZ2214	Centre Mount, swing away to LH, Nucleus
ARZ2241	Centre Mount, swing away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
ARZ2242	Centre Mount, swing away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
ARM2227	10-Way Switch Bracket RH, Mounted for Attendant
ARM2228	10-Way Switch Bracket LH, Mounted for Attendant
ARM1437	Light system with indicators & LED Headlights (included with 10 km/h motor)
ARM1312	Cane holder RH (In combination with Attendant bracket or Dual control the cane holder has to be mounted on the opposite side of the remote)
ARM1363	Cane holder LH (In combination with Attendant bracket or Dual control the cane holder has to be mounted on the opposite side of the remote)
ARM1404	Tray Swing Away RH
ARM1405	Tray Swing Away LH



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):