



Yes, you can:

Tel: 0800 043 0881
Email: salesUK@dhcare.com
Online Spares available at: www.dhcare.co.uk

Power Retail Prescription Form

September 2026
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

TDX SP2 HD Modulite

Standard Features

- Max User Weight up to 180kg (28 stone)
- Crash tested ISO7176 - 19
- Highly adjustable Modulite Seating System
- Telescopic adjustable seat widths from 330mm - 530mm (13" - 21") with an additional 50mm, (2") by arms
- Seat depths 410mm - 560mm stepless adjustment (16" - 22")
- Tension Adjustable backrest
- Standard T-Armrests, height and depth adjustable, detachable, with short armpads
- Manual settable Seat Tilt (0° - 20°)
- Manual settable Backrest Recline (0° - 30° in 3.75° steps)
- Standard 80° hangers with fixed footplate
- LiNX DLX-REM110 with standard bracket
- Selection of 5 shroud and rim colours
- Centre (Driving Wheels) 355mm x 76mm (14" x 3") Black Tyres
- Suspension and Sure Step Stability Technology
- 6 km/h (4mph) 4 Pole Motors
- 8A (200 - 250V) Charger
- 73.5 Ah Batteries
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

Standard Base price: £7,995



Please make your selection with a cross in the box ☐ Free of charge ☒ standard feature POA special request

Customer weight = KG

Code	Description	
BASE MODEL*		
ARJ0000	TDX SP2 Heavy Duty	£7,995
ARZ1300	Special TDX SP2 Heavy Duty (Extended Lead Time). Please provide details in comments section.	POA

USER WEIGHT*			
ARJ1918		User weight less 100 kg	○
ARJ1919		User weight 101 to 125 kg	○
ARJ1921		User weight 126 to 136 kg	○
ARJ1922		User weight 137 to 150 kg (NC 1092, 2052, 2053, 2055, 2059, 2060)	○
ARJ1923		User weight 151 to 160 kg (NC 2192, 2193, 2194, 1092, 2052, 2053, 2055, 2059, 2060)	○
ARJ1924		User weight 161 to 180 kg (NC 2192, 2193, 2194, 1092, 1093, 1094, 1909, 2052, 2053, 2055, 2057, 2059, 2058, 2060)	○
SEATING TYPE*			
ARJ2100		Modulite Seat Unit (NC 2400, 2401, 2402, 2403, 4520, 4521, 2410, 2411, 2412, 2328, 2329, 2330, 2331, 4510, 4530, 4540, 4550)	◆
ARJ2103		Modulite Support+ (NC 2132, 2140, 2142, 2143, 2144, 2147, 1980, 1981, 1982, 1983, 1980, 1981, 1982, 1983, 1985, 2163, 2170, 2171, 2172, 2173, 2174, 2175, 2176, 2177, 2178, 2168, 2169, 2179, 2380, 2381, 2382, 2383, 2385, 2386, 2387, 2389, 2390, 2391, 2393, 2394, 2395, 2151, 2152, 2155, 2195, 2196, 2197, 2199, 2184, 2186, 2188, 2185, 2187, 2189, 1350, 1427, 1432, 1495, 1496, 2324, 2325, 2327)	○
SEAT FRAME*			
ARJ2111		Telescopic Seat Frame + Plate with Side Rails. Width 380 mm to 530 mm. Additional 50 mm width by arm rest adjustment.	◆
SEAT WIDTH*			
ARZ1001		Seat width 330 mm. Seat width achieved on a 380 mm seat frame (2121) with inrigged armrests (2183)	○
ARJ2121		380 mm Seat Plate. Can be reduced to 330 mm in combination with inrigged armrests. (NC 2143, 2701, 2703, 2705)	○
ARJ2122		430 mm Seat Plate. Can be reduced to 380 mm in combination with inrigged armrests. (NC 2703, 2705)	○
ARJ2123		480 mm Seat Plate (NC 2702, 2704, 2706)	○
ARJ2124		530 mm Seat Plate (NC 2174, 2175, 2176, 2177, 2178, 2143, 2702, 2704, 2706)	○
BODY FIT KIT			
UNZ2106		Body Fit Kit - Outrigging Backrest Brackets and Pivot Point (outrig: 1" /25mm per side) (NC 2124, 2192, 2194)	○
SEAT DEPTH RAIL LENGTH*			
ARJ2129		Standard Depth Package (410mm - 510mm. Simple stepless adjustment from the side) (NC 2138)	○
ARJ2130		Extended Depth Package (510mm - 560mm. Simple stepless adjustment from the side) (NC 2132, 2134)	○
SEAT DEPTH*			
ARJ2132		410 mm (NC 2196, 2197, 2701, 2702, 2703, 2704, 2705, 2706)	○
ARJ2134		460 mm	○
ARJ2136		510 mm	○
ARJ2138		560 mm (NC 2143)	○
SEAT OPTION			
ARZSWO		"H" Frame Interface Kit	£495
SEAT CUSHION*			
ARJ2142		Flo-tech Cushion (Lite - 75mm up to 127kg) or (Contour - 100mm up to 191kg) (NC 1981, 1982, 1983, 1980, 1981, 1982, 1983, 1985)	£120
ARJ2143		Flo-tech Contour Visco; (2-way stretch material, 430 mm x 430 mm, up to 125kg) (NC 1981, 1982, 1983, 1980, 1981, 1982, 1983, 1985)	£169
ARJ2147		Matrx Libra (90mm - standard up to 136kg or HD up to 227kg). Extened lead time.	£410
ARJ2144		Matrx Vi Cushion (Mesh fabric) (NC 1981, 1982, 1983, 1980, 1981, 1982, 1983)	£400
ARJ2149		Support+ Cushion	£270
ARJ2140		Without seat cushion (NC 1981, 1982, 1983, 1980, 1981, 1982, 1983, 1985)	○

1st POSITIONING ACCESORIES FOR MATRX LIBRA			
ARJ1980		Tapered Wedges (Pair)	£30
ARJ1981		Abductor	£23
ARJ1982		Lateral Wedge	£20
ARJ1983		Anterior Seat Wedges	£27
2nd POSITIONING ACCESORIES FOR MATRX LIBRA			
ARJ1980		Tapered Wedges (Pair)	£30
ARJ1981		Abductor	£23
ARJ1982		Lateral Wedge	£20
ARJ1983		Anterior Seat Wedges	£27
POSITIONING ACCESSORIES FOR MATRX VI			
ARJ1985		Pelvic Obliquity Kit (Contains 2x positioning accessories)	£54
BACKREST SUPPORT*			
ARJ2163		Backposts short (480 mm). For Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery. Width adjustable cast bracket.	(NC 2173, 2052, 2053) ○
ARJ2162		Backposts tall (540 mm) For Flex3 backrests/Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery. Width adjustable cast bracket.	◆
PUSH HANDLES FOR BACKPOSTS*			
ARJ2107		Backposts WITH push handles	◆
ARJ2108		Backposts WITHOUT push handles, Push handles will be cut down on the vertical section.	(NC 1332, 2220, 2219, 2245, 2227, 2228) ○
BACKREST UPHOLSTERY / CUSHION*			
ARJ2171		Tension Adjustable Upholstery with Nylon overlay. Short = 480 mm / 7 straps. Long = 540 mm / 8 straps. Measured from seat plate to top of upholstery.	(NC 1445, 1496, 2155, 2052, 2053) ◆
ARJ2179		Flex3 Tension Adjustable (540 mm or 480 mm)	(NC 2151, 2192, 2193, 2194, 2199, 2155, 2128) £155
ARJ2172		Flex3 adjustable Backrest Plates (without Backrest Cushion)	(NC 2163, 2128) £250
ARJ2173		Flex3 adjustable Backrest Plates with Laguna Backrest cushion	(NC 2163, 2128) £400
ARJ2400		Support+ Standard Backrest (Spacer fabric)	(NC 2152, 2155) £330
ARJ2401		Support+ Active Strapped Backrest (Spacer fabric)	(NC 2152, 2155) £400
ARJ2402		Support+ Passive Strapped Backrest (Spacer fabric)	(NC 2152, 2155) £400
ARJ2403		Support+ Extra Deep Backrest (Spacer fabric)	(NC 2152, 2155) £450
ARJ2174		Matrx Posture Back, height 400 mm (Mesh fabric) Corresponds to SW setting chosen	(NC 1496, 2124) £353
ARJ2175		Matrx Posture Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen	(NC 1496, 2124) £353
ARJ2176		Matrx Elite E2 Back , height 500 mm (Mesh fabric) Corresponds to SW setting chosen	(NC 1496, 2124) £494
ARJ2177		Matrx Elite TR Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen	(NC 1496, 2124) £496
ARJ2178		Matrx Elite E2 Deep Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen	(NC 1496, 2124) £546
ARJ2168		Matrx MAC Back - Foam Cushion [Mesh fabric]	(NC 1445) £475
ARJ2169		Matrx MAC Back - Hug Stamina Fibre Cushion [Mesh fabric]	(NC 1445) £575
ARJ2170		Without Backrest Upholstery (Back Posts/Back Frame only)	(NC 1496, 2151, 2152, 2155) ○
SUPPORT+ BACKREST HEIGHT (Only applicable for Support+ Back options)			
ARJ4520		Tall - 570 mm (Measured from Seat Plate)	○
ARJ4521		Short - 500 mm (Measured from Seat Plate)	○

MAC BACK BACKREST HEIGHT (Only applicable for MAC Back options)				
ARJ2380		405 mm (16 inch) MAC BH	(NC 1445)	○
ARJ2381		455 mm (18 inch) MAC BH	(NC 1445)	○
ARJ2382		505 mm (20 inch) MAC BH	(NC 1445)	○
ARJ2383		560 mm (22 inch) MAC BH	(NC 1445)	○
MAC BACK LATERAL CONTOUR DEPTH LH (Only applicable for MAC Back options)				
ARJ2385		90 mm (3.5 inch) MAC Contour Depth LH	(NC 1445)	£100
ARJ2386		150 mm (6 inch) MAC Contour Depth LH	(NC 1445)	£100
ARJ2387		230 mm (9 inch) MAC Contour Depth LH	(NC 1445)	£100
MAC BACK LATERAL CONTOUR DEPTH RH (Only applicable for MAC Back options)				
ARJ2389		90 mm (3.5 inch) MAC Contour Depth RH	(NC 1445)	£100
ARJ2390		150 mm (6 inch) MAC Contour Depth RH	(NC 1445)	£100
ARJ2391		230 mm (9 inch) MAC Contour Depth RH	(NC 1445)	£100
MAC BACK BACKREST WIDTH (Only applicable for MAC Back options)				
ARJ2393		Backrest -25 mm (-1 inch) to seat width	(NC 1445)	○
ARJ2394		Backrest matches seat width	(NC 1445)	○
ARJ2395		Backrest +25 mm (+1 inch) to seat width	(NC 1445)	○
UPHOLSTERY COLOUR*				
ARJ2151		Black Nylon Fabric TR18		◆
ARJ2152		Black Multi Stretch PU Coated Wipe Down Fabric TR26	(NC 2171)	○
ARJ2155		Black Spacer Fabric TR36	(NC 2171, 2179)	○
ARX1002		No Upholstery Colour required - No backrest upholstery required	(NC 2179, 2172, 2173, 2174, 2175, 2176, 2177, 2178)	○
ARMRESTS*				
ARJ2191		T-Arm detachable, Standard. Adjustment range from seat plate: 245 - 315 mm		◆
ARJ2181		T-Arm detachable, Extended height. Adjustment range from seat plate: 295 - 365 mm		£40
ARJ2183		Inrigged Armrests. Reduces SW by 25 mm each side	(NC 2185, 2186, 2187, 2188, 2189)	£30
ARJ2193		Flip-up Armrest. Height and depth adjustable (height adjustable from 230mm - 360mm). 150 kg max user weight.	(NC 2176, 2199)	£65
ARJ2192		Following Armrest (height adjustable from 230mm - 300mm). 150 kg max user weight	(NC 1445, 2179)	£185
ARJ2194		Following Armrest (height adjustable from 300mm - 360mm). 150 kg max user weight	(NC 1445, 2179)	£185
ARMPADS*				
ARJ2195		Short Armpad (325mm)		◆
ARJ2196		Long Armpad (375mm)		○
ARJ2197		Long Ergonomic Armpad (397mm)		○
ARJ2410		Support+ Standard Armpad Short (80 mm x 300 mm)		£100
ARJ2411		Support+ Standard Armpad Long (80 mm x 400 mm)		£110
ARJ2412		Support+ P-Shaped Armpad (120 mm x 400 mm)		£130
ARJ2199		Soft Armpad (REA Azalea Wide Pad)		£95
HIP SUPPORTS LH				
ARJ2184		Hip Support LH, pad size 140 x 100 mm	(NC 2702, 2703, 2704, 2705, 2706)	£250
ARJ2186		Hip Support LH, pad size 230 x 100 mm	(NC 2702, 2703, 2704, 2705, 2706)	£250
ARJ2188		Hip Support LH, pad size 330 x 100 mm	(NC 2702, 2703, 2704, 2705, 2706)	£250

HIP SUPPORTS RH			
ARJ2185		Hip Support RH, pad size 140 x 100 mm (NC 2702, 2703, 2704, 2705, 2706)	£250
ARJ2187		Hip Support RH, pad size 230 x 100 mm (NC 2702, 2703, 2704, 2705, 2706)	£250
ARJ2189		Hip Support RH, pad size 330 x 100 mm (NC 2702, 2703, 2704, 2705, 2706)	£250
MOUNTING HARDWARE FOR HIP SUPPORTS			
ARJ2164		Mounting hardware for hip supports (6.1" / 155 mm high)	○
ARJ2165		Mounting hardware for hip supports (8.1" / 206mm high)	○
HEADREST / NECKREST (Adequate headrest option must be selected)*			
ARJ1350		Neckrest depth adjustable	£80
ARJ1432		Headrest sideways adjustable	◆
ARJ1427		Headrest with cheek support	£225
ARJ2328		Support+ Flex Headrest (Medium). Complete with Mounting Hardware & Stem (_2332).	£250
ARJ2329		Support+ Flex Headrest (Large). Complete with Mounting Hardware & Stem (_2332).	£270
ARJ2330		Support+ Standard Headrest (Medium). Complete with Mounting Hardware & Stem (_2332).	£230
ARJ2331		Support+ Standard Headrest (Large). Complete with Mounting Hardware & Stem (_2332).	£250
ARJ2324		Elan Standard Pad 255 mm / 10" wide w/standard cover (ESP26). Complete with Elan Headrest Hardware	£325
ARJ2325		Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (EOP23). Complete with Elan Headrest Hardware	£325
ARJ2327		Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard cover. Complete with Elan Headrest Hardware	£440
ARJ1496		Non Invacare Headrest but profile mount fitting Flex3 hard back	○
ARJ1495		Non Invacare Headrest - no mounting hardware. With Backbrace/Pushbar dia. 22mm in conjunction with telescopic seat frame/backposts	○
SUPPORT+ ACCESSORIES			
ARJ4510		Support+ Wipedown Cushion Cover (PU Coated Type)	£30
ARJ4530		Support+ Wipedown Backrest Cover (PU Coated Type)	£30
ARJ4540		Support+ Wipedown Armpad Covers (PU Coated Type)	£30
ARJ4550		Support+ Wipedown Headrest Cover (PU Coated Type)	£30
LATERAL SUPPORT			
ARJ1445		Lateral Support with extended arm, multi-adjustable pad (NC 2171)	£165

MANUAL / POWERED SEATING OPTIONS*			
ARJ1091		Manually settable Seat Tilt 0° - 20°	◆
ARJ1092		Powered "fixed pivot" Tilt 20° (136kg max user weight)	£400
ARJ2055		Powered "fixed pivot" Tilt 20° & Sit Smarter (136kg max user weight) (NC 2201, 2054)	£540
ARJ1095		Powered "fixed pivot" Tilt 20°, heavy duty (180kg max user weight)	£400
ARJ2056		Powered "fixed pivot" Tilt 20° & Sit Smarter, heavy duty (180kg max user weight) (NC 2201, 2054)	£540
ARJ1093		Powered Tilt (0° - 30°) with CoG shift (160kg max user weight) (NC 2201)	£1,100
ARJ2057		Powered Tilt (0° - 30°) with CoG shift & Sit Smarter (160kg max user weight) (NC 2201)	£1,490
ARJ2052		Powered Tilt (0°-50°) with CoG shift (136kg max user weight) (NC 2201)	£1,450
ARJ2059		Powered Tilt (0° - 50°) with CoG shift & Sit Smarter (136kg max user weight) (NC 2201, 2171)	£2,000
ARJ1094		Seat Lifter with Powered Tilt (300mm, 30°) with CoG shift (160kg max user weight) (NC 2201)	£1,900
ARJ2058		Seat Lifter with Powered Tilt (300 mm, 30°) with CoG shift & Sit Smarter (160kg max user weight) (NC 2201, 2220)	£2,570
ARJ2053		Seat Lifter with Powered Tilt (300mm, 45°) with CoG shift (136kg max user weight) (NC 2201)	£2,500
ARJ2060		Seat Lifter with Powered Tilt (300mm, 45°) with CoG & Sit Smarter (136kg max user weight) (NC 2201, 2220, 2171)	£3,370
SEAT TO GROUND HEIGHT*			
ARJ0354		450mm (17.75") (from floor to front of seat plate at 0° tilt) (NC 1094, 2052, 2053)	○
ARJ0364		475mm (18.75") (from floor to front of seat plate at 0° tilt)	○
SEAT TILT OPTION			
ARJ2054		5° Fixed Anterior Tilt, Assist mount. Will reduce max tilt angle by 5° (NC 1092, 1095)	£400
BACKREST RECLINER*			
ARJ0975		Manual Settable Recliner (0° - 30°, in 3.75° steps) Foldable	◆
ARJ0977		Powered Recliner (0° - 30°) Foldable	£500
ARJ0979		Powered Recliner 30° & Sit Smarter	£670

LEG REST TYPE*				
ARJ0556		Centre Mounted: Manual Centre Mounted Legrest, angle settable. Length 280 - 385 mm	(NC 0517, 0512, 0504, 0513, 0501)	£370
ARJ0589		Centre Mounted: LNX Powered Centre Mounted Legrest, length compensating. Length 340 - 410 mm	(NC 2201, 0504, 0512, 0504, 0513, 0501)	£1,620
ARJ0600		Centre Mounted: LNX Powered Centre Mounted Legrest & Sit Smarter, length compensating. Length 340 - 410 mm	(NC 2201, 0512, 0504, 0513, 0501)	£2,190
UNZ0592		Centre Mounted: LNX Powered Centre Mounted Legrest with Telescoping Foot Platform (lowers to the ground). Length: 380-419mm (Selectable length depends on STF height) Custom-made, requires customised setting.	(NC 2201, 0504, 0512, 0504, 0513, 0501)	POA
ARJ0526		Side Mounted: Standard hanger 80° pair. Lower leg length: 350 - 450 mm	(NC 0501)	◆
ARJ0503		Side Mounted: Vari F pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with tools	(NC 0555, 0501)	£120
ARJ0533		Side Mounted: Vari A pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with handwheel and calf pad	(NC 0554, 0555, 0501)	£235
ARJ0531		Side Mounted: ADE pair. Lower leg length: 290 - 460 mm. Powered angle adjustment: 80 - 0° with calf pad	(NC 2201, 0554, 0555, 0501)	£660
ARJ0530		Side Mounted: ADE pair & Sit Smarter. Lower leg length: 290 - 460 mm. Powered angle adjustment: 80 - 0° with calf pad	(NC 2201, 2365, 0559, 0516, 0517, 0554, 0555, 0501)	£900
ARJ0501		Without Leg Rests	(NC 0504, 0513, 2365, 0559, 0516, 0517, 0559, 0516, 0517)	○
FOOTPLATE TYPE*				
ARJ2365		Angle & width adjustable footplates S (127 W x 190 mm D / 5"x7.5")	(NC 0503, 0533, 0531)	£105
ARJ0559		Angle & width adjustable footplates M (140 W x 240 mm D / 5.5"x9.5")	(NC 0503, 0533, 0531)	○
ARJ0516		Angle & width adjustable footplates L (150 W x 290 mm D / 6"x11.5")	(NC 0503, 0533, 0531)	£135
ARJ0517		Angle adjustable foot platform M (305 W x 275 mm D / 12"x10.75")	(NC 0526, 0503, 0533, 0531)	£270
ARJ0512		Fixed Footplate pair	(NC 0533, 0531)	○
ARJ0504		Angle Adjustable Footplate pair	(NC 0526)	£75
ARJ0513		Angle and depth adjustable Footplate pair		£120
ARJ0501		Without footplates for option 0501 without Leg Rests	(NC 0503, 0533, 0531)	○
AMPUTEE SUPPORT				
ARJ0554		Amputee Support RH	(NC 0556, 0589, 0592)	£195
ARJ0555		Amputee Support LH	(NC 0556, 0589, 0592)	£195
LINX CONTROLS*				
ARJ2201		DLX-REM110 (Drive Only)	(NC 1095, 1093, 1094, 0977, 0589, 0531, 2209, 2210, 2213, 2214, 2270, 2219, 2245, 1909)	◆
ARJ2203		DLX-REM211 (Drive & Actuators)	(NC 2209, 2210, 2213, 2214, 2231, 2232, 2270, 1909)	○
ARJ2204		DLX-REM216 (Drive, Actuators & Lights)	(NC 2210, 2213, 2214, 2231, 2232, 2270)	£200
ARJ2205		DLX-REM400 (Touch Screen Interface. Drive, Act. & Lights. Configurable jack sockets)	(NC 2210, 2213, 2215)	£425
ARJ2206		DLX-REM500 (Touch Screen Interface Display Unit - requires a user input/alternative joystick).	(NC 1208, 1707, 1708, 1709, 1710, 2220, 2270)	£765
REM400 TOGGLE OPTION				
ARJ2270		Assignable Toggle Switch Module for REM400 (DLX-SW400TGL)	(NC 2203, 2204)	£315

CONTROL BRACKETRY*			
ARJ1207		Standard Bracket RH	○
ARJ1208		Standard Bracket LH	○
ARJ1707		Swing Away RH Bracket	£135
ARJ1708		Swing Away LH Bracket	£135
ARJ1709		Matrx Resolve Swing in/out Bracket RH	£245
ARJ1710		Matrx Resolve Swing in/out Bracket LH	£245
ARJ2213		Centre Mount Swing Away to RH Nucleus Bracket (NC 1405)	£495
ARJ2214		Centre Mount Swing Away to LH Nucleus Bracket (NC 1405)	£495
ARJ2209		REM500 Swing Away RH Multi-Adjustable Bracket (NC 2203, 2204, 2205)	○
ARJ2210		REM500 Swing Away LH Multi-Adjustable Bracket (NC 2203, 2204, 2205)	○
ARJ2227		Bracket RH mounted for Attendant	£65
ARJ2228		Bracket LH mounted for Attendant	£65
ARZ2239		Individual mounting option per request (Special). Requires additional information/specification.	POA
DUAL CONTROLS			
ARJ2220		LiNX DLX-REM050 (Dual Control used as secondary remote) (NC 1221, 1222, 2227, 2228, 2205, 2206, 2245, 2128)	£395
ARJ2219		LiNX DLX-ACU200 (Compact Dual Control) (NC 1222, 2227, 2228, 2128, 2201)	£395
ARJ2245		LiNX DLX-IDC (Intuitive Dual Control) (NC 2228, 2128, 2201)	£745
ARJ2420		mo-vis Scoot Control w/ dedicated bracket included. (NC 2108, 1221, 1222)	£3,420
DUAL CONTROL BRACKETRY (Only applicable if Dual Control is selected)			
ARJ1207		Dual Control Standard Bracket RH	○
ARJ1208		Dual Control Standard Bracket LH	○
ARJ1221		IDC Dual Control Detachable Bracket RH (NC 2220, 2227)	○
ARJ1222		IDC Dual Control Detachable Bracket LH (NC 2220, 2228)	○
ARJ9260		Bracket for mo-vis Scoot (NC 2219, 2245)	○
COMPACT REMOTES			
ARJ2231		LiNX DLX-CR400 (Compact Remote) (NC 2203, 2204, 2220)	£825
ARJ2232		LiNX DLX-CR400LF (Compact Remote Low Force) (NC 2203, 2204, 2220)	£1,075
ARJ3000		mo-vis Micro LiNX (NC 2201, 2203, 2204)	£2,960
ARJ3001		mo-vis Multi LiNX (NC 4401, 4402, 2201, 2203, 2204)	£2,640
ARJ3110		mo-vis All Round Low Force Omni (NC 3221, 3222, 4400, 4401, 4402, 2201, 2203, 2204)	£2,200
ARJ3111		mo-vis All Round Omni (NC 3410, 3411, 3412, 3413, 3420, 3421, 3422, 3423, 3200, 3221, 3222, 4400, 4401, 4402, 2201, 2203, 2204)	£1,600
ARJ3112		mo-vis Heavy Duty Omni (NC 3410, 3411, 3412, 3413, 3420, 3421, 3422, 3423, 3200, 3221, 3222, 4300, 4301, 4302, 4400, 4401, 4402, 2201, 2203, 2204)	£3,000
COMPACT REMOTE BRACKETRY (Only applicable if Compact Remote is selected)			
ARJ2207		RH Mount (NC 2201, 2203, 2204)	○
ARJ2208		LH Mount (NC 2201, 2203, 2204)	○
ARJ2213		Centre Mount Swing Away to RH (Nucleus) (NC 1405, 1415, 2201, 2203, 2204)	£495
ARJ2214		Centre Mount Swing Away to LH (Nucleus) (NC 1405, 1339, 2210, 2203, 2204)	£495
ARZ2215		Individual mounting option per request (Special). Requires additional information/specification.	POA
ARJ3300		Maxx Resolve RH (NC 4301, 4302)	£270
ARJ3301		Maxx Resolve LH (NC 4301, 4302)	£270
ARJ3310		BodyPoint Trilock Midline RH (NC 4301, 4302)	£480
ARJ3311		BodyPoint Trilock Midline LH (NC 4301, 4302)	£480
ARJ3320		Tray Mounted (NC 4301, 4302)	○
ARZ3330		Individual mounting option per request (Special). Requires additional information/specification.	POA

CHIN CONTROLS				
ARJ3000		mo-vis Micro LiNX	(NC 2201, 2203, 2204)	£2,960
ARJ3001		mo-vis Multi LiNX	(NC 4401, 4402, 2201, 2203, 2204)	£2,640
ARJ3110		mo-vis All Round Low Force Omni (consider option 3810 for switch access)	(NC 3221, 3222, 4400, 4401, 4402, 2201, 2203, 2204)	£2,200
ARJ3111		mo-vis All Round Omni (consider option 3810 for switch access)	(NC 3410, 3411, 3412, 3413, 3420, 3421, 3422, 3423, 3200, 3221, 3222, 4400, 4401, 4402, 2201, 2203, 2204)	£1,600
CHIN CONTROL BRACKETRY (Only applicable if Chin Control is selected)				
ARJ3400		mo-vis Chin Control Harness (One Size)	(NC 4301, 4302)	£1,000
ARJ3410		mo-vis Manual Lateral Swing Away RH Mounted	(NC 4301, 4302)	£1,520
ARJ3411		mo-vis Manual Lateral Swing Away LH Mounted	(NC 4301, 4302)	£1,520
ARJ3412		mo-vis Manual Vertical Swing Away RH Mounted	(NC 4301, 4302)	£1,520
ARJ3413		mo-vis Manual Vertical Swing Away LH Mounted	(NC 4301, 4302)	£1,520
ARJ3420		mo-vis Powered Lateral Swing Away RH Mounted	(NC 4301, 4302)	£2,400
ARJ3421		mo-vis Powered Lateral Swing Away LH Mounted	(NC 4301, 4302)	£2,400
ARJ3422		mo-vis Powered Vertical Swing Away RH Mounted	(NC 4301, 4302)	£2,400
ARJ3423		mo-vis Powered Vertical Swing Away LH Mounted	(NC 4301, 4302)	£2,400
ARZ3330		Individual mounting option per request (Special). Requires additional information/specification.		POA
HEAD CONTROLS				
ARJ3200		mo-vis Proportional Head Control Omni	(NC 2201, 2203, 2204)	£4,500
ARJ3221		mo-vis Switched Head Control Omni	(NC 2201, 2203, 2204)	£4,500
ARJ3222		mo-vis Switched Head Control w/ Sip & Puff Omni	(NC 2201, 2203, 2204)	£4,750
FOOT CONTROL				
ARJ3010		mo-vis Foot Control LiNX RH Side	(NC 2203, 2204, 2220)	£3,300
ARJ3011		mo-vis Foot Control LiNX LH Side	(NC 2203, 2204, 2220)	£3,300
LiNX 10-WAY SWITCH				
ARJ2251		10-Way Switch DLX-SW10-A for direct Powered Seating Access		£500
LiNX 10-WAY SWITCH BRACKETRY (Only applicable if 10-Way Switch is selected)				
ARJ2217		Bracket RH, mounted for user	(NC 2213)	○
ARJ2218		Bracket LH, mounted for user	(NC 2214)	○
ARJ2227		Bracket RH, mounted for attendant		£65
ARJ2228		Bracket LH, mounted for attendant		£65
LiNX ACCESSORY CABLE				
ARJ2252		Loom providing jack sockets for Drive and Mode (includes connection for 10-way switch)		£100
LiNX USB CHARGER PORT				
ARJ2253		USB Charger Port DLX-USB02-A	(NC 2701, 2703, 2704, 2705, 2706)	£160
LiNX USB CHARGER PORT BRACKETRY (Only applicable if USB Charger is selected)				
ARJ2207		RH Mount		○
ARJ2208		LH Mount		○
FURTHER LiNX / MO-VIS COMPONENTS / ACCESSORIES				
ARJ2302		Not Necessary	(NC 2254, 1711, 1712, 1713, 1714, 1715, 1716)	◆
ARJ2300		To be enclosed only. Standard 300 mm Cable Length for mo-vis Switches		○
ARJ2305		To be enclosed only. Extended 1500 mm Cable Length for mo-vis Switches		○
ARZ2301		To be connected and programmed (Special)		POA
BUDDY BUTTONS				
ARJ1685		Buddy Button. Red / Yellow / Green / Blue / Black provided. Interchangeable colours.		£90

EGG SWITCH			
ARJ2254		Egg Switch, black	£90
PIKO BUTTONS			
ARJ1711		Red Piko Button	£160
ARJ1712		Yellow Piko Button	£160
ARJ1713		Green Piko Button	£160
ARJ1714		Orange Piko Button	£160
ARJ1715		Blue Piko Button	£160
ARJ1716		Black Piko Button	£160
TOGGLES / BUTTONS			
ARZ2261		Stereo Toggle - direct control for Tilt only	£270
ARZ2262		Stereo Button - direct control for Tilt only	£255
ARZ2263		Four-Way Toggle - direct control for Tilt & Recline	£330
ARZ2264		Four-Way Button - direct control for Tilt & Recline	£355
ARZ2207		Toggle / Button RH Mount	○
ARZ2208		Toggle / Button LH Mount	○
ARZ2225		Without Mounting	○
LINX INPUT MODULE			
UNZ2273		DLX-IN200 Module (Jack port for external power button) (NC 2220)	£270
UNZ9249		DLX-IN500 Module (for DB9 Drive input. Only compatible with DLX-REM400 / DLX-REM500)	£785
MO-VIS TWISTER BASIC			
ARJ3500		Black	£140
ARJ3501		Red	£140
ARJ3502		Yellow	£140
ARJ3503		Blue	£140
ARJ3504		Green	£140
MO-VIS TWISTER BASIC MOUNTING (Only applicable for mo-vis Twister Basic)			
ARJ3520		Mounting Set	£50
ARJ3525		Table Mounting Set	£50
MO-VIS TWISTER			
ARJ3600		Black	£160
ARJ3601		Red	£160
ARJ3602		Yellow	£160
ARJ3603		Blue	£160
ARJ3604		Green	£160
MO-VIS TWISTER ON BENDED TUBE INC. D6 CLAMP			
ARJ3700		Black	£220
ARJ3701		Red	£220
ARJ3702		Yellow	£220
ARJ3703		Blue	£220
ARJ3704		Green	£220

MO-VIS SATELLITE TWISTER (Only applicable for mo-vis All Round Remotes)			
ARJ3800		Black	£200
ARJ3801		Red	£200
ARJ3802		Yellow	£200
ARJ3803		Blue	£200
ARJ3804		Green	£200
ARJ3810		All Round Chin Control Set (includes All Round Joystick Ball and 2x Satellite Switches)	£420
MO-VIS TWISTER GOOSENECK			
ARJ3900		Black	£260
ARJ3901		Red	£260
ARJ3902		Green	£260
MO-VIS Q2M TWIST & SWING			
ARJ4000		Black	£500
ARJ4001		Red	£500
ARJ4002		Yellow	£500
ARJ4003		Blue	£500
ARJ4004		Green	£500
MO-VIS TWISTER PRO			
ARJ4100		Black	£160
ARJ4101		Red	£160
ARJ4102		Green	£160
MO-VIS MULTI SWITCH			
ARJ4200		Multi Switch	£620
ARJ4201		Multi Switch Proximity Sensors	£860
ARJ4202		Multi Switch Proximity Set	£840
MO-VIS HANDWARMER			
ARJ4300		Hand Warmer	£800
ARJ4301		Hand Warmer 2 Dome Hand Hood	£400
ARJ4302		Hand Warmer 1 Dome Hand Hood	£300
MO-VIS JOYSTICK KNOBS			
ARJ4400		Micro Joystick Ball	£60
ARJ4401		Micro Joystick Cup 10.1mm	£60
ARJ4402		Micro Joystick Cup 15.5mm	£60
ARJ4403		Multi Joystick Ball 4mm	£60
ARJ4404		Multi Joystick Cup 24mm	£60
ARJ4405		All Round Joystick Ball 6.35mm	£60
ARJ4406		All Round Joystick Cup 37mm	£60
BODYPOINT JOYSTICK KNOBS			
ARJ2721		3" U-Shaped, Joystick Handle, Bodypoint PC101A	£75
ARJ2722		4" U-Shaped, Joystick Handle, Bodypoint PC102A	£75
ARJ2723		Domed Rubber, Joystick Handle, Bodypoint PC107A	£75
ARJ2724		Mushroom Grooved, Joystick Handle, Bodypoint PC110A	£75

SHROUD COLOUR*			
ARJ2640		Obsidian Black Frame Shrouds	○
ARJ2641		Ruby Red Frame Shrouds	○
ARJ2642		Purple Passion Frame Shrouds	○
ARJ2643		Peacock Blue Frame Shrouds	○
ARJ2644		Arctic Blue Frame Shrouds	○
RIM COLOUR*			
ARJ2650		Obsidian Black Frame Shrouds	○
ARJ2651		Ruby Red Frame Shrouds	○
ARJ2652		Purple Passion Frame Shrouds	○
ARJ2653		Peacock Blue Frame Shrouds	○
ARJ2654		Arctic Blue Frame Shrouds	○
FORKS*			
ARJ2631		Single sided fork	£300
ARJ2632		Double sided fork	○
CENTRE WHEELS (DRIVE WHEELS)*			
ARJ0704		DW 14" (355 mm) x 3" (75 mm), Pneumatic	○
ARJ0705		DW 14" (355 mm) x 3" (75 mm), Solid inlay, dual density foam insert	○
ARJ0706		DW 14" (355 mm) x 3" (75 mm), Puncture protected, gel-sealed	○
MOTOR TYPE*			
ARJ1294		3 km/h (2mph)	○
ARJ1282		6 km/h (4mph) (4 pole)	◆
ARJ1909		12 km/h (8mph) (includes light system and indicators) G-Trac recommended	£600
BATTERY*			
ARJ1260		73.5 Ah (Gel)	◆
G TRAC			
ARJ1215		G-Trac (ensures the powerchair takes the intended direction as indicated by the user control, without the need for continual correction by the client on uneven or sloping terrain)	£400
CHARGER*			
ARJ1343		OptiCharge Charger 8A (200 - 250V)	◆
ARJ1313		OptiCharge Charger 10A (200 - 250V) Charging time will be 20% faster	£155
LIGHTS*			
ARJ1437		LED Light System with Indicators (30 LUX) (NC 1909)	£450
ARJ1254		Without light system (NC 1909)	◆
ARX1006		Lights Included with Motor Package 10 & 12kmh (NC 1282)	○
TRANSPORTATION KIT*			
ARJ1450		Preparation for DAHL-Docking station (Adequate headrest is compulsory). Includes a mounting bracket in the battery box. Docking Station for car implementation has to be ordered separately with DAHL Engineering.	£580
ARJ1493		4-Point tie-down (For wheelchair transportation without user)	○
ARJ1494		4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Adequate headrest is compulsory.	◆

FIXATION BELT*				
ARJ2113		Lapbelt (adjustable both sides) adjustable mounting on seat rail extrusion	(NC 2702, 2703, 2704, 2705, 2706)	○
ARJ1332		H-Belt		£75
ARJ2701		Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size L (padded length 620mm)	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 2113)	£95
ARJ2702		Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size M (padded length 460mm)	(NC 2123, 2124, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 2113)	£95
ARJ2703		Bodypoint 2-Point Hip Belt, Centre pull - push button, Cinch mount size L (padded length 620mm)	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 2113)	£90
ARJ2704		Bodypoint 2-Point Hip Belt, Centre pull - push button, Cinch mount size M (padded length 460mm)	(NC 2123, 2124, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 2113)	£90
ARJ2705		Bodypoint 2-Point Hip Belt, Dual pull Side release button, Cinch mount size L (padded length 620mm)	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 2113)	£90
ARJ2706		Bodypoint 2-Point Hip Belt, Dual pull Side release button, Cinch mount size M (padded length 460mm)	(NC 2123, 2124, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 2113)	£90
REAR VIEW MIRROR				
ARJ1415		Rear view Mirror RH	(NC 2213)	£45
ARJ1339		Rear view Mirror LH	(NC 2214)	£45
BAG HOOKS				
ARJ2128		Bag Hooks (Max load up to 6kg = impact on max user weight)		£80
CANE HOLDER				
ARJ1312		Cane holder RH		£45
ARJ1363		Cane holder LH		£45
TRAY				
ARJ1404		Tray, Swing Away RH	(NC 1707, 2213, 2214, 2209)	£200
ARJ1405		Tray, Swing Away LH	(NC 1708, 2213, 2214, 2210)	£200

TOTAL PRICE				
-------------	--	--	--	--

Non Conforming Options:

ARJ2121	380 mm Seat Plate. Can be reduced to 330 mm in combination with inrigged armrests.
ARJ2122	430 mm Seat Plate. Can be reduced to 380 mm in combination with inrigged armrests.
ARJ2123	480 mm Seat Plate
ARJ2124	530 mm Seat Plate
ARJ2132	410 mm
ARJ2134	460 mm
ARJ2138	560 mm
ARJ2142	Flo-tech Cushion (Lite - 75mm up to 127kg) or (Contour - 100mm up to 191kg)
ARJ2143	Flo-tech Contour Visco; (2-way stretch material, 430 mm x 430 mm, up to 125kg)
ARJ2147	Matrx Libra (90mm - standard up to 136kg or HD up to 227kg). Extended lead time.
ARJ2144	Matrx Vi Cushion (Mesh fabric)
ARJ2140	Without seat cushion
ARJ1980	Tapered Wedges (Pair)
ARJ1981	Abductor
ARJ1982	Lateral Wedge
ARJ1983	Anterior Seat Wedges
ARJ1980	Tapered Wedges (Pair)
ARJ1981	Abductor
ARJ1982	Lateral Wedge
ARJ1983	Anterior Seat Wedges
ARJ1985	Pelvic Obliquity Kit (Contains 2x positioning accessories)
ARJ2163	Backposts short (480 mm). For Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery. Width adjustable cast bracket.
ARJ2108	Backposts WITHOUT push handles, Push handles will be cut down on the vertical section.
ARJ2171	Tension Adjustable Upholstery with Nylon overlay. Short = 480 mm / 7 straps. Long = 540 mm / 8 straps. Measured from seat plate to top of upholstery.
ARJ2179	Flex3 Tension Adjustable (540 mm or 480 mm)
ARJ2172	Flex3 adjustable Backrest Plates (without Backrest Cushion)
ARJ2173	Flex3 adjustable Backrest Plates with Laguna Backrest cushion
ARJ2400	Support+ Standard Backrest (Spacer fabric)
ARJ2401	Support+ Active Strapped Backrest (Spacer fabric)
ARJ2402	Support+ Passive Strapped Backrest (Spacer fabric)
ARJ2403	Support+ Extra Deep Backrest (Spacer fabric)
ARJ2174	Matrx Posture Back, height 400 mm (Mesh fabric) Corresponds to SW setting chosen
ARJ2175	Matrx Posture Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen
ARJ2176	Matrx Elite E2 Back , height 500 mm (Mesh fabric) Corresponds to SW setting chosen
ARJ2177	Matrx Elite TR Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen
ARJ2178	Matrx Elite E2 Deep Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen
ARJ2168	Matrx MAC Back - Foam Cushion [Mesh fabric]
ARJ2169	Matrx MAC Back - Hug Stamina Fibre Cushion [Mesh fabric]

ARJ2170	Without Backrest Upholstery (Back Posts/Back Frame only)
ARJ4520	Tall - 570 mm (Measured from Seat Plate)
ARJ4521	Short - 500 mm (Measured from Seat Plate)
ARJ2380	405 mm (16 inch) MAC BH
ARJ2381	455 mm (18 inch) MAC BH
ARJ2382	505 mm (20 inch) MAC BH
ARJ2383	560 mm (22 inch) MAC BH
ARJ2385	90 mm (3.5 inch) MAC Contour Depth LH
ARJ2386	150 mm (6 inch) MAC Contour Depth LH
ARJ2387	230 mm (9 inch) MAC Contour Depth LH
ARJ2389	90 mm (3.5 inch) MAC Contour Depth RH
ARJ2390	150 mm (6 inch) MAC Contour Depth RH
ARJ2391	230 mm (9 inch) MAC Contour Depth RH
ARJ2393	Backrest -25 mm (-1 inch) to seat width
ARJ2394	Backrest matches seat width
ARJ2395	Backrest +25 mm (+1 inch) to seat width
ARJ2151	Black Nylon Fabric TR18
ARJ2152	Black Multi Stretch PU Coated Wipe Down Fabric TR26
ARJ2155	Black Spacer Fabric TR36
ARJ2193	Flip-up Armrest. Height and depth adjustable (height adjustable from 230mm - 360mm). 150 kg max user weight.
ARJ2192	Following Armrest (height adjustable from 230mm - 300mm). 150 kg max user weight
ARJ2194	Following Armrest (height adjustable from 300mm - 360mm). 150 kg max user weight
ARJ2195	Short Armpad (325mm)
ARJ2196	Long Armpad (375mm)
ARJ2197	Long Ergonomic Armpad (397mm)
ARJ2410	Support+ Standard Armpad Short (80 mm x 300 mm)
ARJ2411	Support+ Standard Armpad Long (80 mm x 400 mm)
ARJ2412	Support+ P-Shaped Armpad (120 mm x 400 mm)
ARJ2199	Soft Armpad (REA Azalea Wide Pad)
ARJ2184	Hip Support LH, pad size 140 x 100 mm
ARJ2186	Hip Support LH, pad size 230 x 100 mm
ARJ2188	Hip Support LH, pad size 330 x 100 mm
ARJ2185	Hip Support RH, pad size 140 x 100 mm
ARJ2187	Hip Support RH, pad size 230 x 100 mm
ARJ2189	Hip Support RH, pad size 330 x 100 mm
ARJ1350	Neckrest depth adjustable
ARJ1432	Headrest sideways adjustable
ARJ1427	Headrest with cheek support
ARJ2328	Support+ Flex Headrest (Medium). Complete with Mounting Hardware & Stem (_2332).
ARJ2329	Support+ Flex Headrest (Large). Complete with Mounting Hardware & Stem (_2332).
ARJ2330	Support+ Standard Headrest (Medium). Complete with Mounting Hardware & Stem (_2332).
ARJ2331	Support+ Standard Headrest (Large). Complete with Mounting Hardware & Stem (_2332).
ARJ2324	Elan Standard Pad 255 mm / 10" wide w/standard cover (ESP26). Complete with Elan Headrest Hardware
ARJ2325	Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (EOP23). Complete with Elan Headrest Hardware

ARJ2327	Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard cover. Complete with Elan Headrest Hardware
ARJ1496	Non Invacare Headrest but profile mount fitting Flex3 hard back
ARJ1495	Non Invacare Headrest - no mounting hardware. With Backbrace/Pushbar dia. 22mm in conjunction with telescopic seat frame/backposts
ARJ4510	Support+ Wipedown Cushion Cover (PU Coated Type)
ARJ4530	Support+ Wipedown Backrest Cover (PU Coated Type)
ARJ4540	Support+ Wipedown Armpad Covers (PU Coated Type)
ARJ4550	Support+ Wipedown Headrest Cover (PU Coated Type)
ARJ1445	Lateral Support with extended arm, multi-adjustable pad
ARJ1092	Powered "fixed pivot" Tilt 20° (136kg max user weight)
ARJ2055	Powered "fixed pivot" Tilt 20° & Sit Smarter (136kg max user weight)
ARJ1095	Powered "fixed pivot" Tilt 20°, heavy duty (180kg max user weight)
ARJ1093	Powered Tilt (0° - 30°) with CoG shift (160kg max user weight)
ARJ2057	Powered Tilt (0° - 30°) with CoG shift & Sit Smarter (160kg max user weight)
ARJ2052	Powered Tilt (0°-50°) with CoG shift (136kg max user weight)
ARJ2059	Powered Tilt (0° - 50°) with CoG shift & Sit Smarter (136kg max user weight)
ARJ1094	Seat Lifter with Powered Tilt (300mm, 30°) with CoG shift (160kg max user weight)
ARJ2058	Seat Lifter with Powered Tilt (300 mm, 30°) with CoG shift & Sit Smarter (160kg max user weight)
ARJ2053	Seat Lifter with Powered Tilt (300mm, 45°) with CoG shift (136kg max user weight)
ARJ2060	Seat Lifter with Powered Tilt (300mm, 45°) with CoG & Sit Smarter (136kg max user weight)
ARJ2054	5° Fixed Anterior Tilt, Assist mount. Will reduce max tilt angle by 5°
ARJ0977	Powered Recliner (0° - 30°) Foldable
ARJ0556	Centre Mounted: Manual Centre Mounted Legrest, angle settable. Length 280 - 385 mm
ARJ0589	Centre Mounted: LNX Powered Centre Mounted Legrest, length compensating. Length 340 - 410 mm
UNZ0592	Centre Mounted: LNX Powered Centre Mounted Legrest with Telescoping Foot Platform (lowers to the ground). Length: 380-419mm (Selectable length depends on STF height) Custom-made, requires customised setting.
ARJ0526	Side Mounted: Standard hanger 80° pair. Lower leg length: 350 - 450 mm
ARJ0503	Side Mounted: Vari F pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with tools
ARJ0533	Side Mounted: Vari A pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with handwheel and calf pad
ARJ0531	Side Mounted: ADE pair. Lower leg length: 290 - 460 mm. Powered angle adjustment: 80 - 0° with calf pad
ARJ0501	Without Leg Rests
ARJ2365	Angle & width adjustable footplates S (127 W x 190 mm D / 5"x7.5")
ARJ0559	Angle & width adjustable footplates M (140 W x 240 mm D / 5.5"x9.5")
ARJ0516	Angle & width adjustable footplates L (150 W x 290 mm D / 6"x11.5)
ARJ0517	Angle adjustable foot platform M (305 W x 275 mm D / 12"x10.75")
ARJ0512	Fixed Footplate pair
ARJ0504	Angle Adjustable Footplate pair
ARJ0513	Angle and depth adjustable Footplate pair
ARJ0501	Without footplates for option 0501 without Leg Rests
ARJ0554	Amputee Support RH
ARJ0555	Amputee Support LH
ARJ2201	DLX-REM110 (Drive Only)

ARJ2203	DLX-REM211 (Drive & Actuators)
ARJ2204	DLX-REM216 (Drive, Actuators & Lights)
ARJ2205	DLX-REM400 (Touch Screen Interface. Drive, Act. & Lights. Configurable jack sockets)
ARJ2206	DLX-REM500 (Touch Screen Interface Display Unit - requires a user input/alternative joystick).
ARJ2270	Assignable Toggle Switch Module for REM400 (DLX-SW400TGL)
ARJ1208	Standard Bracket LH
ARJ1707	Swing Away RH Bracket
ARJ1708	Swing Away LH Bracket
ARJ1709	Matrx Resolve Swing in/out Bracket RH
ARJ1710	Matrx Resolve Swing in/out Bracket LH
ARJ2213	Centre Mount Swing Away to RH Nucleus Bracket
ARJ2214	Centre Mount Swing Away to LH Nucleus Bracket
ARJ2209	REM500 Swing Away RH Multi-Adjustable Bracket
ARJ2210	REM500 Swing Away LH Multi-Adjustable Bracket
ARJ2227	Bracket RH mounted for Attendant
ARJ2228	Bracket LH mounted for Attendant
ARJ2220	LiNX DLX-REM050 (Dual Control used as secondary remote)
ARJ2219	LiNX DLX-ACU200 (Compact Dual Control)
ARJ2245	LiNX DLX-IDC (Intuitive Dual Control)
ARJ1208	Dual Control Standard Bracket LH
ARJ1221	IDC Dual Control Detachable Bracket RH
ARJ1222	IDC Dual Control Detachable Bracket LH
ARJ2231	LiNX DLX-CR400 (Compact Remote)
ARJ2232	LiNX DLX-CR400LF (Compact Remote Low Force)
ARJ2213	Centre Mount Swing Away to RH (Nucleus)
ARJ2214	Centre Mount Swing Away to LH (Nucleus)
ARZ2215	Individual mounting option per request (Special). Requires additional information/specification.
ARJ3410	mo-vis Manual Lateral Swing Away RH Mounted
ARJ3411	mo-vis Manual Lateral Swing Away LH Mounted
ARJ3412	mo-vis Manual Vertical Swing Away RH Mounted
ARJ3413	mo-vis Manual Vertical Swing Away LH Mounted
ARJ3420	mo-vis Powered Lateral Swing Away RH Mounted
ARJ3421	mo-vis Powered Lateral Swing Away LH Mounted
ARJ3422	mo-vis Powered Vertical Swing Away RH Mounted
ARJ3423	mo-vis Powered Vertical Swing Away LH Mounted
ARJ3200	mo-vis Proportional Head Control Omni
ARJ3221	mo-vis Switched Head Control Omni
ARJ3222	mo-vis Switched Head Control w/ Sip & Puff Omni
ARJ2227	Bracket RH, mounted for attendant
ARJ2228	Bracket LH, mounted for attendant
ARJ2253	USB Charger Port DLX-USB02-A
ARJ2254	Egg Switch, black
ARJ1711	Red Piko Button
ARJ1712	Yellow Piko Button

ARJ1713	Green Piko Button
ARJ1714	Orange Piko Button
ARJ1715	Blue Piko Button
ARJ1716	Black Piko Button
ARJ4300	Hand Warmer
ARJ4301	Hand Warmer 2 Dome Hand Hood
ARJ4302	Hand Warmer 1 Dome Hand Hood
ARJ4400	Micro Joystick Ball
ARJ4401	Micro Joystick Cup 10.1mm
ARJ4402	Micro Joystick Cup 15.5mm
ARJ1282	6 km/h (4mph) (4 pole)
ARJ1909	12 km/h (8mph) (includes light system and indicators) G-Trac recommended
ARJ2113	Lapbelt (adjustable both sides) adjustable mounting on seat rail extrusion
ARJ1332	H-Belt
ARJ2701	Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size L (padded length 620mm)
ARJ2702	Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size M (padded length 460mm)
ARJ2703	Bodypoint 2-Point Hip Belt, Centre pull - push button, Cinch mount size L (padded length 620mm)
ARJ2704	Bodypoint 2-Point Hip Belt, Centre pull - push button, Cinch mount size M (padded length 460mm)
ARJ2705	Bodypoint 2-Point Hip Belt, Dual pull Side release button, Cinch mount size L (padded length 620mm)
ARJ2706	Bodypoint 2-Point Hip Belt, Dual pull Side release button, Cinch mount size M (padded length 460mm)
ARJ1415	Rear view Mirror RH
ARJ1339	Rear view Mirror LH
ARJ2128	Bag Hooks (Max load up to 6kg = impact on max user weight)
ARJ1405	Tray, Swing Away LH



For a more complex customer requirement,
please email: salesUK@dhcare.com

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):