



**Yes, you can:**

# Power Invacare Ireland Prescription Form

Tel: 01656 776222  
Email: ordersuk@invacare.com  
Online Spares available at: www.invacare.co.uk

07 JANUARY 2026  
Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## TDX SP2 Low Rider

### Standard Features

- Max User Weight up to 160kg (25.2 stone)
- Crash tested ISO7176 - 19
- Highly adjustable Modulite Seating System
- Telescopic adjustable seat widths from 330mm - 530mm (13" - 21") with an additional 50mm, (2") by arms
- Seat depths 410mm - 560mm stepless adjustment (16" - 22")
- Tension Adjustable backrest
- Standard T-Armrests, height and depth adjustable, detachable, with short ampads
- Seat lifter with Powered Tilt (300mm, 30°) CoG shift
- Manual settable Backrest Recline (0° - 30° in 3.75° steps)
- Standard 80° hangers with fixed footplate
- LiNX DLX-REM211 with standard bracket
- Selection of 5 shroud and rim colours
- Centre (Driving Wheels) 355mm x 76mm (14" x 3") Black Tyres
- Suspension and Sure Step Stability Technology
- 6 km/h (4mph) 4 Pole Motors
- 8A (200 - 250V) Charger
- 60 or 73.5 Ah Batteries (STGH dependant)
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

**Standard Base price: €10,409**



**Q'STRAIT**  
**QLK-150**  
Compatible

Please make your selection with a cross in the box    ○ Free of charge    ◆ standard feature    POA special request

Customer weight =                      KG

Code	Description	
<b>BASE MODEL*</b>		
ARQ0071	TDX SP2 Low Rider. 73Ah Batteries. 420mm Seat to Ground Height	○
ARQ0072	TDX SP2 Low Rider. 60Ah Batteries. 403mm Seat to Ground Height	○
ARZ0075	Special TDX SP2 Low Rider. 73Ah Batteries. 420mm Seat to Ground Height (Extended Lead Time). Please provide details in comments section.	POA
ARZ0076	Special TDX SP2 Low Rider. 60Ah Batteries. 403mm Seat to Ground Height (Extended Lead Time). Please provide details in comments section.	POA

<b>USER WEIGHT*</b>			
ARQ1918		User weight less 100 kg	○
ARQ1919		User weight 101 to 125 kg	○
ARQ1921		User weight 126 to 136 kg (NC 1350, 1432, 1427)	○
ARQ1922		User weight 137 to 150 kg (NC 1350, 1432, 1427, 2053)	○
ARQ1923		User weight 151 to 160 kg (NC 2192, 2193, 2194, 1908, 1350, 1432, 1427, 2053)	○
<b>SEATING TYPE*</b>			
ARQ2100		Invacare® Modulite Seat Unit	€10,409
<b>SEAT FRAME*</b>			
ARQ2111		Telescopic Seat Frame + Plate with Side Rails. Width 380 mm to 530 mm. Additional 50 mm width by arm rest adjustment.	◆
<b>SEAT WIDTH*</b>			
ARQ2122		430 mm Seat Plate (NC 2703, 2705)	○
ARQ2123		480 mm Seat Plate (NC 2702, 2704, 2706)	○
ARQ2124		530 mm Seat Plate (NC 2174, 2175, 2176, 2177, 2178, 2143, 2702, 2704, 2706)	○
<b>SEAT DEPTH RAIL LENGTH*</b>			
ARQ2129		Standard Depth Package (410mm - 510mm. Simple stepless adjustment from the side) (NC 2138)	○
ARQ2130		Extended Depth Package (510mm - 560mm. Simple stepless adjustment from the side) (NC 2132, 2134)	○
<b>SEAT DEPTH*</b>			
ARQ2132		410 mm (NC 2196, 2197, 2701, 2702, 2703, 2704, 2705, 2706)	○
ARQ2134		460 mm	○
ARQ2136		510 mm	○
ARQ2138		560 mm (NC 2143)	○
<b>SEAT OPTION</b>			
ARZSWO		"H" Frame Interface Kit (NC 2143)	€507
<b>SEAT CUSHION*</b>			
ARQ2142		Matrx Flo-tech Lite Cushion (NC 1981, 1982, 1983, 1980, 1981, 1982, 1983)	€132
ARQ2143		Matrx Flo-tech Contour Visco; (2-way stretch material, 430 mm x 430 mm only) (NC 1981, 1982, 1983, 1980, 1981, 1982, 1983)	€132
ARQ2147		Matrx Libra Cushion, 90 mm thick (Extended lead time)	€492
ARQ2140		Without seat cushion (NC 1981, 1982, 1983, 1980, 1981, 1982, 1983)	○
<b>POSITIONING ACCESORIES FOR MATRX LIBRA</b>			
ARQ1980		Tapered Wedges (Pair)	€28
ARQ1981		Abductor	€58
ARQ1982		Lateral Wedge	€46
ARQ1983		Anterior Seat Wedges	€43
<b>BACKREST SUPPORT*</b>			
ARQ2162		Backposts tall (540 mm) For Flex3 backrests/Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery. Width adjustable cast bracket.	○
ARQ2163		Backposts short (480 mm). For Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery. Width adjustable cast bracket. (NC 2173, 2053)	○
<b>PUSH HANDLES FOR BACKPOSTS*</b>			
ARQ2107		Backposts WITH push handles	○
ARQ2108		Backposts WITHOUT push handles, Push handles will be cut down on the vertical section. (NC 1332, 2220, 2219, 2245, 2227, 2228)	€50

<b>BACKREST UPHOLSTERY / CUSHION*</b>			
ARQ2171		Tension Adjustable Upholstery with Nylon overlay. Short = 480 mm / 7 straps. Long = 540 mm / 8 straps. Measured from seat plate to top of upholstery.	(NC 1445, 1496, 2155, 2230, 2247) ◆
ARQ2179		Flex3 Tension Adjustable (540 mm or 480 mm)	(NC 2151, 2192, 2193, 2194, 2199, 2155, 2230, 2247, 2128) €138
ARQ2172		Flex3 Adjustable Backrest Plates (without Backrest Cushion)	(NC 2163, 2128) €265
ARQ2173		Flex3 Adjustable Backrest Plates with Laguna Backrest cushion	(NC 2163, 2128) €460
ARQ2174		Matrx Posture Back, height 400 mm (Mesh fabric) Corresponds to SW setting chosen	(NC 1496, 2124) €1,038
ARQ2175		Matrx Posture Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen	(NC 1496, 2124) €1,038
ARQ2176		Matrx Elite E2 Back , height 500 mm (Mesh fabric) Corresponds to SW setting chosen	(NC 1496, 2124) €668
ARQ2177		Matrx Elite TR Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen	(NC 1496, 2124) €572
ARQ2178		Matrx Elite E2 Deep Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen	(NC 1496, 2124) €739
ARQ2168		Matrx MAC Back - Foam Cushion [Mesh fabric]	(NC 1445) €650
ARQ2169		Matrx MAC Back - Hug Stamina Fibre Cushion [Mesh fabric]	(NC 1445) €750
ARQ2170		Without Backrest Upholstery (Back Posts/Back Frame only)	(NC 1496, 2230, 2247, 2151, 2152, 2155) ○
<b>MAC BACK BACKREST HEIGHT (Only applicable for MAC Back options)</b>			
ARQ2380		405 mm (16 inch) MAC BH	(NC 1445) ○
ARQ2381		455 mm (18 inch) MAC BH	(NC 1445) ○
ARQ2382		505 mm (20 inch) MAC BH	(NC 1445) ○
ARQ2383		560 mm (22 inch) MAC BH	(NC 1445) ○
<b>MAC BACK LATERAL CONTOUR DEPTH LH (Only applicable for MAC Back options)</b>			
ARQ2385		90 mm (3.5 inch) MAC Contour Depth LH	(NC 1445) €130
ARQ2386		150 mm (6 inch) MAC Contour Depth LH	(NC 1445) €130
ARQ2387		230 mm (9 inch) MAC Contour Depth LH	(NC 1445) €130
<b>MAC BACK LATERAL CONTOUR DEPTH RH (Only applicable for MAC Back options)</b>			
ARQ2389		90 mm (3.5 inch) MAC Contour Depth RH	(NC 1445) €130
ARQ2390		150 mm (6 inch) MAC Contour Depth RH	(NC 1445) €130
ARQ2391		230 mm (9 inch) MAC Contour Depth RH	(NC 1445) €130
<b>MAC BACK BACKREST WIDTH (Only applicable for MAC Back options)</b>			
ARQ2393		Backrest -25 mm (-1 inch) to seat width	(NC 1445) ○
ARQ2394		Backrest matches seat width	(NC 1445) ○
ARQ2395		Backrest +25 mm (+1 inch) to seat width	(NC 1445) ○
<b>UPHOLSTERY COLOUR*</b>			
ARQ2151		Black Nylon Fabric TR18	○
ARQ2152		Black Multi Stretch PU Coated Wipe Down Fabric TR26	(NC 2171) ○
ARQ2155		Black Spacer Fabric TR36	(NC 2171, 2179) ○
ARX1002		No Upholstery Colour required - No backrest upholstery required	(NC 2179, 2172, 2173, 2174, 2175, 2176, 2177, 2178) ○

<b>ARMRESTS*</b>			
ARQ2191		T-Arm detachable, Standard. Not width adjustable for TDX Low Rider models.	◆
ARQ2181		T-Arm detachable, Extended height. Not width adjustable for TDX Low Rider models.	€50
ARQ2183		Inrigged Armrests. Reduces SW by 25 mm each side (NC 2185, 2186, 2187, 2188, 2189)	€29
ARQ2193		Flip-up Armrest. Height and depth adjustable (height adjustable from 230mm - 360mm). 150 kg max user weight. (NC 2176, 2199)	€57
ARQ2192		Following Armrest (height adjustable from 230mm - 300mm). 150 kg max user weight (NC 1445, 2179)	€210
ARQ2194		Following Armrest (height adjustable from 300mm - 360mm). 150 kg max user weight (NC 1445, 2179)	€210
<b>ARMPADS*</b>			
ARQ2195		Short Armpad (325mm)	◆
ARQ2196		Long Armpad (375mm)	○
ARQ2197		Long Ergonomic Armpad (397mm) (NC 2247)	○
ARQ2199		Soft Arpad (REA Azalea Wide Pad)	€108
<b>HIP SUPPORTS LH</b>			
ARQ2184		Hip Support LH, pad size 140 x 100 mm (NC 2702, 2703, 2704, 2705, 2706, 0072, 0076)	€286
ARQ2186		Hip Support LH, pad size 230 x 100 mm (NC 2702, 2703, 2704, 2705, 2706, 0072, 0076)	€286
ARQ2188		Hip Support LH, pad size 330 x 100 mm (NC 2702, 2703, 2704, 2705, 2706, 0072, 0076)	€286
<b>HIP SUPPORTS RH</b>			
ARQ2185		Hip Support RH, pad size 140 x 100 mm (NC 2702, 2703, 2704, 2705, 2706, 0072, 0076)	€286
ARQ2187		Hip Support RH, pad size 230 x 100 mm (NC 2702, 2703, 2704, 2705, 2706, 0072, 0076)	€286
ARQ2189		Hip Support RH, pad size 330 x 100 mm (NC 2702, 2703, 2704, 2705, 2706, 0072, 0076)	€286
<b>MOUNTING HARDWARE FOR HIP SUPPORTS</b>			
ARQ2164		Mounting hardware for hip supports (6,1" / 155 mm high) (NC 0076)	○
ARQ2165		Mounting hardware for hip supports (8.1" / 206mm high) (NC 0076)	○
<b>HEADREST / NECKREST (Adequate headrest option must be selected)*</b>			
ARQ1350		Neckrest depth adjustable (NC 2247)	€86
ARQ1432		Headrest sideways adjustable (NC 2247)	○
ARQ1427		Headrest with cheek support (NC 2247)	€239
ARQ2324		Elan Standard Pad 255 mm / 10" wide w/standard cover (ESP26). Complete with Elan Headrest Hardware (NC 2247)	€406
ARQ2325		Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (EOP23). Complete with Elan Headrest Hardware (NC 2247)	€406
ARQ2327		Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard cover. Complete with Elan Headrest Hardware (NC 2247)	€548
ARQ1496		Non Invacare Headrest but profile mount fitting Flex3 hard back	○
ARQ1495		Non Invacare Headrest - no mounting hardware. With Backbrace/Pushbar dia. 22mm in conjunction with telescopic seat frame/backposts	○
<b>LATERAL SUPPORT</b>			
ARQ1445		Lateral Support with extended arm, multi-adjustable pad (NC 2171)	€188
<b>MANUAL / POWERED SEATING OPTIONS*</b>			
ARQ1094		Seat Lifter with Powered Tilt (300mm, 30°) with CoG shift (Seat height 460mm - 490mm, 160kg max user weight)	◆
ARQ2053		Seat Lifter with Powered Tilt (300mm, 45°) with CoG shift (Seat height 460mm, 136kg max user weight)	€1,677

<b>SEAT TILT OPTION</b>			
ARQ2054		5° Fixed Anterior Tilt, Assist mount. Will reduce max tilt angle by 5°	€500
<b>BACKREST RECLINER*</b>			
ARQ0975		Manual Settable Recliner (0° - 30°, in 3.75° steps) Foldable	◆
ARQ0977		Powered Recliner (0° - 30°) Foldable	€355
<b>LEG REST TYPE*</b>			
ARQ0556		Centre Mounted: Manual Centre Mounted Legrest, angle settable. Length 280 - 385 mm	(NC 0512, 0504, 0513) €337
ARQ0589		Centre Mounted: LNX Powered Centre Mounted Legrest, length compensating. Length 340 - 410 mm	(NC 0504, 0512, 0504, 0513) €2,013
ARQ0526		Side Mounted: Standard hanger 80° pair. Lower leg length: 350 - 450 mm	◆
ARQ0503		Side Mounted: Vari F pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with tools	○
ARQ0533		Side Mounted: Vari A pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with handwheel and calf pad	€215
ARQ0531		Side Mounted: ADE pair. Lower leg length: 290 - 460 mm. Powered angle adjustment: 80 - 0° with calf pad	€740
ARQ0501		Without Leg Rests	(NC 0504, 0513, 0559, 0559) ○
<b>FOOTPLATE TYPE*</b>			
ARQ0559		Angle & width adjustable footplates M (140 W x 240 mm D / 5.5"x9.5")	(NC 0503, 0533, 0531) ○
ARQ0512		Fixed Footplate pair	(NC 0533, 0531) ◆
ARQ0504		Angle Adjustable Footplate pair	(NC 0526) €79
ARQ0513		Angle and depth adjustable Footplate pair	€128
ARQ0501		Without footplates for option 0501 without Leg Rests	(NC 0503, 0533, 0531) ○
<b>LINX CONTROLS*</b>			
ARQ2203		DLX-REM211 (Drive & Actuators)	(NC 2209, 2210, 2213, 2214, 2231, 2232, 2230, 2247, 2233, 2240, 2270, 2258, 1908, 1909) ◆
ARQ2204		DLX-REM216 (Drive, Actuators & Lights)	(NC 2210, 2213, 2214, 2231, 2232, 2230, 2247, 2233, 2240, 2270, 2258) €162
ARQ2205		DLX-REM400 (Touch Screen Interface. Drive, Act. & Lights. Configurable jack sockets)	(NC 2210, 2213, 2215) €495
ARQ2206		DLX-REM500 (Touch Screen Interface Display Unit - requires a user input/alternative joystick).	(NC 1208, 1707, 1708, 1709, 1710, 2227, 2228, 2270) €909
<b>REM400 TOGGLE OPTION</b>			
ARQ2270		Assignable Toggle Switch Module for REM400 (DLX-SW400TGL)	(NC 2203, 2204) €329

<b>CONTROL BRACKETRY*</b>			
ARQ1207		Standard Bracket RH	○
ARQ1208		Standard Bracket LH	○
ARQ1707		Swing Away RH Bracket	€159
ARQ1708		Swing Away LH Bracket	€159
ARQ1709		Matrx Resolve Swing in/out Bracket RH	€281
ARQ1710		Matrx Resolve Swing in/out Bracket LH	€281
ARQ2213		Centre Mount Swing Away to RH Nucleus Bracket	€739
ARQ2214		Centre Mount Swing Away to LH Nucleus Bracket	€739
ARQ2209		REM500 Swing Away RH Multi-Adjustable Bracket (NC 2203, 2204, 2205)	○
ARQ2210		REM500 Swing Away LH Multi-Adjustable Bracket (NC 2203, 2204, 2205)	○
ARQ2227		Bracket RH mounted for Attendant (NC 2220, 2128)	€71
ARQ2228		Bracket LH mounted for Attendant (NC 2220, 2128)	€71
ARZ2239		Individual mounting option per request (Special). Requires additional information/specification.	POA
<b>DUAL CONTROLS</b>			
ARQ2220		LiNX DLX-REM050 (Dual Control used as secondary remote) (NC 1221, 1222, 2227, 2228, 2205, 2206, 2128, 2230, 2247, 2233, 2240)	€489
ARQ2219		LiNX DLX-ACU200 (Compact Dual Control) (NC 1222, 2227, 2228, 2128)	€489
ARQ2245		LiNX DLX-IDC (Intuitive Dual Control) (NC 2219, 2227, 2228, 2128)	€748
UNZ9259		mo-vis Scoot Control (Only available via Special) (NC 2108, 0071, 0072, 1221, 1222)	€3,020
<b>DUAL CONTROL BRACKETRY (Only applicable if Dual Control is selected)</b>			
ARQ1207		Dual Control Standard Bracket RH	○
ARQ1208		Dual Control Standard Bracket LH	○
ARQ1221		IDC Dual Control Detachable Bracket RH (NC 2219, 2227, 2228)	○
ARQ1222		IDC Dual Control Detachable Bracket LH (NC 2219, 2227, 2228)	○
UNZ9260		Individual mounting option per request (Special). Requires additional information/specification. (NC 0072)	POA
<b>COMPACT REMOTES</b>			
ARQ2231		LiNX DLX-CR400 (Compact Remote) (NC 2203, 2204)	€845
ARQ2232		LiNX DLX-CR400LF (Compact Remote Low Force) (NC 2203, 2204)	€1,199
UNZ9250		mo-vis Micro LiNX (Only available via Special) (NC 2203, 2204, 0071, 0072, 2207, 2208, 2213, 2214, 2235, 2236, 2241, 2242, 2220)	€3,430
UNZ9251		mo-vis Micro Omni (Only available via Special. DLX-IN500 module needed) (NC 2203, 2204, 0071, 0072, 2207, 2208, 2213, 2214, 2235, 2236, 2241, 2242, 2220)	€3,560
UNZ9252		mo-vis Multi LiNX (Only available via Special) (NC 2203, 2204, 0071, 0072, 2207, 2208, 2213, 2214, 2235, 2236, 2241, 2242, 2220)	€3,060
UNZ9253		mo-vis Multi Omni (Only available via Special. DLX-IN500 module needed) (NC 2203, 2204, 0071, 0072, 2207, 2208, 2213, 2214, 2235, 2236, 2241, 2242, 2220)	€3,285
UNZ9254		mo-vis All Round Omni (Only available via Special. DLX-IN500 module needed) (NC 2203, 2204, 0071, 0072, 2207, 2208, 2213, 2214, 2235, 2236, 2241, 2242, 2220)	€1,820
UNZ9255		mo-vis All Round Low Force Omni (Only available via Special. DLX-IN500 module needed) (NC 2203, 2204, 0071, 0072, 2207, 2208, 2213, 2214, 2235, 2236, 2241, 2242, 2220)	€2,630
UNZ9256		mo-vis Heavy Duty Omni (Only available via Special. DLX-IN500 module needed) (NC 2203, 2204, 0071, 0072, 2207, 2208, 2213, 2214, 2235, 2236, 2241, 2242, 2220)	€4,160
ARZ2233		ASL MEC (Proportional MEC Joystick ASL130) (NC 2203, 2204)	€4,230
ARZ2240		ASL Heavy Duty (Proportional Joystick w mode switch ASL133) (NC 2203, 2204)	€2,683

<b>COMPACT REMOTE BRACKETRY (Only applicable if Compact Remote is selected)</b>				
ARQ2207		RH Mount	(NC 2203, 2204)	○
ARQ2208		LH Mount	(NC 2203, 2204)	○
ARQ2213		Centre Mount Swing Away to RH (Nucleus)	(NC 2203, 2204)	€739
ARQ2214		Centre Mount Swing Away to LH (Nucleus)	(NC 2203, 2204)	€739
ARZ2235		Tray Mount RH (ASL603, 100mm x 150mm)	(NC 2203, 2204)	○
ARZ2236		Tray Mount LH (ASL603, 100mm x 150mm)	(NC 2203, 2204)	○
ARZ2241		Centre Mount, Swing Away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 2203, 2204)	€739
ARZ2242		Centre Mount, Swing Away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 2203, 2204)	€739
ARZ2215		Individual mounting option per request (Special). Requires additional information/specification.	(NC 2203, 2204)	POA
<b>CHIN CONTROLS</b>				
UNZ9250		mo-vis Micro LiNX (Only available via Special)	(NC 0072, 2203, 2204, 2220)	€3,430
UNZ9251		mo-vis Micro Omni (Only available via Special. DLX-IN500 module needed)	(NC 0072, 2203, 2204, 2220)	€3,560
UNZ9252		mo-vis Multi LiNX (Only available via Special)	(NC 0072, 2203, 2204, 2220)	€3,060
UNZ9253		mo-vis Multi Omni (Only available via Special. DLX-IN500 module needed)	(NC 0072, 2203, 2204, 2220)	€3,285
UNZ9254		mo-vis All Round Omni (Only available via Special. DLX-IN500 module needed)	(NC 0072, 2203, 2204, 2220)	€1,820
UNZ9255		mo-vis All Round Low Force Omni (Only available via Special. DLX-IN500 module needed)	(NC 0072, 2203, 2204, 2220)	€2,630
<b>CHIN CONTROL BRACKETRY (Only applicable if Chin Control is selected)</b>				
UNZ9261		mo-vis Chin Control Harness (One Size. Special Request Only)	(NC 0072)	€1,255
UNZ9262		mo-vis Manual Swing Away (Special Mounting Only. Please specify what side mounting required)	(NC 0072)	€1,820
UNZ9263		mo-vis Powered Swing Away (Special Mounting Only. Please specify what side mounting required)	(NC 0072)	€3,005
<b>HEAD CONTROLS</b>				
ARQ2230		ASL104 ATOM Head Control with ultrasonic proximity sensors with standard wings	(NC 2179, 2197, 2203, 2204, 2220, 2233, 2324, 2325, 2327)	€4,497
ARQ2247		ASL104P PROTON Head Control with ultrasonic proximity sensors with PROTON wings	(NC 2179, 2197, 2203, 2204, 2220, 2233, 2324, 2325, 2327)	€5,248
UNZ9257		mo-vis Proportional Head Control Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 2233, 2324, 2325, 2327, 2240, 2175, 2171, 2170, 2197, 1350, 1432, 1427, 0071, 0072, 2220)	€9,135
UNZ9258		mo-vis Switched Head Control Omni (Only available via Special. DLX-IN500 module needed)	(NC 2203, 2204, 2233, 2324, 2325, 2327, 2240, 2175, 2171, 2170, 2197, 1350, 1432, 1427, 0071, 0072, 2220)	€5,435
<b>FOOT CONTROL</b>				
UNZ9264		mo-vis Foot Control LiNX (Only available via Special)	(NC 0072, 2203, 2204, 2220)	€3,800
UNZ9265		mo-vis Foot Control Omni (Only available via Special. DLX-IN500 module needed)	(NC 0072, 2203, 2204, 2220)	€3,615
<b>FOOT CONTROL BRACKETRY (Only applicable if Foot Control is selected)</b>				
UNZ9266		Individual mounting option per request (Special). Requires additional information/specification.	(NC 0072)	POA
<b>LINX 10-WAY SWITCH</b>				
ARQ2251		10-Way Switch DLX-SW10-A for direct Powered Seating Access		€539
<b>LINX 10-WAY SWITCH BRACKETRY (Only applicable if 10-Way Switch is selected)</b>				
ARQ2217		Bracket RH, mounted for user	(NC 2213)	○
ARQ2218		Bracket LH, mounted for user	(NC 2214)	○
ARQ2227		Bracket RH, mounted for attendant		○
ARQ2228		Bracket LH, mounted for attendant		○
<b>LINX ACCESSORY CABLE</b>				
ARQ2252		Loom providing jack sockets for Drive and Mode (includes connection for 10-way switch)		€110

<b>LINX USB CHARGER PORT</b>			
ARQ2253	<input type="checkbox"/>	USB Charger Port DLX-USB02-A (NC 2701, 2703, 2704, 2705, 2706)	€162
<b>LINX USB CHARGER PORT BRACKETRY (Only applicable if USB Charger is selected)</b>			
ARQ2207	<input type="checkbox"/>	RH Mount	○
ARQ2208	<input type="checkbox"/>	LH Mount	○
<b>FURTHER LINX COMPONENTS / ACCESSORIES</b>			
ARQ2302	<input type="checkbox"/>	Not Necessary (NC 1686, 1687, 1690, 2254, 1711, 1712, 1713, 1714, 1715, 1716)	◆
ARQ2300	<input type="checkbox"/>	To be enclosed only	○
ARZ2301	<input type="checkbox"/>	To be connected and programmed (Special)	POA
<b>BUDDY BUTTONS</b>			
ARQ1685	<input type="checkbox"/>	Red Buddy Button	€90
ARQ1686	<input type="checkbox"/>	Yellow Buddy Button	€90
ARQ1687	<input type="checkbox"/>	Green Buddy Button	€90
ARQ1690	<input type="checkbox"/>	Blue Buddy Button	€90
<b>EGG SWITCH</b>			
ARQ2254	<input type="checkbox"/>	Egg Switch, black	€90
<b>PIKO BUTTONS</b>			
ARQ1711	<input type="checkbox"/>	Red Piko Button	€161
ARQ1712	<input type="checkbox"/>	Yellow Piko Button	€161
ARQ1713	<input type="checkbox"/>	Green Piko Button	€161
ARQ1714	<input type="checkbox"/>	Orange Piko Button	€161
ARQ1715	<input type="checkbox"/>	Blue Piko Button	€161
ARQ1716	<input type="checkbox"/>	Black Piko Button	€161
<b>TOGGLES / BUTTONS</b>			
ARZ2261	<input type="checkbox"/>	Stereo Toggle - direct control for Tilt only	€273
ARZ2262	<input type="checkbox"/>	Stereo Button - direct control for Tilt only	€263
ARZ2263	<input type="checkbox"/>	Four-Way Toggle - direct control for Tilt & Recline	€342
ARZ2264	<input type="checkbox"/>	Four-Way Button - direct control for Tilt & Recline	€365
ARZ2207	<input type="checkbox"/>	Toggle / Button RH Mount	○
ARZ2208	<input type="checkbox"/>	Toggle / Button LH Mount	○
ARZ2225	<input type="checkbox"/>	Without Mounting	○
<b>BODYPOINT JOYSTICK KNOBS</b>			
ARQ2721	<input type="checkbox"/>	3" U-Shaped, Joystick Handle, Bodypoint PC101A	€75
ARQ2722	<input type="checkbox"/>	4" U-Shaped, Joystick Handle, Bodypoint PC102A	€75
ARQ2723	<input type="checkbox"/>	Domed Rubber, Joystick Handle, Bodypoint PC107A	€75
ARQ2724	<input type="checkbox"/>	Mushroom Grooved, Joystick Handle, Bodypoint PC110A	€75
<b>REMOTE STOP SWITCH</b>			
ARZ2258	<input type="checkbox"/>	Remote Stop Switch (ASL504)	€1,070
<b>LINX INPUT MODULE</b>			
UNZ2273	<input type="checkbox"/>	DLX-IN200 Module (Jack port for external power button) (NC 0072, 2220)	€303
UNZ9249	<input type="checkbox"/>	DLX-IN500 Module (for DB9 Drive input. Only compatible with DLX-REM400 / DLX-REM500) (NC 0072)	€814

<b>SHROUD COLOUR*</b>			
ARQ2640		Obsidian Black Frame Shrouds	○
ARQ2641		Ruby Red Frame Shrouds	○
ARQ2642		Purple Passion Frame Shrouds	○
ARQ2643		Peacock Blue Frame Shrouds	○
ARQ2644		Arctic Blue Frame Shrouds	○
<b>RIM COLOUR*</b>			
ARQ2650		Obsidian Black Rim Inserts	○
ARQ2651		Ruby Red Rim Inserts	○
ARQ2652		Purple Passion Rim Inserts	○
ARQ2653		Peacock Blue Rim Inserts	○
ARQ2654		Arctic Blue Rim Inserts	○
<b>FORKS*</b>			
ARQ2631		Single sided fork (New) includes Black castor	€358
ARQ2632		Double sided fork	○
<b>CENTRE WHEELS (DRIVE WHEELS)*</b>			
ARQ0704		DW 14" (355 mm) x 3" (75 mm), Pneumatic	○
ARQ0705		DW 14" (355 mm) x 3" (75 mm), Solid inlay, dual density foam insert	◆
ARQ0706		DW 14" (355 mm) x 3" (75 mm), Puncture protected, gel-sealed	○
<b>MOTOR TYPE*</b>			
ARQ1294		3 km/h (2mph)	○
ARQ1282		6 km/h (4mph) (4 pole)	◆
ARQ1908		10 km/h (6mph) (includes light system and indicators) G-Trac recommended	○
ARQ1909		12 km/h (8mph) (includes light system and indicators) G-Trac recommended	○
<b>G TRAC</b>			
ARQ1215		G-Trac (ensures the powerchair takes the intended direction as indicated by the user control, without the need for continual correction by the client on uneven or sloping terrain)	€908
<b>CHARGER*</b>			
ARQ1343		OptiCharge Charger 8A (200 - 250V)	○
ARQ1313		OptiCharge Charger 10A (200 - 250V) Charging time will be 20% faster	€174
<b>LIGHTS*</b>			
ARQ1437		LED Light System with Indicators (30 LUX)	(NC 1908, 1909) €519
ARQ1254		Without Light system	(NC 1909) ◆
ARX1006		Lights Included with Motor Package 10 & 12kmh	(NC 1282) ○
<b>TRANSPORTATION KIT</b>			
ARQ1450		TEMPORARILY UNAVAILABLE ON TDX MODELS UNTIL FURTHER NOTICE - Preparation for DAHL-Docking station.	(NC 0072, 0075, 0076) €624
ARQ1494		4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19). Adequate headrest is compulsory.	○

<b>FIXATION BELT*</b>			
ARQ1329		Lapbelt (adjustable both sides) fitted to base of backrest	(NC 2702, 2703, 2704, 2705, 2706) ○
ARQ2113		Lapbelt (adjustable both sides) adjustable mounting on seat rail extrusion	○
ARQ1332		H-Belt	€80
ARQ2701		Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size L (padded length 620mm)	(NC 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113) €117
ARQ2702		Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size M (padded length 460mm)	(NC 2123, 2124, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113) €117
ARQ2703		Bodypoint 2-Point Hip Belt, Centre pull - push button, Cinch mount size L (padded length 620mm)	(NC 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113) €85
ARQ2704		Bodypoint 2-Point Hip Belt, Centre pull - push button, Cinch mount size M (padded length 460mm)	(NC 2123, 2124, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113) €85
ARQ2705		Bodypoint 2-Point Hip Belt, Dual pull Side release button, Cinch mount size L (padded length 620mm)	(NC 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113) €85
ARQ2706		Bodypoint 2-Point Hip Belt, Dual pull Side release button, Cinch mount size M (padded length 460mm)	(NC 2123, 2124, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113) €85
<b>REAR-VIEW MIRROR</b>			
ARQ1415		Rear View Mirror RH	(NC 2213) ○
ARQ1339		Rear View Mirror LH	(NC 2214) ○
<b>BAG HOOKS</b>			
ARQ2128		Bag Hooks (Max load up to 6kg = impact on max user weight)	€95
<b>CANE HOLDER</b>			
ARQ1312		Cane holder RH	€25
ARQ1363		Cane holder LH	€25
<b>TRAY</b>			
ARQ1404		Tray, Swing Away RH	(NC 1707, 2213, 2214, 2209) €373
ARQ1405		Tray, Swing Away LH	(NC 1708, 2213, 2214, 2210) €373

<b>TOTAL PRICE</b>	
--------------------	--

**Non Conforming Options:**

ARQ0071	TDX SP2 Low Rider. 73Ah Batteries. 420mm Seat to Ground Height
ARQ0072	TDX SP2 Low Rider. 60Ah Batteries. 403mm Seat to Ground Height
ARZ0075	Special TDX SP2 Low Rider. 73Ah Batteries. 420mm Seat to Ground Height (Extended Lead Time). Please provide details in comments section.
ARZ0076	Special TDX SP2 Low Rider. 60Ah Batteries. 403mm Seat to Ground Height (Extended Lead Time). Please provide details in comments section.
ARQ2122	430 mm Seat Plate
ARQ2123	480 mm Seat Plate
ARQ2124	530 mm Seat Plate
ARQ2132	410 mm
ARQ2134	460 mm
ARQ2138	560 mm
ARQ2143	Matrx Flo-tech Contour Visco; (2-way stretch material, 430 mm x 430 mm only)
ARQ1980	Tapered Wedges (Pair)
ARQ1981	Abductor
ARQ1982	Lateral Wedge
ARQ1983	Anterior Seat Wedges
ARQ2163	Backposts short (480 mm). For Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery. Width adjustable cast bracket.
ARQ2108	Backposts WITHOUT push handles, Push handles will be cut down on the vertical section.
ARQ2171	Tension Adjustable Upholstery with Nylon overlay. Short = 480 mm / 7 straps. Long = 540 mm / 8 straps. Measured from seat plate to top of upholstery.
ARQ2179	Flex3 Tension Adjustable (540 mm or 480 mm)
ARQ2172	Flex3 Adjustable Backrest Plates (without Backrest Cushion)
ARQ2173	Flex3 Adjustable Backrest Plates with Laguna Backrest cushion
ARQ2174	Matrx Posture Back, height 400 mm (Mesh fabric) Corresponds to SW setting chosen
ARQ2175	Matrx Posture Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen
ARQ2176	Matrx Elite E2 Back , height 500 mm (Mesh fabric) Corresponds to SW setting chosen
ARQ2177	Matrx Elite TR Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen
ARQ2178	Matrx Elite E2 Deep Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen
ARQ2170	Without Backrest Upholstery (Back Posts/Back Frame only)
ARQ2151	Black Nylon Fabric TR18
ARQ2152	Black Multi Stretch PU Coated Wipe Down Fabric TR26
ARQ2155	Black Spacer Fabric TR36
ARQ2193	Flip-up Armrest. Height and depth adjustable (height adjustable from 230mm - 360mm). 150 kg max user weight.
ARQ2192	Following Armrest (height adjustable from 230mm - 300mm). 150 kg max user weight
ARQ2194	Following Armrest (height adjustable from 300mm - 360mm). 150 kg max user weight
ARQ2196	Long Armpad (375mm)
ARQ2197	Long Ergonomic Armpad (397mm)
ARQ2199	Soft Arpad (REA Azalea Wide Pad)
ARQ2184	Hip Support LH, pad size 140 x 100 mm
ARQ2186	Hip Support LH, pad size 230 x 100 mm
ARQ2188	Hip Support LH, pad size 330 x 100 mm

ARQ2185	Hip Support RH, pad size 140 x 100 mm
ARQ2187	Hip Support RH, pad size 230 x 100 mm
ARQ2189	Hip Support RH, pad size 330 x 100 mm
ARQ1350	Neckrest depth adjustable
ARQ1432	Headrest sideways adjustable
ARQ1427	Headrest with cheek support
ARQ2324	Elan Standard Pad 255 mm / 10" wide w/standard cover (ESP26). Complete with Elan Headrest Hardware
ARQ2325	Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (EOP23). Complete with Elan Headrest Hardware
ARQ2327	Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard cover. Complete with Elan Headrest Hardware
ARQ1496	Non Invacare Headrest but profile mount fitting Flex3 hard back
ARQ1445	Lateral Support with extended arm, multi-adjustable pad
ARQ2053	Seat Lifter with Powered Tilt (300mm, 45°) with CoG shift (Seat height 460mm, 136kg max user weight)
ARQ0526	Side Mounted: Standard hanger 80° pair. Lower leg length: 350 - 450 mm
ARQ0503	Side Mounted: Vari F pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with tools
ARQ0533	Side Mounted: Vari A pair. Lower leg length: 290 - 460 mm. Angle adjustment: 70 - 0° with handwheel and calf pad
ARQ0531	Side Mounted: ADE pair. Lower leg length: 290 - 460 mm. Powered angle adjustment: 80 - 0° with calf pad
ARQ0559	Angle & width adjustable footplates M (140 W x 240 mm D / 5.5"x9.5")
ARQ0512	Fixed Footplate pair
ARQ0504	Angle Adjustable Footplate pair
ARQ0513	Angle and depth adjustable Footplate pair
ARQ2203	DLX-REM211 (Drive & Actuators)
ARQ2204	DLX-REM216 (Drive, Actuators & Lights)
ARQ2205	DLX-REM400 (Touch Screen Interface. Drive, Act. & Lights. Configurable jack sockets)
ARQ2206	DLX-REM500 (Touch Screen Interface Display Unit - requires a user input/alternative joystick).
ARQ2270	Assignable Toggle Switch Module for REM400 (DLX-SW400TGL)
ARQ1208	Standard Bracket LH
ARQ1707	Swing Away RH Bracket
ARQ1708	Swing Away LH Bracket
ARQ1709	Matrx Resolve Swing in/out Bracket RH
ARQ1710	Matrx Resolve Swing in/out Bracket LH
ARQ2213	Centre Mount Swing Away to RH Nucleus Bracket
ARQ2214	Centre Mount Swing Away to LH Nucleus Bracket
ARQ2209	REM500 Swing Away RH Multi-Adjustable Bracket
ARQ2210	REM500 Swing Away LH Multi-Adjustable Bracket
ARQ2227	Bracket RH mounted for Attendant
ARQ2228	Bracket LH mounted for Attendant
ARQ2220	LiNX DLX-REM050 (Dual Control used as secondary remote)
ARQ2219	LiNX DLX-ACU200 (Compact Dual Control)
ARQ2245	LiNX DLX-IDC (Intuitive Dual Control)
ARQ1208	Dual Control Standard Bracket LH
ARQ1221	IDC Dual Control Detachable Bracket RH
ARQ1222	IDC Dual Control Detachable Bracket LH

ARQ2231	LiNX DLX-CR400 (Compact Remote)
ARQ2232	LiNX DLX-CR400LF (Compact Remote Low Force)
ARZ2233	ASL MEC (Proportional MEC Joystick ASL130)
ARZ2240	ASL Heavy Duty (Proportional Joystick w mode switch ASL133)
ARQ2207	RH Mount
ARQ2208	LH Mount
ARQ2213	Centre Mount Swing Away to RH (Nucleus)
ARQ2214	Centre Mount Swing Away to LH (Nucleus)
ARZ2235	Tray Mount RH (ASL603, 100mm x 150mm)
ARZ2236	Tray Mount LH (ASL603, 100mm x 150mm)
ARZ2241	Centre Mount, Swing Away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
ARZ2242	Centre Mount, Swing Away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)
ARZ2215	Individual mounting option per request (Special). Requires additional information/specification.
ARQ2230	ASL104 ATOM Head Control with ultrasonic proximity sensors with standard wings
ARQ2247	ASL104P PROTON Head Control with ultrasonic proximity sensors with PROTON wings
ARQ2227	Bracket RH, mounted for attendant
ARQ2228	Bracket LH, mounted for attendant
ARQ2253	USB Charger Port DLX-USB02-A
ARQ2207	RH Mount
ARQ2208	LH Mount
ARQ1686	Yellow Buddy Button
ARQ1687	Green Buddy Button
ARQ1690	Blue Buddy Button
ARQ2254	Egg Switch, black
ARQ1711	Red Piko Button
ARQ1712	Yellow Piko Button
ARQ1713	Green Piko Button
ARQ1714	Orange Piko Button
ARQ1715	Blue Piko Button
ARQ1716	Black Piko Button
ARZ2207	Toggle / Button RH Mount
ARZ2208	Toggle / Button LH Mount
ARZ2258	Remote Stop Switch (ASL504)
ARQ1282	6 km/h (4mph) (4 pole)
ARQ1908	10 km/h (6mph) (includes light system and indicators) G-Trac recommended
ARQ1909	12 km/h (8mph) (includes light system and indicators) G-Trac recommended
ARQ1329	Lapbelt (adjustable both sides) fitted to base of backrest
ARQ2113	Lapbelt (adjustable both sides) adjustable mounting on seat rail extrusion
ARQ1332	H-Belt
ARQ2701	Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size L (padded length 620mm)
ARQ2702	Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size M (padded length 460mm)
ARQ2703	Bodypoint 2-Point Hip Belt, Centre pull - push button, Cinch mount size L (padded length 620mm)
ARQ2704	Bodypoint 2-Point Hip Belt, Centre pull - push button, Cinch mount size M (padded length 460mm)
ARQ2705	Bodypoint 2-Point Hip Belt, Dual pull Side release button, Cinch mount size L (padded length 620mm)

---

ARQ2706	Bodypoint 2-Point Hip Belt, Dual pull Side release button, Cinch mount size M (padded length 460mm)
ARQ2128	Bag Hooks (Max load up to 6kg = impact on max user weight)

---



For a more complex customer requirement,  
please email: [orders@invacare.com.au](mailto:orders@invacare.com.au)

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):