



Yes, you can:

Powerchairs Ireland Prescription Form

NOVEMBER 2023
Customer Ref:

Please refer to www.invacare.ie for Local Distributor

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

TDX SP2 Low Rider

Standard Features

- Max user weight 160kg (25.2 Stone)
- Crash tested ISO7176-Part 19
- Compact, manoeuvrable, CWD TDX SP2
- Telescopic adjustable Seat Widths from 330mm - 530mm (13" - 21")
- Seat Depths 410mm - 560mm Stepless Adjustment (16" - 22")
- Tension Adjustable backrest
- Standard T-Armrests, detachable with Short Armpad
- Seat lifter with powered Tilt (300mm, 30°) CoG Shift
- Choice of seat height either 420mm or 403mm
- Manual settable Recliner 0° - 30° , in 3.75° steps
- Standard 80° Hanger with fixed footplate
- LiNX Control REM211 with standard bracket
- Selection of 10 Stunning Shroud and Rim colours
- Centre (Driving Wheels) 355mm x 76mm (14" x 3") Black Tyres
- Suspension and Sure Step Technology
- 6 km/h 4 Pole Motors
- 8A (200 - 25V) Charge
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

Standard Base price: €10,409



Q'STRAIT
QLK-150
Compatible

Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code	Description	
BASE PRODUCT*		
ARQ0071	TDX SP2 Low Rider 73ah Batteries 420mm seat to ground height	○
ARQ0072	TDX SP2 Low Rider 60ah Batteries 403mm seat to ground height	○
ARZ0071	Special TDX SP2 Low Rider 73ah Batteries 420mm seat to ground height	POA
ARZ0072	Special TDX SP2 Low Rider 60ah Batteries 403mm seat to ground height	POA

User Weight*			
ARQ1918	<input type="checkbox"/>	User weight less 100 kg	○
ARQ1919	<input type="checkbox"/>	User weight 101 to 125 kg	○
ARQ1921	<input type="checkbox"/>	User weight 126 to 136 kg (NC 1350, 1432, 1427)	○
ARQ1922	<input type="checkbox"/>	User weight 137 to 150 kg (NC 1350, 1432, 1427)	○
ARQ1923	<input type="checkbox"/>	User weight 151 to 160 kg (NC 2192, 2193, 2194, 1908, 1350, 1432, 1427)	○
SEATING TYPE*			
ARQ2100	<input type="checkbox"/>	Modulite Seating Unit	€10,409
SEAT FRAME*			
ARQ2111	<input type="checkbox"/>	Telescopic Seat Frame with metal seat plate	◆
SEAT WIDTH*			
ARQ2122	<input type="checkbox"/>	Set at 430mm (17") (adjustable from 430mm to 480mm (17" - 19)) (NC 2703, 2705)	○
ARQ2123	<input type="checkbox"/>	Set at 480mm(19") (adjustable from 480mm to 530mm (19" - 21")) (NC 2702, 2704, 2706)	○
ARQ2124	<input type="checkbox"/>	Set at 530mm (21") (adjustable from 530mm to 580mm (21" - 23")) (NC 2174, 2175, 2176, 2177, 2178, 2143, 2702, 2704, 2706)	○
SEAT DEPTH RAIL LENGTH*			
ARQ2129	<input type="checkbox"/>	Standard Depth Package (16" - 20", 410mm - 510mm) (NC 2138)	○
ARQ2130	<input type="checkbox"/>	Extended Depth Package (20" - 22", 510mm - 560mm) (NC 2132, 2134)	○
SEAT DEPTH*			
ARQ2132	<input type="checkbox"/>	410mm (16") (410mm - 510mm) (NC 2196, 2197, 2701, 2702, 2703, 2704, 2705, 2706)	○
ARQ2134	<input type="checkbox"/>	460mm (18") (410mm - 510mm)	○
ARQ2136	<input type="checkbox"/>	510mm (20") (410mm - 510mm)	○
ARQ2138	<input type="checkbox"/>	560mm (22") (510mm - 560mm) (NC 2143)	○
SEAT OPTION			
ARQSWO	<input type="checkbox"/>	"H" Frame Interface Kit (NC 2143)	€507
SEAT CUSHION*			
ARQ2142	<input type="checkbox"/>	Matrx Flo-tech Lite	€132
ARQ2143	<input type="checkbox"/>	Matrx Flo-tech Contour Visco (NC 2136, 2124, 2138)	€132
ARQ2147	<input type="checkbox"/>	Matrx Libra 90mm thick	€492
ARQ2140	<input type="checkbox"/>	Without seat cushion	○
BACKREST SUPPORT*			
ARQ2162	<input type="checkbox"/>	Backposts 540mm (Tall) must be selected for Flex3 adjustable backrest plates	○
ARQ2163	<input type="checkbox"/>	Backposts 480mm (Short) (NC 2173)	○
PUSH HANDLES FOR BACKPOSTS*			
ARQ2107	<input type="checkbox"/>	Backposts WITH push handles	○
ARQ2108	<input type="checkbox"/>	Backposts WITHOUT push handles (NC 1332, 2220, 2219, 2245, 2227, 2228)	€50

BACKREST UPHOLSTERY / CUSHION*			
ARQ2171		Basic Tension Adjustable Backrest	(NC 1445, 1496, 2155, 2230, 2247) ◆
ARQ2179		Flex3 Tension Adjustable (540mm or 480mm) Black multi stretch PU coated fabric only	(NC 2151, 2192, 2193, 2194, 2199, 2229, 2155, 2230, 2247, 2128) €138
ARQ2172		Flex3 Adjustable Backrest Plates (without Backrest Cushion)	(NC 2163, 2128) €265
ARQ2173		Flex3 Adjustable Backrest Plates with Laguna Backrest cushion	(NC 2163, 2128) €460
ARQ2174		Matrx Posture Back, height 400mm (Mesh fabric)	(NC 1496, 2124) €572
ARQ2175		Matrx Posture Back, height 500mm (Mesh fabric)	(NC 1496, 2124) €572
ARQ2176		Matrx Elite Back, height 500mm (Mesh fabric)	(NC 1496, 2124) €572
ARQ2177		Matrx Elite TR Back, height 500mm (Mesh fabric)	(NC 1496, 2124) €572
ARQ2178		Matrx Elite Deep Back, height 500mm (Mesh fabric)	(NC 1496, 2124) €741
ARQ2170		Without Backrest Upholstery	(NC 1496, 2230, 2247) ○
UPHOLSTERY TYPE*			
ARQ2151		Black fabric	○
ARQ2152		Black multi stretch PU coated fabric	(NC 2171) ○
ARQ2155		Black Spacer fabric	(NC 2171, 2179) ○
ARMREST*			
ARQ2191		T-Armrest Standard, Range of adjustment	◆
ARQ2181		T-Armrest Extended height	€50
ARQ2183		In-rigged Armrests	(NC 2185, 2186, 2187, 2188, 2189) €29
ARQ2193		Flip-up Armrest (from 230mm - 360mm)	(NC 2176, 2199) €57
ARQ2192		Following Armrest (from 230mm - 300mm)	(NC 1445, 2179) €210
ARQ2194		Following Armrest (from 300mm - 360mm)	(NC 1445, 2179) €210
ARMPAD*			
ARQ2195		Short Armpad (325mm)	◆
ARQ2196		Long Armpad (375mm)	○
ARQ2197		Long Ergonomic Armpad (397mm)	(NC 2247) ○
ARQ2199		Soft Arpad (REA Azalea Wide Pad)	€108
HIP SUPPORTS LH			
ARQ2184		SMALL 140mm x 100mm	(NC 2702, 2703, 2704, 2705, 2706) €286
ARQ2186		MEDIUM 230mm x 100mm	(NC 2702, 2703, 2704, 2705, 2706) €286
ARQ2188		LARGE 330mm x 100mm	(NC 2702, 2703, 2704, 2705, 2706) €286
HIP SUPPORTS RH			
ARQ2185		SMALL 140mm x 100mm	(NC 2702, 2703, 2704, 2705, 2706) €286
ARQ2187		MEDIUM 230mm x 100mm	(NC 2702, 2703, 2704, 2705, 2706) €286
ARQ2189		LARGE 330mm x 100mm	(NC 2702, 2703, 2704, 2705, 2706) €286
MOUNTING HARDWARE FOR HIP SUPPORTS			
ARQ2165		Mounting hardware for hip supports (8.1" / 206mm high)	○

HEADREST / NECKREST*			
ARQ1350		Neckrest height and depth adjustable, Max User weight 136kg (NC 2247)	€86
ARQ1432		Headrest, height, depth and sideways adjustable, Max User weight 136kg (NC 2247)	○
ARQ1427		Headrest with cheek support, height and depth adjustable, Max User weight 136kg (NC 2247)	€239
ARQ2324		Elan Standard Pad 355mm (10") wide (NC 2247)	€406
ARQ2325		Elan Occipital Pad 230mm x 180mm (9" x 7") (NC 2247)	€406
ARQ2327		Elan 4 Point Headrest Pad 280mm x 255mm (11" x 10") (NC 2247)	€548
ARQ1496		All Bracketry as above option without Headrest cushion	○
ARQ1495		Push bar dia. 22mm	○
LATERAL SUPPORT			
ARQ1445		Multi-adjustable pad (Azalea type, multi stretch PU coated fabric) (NC 2171)	€188
SEAT LIFTER / SEAT TILT*			
ARQ1094		Seat lifter with Electrical Tilt	◆
SEAT TILT OPTION			
ARQ2054		5° Fixed Anterior Tilt Assist mount	€500
BACKREST RECLINER*			
ARQ0975		Manual Settable Recliner	◆
ARQ0977		Powered Recliner (0° - 30°) Foldable	€355
FOOTREST HANGER TYPE*			
ARQ0556		Centre Mounted: Manual : Hanger, centre mounted (NC 0504, 0513)	€337
ARQ0589		Centre Mounted : LNX Powered: Centre mounted legrest (NC 0504, 0513)	€2,013
ARQ0526		Side Mounted : Standard 80° hangers (NC 0559)	◆
ARQ0503		Side Mounted : Vari F : Angle adjustable	○
ARQ0533		Side Mounted : Vari A : Angle adjustable (NC 0559)	€215
ARQ0531		Side Mounted : ADE Powered (NC 0513, 0559)	€740
ARQ0501		Without Hanger (NC 0504, 0513, 0559)	○
FOOTPLATE TYPE*			
ARQ0512		Fixed Footplate (NC 0533, 0556, 0589, 0501)	◆
ARQ0504		Angle Adjustable Footplate (NC 0589, 0501)	€79
ARQ0513		Angle and depth adjustable Footplate (Adapt) (NC 0589, 0501)	€128
ARQ0559		Angle Adjustable Footplate for Centre mounted legrest (NC 0503, 0533, 0531, 0501)	○
ARQXXXX		Without footplates f (NC 0589, 0526, 0503, 0533, 0531)	○
CONTROLS*			
ARQ2203		REM211 Drive & Actuators (NC 2209, 2210, 2213, 2214, 2229, 2231, 2232, 2230, 2247, 2233, 2240, 2270, 2258, 1908, 1909)	◆
ARQ2204		REM216 Drive, Actuators & Lights (NC 2210, 2213, 2214, 2229, 2231, 2232, 2230, 2247, 2233, 2240, 2270, 2258)	€162
ARQ2205		REM400, touch-screen Interface, built in Bluetooth HID (NC 2210, 2213, 2215)	€495
ARQ2206		REM500 (NC 1208, 1707, 1708, 2227, 2228, 2270)	€909

CONTROL BRACKETRY*			
ARQ1207		Standard Bracket RH	○
ARQ1208		Standard Bracket LH	○
ARQ1707		Swing Away RH Bracket	€159
ARQ1708		Swing Away LH Bracket	€159
ARQ2213		Centre mount swing away to RH Nucleus	€739
ARQ2214		Centre mount swing away to LH Nucleus	€739
ARQ2227		Bracket RH mounted for Attendant (NC 2203, 2204, 2205)	€71
ARQ2228		Bracket LH mounted for Attendant (NC 2203, 2204, 2205)	€71
ARZ2239		Individual mounting option per (special)	POA
ARQ2209		REM500 Swing away RH multi-adjustable (NC 2203, 2204, 2205)	○
ARQ2210		REM500 Swing away LH multi-adjustable (NC 2203, 2204, 2205)	○
REM400 TOGGLE OPTION			
ARQ2270		Assignable Toggle Switch Module for REM400, (NC 2203, 2204)	€329
CHIN CONTROL			
ARQ2229		Proportional Chin Joystick (NC 2203, 2204, 2205)	€3,106
CHIN CONTROL BRACKETRY			
ARQ1289		Chin Control Swing Away Bracket RH Manual	○
ARQ1288		Chin Control Swing Away Bracket LH Manual	○
ARQ1287		Chin Control Swing Away Bracket RH Powered	€2,500
ARQ1286		Chin Control Swing Away Bracket LH Powered	€2,500
DUAL CONTROL			
ARQ2220		REM050 Dual Control (NC 1221, 1222, 2227, 2228, 2205, 2206, 2245, 2128, 2230, 2247, 2233, 2240)	€489
ARQ2219		Compact Dual Control (DLX-ACU200) (NC 1222, 2227, 2228, 2245, 2128)	€489
DUAL CONTROL BRACKETRY			
ARQ1207		Dual Control Standard Bracket RH	○
ARQ1208		Dual Control Standard Bracket LH	○
INVACARE LiNX IDC DUAL CONTROL			
ARQ2245		IDC Dual control (NC 2219, 2227, 2228, 2128)	€748
INVACARE LiNX IDC DUAL CONTROL BRACKETRY			
ARQ1221		IDC Dual Control Detachable Bracket RH (NC 2219, 2227, 2228)	○
ARQ1222		IDC Dual Control Detachable Bracket LH (NC 2219, 2227, 2228)	○
LINX COMPACT REMOTE			
ARQ2231		Proportional Compact Remote DLX-CR400 (NC 2203, 2204)	€845
ARQ2232		Proportional Compact Remote Low Force DLX-CR400LF (NC 2203, 2204)	€1,199
LINX COMPACT REMOTE BRACKET			
ARQ2207		Proportional Compact RH Mount Bracket (NC 2203, 2204)	○
ARQ2208		Proportional Compact LH Mount Bracket (NC 2203, 2204)	○
ARQ2213		Proportional Compact Centre Mount swing away to RH (Nucleus) (NC 2203, 2204)	€739
ARQ2214		Proportional Compact Centre Mount swing away to LH (Nucleus) (NC 2203, 2204)	€739
ARZ2215		Individual mounting option per request (Special) (NC 2203, 2204)	POA
LINX HEAD CONTROL			
ARQ2230		Head Control with ultrasonic proximity sensors (NC 2179, 2197, 2203, 2204, 2220, 2229, 2233, 2324, 2325, 2327)	€4,497
ARQ2247		Head Control with ultrasonic proximity sensors (NC 2179, 2197, 2203, 2204, 2220, 2229, 2233, 2324, 2325, 2327)	€5,248
LINX MICRO EXTREMITY CONTROL			
ARZ2233		MEC Proportional Joystick (NC 2203, 2204)	€4,230

LINX MICRO EXTREMITY CONTROL BRACKET				
ARZ2235		MEC Tray Mount RH	(NC 2203, 2204)	○
ARZ2236		MEC Tray Mount LH	(NC 2203, 2204)	○
ARZ2213		Centre Mount, swing away to RH, Nucleus	(NC 2203, 2204)	€739
ARZ2214		Centre Mount, swing away to LH, Nucleus	(NC 2203, 2204)	€739
ARZ2241		Centre Mount, swing away to RH, Nucleus	(NC 2203, 2204)	€739
ARZ2242		Centre Mount, swing away to LH, Nucleus	(NC 2203, 2204)	€739
ARZ2215		Individual mounting option per request (Special)	(NC 2203, 2204)	POA
LINX HEAVY DUTY REMOTE				
ARZ2240		Proportional Heavy Duty Remote Joystick with mode switch	(NC 2203, 2204)	€2,683
LINX HEAVY DUTY REMOTE BRACKET				
ARZ2235		Tray Mount RH, (ASL603, 100mm x 150mm)	(NC 2203, 2204)	○
ARZ2236		Tray Mount LH, (ASL603, 100mm x 150mm)	(NC 2203, 2204)	○
ARZ2213		Centre Mount, swing away to RH, Nucleus	(NC 2203, 2204)	€739
ARZ2214		Centre Mount, swing away to LH, Nucleus	(NC 2203, 2204)	€739
ARZ2241		Centre Mount, swing away to RH, Nucleus	(NC 2203, 2204)	€739
ARZ2242		Centre Mount, swing away to LH, Nucleus	(NC 2203, 2204)	€739
ARZ2215		Individual mounting option per request (Special)	(NC 2203, 2204)	POA
LINX 10-WAY SWITCH				
ARQ2251		10-Way Switch DLX-SW10-A for direct Powered Seating Access		€539
LINX 10-WAY SWITCH BRACKET				
ARQ2217		10-Way Switch Bracket RH, Mounted for User	(NC 2213)	○
ARQ2218		10-Way Switch Bracket LH, Mounted for User	(NC 2214)	○
ARQ2227		10-Way Switch Bracket RH, Mounted for Attendant		○
ARQ2228		10-Way Switch Bracket LH, Mounted for Attendant		○
LINX ACCESSORY				
ARQ2252		Loom providing jack sockets for Drive and Mode		€110
LINX USB CHARGER PORT				
ARQ2253		USB Charger Port DLX-USB02-A	(NC 2701, 2703, 2704, 2705, 2706)	€162
LINX USB CHARGER PORT BRACKET				
ARQ2207		USB Charger Port RH Mount		○
ARQ2208		USB Charger Port LH Mount		○
FURTHER LINX COMPONENTS / ACCESSORIES				
ARQ2300		To be enclosed only		○
ARZ2301		To be connected and programmed (special)		POA
BUDDY BUTTONS				
ARQ1685		Red Buddy Button		€90
ARQ1686		Yellow Buddy Button		€90
ARQ1687		Green Buddy Button		€90
ARQ1690		Blue Buddy Button		€90
EGG SWITCH				
ARQ2254		Egg Switch, black		€90

PIKO BUTTONS		
ARQ1711	<input type="checkbox"/>	Red Piko Button €161
ARQ1712	<input type="checkbox"/>	Yellow Piko Button €161
ARQ1713	<input type="checkbox"/>	Green Piko Button €161
ARQ1714	<input type="checkbox"/>	Orange Piko Button €161
ARQ1715	<input type="checkbox"/>	Blue Piko Button €161
ARQ1716	<input type="checkbox"/>	Black Piko Button €161
TOGGLES / BUTTONS		
ARZ2261	<input type="checkbox"/>	Stereo Toggle - direct control for Tilt only €273
ARZ2262	<input type="checkbox"/>	Stereo Button - direct control for Tilt only €263
ARZ2263	<input type="checkbox"/>	Four-Way Toggle - direct control for Tilt & Recline €342
ARZ2264	<input type="checkbox"/>	Four-Way Button - direct control for Tilt & Recline €365
ARZ2207	<input type="checkbox"/>	Toggle / Button RH Mount ○
ARZ2208	<input type="checkbox"/>	Toggle / Button LH Mount ○
ARZ2225	<input type="checkbox"/>	Without Mounting ○
BODYPOINT JOYSTICK KNOBS		
ARQ2721	<input type="checkbox"/>	3" U-Shaped, Joystick Handle, Bodypoint PC101A €75
ARQ2722	<input type="checkbox"/>	4" U-Shaped, Joystick Handle, Bodypoint PC102A €75
ARQ2723	<input type="checkbox"/>	Domed Rubber, Joystick Handle, Bodypoint PC107A €75
ARQ2724	<input type="checkbox"/>	Mushroom Grooved, Joystick Handle, Bodypoint PC110A €75
REMOTE STOP SWITCH		
ARZ2258	<input type="checkbox"/>	Remote Stop Switch (ASL504) €1,070
SHROUD COLOUR*		
ARQ0228	<input type="checkbox"/>	Wet Black ○
ARQ2601	<input type="checkbox"/>	SpaceStation Silver ○
ARQ2602	<input type="checkbox"/>	Island Blue ○
ARQ2603	<input type="checkbox"/>	Grape JellyBean (Purple) ○
ARQ2607	<input type="checkbox"/>	Rockstar Red ○
ARQ2608	<input type="checkbox"/>	Lights Out Black - Matt ○
ARQ2609	<input type="checkbox"/>	Invacare Blue ○
RIM COLOUR*		
ARQ2611	<input type="checkbox"/>	SpaceStation Silver ○
ARQ2612	<input type="checkbox"/>	Island Blue ○
ARQ2613	<input type="checkbox"/>	Grape JellyBean (Purple) ○
ARQ2617	<input type="checkbox"/>	Rockstar Red ○
ARQ2618	<input type="checkbox"/>	Lights Out Black - Matt ○
ARQ2619	<input type="checkbox"/>	Invacare Blue ○
ARQ2620	<input type="checkbox"/>	Wet Black ○
ARQ2621	<input type="checkbox"/>	Without Rim insert ○
FORKS*		
ARQ2631	<input type="checkbox"/>	Single sided fork (New) includes Black castor €358
ARQ2632	<input type="checkbox"/>	Double sided fork ○

CENTRE WHEELS (DRIVING WHEELS)*			
ARQ0704		355mm x 76mm (14" x 3") Black Tyres, Pneumatic	○
ARQ0733		355mm x 76mm (14" x 3") Black Tyres with Gel Foam Inserts	○
MOTOR TYPE*			
ARQ1294		3 km/h (2mph)	○
ARQ1282		6 km/h (4mph) (4 pole)	◆
ARQ1908		10 km/h (6mph) (includes light system and indicators) G-Trac recommended	○
ARQ1909		12 km/h (8mph) (includes light system and indicators) G-Trac recommended	○
G TRACK			
ARQ1215		G-Trac	€908
CHARGER*			
ARQ1343		OptiCharge Charger 8A (200 - 250V)	○
ARQ1313		OptiCharge Charger 10A (200 - 250V)	€174
LIGHTS*			
ARQ1437		Light system with Indicators & LED Headlights	(NC 1908, 1909) €519
ARQ1254		Without Light system	(NC 1909) ◆
ARQ1006		Lights Included with Motor Package 10 & 12kmh	(NC 1282) ○
TRANSPORTATION KIT			
ARQ1450		Preparation for DAHL Docking station:	€624
ARQ1494		4-Point tie-down	○
FIXATION BELT*			
ARQ1329		Lapbelt (adjustable both sides)	(NC 2702, 2703, 2704, 2705, 2706) ○
ARQ2113		Lapbelt (adjustable both sides)	○
ARQ1332		H-Belt	€80
ARQ2701		Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size L (padded length 620mm)	(NC 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113) €117
ARQ2702		Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size M (padded length 460mm)	(NC 2123, 2124, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113) €117
ARQ2703		Bodypoint 2-Point Hip Belt, Centre pull - push button	(NC 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113) €85
ARQ2704		Bodypoint 2-Point Hip Belt, Centre pull - push button	(NC 2123, 2124, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113) €85
ARQ2705		Bodypoint 2-Point Hip Belt, Dual pull Side release	(NC 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113) €85
ARQ2706		Bodypoint 2-Point Hip Belt, Dual pull Side release	(NC 2123, 2124, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113) €85
REAR-VIEW MIRROR			
ARQ1415		Rear View Mirror RH	(NC 2213, 1287) ○
ARQ1339		Rear View Mirror LH	(NC 2214, 1286) ○
BAG HOOKS			
ARQ2128		Bag Hooks (Max load up to 6kg = impact on max user weight)	€95
CANE HOLDER			
ARQ1312		Cane holder RH	€25
ARQ1363		Cane holder LH	€25
TRAY			
ARQ1404		Tray Swing Away RH	(NC 1707, 2213, 2214, 2209, 1287) €373
ARQ1405		Tray Swing Away LH	(NC 1708, 2213, 2214, 2210, 1286) €373

TOTAL PRICE	
--------------------	--

Non Conforming Options:

ARQ2122	Set at 430mm (17") (adjustable from 430mm to 480mm (17" - 19))
ARQ2123	Set at 480mm(19") (adjustable from 480mm to 530mm (19" - 21"))
ARQ2124	Set at 530mm (21") (adjustable from 530mm to 580mm (21" - 23"))
ARQ2132	410mm (16") (410mm - 510mm)
ARQ2134	460mm (18") (410mm - 510mm)
ARQ2136	510mm (20") (410mm - 510mm)
ARQ2138	560mm (22") (510mm - 560mm)
ARQ2143	Matrx Flo-tech Contour Visco
ARQ2163	Backposts 480mm (Short)
ARQ2171	Basic Tension Adjustable Backrest
ARQ2179	Flex3 Tension Adjustable (540mm or 480mm) Black multi stretch PU coated fabric only
ARQ2173	Flex3 Adjustable Backrest Plates with Laguna Backrest cushion
ARQ2174	Matrx Posture Back, height 400mm (Mesh fabric)
ARQ2175	Matrx Posture Back, height 500mm (Mesh fabric)
ARQ2176	Matrx Elite Back, height 500mm (Mesh fabric)
ARQ2177	Matrx Elite TR Back, height 500mm (Mesh fabric)
ARQ2178	Matrx Elite Deep Back, height 500mm (Mesh fabric)
ARQ2151	Black fabric
ARQ2155	Black Spacer fabric
ARQ2193	Flip-up Armrest (from 230mm - 360mm)
ARQ2192	Following Armrest (from 230mm - 300mm)
ARQ2194	Following Armrest (from 300mm - 360mm)
ARQ2196	Long Armpad (375mm)
ARQ2197	Long Ergonomic Armpad (397mm)
ARQ2199	Soft Arpad (REA Azalea Wide Pad)
ARQ2184	SMALL 140mm x 100mm
ARQ2186	MEDIUM 230mm x 100mm
ARQ2188	LARGE 330mm x 100mm
ARQ2185	SMALL 140mm x 100mm
ARQ2187	MEDIUM 230mm x 100mm
ARQ2189	LARGE 330mm x 100mm
ARQ1350	Neckrest height and depth adjustable, Max User weight 136kg
ARQ1432	Headrest, height, depth and sideways adjustable, Max User weight 136kg
ARQ1427	Headrest with cheek support, height and depth adjustable, Max User weight 136kg
ARQ2324	Elan Standard Pad 355mm (10") wide
ARQ2325	Elan Occipital Pad 230mm x 180mm (9" x 7")
ARQ2327	Elan 4 Point Headrest Pad 280mm x 255mm (11" x 10")
ARQ1496	All Bracketry as above option without Headrest cushion
ARQ1445	Multi-adjustable pad (Azalea type, multi stretch PU coated fabric)
ARQ0556	Centre Mounted: Manual : Hanger, centre mounted
ARQ0589	Centre Mounted : LNX Powered: Centre mounted legrest
ARQ0526	Side Mounted : Standard 80° hangers
ARQ0503	Side Mounted : Vari F : Angle adjustable

ARQ0533	Side Mounted : Vari A : Angle adjustable
ARQ0531	Side Mounted : ADE Powered
ARQ0501	Without Hanger
ARQ0504	Angle Adjustable Footplate
ARQ0513	Angle and depth adjustable Footplate (Adapt)
ARQ0559	Angle Adjustable Footplate for Centre mounted legrest
ARQ2203	REM211 Drive & Actuators
ARQ2204	REM216 Drive, Actuators & Lights
ARQ2205	REM400, touch-screen Interface, built in Bluetooth HID
ARQ2206	REM500
ARQ1208	Standard Bracket LH
ARQ1707	Swing Away RH Bracket
ARQ1708	Swing Away LH Bracket
ARQ2213	Centre mount swing away to RH Nucleus
ARQ2214	Centre mount swing away to LH Nucleus
ARQ2227	Bracket RH mounted for Attendant
ARQ2228	Bracket LH mounted for Attendant
ARQ2209	REM500 Swing away RH multi-adjustable
ARQ2210	REM500 Swing away LH multi-adjustable
ARQ2270	Assignable Toggle Switch Module for REM400,
ARQ2229	Proportional Chin Joystick
ARQ1287	Chin Control Swing Away Bracket RH Powered
ARQ1286	Chin Control Swing Away Bracket LH Powered
ARQ2220	REM050 Dual Control
ARQ2219	Compact Dual Control (DLX-ACU200)
ARQ1208	Dual Control Standard Bracket LH
ARQ2245	IDC Dual control
ARQ1221	IDC Dual Control Detachable Bracket RH
ARQ1222	IDC Dual Control Detachable Bracket LH
ARQ2231	Proportional Compact Remote DLX-CR400
ARQ2232	Proportional Compact Remote Low Force DLX-CR400LF
ARQ2213	Proportional Compact Centre Mount swing away to RH (Nucleus)
ARQ2214	Proportional Compact Centre Mount swing away to LH (Nucleus)
ARZ2215	Individual mounting option per request (Special)
ARQ2230	Head Control with ultrasonic proximity sensors
ARQ2247	Head Control with ultrasonic proximity sensors
ARZ2233	MEC Proportional Joystick
ARZ2213	Centre Mount, swing away to RH, Nucleus
ARZ2214	Centre Mount, swing away to LH, Nucleus
ARZ2215	Individual mounting option per request (Special)
ARZ2240	Proportional Heavy Duty Remote Joystick with mode switch
ARZ2213	Centre Mount, swing away to RH, Nucleus
ARZ2214	Centre Mount, swing away to LH, Nucleus
ARZ2215	Individual mounting option per request (Special)

ARQ2227	10-Way Switch Bracket RH, Mounted for Attendant
ARQ2228	10-Way Switch Bracket LH, Mounted for Attendant
ARQ2253	USB Charger Port DLX-USB02-A
ARZ2258	Remote Stop Switch (ASL504)
ARQ1282	6 km/h (4mph) (4 pole)
ARQ1908	10 km/h (6mph) (includes light system and indicators) G-Trac recommended
ARQ1909	12 km/h (8mph) (includes light system and indicators) G-Trac recommended
ARQ1329	Lapbelt (adjustable both sides)
ARQ2113	Lapbelt (adjustable both sides)
ARQ1332	H-Belt
ARQ2701	Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size L (padded length 620mm)
ARQ2702	Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size M (padded length 460mm)
ARQ2703	Bodypoint 2-Point Hip Belt, Centre pull - push button
ARQ2704	Bodypoint 2-Point Hip Belt, Centre pull - push button
ARQ2705	Bodypoint 2-Point Hip Belt, Dual pull Side release
ARQ2706	Bodypoint 2-Point Hip Belt, Dual pull Side release
ARQ2128	Bag Hooks (Max load up to 6kg = impact on max user weight)



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):