



**Yes, you can:**

Tel: 01656 776222  
 Email: ordersuk@invacare.com  
 Online Spares available at: www.invacare.co.uk

JANUARY 2025  
 Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## TDX SP2 MAX

### Standard Features



- Max user weight 180kg (28 stone)
- Crash tested ISO 7176 -Part 19
- Compact, manoeuvrable, CWD TDX SP2
- Seat Widths from 490mm - 610mm
- Seat Depths 560mm - 630mm
- Tension Adjustable backrest
- Armrest with reinforced depth adjustable bracket with Short armpad
- Manual settable seat Tilt
- Manual settable backrest Recline
- VARI F hangers with fixed footplate
- LiNX Control REM110 with standard bracket
- Selection of 10 Stunning Shroud and Rim colours
- Centre (Driving Wheels) 355mm x 76mm (14" x 3") Black Tyres
- Suspension and Sure Step Technology
- 6 km/h (4mph) 4 Pole Motors
- 8A (200 - 250V) Charger
- 73.5 Ah Batteries
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

**Standard Base price: £5,382**



**Q'STRAIT QLK-150**  
Compatible

Please make your selection with a cross in the box    ○ Free of charge    ◆ standard feature    POA special request

Customer weight =                      KG

Code	Description	
<b>CHAIR TYPE*</b>		
ARH0000	<input type="checkbox"/> TDX SP2 Max	£5,382
ARZ0000	<input type="checkbox"/> Special TDX SP2 Max (Please provide details in special section)	POA
<b>User Weight*</b>		
ARH1924	<input type="checkbox"/> User weight up to 180 kg	◆
<b>SEATING TYPE*</b>		
ARH1145	<input type="checkbox"/> MAX Seating system	◆

<b>SEAT WIDTH SETTING*</b>			
ARH0378	<input type="checkbox"/>	Seat Width of Seat Frame 550mm (+100mm adjustment in arms)	<input type="radio"/>
ARH0380	<input type="checkbox"/>	Seat Width of Seat Frame 610mm (+100mm adjustment in arms) (NC 0533)	<input type="radio"/>
<b>SEAT DEPTH SETTING SETTABLE IN 3 STEPS*</b>			
ARH0396	<input type="checkbox"/>	Seat Depth 560mm	<input type="radio"/>
ARH0397	<input type="checkbox"/>	Seat Depth 595mm	<input type="radio"/>
ARH0398	<input type="checkbox"/>	Seat Depth 630mm	<input type="radio"/>
<b>SEAT CUSHION*</b>			
ARH1484	<input type="checkbox"/>	MAX Seat Cushion	<input type="radio"/>
ARH1114	<input type="checkbox"/>	Without Seat Cushion (10mm seat board only)	<input type="radio"/>
<b>ADAPTOR TO REDUCE THE SEAT DEPTH BY 130MM</b>			
ARH0399	<input type="checkbox"/>	Adaptor to shorten the seat depth by 130mm (Required if seat width 490mm is selected in combination with seat depth 595mm or 630mm)	<input type="radio"/>
<b>BACKREST UPHOLSTERY*</b>			
ARH1106	<input type="checkbox"/>	Tension Adjustable	<input type="checkbox"/>
<b>UPHOLSTERY COLOUR*</b>			
ARH0278	<input type="checkbox"/>	Black multi stretch PU coated fabric	<input type="checkbox"/>
<b>ARMREST INCLUDING REINFORCED DEPTH ADJUSTABLE BRACKET*</b>			
ARH0635	<input type="checkbox"/>	Fixed Armrest, 250mm - 340mm tool adjustable in height, detachable	<input type="checkbox"/>
<b>ARMPADS*</b>			
ARH0647	<input type="checkbox"/>	Short armpad (325mm)	<input type="checkbox"/>
ARH0649	<input type="checkbox"/>	Long armpad (375mm)	<input type="radio"/>
<b>PUSH BAR*</b>			
ARH1307	<input type="checkbox"/>	Push bar	<input type="checkbox"/>
<b>HEADREST - Adequate headrest option must be selected*</b>			
ARH1432	<input type="checkbox"/>	Headrest: height, depth and sideways adjustable (Azalea type)	<input type="checkbox"/>
ARH1427	<input type="checkbox"/>	Headrest with cheek support: height and depth adjustable (Azalea type)	£199
ARH1495	<input type="checkbox"/>	Push bar dia. 22mm in conjunction with telescopic seat frame / backposts, No Headrest - no mounting hardware.	<input type="radio"/>
<b>SEAT TILT*</b>			
ARH1091	<input type="checkbox"/>	Manually settable Seat Tilt (Max user weight 180kg)	<input type="checkbox"/>
ARH1095	<input type="checkbox"/>	Electric heavy duty basic Tilt (0° - 20°) select either 470mm or 495mm seat to ground height. (Max user weight 180kg)	£384
<b>SEAT HEIGHT SETTING*</b>			
ARH1024	<input type="checkbox"/>	Fixed seat height 470mm (18.5")	<input type="radio"/>
ARH1046	<input type="checkbox"/>	Fixed seat height 495mm (19.5")	<input type="radio"/>
<b>BACKREST RECLINER*</b>			
ARH0958	<input type="checkbox"/>	Manually settable recliner by spindle (0° - 25°)	<input type="checkbox"/>
ARH0947	<input type="checkbox"/>	Electrical recliner (0° - 25°)	£444

FOOTREST HANGER TYPE*				
ARH0556		Centre Mounted: Manual : Hanger, centre mounted, angle settable (0° - 90°) : knee to heel length 390 - 510mm, (including calf pad and angle adjustable footplate)	(NC 0504, 0513, 0599)	£328
ARH0589		Centre Mounted : LNX Powered: Centre mounted legrest, length compensating (including calf pad and angle adjustable footplate)	(NC 0504, 0513, 0599)	£1,442
ARH0564		Reinforced hangers, fixed angle (not detachable without tools) incl Calf pads		○
ARH0503		Side Mounted : Vari F : Angle adjustable with tools 0° - 70°, length 290mm - 460mm		◆
ARH0533		Side Mounted : Vari A : Angle adjustable with hand wheel 0° - 70°, length 290mm - 460mm	(NC 0559)	£206
ARH0531		Side Mounted : ADE Powered : Angle adjustable powered by joystick, 0° - 80° length 290mm - 460mm	(NC 0513, 0599, 0559)	£588
FOOTPLATE TYPE*				
ARH0512		Fixed Footplate	(NC 0533, 0556, 0589)	○
ARH0504		Angle Adjustable Footplate	(NC 0589)	£65
ARH0513		Angle and depth adjustable Footplate (Adapt)	(NC 0589)	£104
ARH0599		Angle adjustable flip-back, heavy-duty footboard (Recommended for SW 550mm)	(NC 0503, 0531, 0533, 0556, 0589)	£294
ARH0559		Angle Adjustable Footplate for Centre mounted legrest	(NC 0503, 0533, 0531)	○
CONTROLS*				
ARH2201		REM110 Drive only	(NC 1437, 1095, 0947, 0589, 0531, 2219, 2245)	◆
ARH2203		REM211 Drive & Actuators	(NC 1437)	○
ARH2204		REM216 Drive, Actuators & Lights		£134
CONTROL BRACKETRY*				
ARH1207		Standard Bracket RH		○
ARH1208		Standard Bracket LH		○
ARH1707		Swing Away RH Bracket		£116
ARH1708		Swing Away LH Bracket		£116
ARH1709		Matrx Resolve Swing in/out bracket RH		£219
ARH1710		Matrx Resolve Swing in/out bracket LH		£219
ARH2227		Bracket RH mounted for Attendant		£56
ARH2228		Bracket LH mounted for Attendant		£56
ARH1289		Tray / Centre control Swing Away Bracket RH	(NC 2203, 2204)	○
ARH1288		Tray / Centre control Swing Away Bracket LH	(NC 2203, 2204)	○
DUAL CONTROL				
ARH2220		REM050 Dual Control (Please select user control and bracket) (Please select dual control bracket below)	(NC 1221, 1222)	£348
ARH2219		Compact Dual Control (DLX-ACU200) (Please select user control and bracket) (Please select dual control bracket below)	(NC 1222, 2201)	£348
DUAL CONTROL BRACKETRY				
ARH1207		Dual Control Standard Bracket RH		○
ARH1208		Dual Control Standard Bracket LH		○
INVACARE LINX IDC DUAL CONTROL				
ARH2245		Invacare IDC Dual control (Please select IDC Dual control bracket below)	(NC 2219, 2227, 2228)	£660
INVACARE LINX IDC DUAL CONTROL BRACKETRY				
ARH1221		IDC Dual control Detachable Bracket RH	(NC 2219, 2227, 2228)	○
ARH1222		IDC Dual control Detachable Bracket LH	(NC 2219, 2227, 2228)	○

<b>LINX 10-WAY SWITCH</b>		
ARH2251		10-Way Switch DLX-SW10-A for direct Powered Seating Access (included loom to connect + additional jack for mode switch) (Please select 10 way switch bracket below) £444
<b>LINX 10-WAY SWITCH BRACKET</b>		
ARH2217		10-Way Switch Bracket RH, Mounted for User ○
ARH2218		10-Way Switch Bracket LH, Mounted for User ○
ARH2227		10-Way Switch Bracket RH, Mounted for Attendant £56
ARH2228		10-Way Switch Bracket LH, Mounted for Attendant £56
<b>LINX ACCESSORY</b>		
ARH2252		Loom providing jack sockets for Drive and Mode (includes connection for 10-way switch) £89
<b>FURTHER LINX COMPONENTS / ACCESSORIES</b>		
ARH2300		To be enclosed only (Accessory loom 2252 required if no REM400 or REM500 is selected) ○
ARZ2301		To be connected and programmed (special) POA
<b>BUDDY BUTTONS</b>		
ARH1685		Red Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £77
ARH1686		Yellow Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £77
ARH1687		Green Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £77
ARH1690		Blue Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £77
<b>EGG SWITCH</b>		
ARH2254		Egg Switch, black (to be enclosed only (Accessory loom 2252 required if no REM400 or REM500 is selected) £77
<b>PIKO BUTTONS</b>		
ARH1711		Red Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £139
ARH1712		Yellow Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £139
ARH1713		Green Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £139
ARH1714		Orange Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £139
ARH1715		Blue Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £139
ARH1716		Black Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected) £139
<b>TOGGLES / BUTTONS</b>		
ARZ2261		Stereo Toggle - direct control for Tilt only £237
ARZ2262		Stereo Button - direct control for Tilt only £226
ARZ2263		Four-Way Toggle - direct control for Tilt & Recline £295
ARZ2264		Four-Way Button - direct control for Tilt & Recline £315
ARZ2207		Toggle / Button RH Mount ○
ARZ2208		Toggle / Button LH Mount ○
ARZ2225		Without Mounting ○

<b>BODYPOINT JOYSTICK KNOBS</b>			
ARH2721		3" U-Shaped, Joystick Handle, Bodypoint PC101A	£64
ARH2722		4" U-Shaped, Joystick Handle, Bodypoint PC102A	£64
ARH2723		Domed Rubber, Joystick Handle, Bodypoint PC107A	£64
ARH2724		Mushroom Grooved, Joystick Handle, Bodypoint PC110A	£64
<b>REMOTE STOP SWITCH</b>			
ARZ2258		Remote Stop Switch (ASL504)	£924
<b>SHROUD COLOUR*</b>			
ARH2640		Obsidian Black Frame Shrouds	○
ARH2641		Ruby Red Frame Shrouds	○
ARH2642		Purple Passion Frame Shrouds	○
ARH2643		Peacock Blue Frame Shrouds	○
ARH2644		Arctic Blue Frame Shrouds	○
<b>RIM COLOUR*</b>			
ARH2650		Obsidian Black Rim Inserts	○
ARH2651		Ruby Red Rim Inserts	○
ARH2652		Purple Passion Rim Inserts	○
ARH2653		Peacock Blue Rim Inserts	○
ARH2654		Arctic Blue Rim Inserts	○
<b>FORKS*</b>			
ARH2631		Single sided fork (New) includes Black castor	£310
ARH2632		Double sided fork	○
<b>CENTRE WHEELS (DRIVING WHEELS)*</b>			
ARH0704		DW 14" (355 mm) x 3" (75 mm), Pneumatic	○
ARH0705		DW 14" (355 mm) x 3" (75 mm), Solid inlay, dual density foam	○
ARH0706		DW 14" (355 mm) x 3" (75 mm), Puncture protected, gel-sealed	○
<b>MOTOR TYPE*</b>			
ARH1294		3 km/h (2mph)	○
ARH1282		6 km/h (4mph) (4 pole)	◆
<b>BATTERY*</b>			
ARH1260		73.5 Ah	◆
<b>G TRAC</b>			
ARH1215		G-Trac (ensures the powerchair takes the intended direction as indicated by the user control, without the need for continual correction by the client on uneven or sloping terrain)	£724
<b>CHARGER*</b>			
ARH1343		OptiCharge Charger 8A (200 - 250V)	◆
ARH1313		OptiCharge Charger 10A (200 - 250V) Charging time will be 20% faster	£134
<b>LIGHTS*</b>			
ARH1437		LED Light System with Indicators (30 LUX)	£422
ARH1254		Without light system	◆

<b>TRANSPORTATION KIT</b>		
ARH1450		Preparation for DAHL Docking station: Docking station for car implementation, please order docking station separately (Adequate headrest is compulsory, see options below) £514
ARH1494		4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory, see options below) ◆
<b>FIXATION BELT*</b>		
ARH1329		MAX Lapbelt (adjustable both sides) fitted to base of backrest ◆
<b>REAR VIEW MIRROR</b>		
ARH1415		Rear-view Mirror RH £37
ARH1339		Rear-view Mirror LH £37
<b>CANE HOLDER</b>		
ARH1312		Cane holder RH (In combination with Attendant bracket or Dual control the cane holder has to be mounted on the opposite side of the remote) £37
ARH1363		Cane holder LH (In combination with Attendant bracket or Dual control the cane holder has to be mounted on the opposite side of the remote) £37

<b>TOTAL PRICE</b>	
--------------------	--

**Non Conforming Options:**

ARH1095	Electric heavy duty basic Tilt (0° - 20°) select either 470mm or 495mm seat to ground height. (Max user weight 180kg)
ARH0947	Electrical recliner (0° - 25°)
ARH0556	Centre Mounted: Manual : Hanger, centre mounted, angle settable (0° - 90°) : knee to heel length 390 - 510mm, (including calf pad and angle adjustable footplate)
ARH0589	Centre Mounted : LNX Powered: Centre mounted legrest, length compensating (including calf pad and angle adjustable footplate)
ARH0503	Side Mounted : Vari F : Angle adjustable with tools 0° - 70°, length 290mm - 460mm
ARH0533	Side Mounted : Vari A : Angle adjustable with hand wheel 0° - 70°, length 290mm - 460mm
ARH0531	Side Mounted : ADE Powered : Angle adjustable powered by joystick, 0° - 80° length 290mm - 460mm
ARH0504	Angle Adjustable Footplate
ARH0513	Angle and depth adjustable Footplate (Adapt)
ARH0599	Angle adjustable flip-back, heavy-duty footboard (Recommended for SW 550mm)
ARH0559	Angle Adjustable Footplate for Centre mounted legrest
ARH2201	REM110 Drive only
ARH2203	REM211 Drive & Actuators
ARH2204	REM216 Drive, Actuators & Lights
ARH2227	Bracket RH mounted for Attendant
ARH2228	Bracket LH mounted for Attendant
ARH2219	Compact Dual Control (DLX-ACU200) (Please select user control and bracket) (Please select dual control bracket below)
ARH2245	Invacare IDC Dual control (Please select IDC Dual control bracket below)
ARH1221	IDC Dual control Detachable Bracket RH
ARH1222	IDC Dual control Detachable Bracket LH
ARH2227	10-Way Switch Bracket RH, Mounted for Attendant
ARH2228	10-Way Switch Bracket LH, Mounted for Attendant
ARH1437	LED Light System with Indicators (30 LUX)



For a more complex customer requirement  
please refer to our Invacare Unique special order form.  
email [specialsuk@invacare.com](mailto:specialsuk@invacare.com)

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):